



# Suggested algorithm for the detection of beer spoilage organisms using GEN-IAL qPCR kits

## 5 potential sources of microbiological spoilage in beer brewing


### Potential sources of primary contamination



### Potential sources of secondary contamination



## Algorithm for in process samples and final product



Other block cycler devices may be suitable as well. Information is available on request. Different kits are available for different real-time thermocyclers:

Tube profile	Real-time device
high profile	ABI 7500, Agilent MX3005P, Quant Studio 5 (0.2 mL)
low profile	Agilent Aria MX, BioRad CFX96, MyGo Pro, ABI 7500 Fast, Quant Studio 5 (0.1 mL)
white strips	BioRad CFX96, LightCycler® 480 II, Analytik Jena qTower <sup>3</sup>



## 1 DNA sample preparation

DNA extraction kit	Intended use	Method	Sample volume
<b>A</b> CSY 0100 GEN-IAL® QuickGEN® Sample preparation in yeast	For beverage samples mainly containing yeast	Centrifugation	1 mL sample
<b>B</b> CSE 0100 GEN-IAL® QuickGEN Sample preparation centrifugation	Enriched or non-filterable samples & colonies	Centrifugation	1 - 50 mL sample
<b>C</b> FSE 0100 GEN-IAL® QuickGEN Sample preparation filtration	All filterable samples with and without enrichment	Filtration	Up to 1 L sample

## 2 Group specific screening + differentiation of bacteria and yeast

Art. No./Test kit	Channel/Target				Comments
	FAM	ROX	HEX	Cy5	
<b>QPP1T0050</b> GEN-IAL® QuickGEN P1 Screening	<i>Lactobacillus</i> spp., <i>Pediococcus</i> spp.	<i>Megasphaera</i> spp., <i>Pectinatus</i> spp.	Yeast	IAC	Possible 1st screening step: Group wise detection of all relevant beer spoilage organisms
<b>QTPYB0096</b> GEN-IAL® QuickGEN First-Beer Yeast and Bacteria Differentiation	Entero- bacteriaceae, wild yeast 1, wild yeast 2	Acetic acid bacteria, <i>Pediococci</i> IAC	<i>Lactobacilli</i> , <i>Pediococci</i> , Bottom-fermented yeast, Top-fermented yeast		For 24 samples; Tubes are pre-coated with the qPCR reagents and contain lyticase – for easier handling and a more convenient lysis step
<b>QTPBF0050</b> GEN-IAL® QuickGEN First-Biofilm	<i>Lactococcus lactis</i> , <i>Leuconostoc mesenteroides</i>	IAC	<i>Pichia anomala</i>		Biofilm forming species
<b>QPP1SD0050</b> GEN-IAL® QuickGEN P1	<i>Lactobacillus</i> spp., <i>Pediococcus</i> spp.	<i>Megasphaera</i> spp., <i>Pectinatus</i> spp.	<i>Saccharomyces cerevisiae var. diastaticus</i>	IAC	Group wise detection of all relevant beer spoilage bacteria and <i>Saccharomyces cerevisiae var. diastaticus</i> . The latter can cause a secondary fermentation

### Screening of bacteria

Art. No./Test kit	Channel/Target				Comments
	FAM	ROX	HEX	Cy5	
<b>QPP1HR0048</b> GEN-IAL® QuickGEN P1 Screening	<i>Lactobacillus</i> spp., <i>Pediococcus</i> spp.	<i>Megasphaera</i> spp., <i>Pectinatus</i> spp.	Hop resistance genes horA/horC	IAC	Additional risk assessment for hop resistance

### Screening of yeast

Art. No./Test kit	Target		Comments
	FAM	HEX	
<b>QTPYD0050</b> GEN-IAL® QuickGEN First-Yeast PCR Kit Dekkera spp.	<i>Dekkera</i> spp.	IAC	Including: <i>D. bruxellensis</i> , <i>D. anomala</i> , <i>D. custersiana</i> , <i>D. naardenensis</i> & <i>D. nanus</i>
<b>TPYUG0050</b> GEN-IAL® Bottom fermented yeast	Bottom fermenting yeast	IAC	Detection of bottom-fermenting yeasts to monitor the purity of yeasts in the quality control process
<b>TPYOG0050</b> GEN-IAL® Top fermented yeast	Top fermenting yeast	IAC	Detection of top-fermenting yeasts to monitor the purity of yeasts in the quality control process

## 3 Identification of bacteria and yeast in one assay

Art. No./Test kit	Target			Comments
	FAM	ROX	HEX	
<b>QTPBD0096</b> GEN-IAL® QuickGEN First-Beer Differentiation PCR Kit	<i>Enterobacteriaceae</i> , <i>P. damnosus</i> , <i>Pectinatus</i> spp., <i>L. brevis</i> / <i>L. parvulus</i> , <i>L. parabrachneri</i> , <i>L. buchneri</i> / <i>L. parabrachneri</i> , <i>L. plantarum</i> / <i>L. paraplantarum</i> , <i>L. acetotolerans</i>	IAC, PTC <i>Saccharomyces cerevisiae var. diastaticus</i> , <i>P. clausenii</i> , <i>L. rossiae</i> , <i>L. paracasei</i> , <i>L. perolens</i> / <i>L. harbinensis</i> ,	<i>Pichia anomala</i> , <i>Pediococcus</i> spp. ( <i>P. acidilactici</i> , <i>P. parvulus</i> , <i>P. inopinatus</i> , <i>P. pentosaceus</i> ), <i>Megasphaera</i> spp., <i>L. lindneri</i> , <i>L. collinoideus</i> / <i>L. paracollinoideus</i> , <i>L. coryniformis</i> , <i>L. backii</i>	Detection of 30 species and identification of 19 species  For 12 samples; Tubes are pre-coated with the qPCR reagents and contain lyticase – for easier handling and a more convenient lysis step

### Identification of bacteria only

Art. No./Test kit	Target			Comments
	FAM	ROX	HEX	
<b>QTPHR0050</b> GEN-IAL® QuickGEN hop resistance genes horA and horC / hitA and orf5 TaqMan™	hor A / hor C	IAC	hit A, ORF 5	Hop resistance genes as genetic markers for beer-spoilage microorganisms
<b>TPPMD0050</b> GEN-IAL® <i>Pectinatus</i> spp./ <i>Megasphaera</i> spp.	<i>Pectinatus</i> spp.	IAC	<i>Megasphaera</i> spp.	Strictly anaerobic bacteria with high beer spoilage potential

### Identification of yeast only

Art. No./Test kit	Target		Comments
	FAM	HEX	
<b>QYDJF0096</b> GEN-IAL® QuickGEN First yeast differentiation PCR Kit	<i>Rhodotorula</i> spp., <i>Candida</i> spp., <i>Saccharomyces ludwigii</i> , <i>Torulasporea delbrückii</i> , <i>Kluyveromyces marxianus</i> , <i>Dekkera</i> spp., <i>Pichia</i> spp.	IAC, PTC <i>Saccharomyces exiguus</i> , <i>Saccharomyces cerevisiae var. diastaticus</i> , <i>Saccharomyces bayanus</i> , <i>Hanseniaspora</i> spp., <i>Debaromyces hansenii</i> , <i>S. pastorianus</i>	12 yeast species  For 12 samples; Tubes are pre-coated with the qPCR reagents and contain lyticase – for easier handling and a more convenient lysis step
<b>TPYDA0050</b> GEN-IAL® <i>Dekkera anomala</i>	<i>Dekkera anomala</i>	IAC	
<b>TPYPA0050</b> GEN-IAL® <i>Pichia anomala</i>	<i>Pichia anomala</i> ( <i>Wickerhamomyces anomalus</i> )	IAC	
<b>TPYSD0050</b> GEN-IAL® <i>Saccharomyces diastaticus</i>	<i>Saccharomyces cerevisiae var. diastaticus</i>	IAC	Can cause a secondary fermentation due to its amylolytic property.
<b>TPYPM0050</b> GEN-IAL® <i>Pichia membranaefaciens</i>	<i>Pichia membranaefaciens</i>	IAC	
<b>QTPYDB0048</b> GEN-IAL® QuickGEN First-Yeast PCR Kit Dekkera bruxellensis	<i>Dekkera bruxellensis</i>	IAC	Precoated tubes containing the qPCR reagents as well as lyticase – for easier handling and a more convenient lysis step
<b>QTPYZB0048</b> GEN-IAL® QuickGEN First-Yeast PCR Kit Zygosaccharomyces bailii	<i>Zygosaccharomyces bailii</i>	IAC	Troublesome spoiling yeast species due to its exceptional tolerance to various stressful conditions Precoated tubes containing the qPCR reagents as well as lyticase – for easier handling and a more convenient lysis step