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QuickGEN PCR Kit

P1 Screening

1. Examples of successful detected microorganisms

Bacteria

Lactobacillus spp.

Lactococcus lactis

Megasphaera spp.

Micrococcus kristinae

Pectinatus spp.

Pediococcus spp.

Yeast

Saccharomyces spp.

Dekkera spp.

Zygosaccharomyces spp.

Rhodotorula spp.

Torulaspora spp.

Kluyveromyces spp.

Debaromyces spp.

Candida spp.

Hanseniaspora spp.

Pichia spp.

2. Intended use

Qualitative detection of bacterial and yeast contaminations in filtrated beer and beer mixes. For products with high amount of yeast use the Screening Kit without yeast detection to find bacterial contaminations more sensitive (Q032).

3. Test principle

The TaqMan[®] real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ; ATTO/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument. To avoid false negative PCR-results an Inhibition Control (ROX/DQ) is amplified together in one reaction vessel with the specific sequence. The system contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit). **The tubes contain lyticase.**

4. Kit contents

The kit contains sufficient reagents for 48 PCR reactions:

1 x Premix	white cap
48 x Dye Strips (freeze-dried, incl. IC-DNA and lyticase)	tube strips
1 x ddH ₂ O	colourless cap
1 x Control-DNA (freeze-dried)	yellow cap

5. Storage conditions

The Control-DNA is freeze-dried and has to be solved in ddH₂O prior to use (see 7.1).

Do **not** freeze the lyophilized Dye Strips and lyophilized Control-DNA.

The PCR-reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected. All reagents are stable for 12 months, if they are stored correctly.

6. Materials required, but not provided

6.1. Instruments

Real-time PCR instrument for low profile PCR tubes

Centrifuge for 1.5 – 2.0 mL reaction tubes

Centrifuge for strips

Pipettes

“Vortex”

6.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

optional: Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction)

7. PCR

7.1. PCR-Setup

When using the kit for the first time, the freeze-dried Control-DNA has to be shortly centrifuged and carefully resolved:

- add 55 μL sterile ddH₂O to the freeze-dried Control-DNA
- after 15 minutes mix well

Before every use thoroughly mix all PCR components and centrifuge briefly.

PCR-reaction Setup:

PCR-Components	volume (μL)
Premix	15.0
Sample-DNA	5.0*
Total volume	20.0

* if using the Simplex easy DNA- or QuickGEN Yeast Sample Preparation Kit:
add 2.5 μL DNA and 2.5 μL sterile ddH₂O

1. Pipette 15 μL of the Premix into each PCR-tube, making sure that, prior to the first filling, the tip has been moistened.
2. Add 5 μL sample DNA, add 5 μL of the Control-DNA for the positive control reaction, add 5 μL of the extraction control and 5 μL of ddH₂O for the negative control* reaction. Use a fresh tip with each DNA filling.
3. Close the tubes immediately and centrifuge them shortly.
4. Place the tubes in the PCR-machine and start run.

Very important: * Please fill up the negative control with 5 μL ddH₂O to avoid unspecific amplification.

Work swiftly to avoid warming up and keep away from light.

7.2 PCR-Program

Note:

For the use of UNG the thermal cycler programm has to be changed according to manufacturers` instructions.

Software 3.3.2

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 30* - 35 x
Cycling Annealing/ Elongation	20 sec	60 °C	

*** for preenriched samples program 30 cycles**

Advanced Settings	Integration time (s)	0.5
	Acquisitions per Cycle	10

Software 3.4.8

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 30* - 35 x
Cycling Annealing/ Elongation	30 sec	60 °C	

*** for preenriched samples program 30 cycles**

Advanced Settings	Integration time (s)	0.5
--------------------------	----------------------	-----

8. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

Lactobacillus/ Pediococcus-DNA: FAM-channel

Megasphaera/ Pectinatus-DNA: HEX-channel

Inhibition Control-DNA: ROX-channel

Yeast-DNA: ATTO-channel


A sample is ***Lactobacillus/ Pediococcus* positive**, if there is a detectable fluorescence increase in the **FAM-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the ROX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is ***Megasphaera/ Pectinatus* positive**, if there is a detectable fluorescence increase in the **HEX-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the ROX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is **yeast positive**, if there is a detectable fluorescence increase in the **ATTO-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the ROX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is **negative**, if there is no detectable fluorescence increase in the different channels and the positive controls have a positive fluorescence signal. The negative controls show no amplification. The Inhibition-Control in the ROX-channel has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

analysis flowchart

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method. GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.