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Advanced Settings

Integration time (s) 0.5
Acquisitions per Cycle 10

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
Dekkera bruxellensis quantitative

1. Intended use

Qualitative and quantitative detection of *Dekkera bruxellensis* in wine and grape must.

2. Test principle

The real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (TaqMan®), which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

The kit contains one specific system for the detection of *Dekkera bruxellensis*. The system emits a maximum fluorescent signal in the FAM/DQ-channel. To avoid false negative PCR-results an Inhibition Control is amplified together in the same reaction vessel with the specific sequence (HEX/DQ-channel). The system contains dUTP. Optional: Use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit). **The tubes contain lyticase.**

3. Kit contents

The kit contains sufficient reagents for 48 reactions:

1 x Premix	white cap
48 x Dye Strips (freeze-dried, incl. Lyticase and IC-DNA)	reaction tubes
1 x ddH ₂ O	colourless cap
1 x Control-DNA (Standard 2, <i>D. bruxellensis</i> , freeze-dried) Standard 2 = 20000 cells/2.5 µL DNA	yellow cap
1 x Solution S (for solving Control-DNA)	green cap

4. Storage conditions

The Control-DNA is freeze-dried, it has to be solved in Solution S prior to use (see 6.1).

Do **not** freeze the lyophilized Dye Strips and lyophilized Control-DNA.

The solved Control-DNA and the Premix should be stored at -20 °C (-4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probe and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR instrument for low profile tubes clear

Centrifuge for 1.5 – 2.0 mL reaction tubes

Centrifuge for strips

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

optional: Uracil N-Glycosylase (0,01 U/μL added to the PCR reaction mix)

6. PCR

6.1. PCR-Setup

When using the kit for the first time, the freeze-dried Control-DNA has to be shortly centrifuged and carefully resolved in Solution S

- add 55 μL Solution S to the freeze-dried Control-DNA
- after 30 minutes mix well

Before every use thoroughly mix all PCR components and centrifuge briefly.

PCR-reaction per sample:

PCR-components	amount (μL)
Premix	17.5
Sample-DNA	2.5
Total volume	20.0

1. Pipette 17.5 μL of the Premix into each PCR-tube, making sure that, prior to the first filling, the tip of the pipette has been moistened.
2. Add 2.5 μL sample DNA, add 2.5 μL of the Control-DNA for the PCR positive control, add 2.5 μL of the extraction control and 2.5 μL of ddH₂O for the negative control* reaction. Use a fresh tip with each DNA filling.
3. Close the tubes immediately and centrifuge them shortly.
4. Place the tubes in the PCR machine and start run.

Very important: * Please fill up the negative control with 2.5 μL ddH₂O to avoid unspecific amplification.

Work swiftly to avoid warming up and keep away from light

6.2 PCR-Program

Note:

For the use of UNG the thermal cycler programm has to be changed according to manufacturers` instructions.

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial Denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	35 sec	63 °C	

For MyGo Pro User

Software 3.3.2

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial Denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	35 sec	63 °C	

Advanced Settings	Integration time (s)	0.5
	Acquisitions per Cycle	10

Software 3.4.8

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial Denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	40 sec	63 °C	

Advanced Settings	Integration time (s)	0.5
--------------------------	----------------------	-----

7. Evaluation

A. qualitative

The evaluation has to be made according to the data analysis programme recommended by the real-time instrument manufacturer.


***Dekkera bruxellensis*-DNA:** FAM-channel

Inhibition Control-DNA: HEX-channel

A sample is ***Dekkera bruxellensis* positive**, if there is a detectable fluorescence increase in the **FAM-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in HEX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). For negative controls it has to be positive.

A sample is **negative**, if there is no detectable fluorescence increase in FAM-channel and the positive controls have a positive fluorescence signal. The negative controls show no amplification in FAM-channel. The Inhibition Control in HEX-channel has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

analysis flowchart

B. quantitative

Quantitative data analysis has to be made according to the data analysis program recommended by the real-time instrument manufacturer. For insufficient results with automatic data analysis, adjust the threshold manually to an optimal efficiency.

In a quantitative real-time PCR assay the ct-values should decrease linearly with ascending DNA concentration per reaction. The software provided with the real-time PCR machine calculates standard curve and slope using the DNA-concentrations stated by the user and the appendant ct-values. The standard row should have an error < 0.1 and an optimal PCR efficiency. The difference between the ct-values of the standard dilutions (at a ratio of 1:10) should be approximately 3.3. Fits the ct-value of the PTC to the ct-value of standard 2 of the standard row? Deviations of 0.5 ct up or down are acceptable. If not, repeat pipetting. The ct-values of the samples with unknown DNA concentrations are now automatically compared to the standard curve and concentrations are assigned.

For accurate quantification, the ct-values of the sample DNA should, as much as possible, lie within the range of the crossing points of the standard series.

For quantification of *Dekkera bruxellensis* concentration per mL sample, the measured values for the samples have to be multiplied with the factor 40. The factor is calculated from the DNA volume of 100 µL after DNA isolation (SEW 0100, CSE 0100) and pipetting of 2.5 µL.

Example:

measured value of $4.61E+03 = 4610 \times 40 = 1.8 \times 10^5$ cfu / mL sample

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.