---

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## Dekkera bruxellensis quantitative

### 1. Intended use

Qualitative and quantitative detection of *Dekkera bruxellensis* in wine and grape must.

### 2. Test principle

The real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (TaqMan<sup>®</sup>), which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

The kit contains one specific system for the detection of *Dekkera bruxellensis*. The system emits a maximum fluorescent signal in the FAM-channel. To avoid false negative PCR-results an inhibition control is amplified together in the same reaction vessel with the specific sequence (HEX/DQ). The system contains dUTP. Optional: Use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit). **The tubes contain lyticase.**

### 3. Kit contents

The kit contains sufficient reagents for 48 reactions:

1 x Premix	white cap
48 x Dye Strips (freeze-dried, tubes incl. lyticase and IC-DNA)	tube strips
6 x Cap Strips	
1 x ddH <sub>2</sub> O	colourless cap
1 x Control-DNA (Standard 2, <i>D. bruxellensis</i> , freeze-dried) Standard 2 =20000 cells/2,5 µL DNA	yellow cap
1 x Solution S (for solving Control-DNA)	green cap

#### **4. Storage conditions**

**The Control-DNA is freeze-dried, it has to be solved in Solution S prior to use (see 6.1).**

Do **not** freeze the lyophilized Dye Strips and lyophilized Control-DNA.

The PCR reagents should be stored at 2 - 8 °C (35 – 46 °F). Keep Premix for storage at - 20 °C (- 4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected. All reagents are stable for 12 months, if they are stored correctly.

#### **5. Materials required, but not provided**

##### **5.1. Instruments**

Real-time PCR machine for low profile tubes

Centrifuge for 1.5 – 2.0 mL reaction tubes

Centrifuge for strips

Pipettes

“Vortex”

##### **5.2. Reagents and plastic ware**

sterile ddH<sub>2</sub>O

sterile filter tips

**optional:** Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction mix)

## 6. PCR

### 6.1. PCR-Setup

*When using the kit for the first time, the freeze-dried Control-DNA has to be shortly centrifuged and carefully resolved in Solution S*

- add 55 µL Solution S to the freeze-dried Control-DNA
- after 30 minutes mix well

**Before every use thoroughly mix all PCR-components and centrifuge briefly.**

Remove the film from the needed tube strips and pipette the PCR-Components. After pipetting close the tube strips with the provided cap strips.

PCR-reaction per sample:

PCR-Components	amount (µL)
Premix	17.5
Sample-DNA	2.5
Total volume	<b>20.0</b>

1. Pipette 17.5 µL of the Premix into each PCR-tube, making sure that, prior to the first filling, the tip of the pipette has been moistened.
2. Add 2.5 µL sample DNA, add 2.5 µL of the Control-DNA for the PCR positive control, add 2.5 µL of the extraction control and 2.5 µL of ddH<sub>2</sub>O for the negative control\* reaction. Use a fresh tip with each DNA filling.
3. Close the tubes immediately and centrifuge them shortly.
4. Place the tubes in the PCR machine and start run.

**Very important: \* Please fill up the negative control with 2.5 µL ddH<sub>2</sub>O to avoid unspecific amplification.**

**Work swiftly to avoid warming up and keep away from light**

## 6.2 PCR-Program

### 6.2.1 PCR-Program LC480

1. Click the button **tool** on the right side in the window **LightCycler 480 Software release 1.5.0. SP1**
2. Click the button **Detection formats** at the left side of the menu bar
3. Click **New** in the window **Detection formats** and name the experiment
4. Open the window **Filter Combination Selection** and choose the following filter combinations: 465-510 / 533-580
5. Open the window **Selected Filter Combination** list and add the following amounts

Excitation Filter	Emission Filter	Name	Melt factor	Quant factor	Max. Integ. Time
465	510	465-510 FAM	1	10	2
533	580	533-580 HEX	1	10	2

6. **Close** the window
7. Click the button **New Experiment** on the right side of the menu bar
8. From the pull-down menu **Detection formats** choose the defined experiment, click the button **Customize** and check the detection formats. All of them have to be activated
9. Click the button **ok**
10. Define the following program:

2. Programm Name: **Heat**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
37	None	00:15:00	4.40		0	0	0
95	None	00:15:00	4.40		0	0	0

4. Programm Name: **Ampli**

Cycles **40** Analysis Mode **Quantification**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
95	None	00:00:10	4.40		0	0	0
63	Single	00:00:35	2.20		0	0	0

5. Programm Name : **Cool**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
40	None	00:00:20	2.20		0	0	0

**Optional:** Saving the program as **run template**. The template button allows to select and apply a template to the currently open object and to save the currently open object as a template. Click the clamp beside the button **Apply Template** . Click **Save as template** and save the file in **templates**.

11. Click the button **Subset editor** on the left side
12. Click the button **+** and **New Subset 1** appears
13. Mark the wells in the **New Subset 1 Settings** window
14. Click the button **Apply**
15. Click the button **Sample editor** on the left side
16. **Very important:** Activate In the window **Step 1 Select Workflow: Abs.Quant**
17. Choose the subset **New Subset** in the window **Step 2 Select Samples**
18. Define your probes in the **Sample table**
19. Switch to the button **Experiment** and start the run with the button **start run**

## 7. Evaluation

### A. qualitative

The evaluation has to be made according to the data analysis programme recommended by the real-time instrument manufacturer.

For LC480: For analysis activate Colour Compensation

***Dekkera bruxellensis*-DNA:** FAM-channel


**Inhibition Control-DNA:** HEX-channel

A sample is ***Dekkera bruxellensis* positive**, if there is a detectable fluorescence increase in the **FAM-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in HEX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). For negative controls it has to be positive.

A sample is negative, if there is no detectable fluorescence increase in FAM-channel and the positive controls have a positive fluorescence signal. The negative controls show no amplification in FAM-channel. The Inhibition Control in HEX-channel has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

# analysis flowchart

LC480: analysis after activating Colour Compensation

## B. quantitative

Quantitative data analysis has to be made according to the data analysis program recommended by the real-time instrument manufacturer. For insufficient results with automatic data analysis, adjust the threshold manually to an optimal efficiency.

In a quantitative real-time PCR assay the ct-values should decrease linearly with ascending DNA concentration per reaction. The software provided with the real-time PCR machine calculates standard curve and slope using the DNA-concentrations stated by the user and the appendant ct-values. The standard row should have an error < 0.1 and an optimal PCR efficiency. The difference between the ct-values of the standard dilutions (at a ratio of 1:10) should be approximately 3.3. Fits the ct-value of the PTC to the ct-value of standard 2 of the standard row? Deviations of 0.5 ct up or down are acceptable. If not, repeat pipetting. The ct-values of the samples with unknown DNA concentrations are now automatically compared to the standard curve and concentrations are assigned.

For accurate quantification, the ct-values of the sample DNA should, as much as possible, lie within the range of the crossing points of the standard series.

For quantification of *Dekkera bruxellensis* concentration per mL sample, the measured values for the samples have to be multiplied with the factor 40. The factor is calculated from the DNA volume of 100 µL after DNA isolation (SEW 0100, CSE 0100) and pipetting of 2.5 µL.

### Example:

measured value of  $4.61E+03 = 4610 \times 40 = 1.8 \times 10^5$  cfu / mL sam

### Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.