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Dekkera bruxellensis quantitative

1. Intended use

Detection of *Dekkera bruxellensis* in wine and grape must by real-time PCR.

2. Test principle

The real-time PCR is based on Hotstart-PCR and sequence-specific dual labelled probes (TaqMan®), which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

The kit contains one specific system for the detection of *Dekkera bruxellensis*. The system emits a maximum fluorescent signal at 518 nm (FAM/DQ, respectively NFQ). To avoid false negative PCR-results an inhibition control is amplified together in the same reaction vessel with the specific sequence (ROX/DQ, respectively NFQ). The system contains dUTP. Optional: Use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit).

3. Kit contents

The kit contains sufficient reagents for 50 PCR reactions.

1 x Premix	white cap
1 x Dye Mix (freeze-dried, incl. IC-DNA)	dark vial, red cap
1 x ddH ₂ O	colourless cap
1 x Control-DNA (Standard 2, <i>D. bruxellensis</i> , freeze dried) Standard 2 = 20000 cells/2.5 µL DNA	yellow cap
1 x Solution S (for solving Control-DNA)	green cap

4. Storage conditions

The Dye Mix is freeze-dried, it has to be solved in ddH₂O prior to use.

The freeze-dried Control-DNA has to be solved in Solution S (see 6.1).

Do **not** freeze the lyophilized Dye Mix and lyophilized Control-DNA.

The PCR reagents should be stored at 2 – 8 °C (35 – 46 °F).

The solved Control-DNA and the Premix should be stored at - 20 °C (- 4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Mix contains the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 6 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR instrument with FAM- and ROX-channel

Centrifuge for 1.5 – 2.0 mL reaction tubes

Centrifuge for multiwell plates or stripes

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile reaction vessels 1.5 – 2.0 mL

sterile optical tubes, plates

sterile filter tips

optional: Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction mix)

Colour Compensation Set for LC480 (GEN-IAL[®], Q800)

LC480: Use a Colour compensation prior to first use of the kit.

6. PCR

6.1. PCR-Setup

When using the kit for the first time, the freeze-dried kit components have to be shortly centrifuged and carefully resolved in ddH₂O or Solution S:

- add 80 µL sterile ddH₂O to freeze-dried Dye Mix
- add 55 µL Solution S to freeze-dried Control-DNA
- after 30 minutes mix well

Before every use thoroughly mix all PCR components and centrifuge briefly.

PCR-reaction per sample:

PCR-components	amount (µL)
Premix	16.0
Dye Mix	1.5
Sample or Control-DNA	2.5
Total volume	20.0

1. Prepare a mastermix by mixing Premix and Dye Mix
2. Multiply said volumes with the number of PCR preparations including controls (positive control, negative control, extraction control), taking into account pipette reserves of approximately 5-10 %.
3. Divide 17.5 µL of the PCR-mastermix among the individual reaction vessels, making sure that, prior to the first filling, the tip of the pipette has been moistened.
4. Add 2.5 µL sample DNA, add 2.5 µL of the Control-DNA for the PCR positive control, add 2.5 µL of the extraction control and 2.5 µL of ddH₂O for the negative control* reaction. Use a fresh tip with each DNA filling.
5. Close the tubes immediately and centrifuge them shortly.
6. Place the tubes in the PCR-machine and start run.

Very important: * Please fill up the negative control with 2.5 µL ddH₂O to avoid unspecific amplification.

Work swiftly to avoid warming up and keep away from light

6.2 PCR-Program

6.2.1 PCR-Program LC480

1. Click the button **tool** on the right side in the window **LightCycler 480 Software release 1.5.0. SP1**.
2. Click the button **Detection formats** at the left side of the menu bar.
3. Click **New** in the window **Detection formats** and name the experiment
4. Open the window **Filter Combination Selection** and choose the following filter combinations: 465-510 and 533-610.
5. Open the window **Selected Filter Combination** list and add the following amounts.

Excitation Filter	Emission Filter	Name	Melt factor	Quant factor	Max. Integ. Time
465	510	465-510 FAM	1	10	2
533	610	533-610 ROX	1	10	2

6. **Close** the window.
7. Click the button **New Experiment** on the right side of the menu bar.
8. From the pull-down menu **Detection formats** choose the defined experiment, click the button **Customize** and check the detection formats. All of them have to be activated.
9. Click the button **ok**
10. Define the following program:

1. Programm Name: **Heat**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
95	None	00:10:00	4.40		0	0	0

2. Programm Name: **Ampli**

Cycles **40** Analysis Mode **Quantification**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
95	None	00:00:15	4.40		0	0	0
63	Single	00:00:35	2.20		0	0	0

3. Programm Name : **Cool**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
40	None	00:00:20	2.20		0	0	0

Optional: Saving the program as **run template**. The template button allows to select and apply a template to the currently open object and to save the currently open object as a template. Click the clamp beside the button **Apply Template**. Click **Save as template** and save the file in **templates**.

11. Click the button **Subset editor** on the left side
12. Click the button **+** and **New Subset 1** appears
13. Mark the wells in the **New Subset 1 Settings** window
14. Click the button **Apply**
15. Click the button **Sample editor** on the left side
16. **Very important:** Activate In the window **Step 1 Select Workflow: Abs.Quant**
17. Choose the subset **New Subset** in the window **Step 2 Select Samples**
18. Define your probes in the **Sample table**
19. Click the button **Experiment** and start the run with **Start run**

6.2.2 PCR-Program for other real-time machines

The Premix contains no ROX. This must be considered according to the settings of the real-time machine. **ABI 7500:** Under “Assign Targets and Samples” in ‘Select the dye to use as the passive reference’ choose „none“.

Mx3500P: Under Instruments in ‘Filter Set Gain Settings’ change the Filter set gain multiplier to CY5 x1, ROX x1, HEX-JOE x4 und FAM x8.

For the use of UNG the thermal cycler program has to be changed according to manufacturers instructions.

Step	Time	Temp.	
Initial denaturation	10 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	35 sec	63 °C	

7. Evaluation

A. qualitative

The evaluation has to be made according to the data analysis programme recommended by the real-time instrument manufacturer.

For LC480: Activate Colour Compensation


***Dekkera bruxellensis*-DNA:** FAM-channel (LC480: 465-510 nm)

Inhibition Control-DNA: ROX-channel (LC480: 533-610 nm)

A sample is ***Dekkera bruxellensis* positive**, if there is a detectable fluorescence increase in the **FAM-channel (LC480: 465-510 nm)** and the negative controls show no amplification. The Inhibition Control in ROX-channel (LC480: 533-610 nm) may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). For negative controls it has to be positive.

A sample is negative, if there is no detectable fluorescence increase in FAM-channel and the positive controls have a positive fluorescence signal. The negative controls show no amplification in FAM-channel. The Inhibition Control in ROX-channel (LC480: 533-610 nm) has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

analysis flowchart

B. quantitative

Quantitative data analysis has to be made according to the data analysis program recommended by the real-time instrument manufacturer. For insufficient results with automatic data analysis, adjust the threshold manually to an optimal efficiency.

In a quantitative real-time PCR assay the ct-values should decrease linearly with ascending DNA concentration per reaction. The software provided with the real-time PCR machine calculates standard curve and slope using the DNA-concentrations stated by the user and the appendant ct-values. The standard row should have an error < 0.1 and an optimal PCR efficiency. The difference between the ct-values of the standard dilutions (at a ratio of 1:10) should be approximately 3.3. Fits the ct-value of the PTC to the ct-value of standard 2 of the standard row? Deviations of 0.5 ct up or down are acceptable. If not, repeat pipetting. The ct-values of the samples with unknown DNA concentrations are now automatically compared to the standard curve and concentrations are assigned.

For accurate quantification, the ct-values of the sample DNA should, as much as possible, lie within the range of the crossing points of the standard series.

For quantification of *Dekkera bruxellensis* concentration per mL sample, the measured values for the samples have to be multiplied with the factor 40. The factor is calculated from the DNA volume of 100 µL after DNA isolation (SEW 0100, CSE 0100) and pipetting of 2.5 µL.

Example:

measured value of $4.61E+03 = 4610 \times 40 = 1.8 \times 10^5$ cfu / mL sample

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.