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Yeast Differentiation

1. Intended use

Differentiation of beer spoilage yeasts in beer and beer mixing drinks.
The following yeast are differentiated:

Saccharomyces exiguus

Saccharomyces cerevisiae var. *diastaticus*

Saccharomyces bayanus / *pastorianus*

Saccharomycodes ludwigii

Torulasporea delbrückii

Klyveromyces marxianus

Debaromyces hansenii

<i>Dekkera</i> spp. :	<i>D. anomala</i>
	<i>D. custersiana</i>
	<i>D. bruxellensis</i>
	<i>D. naardenensis</i>
	<i>D. nanus</i>
<i>Rhodotorula</i> spp. :	<i>R. mucilaginoso</i>
	<i>R. glutinis</i>
	<i>R. graminis</i>
	<i>R. toruloides</i>
	<i>R. bacarum</i>
<i>Candida</i> spp. :	<i>C. albicans</i>
	<i>C. parapsilosis</i>
	<i>C. glabrata</i>
	<i>C. tropicalis</i>
	<i>C. sake</i>
	<i>C. intermedia</i>
	<i>C. pulcherrima</i>
<i>Pichia</i> spp. :	<i>P. anomala</i>
	<i>P. fermentans</i>
	<i>P. membranaefaciens</i>
<i>Hanseniaspora</i> spp. :	<i>H. guillermondii</i>
	<i>H. uvarum</i>
	<i>H. osmophila</i>

2. Test principle

The TaqMan[®] real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument. To avoid false negative PCR-results an Inhibition Control is amplified together in one tube with the specific sequence. The system contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit). **The tubes contain lyticase.**

3. Kit contents

The kit contains sufficient reagents for 96 PCR reactions (12 samples):

1 x Premix	white cap
12 x Dye Strips (freeze-dried, tubes contain lyticase)	tube strips
1 x ddH ₂ O (for PCR negative controls)	colourless cap

4. Storage conditions

Do **not** freeze the lyophilized Dye Strips.

The PCR reagents should be stored at 2 - 8 °C (35 – 46 °F). Keep Premix for storage at - 20 °C (- 4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR instrument for low profile tubes

Centrifuge for tube strips

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

optional: Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction mix)

6. PCR

6.1. PCR-Setup

Before every use thoroughly mix Premix and centrifuge briefly.

PCR-reaction setup:

PCR-Components	Volume (μL)
Premix	17.5
Sample-DNA	2.5
Total volume	20.0

Pipetting schema:

Tube 1*	Premix 17.5 μL	ddH₂O 2.5 μL <u>Attention:</u> pipette no sample-DNA
Tube 2*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 3*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 4*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 5*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 6*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 7*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 8*	Premix 17.5 μL	sample-DNA 2.5 μL

* please notice the orientation of the tubes: they have to be opened to the right side

6.2 PCR-Program

Note: For the use of UNG the thermal cycler program has to be changed according to manufacturers` instructions.

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30*-35x
Cycling Annealing/ Elongation	20 sec	63 °C	

* for preenriched samples program 30 cycles

For MyGo Pro User:

Software 3.3.2

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30*-35x
Cycling Annealing/ Elongation	20 sec	63 °C	

* for preenriched samples program 30 cycles

Advanced Settings	Integration time (s)	0.5
	Acquisitions per Cycle	10

Software 3.4.8

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30*-35x
Cycling Annealing/ Elongation	30 sec	63 °C	

* for preenriched samples program 30 cycles

Advanced Settings	Integration time (s)	0.5
--------------------------	----------------------	-----

7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

The negative control (NTC, Tube 1) has to be negative in the FAM channel and positive in the HEX-channel (Inhibition Control). Is the Inhibition Control in the HEX-channel negative and/or the results in the FAM-channel positive, PCR has to be repeated.

The Inhibition Control (Tube 7) has to be positive in the HEX-channel. Is the Inhibition Control in the HEX-channel negative inhibitors or high amount of DNA are in the sample reaction.

The Positive Control (Tube 8) has to be positive in the HEX-channel. If it is negative, PCR has to be repeated.

PCR-Analysis:

please notice the orientation of the tubes: they have to be opened to the right side

	Channels	Species	Signal
Tube 1	FAM HEX	NTC Inhibition Control	negative positive
Tube 2	FAM HEX	Rhodotorula spp. Saccharomyces exiguus	
Tube 3	FAM HEX	Candida spp. Saccharomyces cerevisiae var. diastaticus	
Tube 4	FAM HEX	Saccharomycodes ludwigii Debaromyces hansenii	
Tube 5	FAM HEX	Torulaspota delbrückii Saccharomyces bayanus / pastorianus	
Tube 6	FAM HEX	Kluyveromyces marxianus Hanseniaspora spp.	
Tube 7	FAM HEX	Dekkera spp. Inhibition Control	positive
Tube 8	FAM HEX	Pichia spp. PTC	positive

Example for interpretation:

If sample shows a positive FAM-signal in Tube 3, *Candida spp.* is identified. If the sample contains mixtures, you have positive signals in different channels.

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.