---

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### Advanced Settings

Integration time (s)	0.5
Acquisitions per Cycle	10

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
## Dekkera spp.

### 1. Intended use

Detection of *Dekkera* spp. in beverages (e.g. beer, wine, grape must) by real-time PCR. The following species are detected: *D. bruxellensis*, *D. anomala*, *D. custersiana*, *D. naardenensis* and *D. nanus*

### 2. Test principle

The real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (TaqMan<sup>®</sup>), which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

The kit contains a specific system for the detection of *Dekkera spp.* The system emits a maximum fluorescent signal in the FAM/DQ channel. To avoid false negative PCR-results an Inhibition Control is amplified together in the same reaction vessel with the specific sequence (HEX/DQ channel). The system contains dUTP.

Optional: Use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit). **The tubes contain lyticase.**

### 3. Kit contents

The kit contains sufficient reagents for 48 reactions:

1 x Premix	white cap
48 x Dye Strips (freeze-dried, incl. IC-DNA and lyticase)	tube strips
1 x ddH <sub>2</sub> O	colourless cap
1 x Control-DNA (freeze-dried)	yellow cap

#### **4. Storage conditions**

**The Control-DNA is freeze-dried and has to be solved in ddH<sub>2</sub>O prior to use (see 6.1).**

Do **not** freeze the lyophilized Dye Strips and lyophilized Control-DNA.

The PCR reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

#### **5. Materials required, but not provided**

##### **5.1. Instruments**

Real-time PCR machine for low profile tubes

Centrifuge for 1.5 – 2.0 mL reaction tubes

Centrifuge for strips

Pipettes

“Vortex”

##### **5.2. Reagents and plastic ware**

sterile ddH<sub>2</sub>O

sterile filter tips

**optional:** Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction mix)

## 6. PCR

### 6.1. PCR-setup

*When using the kit for the first time, the freeze-dried Control-DNA has to be shortly centrifuged and carefully resolved:*

- add 55  $\mu\text{L}$  sterile ddH<sub>2</sub>O to the freeze-dried Control-DNA
- after 15 minutes mix well

**Before every use thoroughly mix all PCR-components and centrifuge briefly.**

PCR-reaction per sample:

PCR-Components	amount ( $\mu\text{L}$ )
Premix	15,0
Sample-DNA	5,0*
Total volume	<b>20,0</b>

\* if using the Simplex Easy DNA- or QuickGEN Yeast Sample Preparation kit:  
add 2.5  $\mu\text{L}$  DNA and 2.5  $\mu\text{L}$  sterile ddH<sub>2</sub>O

1. Pipette 15  $\mu\text{L}$  of the Premix into each PCR-tube, making sure that, prior to the first filling, the tip has been moistened.
2. Add 5  $\mu\text{L}$  sample DNA, add 5  $\mu\text{L}$  of the Control-DNA for the PCR positive control, add 5  $\mu\text{L}$  of the extraction control and 5  $\mu\text{L}$  of ddH<sub>2</sub>O for the negative control\* reaction. Use a fresh tip with each DNA filling.
3. Close the tubes immediately and centrifuge them shortly.
4. Place the tubes in the PCR-machine and start run.

**Very important: \* Please fill up the negative control with 5  $\mu\text{L}$  ddH<sub>2</sub>O to avoid unspecific amplification.**

**Work swiftly to avoid warming up and keep away from light.**

## 6.2 PCR-Program

**Note:** For the use of UNG the thermal cycler programmes have to be changed according to manufacturer's instructions.

Step	Time	Temp.	
Lyticase treatment	15 min	37°C	
Initial Denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	20 sec	63 °C	

### For MyGo Pro User

#### Software 3.3.2

Step	Time	Temp.	
Lyticase treatment	15 min	37°C	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	20 sec	63 °C	

<b>Advanced Settings</b>	Integration time (s)	0.5
	Acquisitions per Cycle	10

#### Software 3.4.8

Step	Time	Temp.	
Lyticase treatment	15 min	37°C	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	30 sec	63 °C	

<b>Advanced Settings</b>	Integration time (s)	0.5
--------------------------	----------------------	-----

## 7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.


**Dekkera spp.-DNA:** FAM-channel

**Inhibition Control-DNA:** HEX-channel

A sample is **Dekkera spp.** positive, if there is a detectable fluorescence increase in the **FAM-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the HEX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). For negative controls it has to be positive.

A sample is negative, if there is no detectable fluorescence increase in the FAM-channel and the positive controls have a positive fluorescence signal. The negative controls have to show no amplification in the FAM-channel. The Inhibition Control in the HEX-channel has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

# analysis flowchart

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method. GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.