# **RIDA<sup>®</sup> QUICK DON RQS ECO**

**Art. No. R5911**

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## Brief information

RIDA<sup>®</sup>QUICK DON RQS ECO (Art. No. R5911) is a quantitative immunochromatographic test with a water extraction for the determination of deoxynivalenol (DON) in grain (wheat, corn, barley, oats).

The test kit contains 20 test strips, each of which can be used for one analysis.

Results are evaluated with the RIDA<sup>®</sup>QUICK SCAN reader (Art. No. ZG5005) or the RIDA<sup>®</sup>SMART APP software (Art. No. ZRSAM1000).

Preparation of samples: homogenize and extract

Time requirement:            Sample preparation (for 10 samples) .... approx. 10 min  
   Test procedure (incubation period) ..... 3 min

Limit of detection:            0.15 mg/kg (ppm)

Limit of quantification:      0.25 mg/kg

Measuring range:            0.25 - 7.5 mg/kg and 7.5 - 50 mg/kg

Specificity:                    The RIDA<sup>®</sup>QUICK DON RQS ECO test detects DON in wheat, corn, barley and oats samples.

The specificity of the RIDA<sup>®</sup>QUICK DON RQS ECO test was established by determining the cross reactivity to the corresponding substances in the buffer system. The specificity in the sample can differ due to the matrix effects of the values determined in the buffer system. Before cross-reactive substances are analyzed, their detection limit and recovery rate in the relevant matrix must be determined by the user. The test cannot distinguish between analytes and cross-reactive substances.

## 1. Intended use

RIDA®QUICK DON RQS ECO is a quantitative immunochromatographic test with a water extraction in a test strip format to determine DON in grains (wheat, corn, barley, oats).

## 2. General information

DON is a trichothecene mycotoxin. It is produced by fungi in the *Fusarium* genus. DON can be found in plant products and grains in particular. DON is the most prevalent toxin of the more than 100 known trichothecenes in Europe and North America.

## 3. Test principle

The immunochromatographic test in the form of a test strip is based on an antigen-antibody reaction. A specific anti-DON antibody detects the DON in the sample and results in a line pattern (test line / control line).

The intensity of the test line depends on the DON concentration of the sample. It increases with increasing DON contamination.

The control band (control line) should be present in all cases in order to prove that the test strip is valid. The control line fades as the DON concentration in the sample increases.

The test strip is evaluated using the RIDA®SMART APP software or the RIDA®QUICK SCAN reader.

## 4. Reagents provided

The reagents in the kit are sufficient for 20 analyses. Every test kit contains:

Component	Lid color	Format	Content
20 x Test strips	-	Ready for use	One analysis per test strip (individually packaged)
1 x Mobile solvent	White	Ready for use	25 ml
1 x RIDA®SMART APP cover	-	Ready for use	Lot-specific evaluation cover for use with RIDA®SMART APP

## 5. Reagents required but not provided

### 5.1. Equipment

- Grinder
- Scale
- Graduated cylinder
- 100 µl - 1000 µl pipette
- Latest version of the RIDA<sup>®</sup>SMART APP software (Art. Nr. ZRSAM1000) activated and installed on a compatible smartphone. Alternatively RIDA<sup>®</sup>QUICK SCAN (Art. Nr. ZG5005) with software version V1.0.1.7 2.0 or higher.

#### Optional

- Shaker (e.g. Rock-it 360, Trilogy<sup>®</sup> Analytical Laboratory Inc. or equivalent)
- Centrifuge (≥ 2000 g)
- Filter paper Munktell Ahlstrom 3hw or equivalent

### 5.2. Reagents

- Distilled or deionized water

## 6. Warnings and precautions for the users

This test should be carried out only by trained personnel. Always strictly adhere to the instructions for use of this test.

The test strips are sensitive to moisture. Moist test strips can have a negative impact on the test result. It is important to protect the strips against moisture, this applies particularly if the test strip packaging has been opened.

This kit may contain harmful substances. Please refer to the component safety information in the material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.de](http://www.r-biopharm.de).

## 7. Storage instructions

Store the test at 2 - 8 °C (35 - 46 °F). Do not freeze components of the test kit.

After the expiration date (see expiration on outer label of test kit), the quality guarantee is no longer valid.

The exchange of individual reagents between kits with different lot numbers is not permitted.

## 8. Sample preparation

**All test kit components must be brought to room temperature (20 - 25 °C / 68 - 77 °F) prior to use. The samples must be prepared at room temperature.**

The samples should be stored in a cool place, protected from light.

Before extraction, grind and homogenize a representative sample (a sample taken in accordance with official sampling regulations).

### 8.1. Measurement range 0.25 - 7.5 mg/kg

- Weigh 5 g of the homogenized sample into a sealable tube and add 25 ml of distilled or deionized water
- Close the tube and shake vigorously for 30 sec (e.g. on Rock-it 360)
- Centrifuge (e.g. 1 min, 2000 g), filter or let settle to get a particle-free supernatant
- Dilution I: Dilute 100 µl of the particle-free supernatant with 500 µl of the mobile solvent and mix thoroughly (e.g. invert tube five times)
- Use 100 µl of Dilution I in the test

### 8.2. Measurement range 7.5 - 50 mg/kg

- Dilution II: Dilute 100 µl of Dilution I (see 8.1) with 700 µl of the mobile solvent and thoroughly mix (e.g. invert tube five times)
- Use 100 µl of Dilution II in the test

#### **Note:**

**The weight of the sample can be increased if necessary. The amount of water used should be adjusted in such cases, e.g., 50 g sample in 250 ml of distilled or deionized water.**

## 9. Test procedure

**All used components must be brought to room temperature (20 - 25 °C / 68 - 77 °F) prior to use. The test must also be performed at room temperature.**

- Remove a RIDA<sup>®</sup>QUICK DON RQS ECO test strip from its foil pouch
- Apply 100 µl of Dilution I or Dilution II (see 8. Sample preparation) to the application area (sample) of the test strip
- Incubate the test strip at room temperature
- Evaluate with the RIDA<sup>®</sup>QUICK SCAN or the RIDA<sup>®</sup>SMART APP after an incubation period of 3 minutes

## 10. Evaluation


The left-hand line in the reaction area is the test line (see fig. 1). Its intensity depends on the DON concentration of the sample.

The right-hand line in the reaction area is a control line and must appear after each test procedure. If this line is not visible, the test result is not valid because of an improperly performed test or deterioration of reagents. In this case, the test should be repeated with a new test strip. Please notify your local distributor if the control line is still missing after repeated tests.

Please note: The control line will only be faintly visible if the sample is highly contaminated with DON.

Test strips must always be evaluated after 3 minutes in order to achieve comparable results.

Fig. 1: RIDA<sup>®</sup>QUICK DON RQS ECO test strip

## 10.1. Evaluation with RIDA<sup>®</sup>SMART APP (Art. No. ZRSAM1000)

First, carefully read the user guide for the RIDA<sup>®</sup>SMART APP. The evaluation must be performed using these specifications.

Brief instructions on how to use the software application are available online at [www.r-biopharm.com](http://www.r-biopharm.com).

The RIDA<sup>®</sup>SMART APP cover enclosed in the test kit is required for test strip evaluation. The RIDA<sup>®</sup>SMART APP cover contains lot-specific information and must be used only with the intended test-kit lot. More information can be found on the quality assurance certificates enclosed in the test kit.

Please make sure the correct RIDA<sup>®</sup>SMART APP cover is present prior to starting the test procedure and test strip evaluation.

RIDA<sup>®</sup>SMART APP software applications for evaluation:

Matrix	Measuring range	RIDA <sup>®</sup> SMART APP application
Wheat	0.25 - 7.5 mg/kg	Wheat 0,25 - 7,5 ppm
	7.5 - 50 mg/kg	Wheat 7,5 - 50 ppm
Corn	0.25 - 7.5 mg/kg	Corn 0,25 - 7,5 ppm
	7.5 - 50 mg/kg	Corn 7,5 - 50 ppm
Barley/oats	0.25 - 7.5 mg/kg	Barley_Oats 0,25 - 7,5 ppm
	7.5 - 50 mg/kg	Barley_Oats 7,5 - 50 ppm

## 10.2. Evaluation with RIDA<sup>®</sup>QUICK SCAN

First, carefully read the user guide for the RIDA<sup>®</sup>QUICK SCAN. The procedure for evaluating test strips is provided in Section 3.

### **Please note:**

The quality assurance certificate included in the test kit contains information for several standard curves. Please choose the appropriate curve information for the matrix and measurement range of interest. To enter the curve parameters manually into the RIDA<sup>®</sup>QUICK SCAN reader use the plain text information on the quality assurance certificate; to enter the data with the barcode scanner use the 2D barcodes. In both cases, the checksum shown on the RIDA<sup>®</sup>QUICK SCAN screen needs to be identical to the respective checksum on the certificate.

RIDA®QUICK SCAN evaluation methods:

Measuring range	RIDA®QUICK SCAN method
0.25 - 7.5 mg/kg	DON Low
7.5 - 50 mg/kg	DON High

- Press the middle button (*Enter*) to activate the RIDA®QUICK SCAN reader
- Make sure that you have selected the correct method (*DON Low* or *DON High*)
- Go to the *Start Scan* menu item and press *Enter*
- Make sure that the checksum on the display matches the respective checksum on the quality assurance certificate included in the test kit
- Select the *Start Scan* menu item again and press *Enter*
- The reader will now evaluate the test strip and display the result
- The test is now complete

## 11. Sensitivity

The RIDA®QUICK DON RQS ECO test is capable of detecting DON contaminations quantitatively at levels between 0.25 - 50 mg/kg (0.25 - 7.5 mg/kg and 7.5 - 50 mg/kg).

**For further product information, please contact [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### R-Biopharm AG

Postanschrift / Postal Address:

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64297 Darmstadt, Germany

Sitz / Corporate Seat: Pfungstadt

Tel.: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-Mail: [info@r-biopharm.de](mailto:info@r-biopharm.de)

[www.r-biopharm.com](http://www.r-biopharm.com)

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are registered trademarks of R-Biopharm AG