# Simplex<sup>®</sup> Easy Spin

DNA-Extraktion aus *Alicyclobacillus spp.*

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# Simplex<sup>®</sup> Easy Spin

## 1. Intended use

Fast *Alicyclobacillus spp.* DNA extraction from e.g. juices, concentrates and tomato products.

## 2. Content

The kit contains sufficient reagents for 50 preparations:

- 1 x Lysis buffer SES
- 1 x Binding buffer
- 3 x Elution buffer
- 50 x Columns incl. Collection tubes

## 3. Materials required, but not provided

### 3.1. Instruments

- Heating block or water bath 37 °C – 95 °C (98.6 °F – 203 °F)
- Microcentrifuge for 1.5 – 2.0 mL reaction tubes
- Pipettes: 10 – 100 µL, 100 – 1000 µL, 0.5 – 10 µL

### 3.2. Reagents and plastic ware

- sterile safe-seal reaction tubes 1.5 – 2.0 mL
- suitable filter tips
- single use gloves
- optional: lysozyme
- Ethanol (96-100 %)

## 4. Storage

All reagents can be stored at room temperature (18-25 °C, 64.4 °F – 77 °F) and are stable up to one year. If there is any precipitate present in the buffers, warm the buffer to 37 °C (98.6 °F) to dissolve the precipitate before use.

**Attention:** *The binding buffer contains Guanidin hydrochloride.*

## 5. Risk and safety phrases (R & S)

<b>Component</b>	<b>Hazard components</b>	<b>Hazard symbol</b>	<b>Risk phrases</b>	<b>Safety phrases</b>
Binding buffer	Guanidin hydrochloride 50 – 66%	 Warnung	R 22-36/38	S26-37/39

### Risk phrases

R 22                    Harmful if swallowed  
R 36/38                Irritating to eyes and skin

### Safety phrases

S 26                    In case of contact with eyes, rinse immediately with plenty of water and seek medical advice  
S 37/39                Wear suitable gloves and eye/face protection

## 6. Indications of deterioration of reagents

In case of accurate handling deterioration unknown.

## 7. DNA-Extraction from *Alicyclobacillus spp.*

### 7.1 Sample preparation according to IFU-method 12

*For non filterable solutions:*

- Add 90 mL BAT-Bouillon to 10 mL product
- For inactivation of vegetative cells and activation of spores incubate the sample for 10 min at 80 °C and cool down to 45 °C (113 °F)
- Incubate the sample (3-7 days, aerobe at 45 °C +/- 1 °, 113 °F)

*For filterable solutions:*

- For inactivation of vegetative cells and activation of spores incubate 100 mL sample for 10 min at 80 °C (176 °F) and cool down to 45 °C (113 °F)
- Membrane filtration of the sample (0.45 µm filter)
- Incubate the filter in BAT-Bouillon (3-7 days, aerobe at 45 °C +/- 1 °, 113 °F)

### 7.2 DNA-Extraction protocol for non filterable solutions (e.g. concentrates)

- add 2 mL of the pre-enriched sample in a 2.0 mL safe-seal reaction tube
- centrifuge 5 min at 1500 x g
- transfer 1 mL supernatant in a 1.5 mL safe-seal reaction tube
- centrifuge 3 min at 18000 x g
- discard the supernatant
- resuspend the pellet in 250 µL Lysis buffer SES
- **Optional:** add 100 U Lysozyme and incubate 30 min at 65 °C (149 °F)
- incubate 30 min at 95 °C (203 °F)
- add 240 µL Binding buffer and 160 µL Ethanol
- vortex 10 sec
- transfer the lysate to the column placed in a collection tube
- centrifuge 2 min at 11000 x g
- transfer the column in a new 1.5 mL safe-seal reaction tube
- add 100 µL Elution buffer (preheated to 95 °C, 203 °F) onto the column
- incubate for 15 min at room temperature
- centrifuge 1 min at 11000 x g
- centrifuge 1 min at 11000 x g
- use 5 µL of the eluted DNA for PCR

### 7.3 DNA-Extraction protocol for filterable solutions

- add 2 mL of the pre-enriched sample in a 2.0 mL safe-seal reaction tube
- centrifuge 5 min at 18000 x g
- discard the supernatant
- resuspend the pellet in 250 µL Lysis buffer SES
  - Optional:** add 100 U Lysozyme and incubate 30 min at 65 °C (149 °F)
- incubate 30 min at 95 °C (203 °F)
- add 240 µL Binding buffer and 160 µL Ethanol
- vortex 10 sec
- transfer the lysate to the column placed in a collection tube
- centrifuge 2 min at 11000 x g
- transfer the column in a new 1.5 mL safe-seal reaction tube
- add 100 µL Elution buffer (preheated to 95 °C, 203 °F) onto the column
- incubate for 5 min at room temperature
- centrifuge 1 min at 11000 x g
- use 5 µL of the eluted DNA for PCR

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.