
**CONGEN**

**SureFast<sup>®</sup>**  
**SARS-CoV-2 PLUS**

Art. No. F7110  
100 rxn

**User Manual**

**March 2021**

 **Inhalt**

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## **Content**

1	General Information .....	8
1.1	Description .....	8
1.2	Limit of Detection .....	8
1.3	RNA Preparation .....	8
1.4	Kit components and storage .....	8
1.5	Additionally required equipment and materials .....	9
1.6	Setup.....	9
1.7	Detection channel Set-up .....	9
2	Qualitative Analysis .....	10
2.1	Protocol .....	10
2.1.1	RNA Preparation .....	10
2.1.2	Preparation of the master-mix .....	10
2.1.3	Preparation of the real-time PCR-mix .....	11
2.2	Interpretation of results .....	11
3	Validation .....	11
4	Further Information .....	12
4.1	Product Information .....	12
4.2	Technical Support .....	12

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## 1 General Information

### 1.1 Description

The SureFast® SARS-CoV-2 PLUS is a real-time RT-PCR for the direct, qualitative detection of novel coronavirus (SARS-CoV-2) RNA. The SureFast® SARS-CoV-2 PLUS real-time RT-PCR in combination with the SureFast® PREP DNA/RNA Virus (Art. No. F1051) has been validated and certified by AOAC (license-no. 022102).

Do not use for in-vitro diagnostics. Also do not use diagnostic procedures according to directive 98/79/EC or regulation (EU) 2017/746.

Each reaction contains an Internal Control RNA (ICR, consisting of MS2-bacteriophage) as an internal control of sample preparation procedure and to monitor possible PCR-inhibition. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride).

The real-time RT-PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time. The technical validation of instruments was performed on Roche LightCycler® 480 II, Qiagen Rotor-Gene Q, Applied Biosystems 7500 Fast Dx, Bio-Rad CFX96, R-Biopharm RIDA®CYCLER, Agilent AriaDx and Agilent Mx3005P.

The external validation according to AOAC was performed on the Bio-Rad CFX96.

#### Note:

Samples tested negative with SureFast®SARS-CoV-2 PLUS still could contain a SARS-CoV-2 contamination below the limit of detection of the assay.

### 1.2 Limit of Detection

The SureFast® SARS-CoV-2 PLUS real-time RT-PCR has a limit of detection of  $\leq 25$  RNA copies.

The assay limit of detection depends on sample matrix, processing grade, RNA preparation and RNA content. The SureFast® PCR systems are very sensitive and therefore even a small amount of target RNA is sufficient for a successful analysis. The concentration of total RNA in the sample does not allow a conclusion on the quantity and quality of the target RNA.

### 1.3 RNA Preparation

For RNA-preparation the use of SureFast® PREP DNA/RNA Virus (Art. No. F1051) is recommended.

### 1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 $\mu$ l	Yellow
2	Enzyme Mix	1 x 80 $\mu$ l	Red
R	Internal Control RNA	2 x 1800 $\mu$ l	Brown
N	PCR Water	1 x 500 $\mu$ l	White
P	Positive Control	1 x 250 $\mu$ l	Light Blue

**Store all reagents at  $-20^{\circ}\text{C}$  and protected from light.**

**1.5 Additionally required equipment and materials**

- RNA-Extraction kit (e.g. SureFast® PREP DNA/RNA Virus Art. No. F1051)
- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumables (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

**1.6 Setup**

	<b>Blockcycler &amp; Rotor-Gene &amp; R-Biopharm RIDA®CYCLER</b>	<b>LightCycler® 480 II</b>
Reverse Transcription	10 min, 58 °C	10 min, 58 °C
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

**1.7 Detection channel Set-up**

<b>Real-time PCR device</b>	<b>Detection</b>	<b>Detection channel</b>	<b>Quencher</b>	<b>Note</b>
<b>Agilent Mx3005P</b>	SARS-CoV-2	FAM	+	
	ICR	HEX	+	
<b>Agilent AriaMx</b>	SARS-CoV-2	FAM	+	
	ICR	HEX	+	
<b>Applied Biosystems 7500 Fast</b>	SARS-CoV-2	FAM	None	Check the passive reference option ROX is none.
	ICR	VIC	None	
<b>Bio-Rad CFX96</b>	SARS-CoV-2	FAM	+	
	ICR	VIC/HEX	+	
<b>R-Biopharm RIDA®CYCLER</b>	SARS-CoV-2	green	+	
	ICR	yellow	+	
<b>Qiagen Rotor-Gene Q</b>	SARS-CoV-2	green	+	<b>Note:</b> Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	ICR	yellow	+	
<b>Roche LightCycler® 480 II</b>	SARS-CoV-2	465-510	+	
	ICR	533-580	+	

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 RNA Preparation

For RNA-preparation the use of SureFast® PREP DNA/RNA Virus (Art. No. F1051) is recommended.

The test assay contains an Internal Control RNA (ICR), which can either be used as PCR inhibition control or as positive extraction control for the sample preparation procedure and as a PCR inhibition control.

If the ICR is used only as a PCR inhibition control, 1 µl per reaction of the ICR should be added to the master-mix (see point 2.1.2).

If the ICR is used as an extraction control for the sample preparation procedure and as PCR inhibition control, 20 µl of the ICR should be added during extraction procedure. The ICR should always be pipetted to the sample with added Lysis Buffer and must not be added directly to the raw sample.

#### 2.1.2 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, negative extraction control, positive control. The test assay contains an Internal Control RNA (ICR), which can either be used as PCR inhibition control or as positive extraction control.

#### Reactions needed for the qualitative SARS-CoV-2 detection:

3 reactions for controls (1x negative control, 1x extraction control, 1x positive control)

For each sample: at least 1 reaction for each sample RNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

#### Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212,3 µl
Enzyme Mix	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20.0 µl</b>	<b>220.0 µl</b>

Mix each master-mix well and centrifuge shortly before use.

#### Example for the calculation and preparation of 10 reactions for ICR only as PCR inhibition control:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
ICR	1.0 µl	11.0 µl
Enzyme Mix	0.7 µl	7.7 µl
<b>Total volume</b>	<b>21.0 µl</b>	<b>231.0 µl</b>

Mix each master-mix well and centrifuge shortly before use.

### **2.1.3 Preparation of the real-time PCR-mix**

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Negative control: pipette 5 µl of PCR Water into the designated tubes/wells and close them.
- Pipette 5 µl of sample RNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

### **2.2 Interpretation of results**

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

SARS-CoV-2 RNA is detected in the FAM-channel. In the VIC/HEX-channel the internal amplification or positive extraction control (ICR) is detected. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the internal amplification control (IAC) is not needed to obtain a positive result of the positive control.

A sample is stated **positive** for SARS-CoV-2, if the sample RNA shows amplification in the FAM-channel.

A sample is stated **negative** for SARS-CoV-2, if the sample RNA shows no amplification in the FAM-channel and if the internal amplification or positive extraction control (VIC/HEX-channel) of the sample is **positive** with a shift in Cp-Value  $\leq 2$  compared to the negative control.

If the sample RNA in the VIC/HEX-Channel shows **no amplification** or a shift in Cp-value  $> 2$  compared to the negative control, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances RNA isolation and purification of the sample need to be improved. Alternatively the RNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific SARS-CoV-2 PCR assay.

## **3 Validation**

The external validation of the SureFast® SARS-CoV-2 PLUS real-time RT-PCR in combination with the SureFast® PREP DNA/RNA Virus was performed by the AOAC Performance Tested Methods (PTM) Program.

The results of the validation study demonstrated no differences between the SureFast® SARS-CoV-2 PLUS test method and the CDC (Center for Disease Control and Prevention) 2019-Novel Coronavirus Real-Time RT-PCR Diagnostic Panel method. The SureFast® SARS-CoV-2 PLUS real-time RT-PCR was approved for stainless steel environmental surface swabs.

## **4 Further Information**

### **4.1 Product Information**

- Detailed information about setup of several real-time PCR devices  
(Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Validation Report upon request
- AOAC certificate (license no. 022102)

### **4.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).