RIDA[®] Clean Extract

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RIDA[®]Clean Extract

1. Intended use

The RIDA[®]Clean Extract is a sample preparation kit for oil-containing emulsions, lubricants, paints and comparable solutions to make them suitable for measurement of ATP with LuciPac[®]Pen AQUA/Lumitester PD-20/PD-30. High concentrations of microorganisms which can have blight on the properties of these emulsions can be detected fast and efficiently using this system of ATP measurement subsequently to sample preparation with RIDA[®]Clean Extract.

2. General

Special microorganisms (bacteria or filamentous fungi) can grow in oil-containing liquids or lubricants or paints and change the nature of the liquid by their metabolism. This may lead to damages of e.g. motors, pumps or gears. The detection of the bacteria has to be done by classical cultivation on agar plates, yet. Results are very often available not earlier then after 2-5 days. Fast and appropriate counter-measures on contaminations are often not possible, anymore. Using RIDA[®]Clean Extract the detection of microbial contaminations can be done within a few minutes.

3. Test principle

As all living organisms bacteria and fungi contain the energy rich compound ATP (adenosine-triphosphate). This substance can be measured using a bioluminescence method. In order to do this ATP has to be released from microbial cells first by a lysis reagent. Released ATP will then be extracted from the oily or turbid liquid phase and transferred to the aqueous phase. The aqueous phase will be taken up by using the sampler stick of LuciPac[®]Pen AQUA and transferred into the test tube. Within the subsequent initiating enzymatic reaction (Luciferase-reaction) with the reagent in the tube, light will be emitted which can be measured with Lumitester PD-20/PD-30. The values will be given as “relative light units” (RLU).

4. Reagents provided

Each sample preparation kit contains sufficient reagents for preparation of 20 samples.

Component	Cap color	Format		Content
Reagent L1	Blue	Ready-to-use	Lysis buffer	20 pcs
Reagent L2	Red	Ready-to-use	Extraction solution	20 pcs
Pipettes		Disposable		20 pcs

5. Materials required but not provided

5.1. Equipment

- Lumitester PD-20 (will not be produced anymore!)
- Alternative: Lumitester PD-30 (Art. No. ZLT1402653)

5.2. Reagents

- LuciPac[®]Pen AQUA (Art. Nr. ZLA1002673)

6. Warnings and precautions for the users

Reagent tubes (L1 and L2) which have been used for extraction procedure will contain residues of treated sample liquids. They should be disposed according to current regulations for waste disposal of the respective liquid. Used LuciPac[®]Pen AQUA tubes can be disposed with the normal domestic waste.

7. Storage instructions

Store RIDA[®]Clean Extract at room temperature. The shelf life of the kit is mentioned as “Expiry” on the outer label. No quality guarantee is accepted after expiry of the kit.


8. Indication of instability or deterioration of reagents

The reagents don't show any visually recognizable advice on a possible decrease of reactivity.


9. Preparation of Samples

1. Take a sample portion using a Pasteur pipette and transfer it to reagent tube L1 (blue closure). Fill up the tube until the liquid is directly underneath the vertical grooves.
2. Close tube L1 tightly by screwing up the closure and mix intensively (shake vigorously or vortex approx. 45 sec). Keep the tube closed and allow the liquid to stand for 1 min to reach a complete lysis.
3. Transfer the total content of tube L1 into tube L2 (red closure) and close tightly by screwing up the closure.
4. Mix intensively for approx. 30 sec (shake vigorously or vortex).
5. Centrifuge tube L2 for 1 min at 6.000 to 13.000 rpm. Alternatively allow the liquid to stand until the phases have separated completely.
6. In case of preparation of oils, cooling lubricants/lubricants, fats and aqueous solutions discard the upper phase completely by decantation. As the lower, clear phase is strongly viscous it will remain in the tube.
In case of preparation of paints and other turbid liquids the clear phase will go up to the top so that decantation is not necessary.
7. Perform sampling from the prepared solution using the sampler stick of LuciPac[®] Pen AQUA as described under “10. Test implementation”.

Flow chart in case of oils, lubricants, fats and aqueous solutions:

Flow chart in case of paints and other turbid liquids:

10. Test implementation

1. In case of oils, cooling lubricants/lubricants, fats and aqueous solutions dip the sampler stick of LuciPac[®]Pen AQUA into the remaining strongly viscous phase until the grooves of the sampling head are wetted completely.

In case of paints and other turbid liquids dip the sampler stick into the upper phase directly until the grooves of the sampling head are wetted completely.

2. Pull out sampler stick and transfer it back to the test tube of LuciPac[®]Pen AQUA. Press it in to the stop. Shake down the reaction liquid into the reaction chamber and rotate tube slightly to solve the reagent (s. instructions for LuciPac[®]Pen AQUA). Perform measurement as described in the instructions for LuciPac[®]Pen AQUA or in the Lumitester PD-20/PD-30 manual, respectively.

11. Results

Warning color code for lubricants, water based emulsions, oily emulsions and paints

Approx. limit value (CFU/ml)	Microbial contamination	Water quality/actions	Color code	Approx. RLU values (undiluted)
$<10^2$	very low	normal / no actions	green	20 - 1.000 level A: passed
$10^3 - 10^4$	medium - strong	warning range / increased control / cleaning actions	orange	1.000 - 10.000 level B: caution - failed
$>10^5$	high - very high	critical / urgent need for immediate action	red	> 10.000 level C: critical

RLU values according to the warning color code come along as the fastest possible guidance to microbial loads and corresponding problems. The shown correlation between RLU values and classical cfu values are representative average values which derive from several measurements and analyses with real life samples.

As RLU values may differ due to type and growth activity of microorganisms the border areas in the table (e.g. 1000 – 10000) are not exactly defined. Values which are more or less „between“ two different warning colors should be taken skeptically as the real amount of microorganisms in the sample may be different somehow. In such cases additional samples should be analyzed verify in which true condition the emulsion currently is.

Additionally every plant that performs ATP measurements in lubricants, paints or comparable liquids should carry out individual calibrations for his specific systems. This makes it possible to correlate the measured RLU values with individual loads of microorganisms.

General advices:

As long as the RLU values are low there is no serious load present. In the case of increased values appearing the reasons for that have to be checked out and it should be clarified if mainly bacteria or fungi are responsible. In cases of critical values it makes sense to perform classical cultivation based methods to determine the cfu value and to identify the microbial species to estimate their hazardous potential.

Additional plant internal experience and measured values as pH, frothing, odor development as well as turbidity, are criteria that may help to decide if urgent action is required or not.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG