


# Lateral flow tests from bioavid with hook line

Reliable identification also of high positive samples

## Field of application

- Hygiene control
- CIP water analysis
- Food analysis



## Benefits



**Extended functionality: included hook line for safe identification of high positive samples**



**Quick and reliable results in 10 minutes**



**All included: all further reagents required for the assay – including a positive control – are contained in the test kit**



**Specifically developed accessories available for hygiene control (swabbing kits) (Art. No. BS800-25 and BS801-25)**




**Also suitable for non-laboratory staff**

## New functionality – tracking the hook effect

The hook effect (or “overload”-effect; Standard EN 15633-1. 2019) is observed when a very high amount of an analyte is present in the sample. In this case, the analyte’s amount exceeds the amount of color-labelled antibodies and the ratio between both needed for the formation of the test band becomes imbalanced. Although the analyte is present, only a faint or even no test band is visible. Hence, there is risk for a falsely low or negative interpretation of a test result. This effect becomes visible by the hook line. A missing hook line indicates high allergen content in a sample. In this case, test users must not interpret a missing test line as negative.

1 to 3 purple lines may appear in the reaction field on the strip. The upper line is the control line **C**. It indicates that the test has been performed correctly. The middle line is the hook line **H**. It is visible in a valid test if the allergen concentration does not exceed ca. 1,000 or 10,000 ppm (depending on test). The lowest line is the test line **T**. It indicates the presence of the target residue in the sample in the range of approx. 1 to 1,000 ppm.



**C and H line present, no test line:**  
negative result (A)


**3 lines present:** positive result (B)

**C line and T line present,**  
**Hook line faint or absent:** high positive result  
(≥ 1,000 ppm or 10,000 ppm) (C)


**Only C line present:** suspected high positive,  
repeat test with higher sample dilution (D)

**No C line:** Test is invalid

## 1 Sample preparation



## 2 Test procedure



## Food allergies

In the case of a food allergy, certain parts of a food (allergens) trigger an immune reaction in humans. This immune reaction causes a wide range of symptoms (e.g. slight prickling in the mouth, rash, anaphylactic shock). Worldwide, about 2 - 8 % of people suffer from a food allergy. Since there is currently no effective therapy for food allergies,

sensitized persons must follow an allergen-free diet. Therefore, labelling allergens on food is essential. For prevention of unintended contamination, food industry must consider a respective allergen management (e.g. testing of raw materials, checking production lines). Especially for the last aspect, lateral flows offer a fast, easy and reliable method.

## Fast replies - Lateral Flow Devices (LFDs)

LFDs are a simple and fast technique combining chromatography and an immunological method to detect the presence or absence of an analyte (here:

allergen). Within 10 minutes, a reliable result can be observed without the need for elaborate equipment. Thus, they are applicable everywhere.

## Lateral flow tests including hook line

Product	Description	No. of tests/amount	Art. No.
<b>Almond</b>			
bioavid Lateral Flow Almond	Qualitative detection of Almond; Hook line upper limit: 1,000 mg/kg (ppm)	20 strips	BLH701-20
<b>Coconut</b>			
bioavid Lateral Flow Coconut	Qualitative detection of Coconut; Hook line upper limit: 10,000 mg/kg (ppm)	20 strips	BLH700-20
<b>Hazelnut</b>			
bioavid Lateral Flow Hazelnut	Qualitative detection of Hazelnut; Hook line upper limit: 1,000 mg/kg (ppm)	20 strips	BLH704-20
<b>Mustard</b>			
bioavid Lateral Flow Mustard	Qualitative detection of Mustard; Hook line upper limit: 1,000 mg/kg (ppm)	20 strips	BLH703-20
<b>Peanut</b>			
bioavid Lateral Flow Peanut	Qualitative detection of Peanut; Hook line upper limit: 1,000 mg/kg (ppm)	20 strips	BLH706-20

Detection limit for all tests is approx. 1 mg/kg (ppm) allergen; result is available after 10 minutes.