# RIDASCREEN® Enro-/Ciprofloxacin

**REF R3111**

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# RIDASCREEN® Enro-/Ciprofloxacin

## Brief information

RIDASCREEN® Enro-/Ciprofloxacin (R3111) is a competitive enzyme immunoassay for the quantitative analysis of enrofloxacin and its active metabolite ciprofloxacin in milk (raw milk, fresh milk and UHT milk) and milk powder (skim milk powder and whole milk powder).

All reagents required for performing the enzyme immunoassay, including the standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 duplicate determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	Milk: homogenize and use directly in the test. Milk powder: dissolve, homogenize, and use in the test.
Time requirement:	Test procedure (incubation time)..... 45 min
Limit of detection:	Fresh milk, UHT milk, raw milk, milk powder: 0.6 µg/l
Detection capability (CC $\beta$ ):	Fresh milk, UHT milk, raw milk, milk powder: 0.8 µg/l
Recovery rate:	Of Ciprofloxacin..... 85 % - 105 % Of 50% Ciprofloxacin + 50% Enrofloxacin ..... 70 % - 80 %
Specificity:	Ciprofloxacin .....100 % Enrofloxacin .....91 % Sarafloxacin .....17 % Danofloxacin .....15 % Norfloxacin .....15 % Difloxacin .....12 % Marbofloxacin.....12 % Ofloxacin, Flumequine, Oxolinic acid..... < 10 %

The specificity of the RIDASCREEN® Enro-/Ciprofloxacin assay was determined by analyzing the cross-reactivity to the respective substances in the buffer system. In samples, the specificity may differ from the values determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection and the recovery rate in the relevant matrix. The test cannot distinguish between analytes and cross-reactive substances.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice. Such lists minimum standards concerning the framework conditions when using test kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website [www.r-biopharm.com/products/food-feed-analysis](http://www.r-biopharm.com/products/food-feed-analysis).

## **Product offer**

RIDA® Ciprofloxacin Spiking Solution (Art. No. R3198)

### **1. Intended use**

RIDASCREEN Enro-/Ciprofloxacin is a competitive enzyme immunoassay for the quantitative analysis of enrofloxacin and its active metabolite ciprofloxacin in milk (raw milk, fresh milk and UHT milk) and milk powder (skim milk powder and whole milk powder).

### **2. General**

The mechanism of action of fluoroquinolones is based on the inhibition of gyrase, an enzyme essential to bacterial DNA synthesis. Thus, fluoroquinolones have very high efficacy against most gram-negative and gram-positive bacteria and are widely used in veterinary medicine to treat cattle, pigs, and poultry.

The use of fluoroquinolones has increased in recent years to prevent infectious diseases from occurring during livestock production. Particularly in the group of fluoroquinolones, the development of resistance was found to be associated mostly with *E. coli*.

Therefore, in the European Union, maximum residue limits (MRL) apply to some substances of this group, as described in COMMISSION REGULATION (EU) No. 37/2010. One representative of this group approved in human and veterinary medicine is enrofloxacin (metabolite: ciprofloxacin). The MRL for the sum of enrofloxacin and ciprofloxacin in milk (cow, sheep, or goat) is 100 µg/kg.

### 3. Test principle

The ELISA test system is based on the antigen-antibody reaction. The wells of the microtiter strips are coated with capture antibodies and anti-ciprofloxacin antibodies. Standards, the sample solutions and enzyme conjugate are added. Free antigen (enrofloxacin/ciprofloxacin) and enzyme conjugate compete for the antibody binding sites (competitive enzyme immunoassay). Any unbound enzyme conjugate is then removed in a washing step. For detection, substrate/chromogen is added to the wells and incubated. Bound enzyme conjugate converts the chromogen into a blue end product. Adding the stop solution causes a color change from blue to yellow. The measurement is made photometrically at 450 nm; the extinction is inversely proportional to the antigen concentration (enrofloxacin / ciprofloxacin) in the sample.

### 4. Reagents provided

Each kit contains sufficient materials for a maximum of 48 duplicate determinations (including standard analyses). Each test kit contains:

Component	Cap color	Format	Volume
<b>Microtiter plate</b>	-	Ready-to-use	96 wells
<b>Dilution buffer</b>	Transparent	Ready-to-use	60 ml
<b>Standard 1</b>	White	Ready-to-use	0 µg/l 1.3 ml
<b>Standard 2</b>	White	Ready-to-use	0.25 µg/l 1.3 ml
<b>Standard 3</b>	White	Ready-to-use	0.75 µg/l 1.3 ml
<b>Standard 4</b>	White	Ready-to-use	2.25 µg/l 1.3 ml
<b>Standard 5</b>	White	Ready-to-use	6.75 µg/l 1.3 ml
<b>Standard 6</b>	White	Ready-to-use	20.25 µg/l 1.3 ml
<b>Wash buffer salt Tween</b>	-	Salt for dissolving	
<b>Conjugate</b>	Red	Ready-to-use	7.5 ml
<b>Substrate/ Chromogen</b> Red Chromogen Pro	Brown	Ready-to-use	13 ml
<b>Stop solution</b>	Yellow	Ready-to-use	14 ml

## **5. Reagents required but not provided**

### 5.1 Equipment

- Microtiter plate spectrophotometer (450 nm)
- Shaker
- Graduated pipettes
- Variable 20 - 200 µl and 200 - 1000 µl micropipettes

### 5.2 Reagents

- Distilled or deionized water

## **6. Warnings and precautions for the users**

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

Ensure the proper and responsible disposal of all reagents and materials after their use. For disposal, please adhere to national regulations.

## **7. Storage instructions**

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

The red-colored substrate/ chromogen is photosensitive, therefore, avoid direct exposure to light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

- Bluish coloration of the red-stained substrate/chromogen prior to test implementation
- A value of less than 0.8 absorbance units ( $A_{450\text{ nm}} < 0.8$ ) for standard 1.

## 9. Preparation of Samples

The samples should be stored in a cool place, protected against light.

### 9.1 Milk (raw milk, fresh milk, UHT milk)

- Use 50 µl of the sample in each well directly in the test
- Dilution factor 1

### 9.2 Milk powder (skim milk powder, whole milk powder)

- Dissolve the milk powder 1:10 in demineralized water, e.g. 10 g milk powder + 90 ml demineralized water
- Shake / vortex until powder is completely dissolved
- Use 50 µl of dissolved milk powder per well in the test
- Dilution factor 1

## 10. Test implementation

### 10.1 Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

As wash buffer, a PBS Tween buffer is needed; please use the enclosed buffer salt (see section 4. Reagents provided). To prepare the buffer, dissolve the entire contents of the pouch in 1 l distilled water. The dissolved wash buffer can be stored for approximately 4 to 6 weeks at 2 - 8 °C.

Alternative: Dissolve the contents of the pouch in 100 ml distilled water (10-fold concentrate). The solution can be stored for approximately 8 - 12 weeks at room temperature (20 - 25 °C). To prepare the ready-to-use solution, mix 1 part of the 10-fold concentrate with 9 parts of distilled water.

Return all reagents to 2 - 8 °C immediately after use.

## 10.2 Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Pipette 50 µl each of the standard and the samples (see section 9. Preparation of Samples) as duplicates into the respective wells. Use new pipette tips for each standard or sample.
3. Pipette 50 µl each of conjugate into the respective wells, mix carefully by shaking the plate per hand, and incubate for 30 minutes at room temperature (20 - 25 °C).
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times) against absorbent paper to ensure complete removal of liquid from the wells. Wash each well using 250 µl diluted wash buffer (see 10.1 Test preparation). Repeat the washing step another two times.
5. Pipette 100 µl each of substrate/ chromogen into the wells. Mix carefully by hand and incubate in the dark for 15 minutes at room temperature (20 - 25 °C).
6. Pipette 100 µl stop solution into each well and mix carefully by hand. Measure the absorbance at 450 nm within 15 minutes after adding the stop solution.

## 11. Results

A special software, the **RIDASOFT® Win.NET** (Art. No. Z9996FF), is optionally available for evaluation of the RIDASCREEN® enzyme immunoassays. The course of the standard curve is shown in the Quality Assurance Certificate (certificate of analysis) enclosed in the test kit.

### Interpretation of results

In order to obtain the enrofloxacin / ciprofloxacin concentration in µg/l (µg/kg) actually contained in a sample, the concentration read from the standard curve must be further multiplied by the corresponding dilution factor. When working according to the instructions, the dilution factors are as follows:

Milk .....	1
Milk powder (referring to dissolved milk powder) .....	1

## Recommendation

In order to ensure a high analytical performance we recommend:

- To analyze each sample material in duplicates.
- To do spike experiments to ensure an accurate and correct test procedure.

Samples measured outside the measuring range at  $> 20 \mu\text{g/l}$  should be diluted 1:10 up to maximal 1:30 with the enclosed sample dilution buffer.

## Further application notes

A summary of the performance data for the ThunderBolt<sup>®</sup> analyzer is available on request.

**Further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## Version overview

Version number	Chapter and title
2020-09-15	Release version

## Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.