r-biopharm®

RIDASCREEN® Total Gluten

Art. No. R7041

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RIDASCREEN[®] Total Gluten

Brief information

RIDASCREEN[®] Total Gluten (R7041) is a sandwich enzyme immunoassay for the quantitative determination of intact (non-hydrolyzed) gluten from gluten containing cereals (wheat, rye and barley) in oat and oat products.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization and extraction
Time requirement:	sample preparation with Cocktail (patented) for 10 samples..... approx. 2 h test implementation (incubation time)50 min
Standard material:	The RIDASCREEN [®] standard material is a total gluten extract from four wheat varieties. The results with this standard material are traceable to the oat samples described in the AOAC SMPR [®] 2017.021.
Limit of detection:	4 mg/kg (ppm) gluten depending on the tested sample matrix
Limit of quantification:	5 mg/kg (ppm) gluten
Specificity:	The used monoclonal antibodies detect gliadins from wheat and corresponding prolamins from rye and barley, high-molecular-weight (HMW) glutenin subunits (GS) from wheat and HMW-secalins from rye as well as low-molecular-weight (LMW)-GS from wheat (see also 2.).

There are not any known cross reactivities.

Cross reactivities of the antibodies used have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP)-Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products for gliadin/gluten determination

RIDASCREEN® Gliadin (R7001)
RIDASCREEN® FAST Gliadin (R7002)
RIDASCREEN® Gliadin competitive (R7021)
RIDASCREEN® FAST Gliadin Sensitive (R7051)
Cocktail (patented) (R7006 or R7016)
Cocktail ECO (R7080)
RIDA® Extraction Solution (colorless) (R7098)
Set of 3 processed Gliadin Assay Controls (R7012)
RIDA® QUICK Gliadin (R7003/R7004/R7005)
SureFood® Allergen Gluten (S3606)

1. Intended use

RIDASCREEN® Total Gluten is a sandwich enzyme immunoassay for the quantitative determination of intact (non-hydrolyzed) gluten from gluten containing cereals (wheat, rye and barley) in oat and oat products.

2. General

Traditionally, cereal proteins are subdivided according to their solubility (Osborne fractionation):

Soluble in	Fraction
Water	Albumins (non gluten proteins)
0.5 M NaCl	Globulins (non gluten proteins)
60% Ethanol	Prolamins (gluten proteins)
Insoluble in 60% Ethanol	Glutelins (gluten proteins)

Unfortunately, the subdivision according to solubility is not an accurate criterion, since co-solubilization and co-precipitation occur frequently. Furthermore, additional factors such as temperature can have a high influence on the result. Especially the subdivision into prolamins and glutelins (which form together the gluten proteins) according to their solubility is very difficult. Therefore, the modern subdivision of the gluten proteins into protein families according to their electrophoretic mobility and their amino acid sequences is much more accurate.

Due to high sequence homologies between gluten proteins within one cereal species and between different species and due to the combination of four different monoclonal antibodies, the RIDASCREEN® Total Gluten enables the detection of almost all gluten proteins:

Wheat			
Protein family	Fraction	Mean content of gluten protein*	Primarily reactive with
α/β-gliadins	Mainly prolamins	33%	R5 antibody
γ-gliadins	Mainly prolamins	27%	R5 antibody
ω1,2-gliadins	Mainly prolamins	4%	R5 antibody
ω5-gliadins	Mainly prolamins	3%	R5 antibody (weak)
LMW-glutenin-subunits	Mainly glutelins	22%	LMW 1 & 2 antibody
HMW-glutenin-subunits	Mainly glutelins	11%	HMW antibody

Rye			
Protein family	Fraction	Mean content of gluten protein*	Primarily reactive with
ω-secalins	Mainly prolamins	18%	R5 antibody
γ40k-secalins	Prolamins and glutelins	25%	R5 antibody
γ75k-secalins	Prolamins and glutelins	48%	R5 antibody
HMW secalins	Mainly glutelins	9%	HMW antibody

Barley			
Protein family	Fraction	Mean content of gluten protein*	Primarily reactive with
B-hordeins	Mainly prolamins	27%	R5 antibody
C-hordeins	Mainly prolamins	36%	R5 antibody
γ-hordeins	Mainly prolamins	32%	R5 antibody
D-hordeins	Mainly glutelins	5%	No detection

* Significant differences between different cultivars of the respective cereal and harvest years possible. Data from Wieser et al. Celiac Disease and Gluten (2014) Elsevier Inc. Amsterdam, ISBN 978-0-12-420220-7, page 107

The use of wheat flour and gluten in foodstuffs is extremely common because of their heat stability and useful effects on e.g. texture, moisture retention and flavor. Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

According to the Codex Alimentarius (CODEX STAN 118/1979), two categories for labeling of food according to the gluten content exist:

- 1) Food products which contain less than 20 mg/kg gluten can be labeled as "gluten-free".
- 2) Food products labeled as "very low gluten" can have a gluten content above 20 and up to 100 mg/kg gluten.

The threshold of 20 mg/kg gluten has been adopted by national legislations in many countries.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies against gluten proteins. After adding the standard solutions or samples to the wells, present gluten will bind to the specific antibodies. The result is an antibody-antigen-complex. Components not bound by the antibodies are removed in a washing step. Then, antibodies conjugated to peroxidase are added. These conjugates are bound to the antibody-antigen-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound conjugates are then removed in a washing step. Enzyme substrate and chromogen are added to the wells and incubated. The bound conjugate converts the substrate/chromogen into a blue product. Addition of stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the gluten concentration of the sample.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format	Concentration	Volume
Microtiter plate		Ready to use		96 wells
Buffer	Transparent	Ready to use		110 ml
Standard 1	Transparent	Ready to use	0 mg/kg gluten*	1.3 ml
Standard 2	Transparent	Ready to use	5 mg/kg gluten*	1.3 ml
Standard 3	Transparent	Ready to use	10 mg/kg gluten*	1.3 ml
Standard 4	Transparent	Ready to use	20 mg/kg gluten*	1.3 ml
Standard 5	Transparent	Ready to use	40 mg/kg gluten*	1.3 ml
Standard 6	Transparent	Ready to use	80 mg/kg gluten*	1.3 ml
Wash buffer	Brown	Concentrate	10x	100 ml
Conjugate	Red	Ready to use		11 ml
Substrate/Chromogen Red Chromogen Pro	Brown	Ready to use		13 ml
Stop Solution	Yellow	Ready to use		14 ml

* The dilution factor 1000 from the sample preparation has already been considered when labeling the standards. Therefore, the gluten concentration of samples can directly be read from the standard curve.

5. Reagents required but not provided

5.1. Equipment:

- Microtiter plate spectrophotometer (450 nm)
- Centrifuge, centrifugal vials (e.g. Brand 10742512)
- Shaker
- Water bath (50 °C / 122 °F)
- Laboratory mincer / grinder, pestle and mortar, ultra-turrax or homogenizer
- Graduated pipettes
- Variable 20 µl - 200 µl and 200 - 1000 µl micropipettes
- Gloves
- Optional: RIDASOFT® Win.NET (Z9996)

5.2. Reagents:

- Distilled or deionized water
- Gluten-free skim milk powder (food quality)
- Cocktail (patented); Art. No.: R7006 (105 ml) or R7016 (1000ml)
- Ethanol solution (80 %): i.e. add 120 ml ethanol p.a. to 30 ml distilled water and shake well

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instructions for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Regarding the handling of ethanol, please refer to the safety instructions of the respective manufacturer.

The Cocktail (patented) contains β-mercaptoethanol. It is recommended to work under a chemical hood

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

To avoid moisture inside the wells, open the foil bag for withdrawal of microwells only after having reached room temperature (20 - 25 °C / 68 - 77 °F).

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided, and further store at 2 - 8 °C (35 - 46 °F).

The reddish substrate/chromogen is light sensitive. Therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- Any bluish coloration of the reddish substrate/chromogen before pipetting into microplate wells
- A value of less than 1.2 absorbance units ($A_{450\text{ nm}} < 1.2$) for standard 6

9. Preparation of samples

Airborne cereal dust and dirty laboratory equipment may lead to gluten contamination of the assay. Hence it is recommended to take the following precautionary measures. Before starting and during the assay wear gloves.

- Clean surfaces, glass vials, mincers and other equipment with 40 - 60 % ethanol or 40 - 60 % 2-propanol
- Carry out the sample preparation in a room isolated from the ELISA procedure
- Check for gliadin / gluten contamination of reagents and equipment with the test strips RIDA[®]QUICK Gliadin (R7003, R7004, or R7005)
- It is recommended to **work under a chemical hood**, because of β-mercaptoethanol in the Cocktail (patented)
- β-mercaptoethanol can disturb the ELISA; therefore, dilute the samples at least 1:1000
- Try to homogenize with a clean spatula if samples gelatinizing during preparation

9.1. General extraction with the Cocktail (patented) (R7006 / R7016) for oat and oat products

- The gluten distribution in oat samples can be very inhomogeneous; furthermore, the samples are difficult to homogenize; therefore, grind and homogenize well at least 200 g
- Weigh in 1 g hereof to a sufficient vial (e.g. centrifugal vials with lid) and add 10 ml of the Cocktail (patented), close the vial and mix well; pay attention to obtain a homogeneous suspension
- Add 30 ml 80 % ethanol (see 5.2.), close the vial, and mix well; pay attention to obtain a homogeneous suspension

The weighted sample can be further increased (>1 g) in case of very inhomogeneous samples. The volume of Cocktail (patented) and 80% ethanol must be increased accordingly. Choose a properly sized vial (e.g. Nalgene bottle 250 ml).

- Incubate for 40 min at 50 °C (122 °F) in a water bath
- Shake for 1 h upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F)
- Centrifugation of the extract: 10 min, at least 2500 g, at room temperature (20 - 25 °C / 68 - 77 °F); alternatively, the extract can be filtered
- Put the particle-free supernatant into a screw top vial for further processing or storage (see remark below)

If it is not possible to achieve a particle-free supernatant with the described centrifugation, it will be necessary to filter the supernatant additionally; alternatively, 2 ml of the extract can be transferred into a reaction tube and centrifuged by using a microcentrifuge (10 min, 10000 g).

- Dilute the sample 1:25 (40 µl + 960 µl) for ELISA testing with the buffer from the test kit: the final dilution factor is 1000
- Use 100 µl per well of the diluted extract within 30 min for ELISA testing

9.2. Extraction with the Cocktail (patented) (R7006 / R7016) for tannin- and polyphenol-containing oat products (e.g. oat products with high content of chocolate, coffee, cocoa, chestnut flour, teff flour, buckwheat, millet or spices)

- The gluten distribution in oat samples can be very inhomogeneous; furthermore, the samples are difficult to homogenize; therefore, grind and homogenize well at least 200 g

- Weigh in 1 g thereof to a sufficient vial (e.g. centrifugal vials with lid), add 1 g skim milk powder and additionally 10 ml of the Cocktail (patented); close the vial and mix well; pay attention to obtain a homogeneous suspension
- Add 30 ml 80 % ethanol (see 5.2.), close the vial and mix well. Please pay attention to obtain a homogeneous suspension

The weighted sample can be further increased (>1 g) in case of very inhomogeneous samples. The volume of Cocktail (patented) and 80% ethanol as well as the amount of skim milk powder must be increased accordingly. Choose a properly sized vial (e.g. Nalgene bottle 250 ml).

- Incubate for 40 min at 50 °C (122 °F) in a water bath
- Shake for 1 h upside down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F)
- Centrifugation of the extract: 10 min, at least 2500 g, at room temperature (20 - 25 °C / 68 - 77 °F); alternatively, the extract can be filtered
- Put the particle-free supernatant into a screw top vial for further processing or storage (see remark below)

If it is not possible to achieve a particle-free supernatant with the described centrifugation, it will be necessary to filter the supernatant additionally; alternatively, 2 ml of the extract can be transferred into a reaction tube and centrifuged by using a microcentrifuge (10 min, 10000 g).

- Dilute the sample 1:25 (40 µl + 960 µl) for ELISA testing with the buffer from the test kit: the final dilution factor is 1000
- Use 100 µl per well of the diluted extract within 30 min for ELISA testing

Remark for 9.1. and 9.2.:

The undiluted supernatant obtained after centrifugation or the filtrate can be stored in a tightly closed vial in the dark at room temperature (20 - 25 °C / 68 - 77 °F) for up to two weeks.

If samples have a higher absorbance than standard 6, they should be further diluted and retested in the ELISA. For this, dilute the extract again 1:25 with buffer as described in 9.1 and 9.2. Then perform the additional dilution of the diluted sample extract with the following buffer:

- 1% Cocktail (patented)
- 3% of an 80% ethanol solution
- 96% buffer

This equals to e.g. 50 µl Cocktail (patented), 150 µl 80% ethanol and 4800 µl buffer.

Using this buffer, the composition of the sample medium is identical to that used for extraction as described in 9.1 or 9.2.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The washing buffer is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml dist. water). Prior to dilution, dissolve any crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 20 - 25 °C (68 - 77 °F) for four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than 3 strips (24 wells) at a time. This avoids a shift by pipetting and guarantees consistent incubation times. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one Cat.-No. 655101 or Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific) should be used as a pre-plate. All standards and samples are pipetted into the uncoated plate first (at least 150 µl per well) and then quickly transferred to the coated microtiter plate with an 8-channel pipette.

It is recommended to pipette the conjugate, the substrate/chromogen and the stop solution with a multi-channel or stepper pipette to avoid a time shift over the plate.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard or sample to separate duplicate wells and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see 10.1.) and pour out the liquid again. Repeat two additional times.

4. Add 100 µl of the conjugate to each well and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl diluted wash buffer (see 10.1.) and pour out the liquid again. Repeat two additional times.
6. Add 100 µl of substrate/chromogen to each well and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 min after addition of stop solution.

11. Results

Special software RIDASOFT® Win.NET (Art. No. Z9996) is available for evaluation of the RIDASCREEN® enzyme immunoassays. The calculation should be done by use of a 4-parameter function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. The document 'Compliance Criteria' which is available upon request provides criteria for evaluation of standard curves. For quality assurance, assay controls should be used.

Compared to the certificate, higher absorbance values ($A_{450\text{ nm}}$) for the standard curve, especially for the zero standard, may be a result of insufficient washing or gluten contamination.

A further dilution and new detection of the samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 6. Please use the buffer mentioned for further sample dilution (see remark for 9.1. and 9.2.).

Please note:

When working in accordance with the regulation stated, the dilution factor is 1000. The gluten concentration can directly be read from the standard curve (the sample dilution factor of 1000 is already taken into account (see 4. Reagents provided).

For sample dilutions of more than 1:1000, the additional dilution factor must be considered for calculation of the gluten concentration.

In General:

Samples tested negative may still contain a gluten contamination below the limit of detection of the assay, or they might contain other components such as starch for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity. When analyzing a non-validated matrix (see validation report), it is recommended to verify the results by spiking experiments.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the Validation Report.

Recommendations for guarantee of high analytical performance:

Adhere to the specified temperature of 20 - 25 °C (68 - 77 °F) when running the assay. Deviations hereof may influence the analyte recovery.

Each sample extract should be analyzed in duplicate.

Use also gluten free and gluten containing (spiked) samples as test controls.

To ensure an accurate result, spike experiments are recommended.

Strong acidic or alkaline samples have to be neutralized prior to analysis.

Confirm results with PCR (e.g. SureFood[®] Allergen Gluten, Art. No. S3606).

In the production of foods such as beer or sourdough, proteins are fragmented. In sandwich ELISAs protein fragments lead to a reduced recovery. Such samples must be analyzed with a competitive ELISA test systems like the RIDASCREEN[®] Gliadin competitive (R7021).

Further applications:

–None

Details about automated assay procedures and further product information are available on request from your local distributor or R-Biopharm AG.

The data correspond to our present state of technology and provide information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

R-Biopharm AG

Postanschrift / Postal Address:

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64297 Darmstadt, Germany

Sitz / Corporate Seat: Pfungstadt

Tel.: +49 (0) 61 51 - 81 02-0

Fax: +49 (0) 61 51 - 81 02-40

E-mail: info@r-biopharm.de

www.r-biopharm.com

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are registered trademarks of R-Biopharm AG