

AFLATOXIN B1 ELISA

5121AFB[20]04.20

**A competitive enzyme immunoassay
for quantitative analysis of Aflatoxin B1
in food and feed samples**

EUROPROXIMA AFLATOXIN B1 ELISA

A competitive enzyme immunoassay for quantitative analysis of Aflatoxin B1 in food and feed samples

TABLE OF CONTENTS


	PAGE:
Brief Information.....	2
1. Introduction.....	2
2. Principle of the Aflatoxin B1 ELISA.....	3
3. Specificity and Sensitivity.....	3
4. Handling and Storage.....	4
5. Kit contents.....	5
6. Equipment and materials required but not provided.....	6
7. Precautions.....	6
8. Sample preparation.....	7
9. Preparation of reagents.....	9
10. Assay Procedure.....	10
11. Interpretation of results.....	11
12. Literature.....	13
13. Ordering information.....	13
14. Revision history.....	13

The quality management system of R-Biopharm Nederland B.V. is
ISO 9001:2015 certified

BRIEF INFORMATION

The Aflatoxin B1 ELISA is a competitive enzyme immunoassay for the screening of food and feed samples on the presence of this mycotoxin. The test is based on antibodies directed against Aflatoxin B1. The ELISA kit contains a 96 microtiter plate as well as all essential reagents including ready-to-use standards to perform the test. Methods for a fast and efficient extraction of aflatoxin B1 from different matrices are included in the kit manual.

1. INTRODUCTION



Chemical structure of Aflatoxin B1

The Aflatoxins (Aflatoxins B₁, B₂, G₁, G₂) are a group of structurally related, extremely toxic compounds produced by the moulds *Aspergillus flavus*, *A. paraciticus* and *A. nomius*. Aflatoxins cause cancer, mainly of the liver but also of the gut, lungs and breast. These moulds occur on food, feed and their ingredients, derived from tropical and sub-tropical areas. The most pronounced contamination has been found in cereals, rice, maize, soy, tree nuts and peanuts.

Within the EU, maximum levels (MLs) for a number of mycotoxins have been set by Commission Regulation (EC) 1881/2006 (which has replaced earlier Regulations (EC) No 1525/98, 194/97 and 466/2001). Concerning the Aflatoxins, this Regulation has been amended by Commission Regulation (EC) No 165/2010. In this amending Regulation, changes in the MLs were needed for alignment with Codex Alimentarius FAO/WHO standards, and new groups of foodstuffs (oilseeds and apricot kernels) are specified. The MLs for Aflatoxin B1 in most of the foodstuffs for direct human consumption vary from 2 to 12 µg/kg (ppb). However, for processed cereal-based foods and baby foods for infants and young children as well as for dietary foods for special medical purposes intended specifically for infants, the ML for Aflatoxin B1 is set at 0.1 µg/kg (0.1 ppb). Especially for this low measuring range EuroProxima can provide the Aflatoxin B1 sensitive ELISA (5121AFBS).

2. PRINCIPLE OF THE AFLATOXIN B1 ELISA

The kit is based on a microtiter plate (12 strips, each 8 wells), precoated with rabbit antibodies to mouse IgG. In one incubation step, specific antibodies (mouse anti-Aflatoxin), enzyme labelled Aflatoxin (enzyme conjugate) and Aflatoxin B1 standards or sample are added to the precoated wells. The specific antibodies are bound by the immobilised rabbit antibodies and at the same time free Aflatoxins (in the standard solution or in the sample) and enzyme labelled Aflatoxin compete for the specific antibody binding sites (competitive enzyme immunoassay).

After an incubation time of one hour, the non-bound (enzyme labelled) reagents are removed in a washing step. The amount of bound enzyme conjugate is visualised by the addition of chromogen substrate (tetramethylbenzidine, TMB). Bound enzyme transforms the chromogen into a coloured product.

The substrate reaction is stopped by the addition of sulfuric acid. The colour intensity is measured photometrically at 450 nm and is inversely proportional to the Aflatoxin B1 concentration in the sample.

3. SPECIFICITY AND SENSITIVITY

The Aflatoxin B1 ELISA utilizes antibodies raised in mouse against protein conjugated Aflatoxin. The reactivity pattern of the antibody is:

Cross- reactions:	Aflatoxin B1	100%
	Aflatoxin B2	20%
	Aflatoxin G1	17%
	Aflatoxin G2	4%

The cross-reactivities are determined in a buffer system. The reported values may be different in samples due to matrix effects.

The test cannot discriminate between analytes and cross-reactive substances

The Limit of detection (LOD) are determined under optimal conditions. Cut-off criteria need critical consideration.

Matrix	Procedure	LOD ppb
cereals	8.1	0.5
rice	8.1	0.4
egg	8.1	0.2
nuts	8.1	0.8
honey	8.1	0.2
mashed fruits	8.1	0.6
edible oils	8.3.2	0.2
feed	8.1	1.0

If the sample is found to be non-compliant, the results shall be verified by re-analysis of the sample using a confirmatory method.

4. HANDLING AND STORAGE

- Kit and kit components should be stored at 2°C to 8°C in a dark place. For repeated use store kit components as specified under chapter 9.
- After the expiry date of the kit and/or components has passed, no further quality guarantee is valid.
- Bring all kit components including the microtiter plate to ambient (room) temperature before use.
- Dilute the kit components immediately before use, but after the components are brought to ambient temperature.
- Avoid condensation in the wells of the plate. Bring the sealed plate to ambient temperature before opening the plate sealing.
- The substrate chromogen solution can be stored in a refrigerator (2°C to 8°C) until the expiry date stated on the label.
- Exposure of the chromogen solution to light should be avoided.

Degeneration of the reagents may have occurred when the following phenomena are observed:


- A blue colouring of the chromogen solution before transferring it into the wells.
- A weak or no colour reaction in the zero standard wells ($E_{450nm} < 0.8$).

5. KIT CONTENTS

Manual

One sealed microtiter plate (12 strips, 8 wells each), coated with antibodies directed against mouse-IgG. Ready-to-use.

Position of the reagents in the kit. For preparation of the reagents see Chapter 9.



1. **Dilution buffer** (20 ml, 4x concentrated)
2. **Rinsing buffer** (30 ml, 20x concentrated)
3. **Substrate solution** (12 ml, Ready-to-use)
4. **Stop solution** (15 ml, Ready-to-use)
5. **Conjugate** (lyophilised , blue cap)
6. **Antibody** (lyophilised, yellow cap)
7. not in use
8. not in use
9. not in use
10. not in use
11. not in use
12. not in use
13. **Zero Standard solution** (2 ml, Ready-to-use)
14. **Standard solution 1** (1ml, Ready-to-use) **0.0157 ng/ml**
15. **Standard solution 2** (1ml, Ready-to-use) **0.0313 ng/ml**
16. **Standard solution 3** (1ml, Ready-to-use) **0.0625 ng/ml**
17. **Standard solution 4** (1ml, Ready-to-use) **0.125 ng/ml**
18. **Standard solution 5** (1ml, Ready-to-use) **0.25 ng/ml**
19. **Standard solution 6** (1ml, Ready-to-use) **0.5 ng/ml**

6. EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- Scales and weighing vessels
- Gloves
- Glass tubes
- Fume hood
- Homogeniser (vortex, mixer)
- Centrifuge (2000 x g)
- Automated microtiter plate washer or 8-channel micropipette 100 – 300 µl
- Microtiter plate shaker
- Microtiter plate reader with 450 nm filter
- Test tubes (10 – 15 ml)
- Micropipettes, 25 – 1000 µl
- Multipipette with 2.5 ml combitips
- Aluminium foil or parafilm
- Filter
- Methanol
- n-Hexane
- Dichloromethane
- Distilled water

7. PRECAUTIONS

- This kit may contain hazardous substances. For hazard notes please refer to the appropriate safety data sheets (SDS).
- Avoid contact of all biological materials with skin and mucous membranes.
- Do not pipette by mouth.
- Do not eat, drink, smoke, store or prepare foods, or apply cosmetics within the designated work area.
- Do not use components past expiration date and do not use components from different lots.
- Each well is ultimately used as an optical cuvette. Therefore, do not touch the under surface of the wells, prevent damage and dirt.
- All components should be completely dissolved before use. Take special attention to the substrate and rinsing buffer, which crystallize at +4°C.
- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting and washing throughout this procedure are necessary to maintain good precision and accuracy.

R-Biopharm Nederland makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm Nederland shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

8. SAMPLE PREPARATION

8.1 Feed and Food

- The first step in the sample preparation process depends on the nature of the compound to be tested. In general a homogenous sample has to be obtained from a representative part of the compound. Alternative techniques such as grinding, pulverizing, etc can be used.
- Approximately 50-100 gram of sample is ground and pulverised into a fine homogenous compound
- To 3 gram of ground sample add 9 ml of 80% methanol
- Shake thoroughly for 10 minutes at 20°C - 25°C
- Centrifuge the sample, 10 minutes, 2000 x g or filtrate (recommended for larger sample volumes).
- An aliquot of 50 µl of the supernatant obtained after centrifugation, or an aliquot of 50 µl of filtrate obtained after filtration, is diluted with 150 µl of dilution buffer to obtain a solution containing 20% methanol.
- Pipette 50 µl of the diluted supernatant in the respective wells of the ELISA plate.

8.2 Alternative method feed and food

- The first step in the sample preparation process depends on the nature of the compound to be tested. In general a homogenous sample has to be obtained from a representative part of the compound. Alternative techniques such as grinding, pulverizing, etc can be used.
- Approximately 50-100 gram of sample is ground and pulverised into a fine homogenous compound
- Add 10 ml distilled water to 1 gram of homogenised ground sample and vortex
- Add 10 ml dichloromethane (CH₂Cl₂) to the mixture
- Shake thoroughly, head over head for 5-10 minutes (Rotor)
- Centrifuge the sample mixture 10 minutes, 2000 x g, 20°C - 25°C.
- Filtrate the layer underneath over a filter
- Evaporate 1 ml of the filtrate under a mild stream of nitrogen at 50°C.
- Dissolve the residue in 1 ml sample dilution buffer (see chapter 9)
- Add 1 ml of n-hexane
- Vortex for 1 minute
- Centrifuge 10 minutes, 2000 x g, 20°C - 25°C.
- Remove the upper layer (n-hexane)
- Pipette 50 µl of the layer underneath in the respective wells of the ELISA plate.

Some matrices need a second defatting step. The liquid appears to be turbid after the first defatting step. In that case repeat the defatting procedure.

Repeat the defatting step till the layer underneath is clear. Use 50 µl of the layer underneath in the respective wells of the ELISA plate.

8.3.1 Procedure to extract fatty samples, like oil, seeds, butter etc.

- Weigh 2 g of finely cut and subsequently homogenized butter in the glass tube.
- Add 8 ml of dichloromethane.
- Mix 30 minutes head over head.
- Filter the butter extract through
- Take 4 ml of the filtrate and evaporate to dryness at 50°C, under a mild stream of nitrogen.
- Dissolve the residue in 1 ml sample dilution buffer (see chapter 9)
- Add 1 ml of n-hexane (for defatting) and vortex.
- Centrifuge 5 minutes, 2000 x g, 20°C - 25°C
- Pipette the layer underneath into a clean glass tube.
- Add 1 ml n-hexane and vortex.
- Centrifuge 5 minutes, 2000 x g, 20°C - 25°C
- Remove the upper layer (n-hexane)
- Pipette 50 µl of the layer underneath in the respective wells of the ELISA plate.

8.3.2 Alternative method: Edible oils

- Add to 1 ml of 1 g oil, 3 ml 80% methanol and 2 ml n-hexane
- Vortex and shake head over head for 10 minutes
- centrifuge 5 minutes, 2000 x g
- Remove the upper layer (n-hexane)
- Dilute 100 µl extract with 400 µl sample dilution buffer (see chapter 9), vortex
- Pipette 50 µl of the dilution in the respective wells of the ELISA plate.

9. PREPARATION OF REAGENTS

Before starting the test, allow the reagents to come to ambient temperature. Any reagents not used should be put back into storage immediately at 2°C - 8°C. Prepare reagents fresh before use.

Microtiter plate

Return unused strips into the zip lock bag with desiccant and store at 2°C - 8°C for use in subsequent assays. Retain also the strip holder.

Rinsing buffer

The rinsing buffer is delivered 20 times concentrated. Prepare dilutions freshly before use. For each strip 20 ml of diluted rinsing buffer is used (1 ml concentrated rinsing buffer + 19 ml distilled water).

Substrate/chromogen solution

The substrate/chromogen solution (ready to use) tends to precipitate at 4°C. Take care that this vial is at room temperature (keep in the dark) and mix the content before pipetting into the wells.

Dilution buffer

This buffer is for dissolving conjugate, antibody and to prepare the sample dilution buffer.

The dilution buffer is 4 times concentrated. Before dilution (20 ml buffer + 60 ml distilled water) the concentrated buffer should be at room temperature (20°C - 25°C) and thoroughly mixed. Concentrated buffer can show precipitates of the contents. Mix well before dilution with distilled water.

Sample dilution buffer

Prepare freshly before use. Prepare the sample dilution buffer by adding 2 ml of 100% methanol to 8 ml of dilution buffer.

Conjugate solution

Reconstitute the vial of lyophilised conjugate (Aflatoxin-HRP) with 4 ml of dilution buffer, mix thoroughly and keep in the dark until use.

Antibody solution

Reconstitute the vial of lyophilised antibodies with 4 ml of dilution buffer, mix thoroughly and keep in the dark until use.

10. ASSAY PROCEDURE

Rinsing protocol

In ELISA's, between each immunological incubation step, unbound components have to be removed efficiently. This is reached by appropriate rinsing. It should be clear that each rinsing procedure must be carried out with care to guarantee good inter- and intra-assay results.

Basically, manual rinsing or rinsing with automatic plate wash equipment can be performed as follows:

Manual rinsing

1. Empty the contents of each well by turning the microtiter plate upside down and remove residual liquid by striking the plate against a paper towel.
2. Fill all the wells to the rims (300 μ l) with rinsing solution.
3. This rinsing cycle (1 and 2) should be carried out 3 times.
4. Turn the plate upside down and empty the wells by a firm short vertical movement.
5. Place the inverted plate on absorbent paper towels and tap the plate firmly to remove residual washing solution in the wells.
6. Take care that none of the wells dry out before the next reagent is dispensed.

Rinsing with automatic microtiter plate wash equipment

When using automatic plate wash equipment, check that all wells can be aspirated completely, that the rinsing solution is nicely dispensed reaching the rim of each well during each rinsing cycle. The washer should be programmed to execute three rinsing cycles.

Assay Protocol

1. Prepare samples according to Chapter 8 (Sample treatment) and prepare reagents according to Chapter 9 (Preparation of reagents).
2. Pipette 100 μ l of the zero standard in duplicate (wells H1, H2, blank).
Pipette 50 μ l of the zero standard in duplicate (wells A1, A2; Bmax).
Pipette 50 μ l of each of the Aflatoxin B1 standard solutions in duplicate (wells B1,2 to I1,2 i.e. 0.0157, 0.0313, 0.0625, 0.125, 0.25 and 0.5 ng/ml).
3. Pipette 50 μ l of each sample solution in duplicate into the remaining wells of the microtiter plate (40 samples; 80 wells).
4. Add 25 μ l of conjugate (Aflatoxin-HRP) to all wells, except wells H1 and H2.
5. Add 25 μ l of antibody solution to all wells, except wells H1 and H2.
6. Seal the microtiter plate and shake the plate for a few seconds.

7. Incubate for 1 hour in the dark at 37°C.
8. Discard the solution from the microtiter plate and wash 3 times with rinsing buffer.
9. Pipette 100 µl of substrate solution into each well.
10. Incubate 30 minutes at 20°C – 25°C.
11. Pipette 100 µl of stop solution to each well.
12. Read the absorbance values immediately at 450 nm.

11. INTERPRETATION OF RESULTS

Subtract the mean optical density (O.D.) of the wells H1 and H2 (Blank) from the individual O.D. of the wells containing the standards and the samples.

The O.D. values of the six standards and the samples (mean values of the duplicates) are divided by the mean O.D. value of the zero standard (Bmax, wells A1 and A2) and multiplied by 100. The zero standard (Bmax) is thus made equal to 100% (maximal absorbance) and the other O.D. values are quoted in percentages of the maximal absorbance.

$$\frac{\text{O.D. standard (or sample)}}{\text{O.D. zero standard (Bmax)}} \times 100 = \text{Percentage maximal absorbance}$$

Calibration curve:

The values (% maximal absorbance) calculated for the standards are plotted on the Y-axis versus the analyte equivalent concentration (ng/ml) on a logarithmic X-axis.

Alternative for calibration curve:

The absorption value of the standards is plotted on the Y-axis versus the concentration on the X-axis. The Y-axis is in logit the Y-axis is logarithmic.




Figure 1 : Example of a calibration curve

The amount of Aflatoxin in the samples is expressed as Aflatoxin B1 equivalents. The Aflatoxin B1 equivalents in the extracts (ng/ml) corresponding to the % maximal absorbance of each extract can be read from the calibration curve.

8.1 Food and feed samples

The aflatoxin B1 equivalents, as read from the standard curve, have to be multiplied by a factor 16 to obtain the aflatoxin B1 content in food and feed samples.

8.2 Alternative method food and feed

The aflatoxin B1 equivalents, as read from the standard curve, have to be multiplied by a factor 10 to obtain the aflatoxin B1 content in food and feed samples.

8.3.1 Procedure to extract fatty samples like, seeds and oil

The concentration of Aflatoxin B1 can directly be read from the calibration curve.

8.3.2 Alternative method edible oils

The aflatoxin B1 equivalents, as read from the standard curve, have to be multiplied by a factor 20 to obtain the aflatoxin B1 content in edible oils.

12. LITERATURE

J.E. Smith, C.W. Lewis and J.G. Anderson. Mycotoxins in human nutrition and health. EU Directorate-General XII, Science, Research and Development, 1990.

J.L. Richard, G.A. Bennett, P.F. Ross and P.E. Nelson. Analysis of naturally occurring mycotoxins in feedstuffs and food. J. Anim. Sci. 71, 2563-2574, 1993.

Commission Regulation (EC) No 1881/2006 of 19 December 2006 setting maximum levels for certain contaminants in foodstuffs. Off. J. European Union. L364 (2006) 5-24.

Commission Regulation (EU) No 165/2010 of 26 February 2010 amending Regulation (EC) No 1881/2006 setting maximum levels for certain contaminants in foodstuffs as regards aflatoxins. Off. J. European Union. L50 (2010) 8-12.

13. ORDERING INFORMATION

For ordering the Aflatoxin B1 ELISA kit, please use cat. code 5121AFB.

14. REVISION HISTORY

The manual is adapted to a new layout of the test kit. Several textual changes are added.