bioavid Lateral Flow Tests inkl. Hook-Linie BLH-Serie

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

bioavid Lateral Flow tests incl. Hook Line BLH-Line

INTENDED USE

The bioavid lateral flow tests with included hook line are immunochromatographic tests for the sensitive and qualitative detection of allergen residues on surfaces (e.g. swab test for the hygiene control in food production lines), in cleansing and process water (CIP water) and food samples. The integrated hook line prevents false negative interpretation of high positive samples.

Time requirement:

Sample preparation (swab, CIP water) approx. 1 min
 Sample preparation (liquid samples)..... approx. 3 min
 Sample preparation (solid samples)..... approx. 5 min
 Test implementation (incubation time) 10 min

Art. No.	Parameter	Limit of Detection (LoD)	
		CIP water (ppm)	Surfaces ($\mu\text{g}/\text{cm}^2$)
BLH700-15	Coconut	1	≤ 0.03
BLH701-15	Almond	1	≤ 0.02
BLH703-15	Mustard	1	≤ 0.08
BLH704-15	Hazelnut	1	≤ 0.02
BLH705-15	Macadamia	0.5	0.1
BLH706-15	Peanut	1	≤ 0.02
BLH707-15	Walnut	1	0.1
BLH708-15	Egg	1	0.015 - 0.08
BLH709-15	Sesame	1	≤ 0.01
BLH710-15	Cashew	1	≤ 0.08
BLH711-15	Pistachio	0.1	0.05
BLH714-15	Casein	1	0.015
BLH716-15	Crustacean	10 - 100	0.15 - 1.5

The presence of cleaners and sanitizers in the swab or CIP water sample, as well the type of food matrix being tested, can affect limit of detection. Detailed information can be found in the validation report.

BACKGROUND INFORMATION

In the case of a food allergy, certain parts of a food (allergens) trigger an immune reaction in humans. This immune reaction is the cause of a variety of different symptoms, which can range from mild tingling in the mouth, to rash, to a life-threatening anaphylactic shock. Worldwide, about 2 - 8 % of people suffer from a food allergy. Since there is currently no effective therapy for food allergies, sensitized persons must follow an allergen-free diet. Therefore, labelling allergens on food is essential. According to Regulation (EU) No. 1169/2011, the 14 most important allergens and their products must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand. For prevention of unintended contamination, food industry must consider a respective allergen management (e.g. testing of raw materials, checking production lines). Lateral flows offer a fast, easy and reliable method for the detection of allergen residues. Summarized information about the usage of LFDs within an allergen management can be found on our infographic, which can be downloaded on www.r-biopharm.com.

HOOK EFFECT

Most important challenges with the safe and fast LFD analysis are false negative results due to very high allergen concentrations in the test sample. The so-called “high dose hook effect” (or “overload effect”; Standard EN 15633-1. 2019) is observed when a very high amount of an analyte is present in the sample. In this case, the development of the test line can be inhibited and there is a risk that the test result will be falsely evaluated as negative (or analyte undetectable). The hook line on the test strip makes this effect visible. A missing hook line indicates a high allergen content in the sample and can no longer be evaluated as negative (see chapter 8).

1. TEST PRINCIPLE

The principle of the test is an antigen-antibody reaction. The reaction vial contains allergen-specific antibodies. After addition of running buffer, the sample (e.g. supernatant, swab sample, CIP water) is added to the reaction vial. If the sample contains the target antigen, an antigen-antibody complex will form. After incubation, a dip stick is placed into the reaction vial. Capillary forces cause the liquid to flow up the dip stick. The

dip stick contains three areas, where reactions take place and lines for visual evaluation can form. First, the control line in order to verify the test run, which should appear every time. Second, the test line, where antigen-specific antibodies are immobilized. If target antigen is present, the previously formed antigen-antibody complexes react with the immobilized antibodies on the test line resulting in the formation of an antibody-antigen-antibody complex. Third, the hook line, where target analyte is immobilized and will appear, as long as target-analyte-specific antibodies are present unbound. The result is read visually. The test is very rapid, reliable and does not require laboratory equipment for most applications.

2. STORAGE

Store the test kit at 2 - 25 °C (35 - 77 °F).

3. REAGENTS PROVIDED

Each kit contains sufficient materials for 15 measurements.

Content	Description	Amount
Reaction vials	Containing labeled antibodies in stabilized dried form, needs to be dissolved with running buffer	15 x
Allergen test strips	1 tube, reclosable	15 x
Dropper bottle with running buffer	8 mL	1 x
Positive control	Needs to be dissolved with 1 mL distilled or deionized water	1 x
Plastic pipettes	With 100 µL mark	15 x
Plastic swab	Used to swab surfaces	16 x
Pre-filled PBS vials	Containing 1.2 mL PBS; used for swabbing (moisten the swab, swab and flush swab)	15 x
Evaluation card		1 x

4. WARNINGS AND PRECAUTIONS FOR THE USERS

Bring the test kit to room temperature (20 - 25 °C / 68 - 77 °F) before use. Do not open the tube with test strips when it is colder than the environment, since strip deterioration may result from moisture.

Once removed, test strips should not be returned to the test strip tube.

The test strips are sensitive to moisture. Therefore, remove the test strips from the test strip tube immediately before use only.

The product / test is only suitable within the scope of its intended use.

Do not use the test kit after the expiration date (see test kit label).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).


5. PREPARATION OF SAMPLES

Positive control (if needed):

Add 1 mL distilled or deionized water to the positive control and leave for 5 minutes. Then apply 0.1 mL as a sample into the assay. The positive control can be stored at $-20\text{ }^{\circ}\text{C}$ ($-4\text{ }^{\circ}\text{F}$) for at least one year. Do not repeatedly freeze and thaw. Dissolved or thawed controls are stable for 2 days at room temperature ($20 - 25\text{ }^{\circ}\text{C}$ / $68 - 77\text{ }^{\circ}\text{F}$).

5.1. Solid food samples:

1. Grind/mix a representative amount of sample (e.g. 5 - 50 g) thoroughly and dissolve the sample 1:10 (1+9) in distilled or deionized water (e.g. 1 g sample + 9 mL distilled or deionized water). Mix until a uniform suspension is obtained.
2. Centrifuge 2 min at $2,000 \times g$ or filter. The supernatant or extract should be as particle-free as possible.

5.2. Liquid food samples:


1. Mix a representative amount of sample (e.g. 5 - 50 mL) thoroughly and dissolve the sample 1:10 (1+9) in distilled or deionized water (e.g. 1 mL sample + 9 mL distilled or deionized water). Mix until a uniform suspension is obtained.
2. Particle-containing extracts should be centrifuge for 2 min at $2,000 \times g$ or filtered. The supernatant or extract should be free of suspended solids.

5.3. CIP water:

1. 0.1 mL of cleansing or process water (CIP water) can directly be used in the assay.

5.4. Surfaces (swab samples):






1. Use one vial (pre-filled with PBS) per sample. Moisten a clean cotton swab with the contained PBS.
2. Swab the area of interest with a crosshatch technique in all directions (see graphic on the right).
3. Release the sample from the swab into the vial from step 1: squeeze the swab against the vial, carefully rotate and move up and down.

6. ASSAY


Note the sample ID on a reaction vial and open it.

1. Add 0.1 mL (3 drops) of running buffer into the reaction vial.
2. Add 0.1 mL of sample (supernatant, CIP water or swab sample) into the reaction vial (single-use plastic pipettes with 0.1 mL scale are contained in the kit) and mix the content by carefully pipetting up and down (approx. 3-times) or close the reaction vial and mix intensively by carefully tapping the vial on the desk a few times. The liquid must turn red.
3. Incubate for **5 minutes**.
4. Immediately before use, take an allergen test strip out of the test strip tube. Immediately close the test strip tube again. Insert the test strip with the arrows pointing down into the reaction vial.
5. Read the result at **5 minutes**.

7. INTERPRETATION OF THE RESULT

The allergen test strip includes three lines in the reaction field: the control line (C, top line), the hook line (H, middle line) and the test line (T, lowest line). In total, 1 to 3 purple lines may appear in the reaction field on the test strip. The control line (C) indicates that the test has been performed correctly and consequently must be present on every valid test strip. The test line (T) indicates the presence of the target residue in the sample. Very high concentrations of allergen residues may reduce the intensity of the test line or suppress its formation completely. The hook line (H) serves as an “alarm line”. If it is missing, a high content of allergen is present in the sample. In this case, a missing test line must not be interpreted as a negative result.

8. LIMITATION

Due to the large amount of different foods, it is not possible for the manufacturer of a test kit to validate all foods for the test. The suitability of the respective food for use in the test must be verified by the user himself.

Some solid samples containing cereal flour and/or thickeners may interfere with a correct development of the hook line if applied in a 1:10 dilution. Thus, very viscous samples should be tested in a 1:20 dilution. In this case, lower sensitivity should be noted due to higher sample dilution.

For testing samples containing polyphenols (e.g. coffee, cocoa, spices, wine) or samples containing high amounts of oil the use of the Absorbent Buffer from bioavid (Art. No. BS810-15) is recommended.

Pure nuts interfere with the development of the hook line and must be diluted 1:20 (instead of 1:10) with distilled or deionized water. In this case, lower sensitivity should be noted due to higher sample dilution.

Bioavid lateral flow tests are qualitative tests. They do not allow quantitative evaluation even in the case of variably developed test lines (T).

9. NOTES

Product-specific information can be found in the validation report.

For verification of unknown food samples, an SOP (Standard Operating Procedure) is available upon request.

Version overview

Version number	Chapter and title
2023-02-01	Release version
2023-07-18	Revision
2023-10-20	Revision
2024-02-27	Current version Changes made: <ul style="list-style-type: none">– Added testsystems BLH705-15 Macadamia and BLH711-15 Pistachio to Intended Use– Adaption in chapter “Limitations” (Addition of suitability of the Absorbent buffer (Art. No. BS810-15) for oily samples)

Disclaimer


The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

1 Sample preparation

2 Test procedure