# bioavid Lateral Flow Tests inkl. Hook-Linie BLH-Serie

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# bioavid Lateral Flow tests incl. Hook Line BLH-Line

## INTENDED USE

The bioavid lateral flow tests with included hook line are immunochromatographic tests for the sensitive and qualitative detection of allergen residues on surfaces (e.g. swab test for the hygiene control in food production lines), in cleansing and process water (CIP water) and food samples. The integrated hook line prevents false negative interpretation of high positive samples.

### Time requirement:

Sample preparation swab test .....	approx. 1 min
Sample preparation CIP water.....	approx. 1 min
Sample preparation (liquid samples) .....	approx. 3 min
Sample preparation (solid samples) .....	approx. 5 min
Test implementation (incubation time).....	10 min

CIP = Clean-In-Place (rinse water from cleaning process)

## LIMIT OF DETECTION:

Art. No.	Parameter	Limit of Detection (LoD)	
		CIP water ( $\mu\text{g/mL}^*$ )	Surfaces ( $\mu\text{g/cm}^2$ / $\mu\text{g/mL}^{**}$ )
BLH700-15	Coconut	0.1	0.1 / 10
BLH701-15	Almond	0.1	0.05 / 5
BLH702-15	Brazil nut	0.05	0.025 / 25
BLH703-15	Mustard	1.0	0.5 / 50
BLH704-15	Hazelnut	0.1	0.05 / 5
BLH705-15	Macadamia	0.5	0.1 / 10
BLH706-15	Peanut	0.5	0.02 / 2
BLH707-15	Walnut	2.0	0.4 / 40
BLH708-15	Egg	0.05	0.05 / 5
BLH709-15	Sesame	0.2	0.02 / 2
BLH710-15	Cashew	0.1	0.08
BLH711-15	Pistachio	0.025	0.05 / 5
BLH712-15	Soy	0.1	0.05 / 5
BLH714-15	Casein	0.01	0.01 / 1
BLH716-15	Crustacean	2.0	0.4 / 40

The presence of cleaners and sanitizers in the swab or CIP water sample, as well the type of food matrix being tested, can affect limit of detection. Detailed information can be found in the validation report.

\* $\mu\text{g/mL}$  = ppm = mg/kg

\*\*Converted according to spiked surface area of 100  $\text{cm}^2$

## BACKGROUND INFORMATION

In the case of a food allergy, certain parts of a food (allergens) trigger an immune reaction in humans. This immune reaction is the cause of a variety of different symptoms, which can range from mild tingling in the mouth, to rash, to a life-threatening anaphylactic shock. Worldwide, about 2 - 8 % of people suffer from a food allergy. Since there is currently no effective therapy for food allergies, sensitized persons must follow an allergen-free diet. Therefore, labelling allergens on food is essential. According to Regulation (EU) No. 1169/2011, the 14 most important allergens and their products must be declared on food labels. Similar regulations exist e.g. in the USA, Canada, Australia and New Zealand. For prevention of unintended contamination, food industry must consider a respective allergen management (e.g. testing of raw materials, checking production lines). Lateral flows offer a fast, easy and reliable method for the detection of allergen residues. Summarized information about the usage of LFDs within an allergen management can be found on our infographic, which can be downloaded on [www.r-biopharm.com](http://www.r-biopharm.com).

## HOOK EFFECT

One of the most important challenges with the safe and fast LFD analysis are false negative results due to very high allergen concentrations in the test sample. The so-called “high dose hook effect” (or “overload effect”; Standard EN 15633-1. 2019) is observed when a very high amount of an analyte is present in the sample. In this case, the development of the test line can be inhibited and there is a risk that the test result will be falsely evaluated as negative (or analyte undetectable). The hook line on the test strip makes this effect visible. A missing hook line indicates a high allergen content in the sample and can no longer be evaluated as negative (see chapter 7).

## 1. TEST PRINCIPLE

The principle of the test is an antigen-antibody reaction. The reaction vial contains allergen-specific antibodies. After addition of running buffer, the sample (e.g. supernatant, swab sample, CIP water) is added to the reaction vial. If the sample contains the target antigen, an antigen-antibody complex will form. After incubation, a dip stick is placed into the reaction vial. Capillary forces cause the liquid to flow up the dip stick. The dip stick contains three areas, where reactions take place and lines for visual evaluation can form. First, the control line in order to verify the test run, which should appear every time.

Second, the test line, where antigen-specific antibodies are immobilized. If target antigen is present, the previously formed antigen-antibody complexes react with the immobilized antibodies on the test line resulting in the formation of an antibody-antigen-antibody complex. Third, the hook line, where target analyte is immobilized and will appear, as long as target-analyte-specific antibodies are present unbound. The result is read visually. The test is very rapid, reliable and does not require laboratory equipment for most applications.

## 2. STORAGE

Store the test kit at 2 - 25 °C (35 - 77 °F).

## 3. REAGENTS PROVIDED

Each kit contains sufficient materials for 15 determinations. Each test kit contains:

Component	Format	Content
<b>Test strip</b>	Ready to use in a container, resealable	15 pieces
<b>Reaction vials</b>	Contains labelled antibodies in stabilized dried form; needs to be dissolved with running buffer	15 pieces
<b>Dropper bottle with running buffer</b>		8 mL
<b>Positive control</b>	Needs to be dissolved with 1 mL distilled or deionized water	1 pieces
<b>Single-use plastic pipettes</b>	With 100 µL mark	15 pieces
<b>Plastic swabs</b>	For sampling surfaces	16 pieces
<b>Pre-filled PBS vials</b>	Contains 1 mL PBS; required for swab test	15 pieces
<b>Evaluation card</b>		1 piece

## 4. WARNINGS AND PRECAUTIONS FOR THE USERS

The product / test is only suitable within the scope of its intended use.

If stored in a refrigerator, bring the test kit to room temperature (20 - 25 °C / 68 - 77 °F) before use. Do not open the tube with test strips when it is colder than the environment, since strip deterioration may result from moisture.

Once removed, test strips should not be returned to the test strip tube.

The test strips are sensitive to moisture and must be absolutely protected against it. Moist test strips can have a negative impact on the result. For this reason, remove the test strips from packaging only after having reached room temperature and immediately prior to use in the test.

Do not use the test kit after the expiration date (see test kit label).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 5. PREPARATION OF SAMPLES


Airborne allergens and dirty laboratory equipment lead to contamination of the assay. Therefore, clean surfaces and equipment (e.g. mincers, vessels) thoroughly before and after each sample.

### **Positive control (if needed):**

Add 1 mL distilled or deionized water to the positive control and leave for 5 minutes. Mix well. Then apply 0.1 mL as a sample into the assay. The positive control can be stored at -20 °C (-4 °F) for at least one year. Do not repeatedly freeze and thaw. Dissolved or thawed controls are stable for 2 days at room temperature (20 - 25 °C / 68 - 77 °F).

### 5.1. Solid food samples:

1. Grind/mix a representative amount of sample (e.g. 5 - 50 g) thoroughly and dissolve the sample 1:10 (1+9) in distilled or deionized water (e.g. 1 g sample + 9 mL distilled or deionized water). Mix until a uniform suspension is obtained.
2. Centrifuge 2 min at 2,000 x g or filter. The supernatant or extract should be as particle-free as possible.

### 5.2. Liquid food samples:


1. Mix a representative amount of sample (e.g. 5 - 50 mL) thoroughly and dissolve the sample 1:10 (1+9) in distilled or deionized water (e.g. 1 mL sample + 9 mL distilled or deionized water). Mix until a uniform suspension is obtained.
2. Particle-containing extracts should be centrifuged for 2 min at 2,000 x g or filtered. The supernatant or extract should be free of suspended solids.

### 5.3. CIP water:

1. 0.1 mL of cleansing or process water (CIP water) can directly be used in the assay.

## 5.4. Surfaces (swab samples):

1. Use one vial (pre-filled with PBS) per sample. Moisten a clean cotton swab with the contained PBS.
2. Swab the area of interest with a crosshatch technique in all directions (see graphic on the right). An area of about 10 x 10 cm<sup>2</sup> is recommended, but it is also possible to take samples from trouble spots (e.g. equipment).
3. Release the sample from the swab into the vial from step 1: squeeze the swab against the vial, carefully rotate and move up and down.






## 6. ASSAY

Note the sample ID on a reaction vial and open it.


1. Add 0.1 mL (3 drops) of running buffer into the reaction vial.
2. Add 0.1 mL of sample (supernatant, CIP water or swab sample) into the reaction vial (single-use plastic pipettes with 0.1 mL scale are contained in the kit) and mix the content by carefully pipetting up and down (approx. 3-times) or close the reaction vial and mix intensively by carefully tapping the vial on the desk a few times. Due to the coloration of the antibodies in the reaction vial, the liquid should turn red when mixing only with the running buffer. However, after adding the sample, depending on the coloration of the sample, the red color may not be visible.
3. Incubate for **5 minutes**.

4. Immediately before use, take an allergen test strip out of the test strip tube. Immediately close the test strip tube again. Insert the test strip with the arrows pointing down into the reaction vial.
5. Read the result at **5 minutes**.

## 7. INTERPRETATION OF THE RESULT

The allergen test strip includes three lines in the reaction field: the control line (C, top line), the hook line (H, middle line) and the test line (T, lowest line). In total, 1 to 3 purple lines may appear in the reaction field on the test strip. The control line (C) indicates that the test has been performed correctly and consequently must be present on every valid test strip. The test line (T) indicates the presence of the target residue in the sample. Very high concentrations of allergen residues may reduce the intensity of the test line or suppress its formation completely. The hook line (H) serves as an “alarm line”. If it is missing, a high content of allergen is present in the sample. In this case, a missing test line must not be interpreted as a negative result.

## 8. LIMITATION

Test results may vary depending on the actual test procedure and the laboratory environment.

The limit of detection depends on sample type and extraction efficiency or the properties of the swabbed surface and the kind of contamination, respectively.

Due to the large amount of different foods, it is not possible for the manufacturer of a test kit to validate all foods for the test. The suitability of the respective food for use in the test must be verified by the user himself.

Some solid samples containing cereal flour and/or thickeners may impact the running behavior if applied in a 1:10 dilution. Thus, very viscous samples should be tested in a 1:20 dilution. In this case, lower sensitivity should be noted due to higher sample dilution.

For testing samples containing polyphenols (e.g. coffee, cocoa, spices, wine) or samples containing high amounts of oil the use of the Absorbent Buffer from bioavid (Art. No. BS810-15) is recommended.

Strongly acidic or basic samples can impair the test system and should be neutralized before extraction.

Pure nuts interfere with the development of the hook line and must be diluted 1:20 (instead of 1:10) with distilled or deionized water. In this case, lower sensitivity should be noted due to higher sample dilution.

Bioavid lateral flow tests are qualitative tests. They do not allow quantitation of the test results based on test or hook line intensity.

## 9. NOTES

For detailed results and further information please refer to the current validation report. In addition, data may be available from comparative laboratory tests and interlaboratory comparisons.

For verification of unknown food samples, an SOP (Standard Operating Procedure) is available upon request.

## 10. RECOMMENDATION

In order to ensure a high analytical performance we recommend:

- To comply with the general quality assurance requirements for laboratories as listed in the standards EN 15633-1 and 15842.
- To adjust the sample's pH value to neutral (pH 6.5 to 7.5) prior to extraction in case of extremely acidic or basic samples.
- To verify and adjust the sample's pH value to neutral after extraction.
- To do spike experiments to ensure an accurate and correct test procedure.
- To use a RIDASCREEN® ELISA for quantification of a positive result.


- To clear unexpected results, use the ELISA RIDASCREEN® or run a PCR (e.g. SureFood®).
- For quick and easy documentation of the test strips, use the “RIDA®SMART APP Allergen” documentation app available free of charge in the Google Play Store (only for Android devices).


## Version overview

Version number	Chapter and title
2023-02-01	Release version
2023-07-18	Revision
2023-10-20	Revision
2024-02-27	Revision
2024-07-30	Revision Changes made: <ul style="list-style-type: none"> <li>– Addition of the new test system bioavid Lateral Flow Soy incl. Hook Line (Art. No. BLH712-15)</li> <li>– Adaptation of the specifications based on new data</li> <li>– Chapter 3: Reformatting table</li> <li>– Chapter 5: Added recommended sampling for the swab test</li> <li>– Chapter 6: Note of red coloration of the liquid depending on sample color after adding the sample to the reaction vial</li> <li>– Additions to chapter 8 and 9</li> <li>– Chapter 10 „Recommendation“ added</li> <li>– Revision of the disclaimer</li> <li>– Insert of bioavid logo</li> </ul>
2024-10-14	Current version Changes made: <ul style="list-style-type: none"> <li>– Addition of the new test system bioavid Lateral Flow Brazil nut incl. Hook Line (Art. No. BLH702-15)</li> <li>– Correction of version overview; version date of the previous version (from 2024-07-22 to 2024-07-30)</li> </ul>

## 1 Sample preparation

## 2 Test procedure

## Disclaimer

1. In conformance with the German Civil Code (“BGB”) R-Biopharm AG provides a limited warranty (“Gewährleistung”) against defects in design and manufacture present at delivery that make the Product unsuitable or unsafe for the contractually specified Product use and location through properly qualified personnel. Such limited warranty warrants Product title and against such material Product defects for 12 months, or if the Product is provided with a shelf life, the length of the declared shelf life, calculated from the date of risk transfer pursuant to delivery terms, and further provided customer provides timely and proper notice of the defect.  
ALL OTHER WARRANTIES OR GUARANTIES, EXPRESS OR IMPLIED, OF ANY KIND ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, THE COURSE OF DEALING BETWEEN THE PARTIES OR OTHERWISE. The responsibility for any express or implied warranties or extensions of R-Biopharm’s own warranty made by dealers, distributors, or resellers is solely with such dealer, distributor, or reseller, and R-Biopharm AG expressly disclaims liability for such warranties.
2. R-Biopharm Products are designed for sole use by properly trained and qualified personnel applying appropriate industry standards and practices. Defects are covered only to the extent that such defects are demonstrably the result of defective material, design, parts, or faulty workmanship, which affect safety or result in substandard performance below specifications. R-Biopharm AG is not liable for any improper use nor the consequences of
  - a. the failure to read, understand and follow use and safety instructions, or use proper preparation and storage;
  - b. the failure to utilize trained and unqualified personnel, suitable samples or sampling techniques;
  - c. chemical, electromagnetic, mechanical, or electrolytic influences outside R-Biopharm AG provided standard perimeters through product specific data sheets, written product descriptions, or as otherwise expressly agreed upon in writing; or
  - d. any combination thereof.
3. R-Biopharm AG is also not liable for any changes or modifications, not carried out by R-Biopharm AG, nor consequences of use, applications, or processing which are unsafe or otherwise inconsistent with intended Product purpose as limited by written Product descriptions and specifications of R-Biopharm AG.
4. R-Biopharm AG’s liability for ordinary breach of contract is limited to repair, replacement, other substitute performance, or refund. The choice of remedies is within R-Biopharm AG’s sole discretion. R-Biopharm AG is not contractually liable for any incidental, or consequential damages, including but not limited to purchaser’s expenses, losses, or damages from loss of good will, frustrated business purposes, sales expectations, investment loss, actual or anticipated lost profits, End User indemnity, or any other business expenditures.
5. The foregoing limited warranty is solely intended to fulfill the warranty requirements (“Gewährleistung”) implied by German Civil Code. It is not intended to extend the scope or period of such Gewährleistung or provide additional warranties. Nor is this Disclaimer otherwise intended to limit or extent mandatory liability for damages in tort resulting from injury to life, body, or health, for intentional or grossly negligent acts, fraud, or due to strict liability under applicable product liability laws. A reversal of the burden of proof or rules of contract interpretation is not intended.