bioavid Lateral Flow Tests incl. Hook-Line

REF **BLH7xx-15**

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Bioavid Lateral Flow Tests incl. Hook Line

REF BLH7xx-15

Brief information and Intended Use

The bioavid lateral flow tests with included hook line (Art. No. BLH7xx-15) are immunochromatographic rapid tests for the qualitative detection of allergen residues on surfaces (e.g. in food production lines), in cleaning and process water (CIP water) and in food. The integrated hook line prevents false negative interpretation of high positive samples.

A test kit contains 15 test strips, for one analysis each.

Time requirement:

Swab sampling:	approx. 1 min
Sample preparation (CIP) water	direct use
Sample preparation solid food samples	approx. 5 min
Sample preparation liquid food samples	approx. 3 min
Test implementation (incubation times)	10 min

CIP = Clean-In-Place (rinse water from cleaning process)

Limit of detection

Results are reported as protein values. The presence of cleaning agents or disinfectants in swab or CIP water samples, as well the food matrix type being tested, affect limits of detection.

Art. No.	Parameter	Limit of detection (LoD)		
		(CIP) water (µg/ml)	Swab test (ng/cm ² * / µg/ml**)	Food samples*** (mg/kg protein)
BLH700-15	Coconut	0.01	1 / 0.1	0.1
BLH701-15	Almond	0.01	2 / 0.2	0.2
BLH702-15	Brazil nut	0.01	4 / 0.4	0.1
BLH703-15	Mustard	0.1	20 / 2	1
BLH704-15	Hazelnut	0.01	1 / 0.1	0.1
BLH705-15	Macadamia	0.005	2 / 0.2	0.1
BLH706-15	Peanut	0.05	5 / 0.5	1
BLH707-15	Walnut	0.25	6 / 0.6	3.6
BLH708-15	Egg	0.05	2,5 / 0.25	0.5
BLH709-15	Sesame	0.05	2 / 0.2	0.5
BLH710-15	Cashew	0.01	1 / 0.1	0.2
BLH711-15	Pistachio	0.1	5 / 0.5	1
BLH712-15	Soy	0.01	2,5 / 0.25	0.25
BLH714-15	Casein	0.2	20 / 2	2
BLH716-15	Crustacean	0.25	50 / 5	2.5

* Calculated from a spiked surface area of 100 cm²

** Refers to the theoretical concentration in the swab buffer (1 mL) after extraction of the swab.

*** Matrix-dependent

1. General information

In the case of a food allergy, certain parts of a food (allergens) trigger an immune reaction in humans. This immune reaction is the cause of a variety of different symptoms, which can range from mild tingling in the mouth, to rash, to a life-threatening anaphylactic shock. Worldwide, about 2 - 8 % of people suffer from a food allergy. Since there is currently no effective therapy for food allergies, sensitized persons must follow an allergen-free diet. Therefore, labelling allergens on food is essential. According to Regulation (EU) No. 1169/2011, the 14 most important allergens and their products must be declared on food labels. Similar regulations apply, for example, in the USA, Canada, Australia, and New Zealand. For prevention of unintended contamination, food industry must consider a respective allergen management (e.g. testing of raw materials, checking production lines). Lateral flows offer a fast, easy and reliable method for the detection of allergen residues. Brief information on the utility of

LFDs within an allergen management is available at [Rapid allergen testing – Food & Feed Analysis](#).

2. Hook effect

One of the most important challenges with the safe and fast LFD analysis are false negative results due to very high allergen concentrations in the test sample. The so-called “high dose hook effect” (or “overload effect”, Standard EN 15633-1. 2019) is observed when a very high amount of analyte is present in the sample. In this case, the development of the test line can be inhibited and there is a risk that the test result will be falsely evaluated as negative (or analyte undetectable). The hook line on the test strip makes this effect visible. A missing hook line indicates a high allergen content in the sample and can no longer be evaluated as negative (see chapter 11).

3. Test principle

Use the test according to the procedure described below. The explanation that follows describes how the test works.

The test is based upon an antigen-antibody reaction. The reaction vial contains allergen-specific antibodies. After addition of the running buffer, the sample (e.g. supernatant, swab sample, (CIP) water) is added to the reaction vial. If the sample contains the target antigen, an antigen-antibody complex forms. After incubation, a strip is placed into the reaction vial. Capillary forces cause the liquid to migrate up the dip stick. Reactions take place at three distinct lines, which allow visual evaluation of the test results.

In the direction of the liquid flow these are in order:

1. The Test Line: Antigen-specific antibodies are immobilized at this line. If the target antigen is present, the previously formed antigen-antibody-complexes are captured by the immobilized antibodies, resulting in a visible line that indicates a positive test result.
2. The Hook Line: Target antigen is immobilized, leading to a visible line as long as antigen-specific antibodies are present unbound. A positive test result at very high antigen concentrations therefore is accompanied by a faint or absent hook line to indicate a high-dose hook effect.

3. The Control Line: Detection-particles bind at this line regardless of whether antigen is present for control purposes. Test results without the presence of the Control Line are always invalid.

Test strips are evaluated visually. The test procedure is rapid and does not require laboratory equipment for useful results in most applications. A video explaining and illustrating the principle of the hook effect can be found on our website ([Rapid allergen testing – Food & Feed Analysis](#)) or on our YouTube channel.

4. Reagents provided

One kit contains sufficient reagents for 15 determinations. Each test kit contains:

Component	Format	Volume
Test strip	Ready-to-use in a resealable container	15 pieces
Reaction vials	Contains labelled antibodies in a stabilized and dried form that must be reconstituted with running buffer	15 pieces
Running buffer	Ready-to-use in a dropper bottle	8 mL
Positive control	Needs to be dissolved with 1 mL distilled or deionized water	1 pieces
Single-use Pasteur pipettes	With 100 µL calibration mark	15 pieces
Cotton swabs	To sample surface	16 pieces
Pre-filled swabbing buffer	Contains 1 mL PBS	15 pieces
Evaluation card		1 piece

5. Additional materials and equipment required, but not provided

- Laboratory gloves

For testing food samples:

- Analytical Scale with a weighing range of not less than 50 g and accuracy of ± 0.01 g)
- Mixer, mortar, ultra-turrax or homogenizer
- Distilled or deionized water
- Centrifuge capable of at least 2,000 x g plus centrifuge vials with lids
- As an alternative: folded filters with a pore size of 8 to 12 µm
- If required: vortexer

- If required: bioavid Absorbent buffer (Art. No. BS810-15); Recommended in particular for samples containing polyphenol (e.g. coffee, cocoa, red wine) or oils. This buffer is also suitable for neutralizing slightly acidic or alkaline samples (Note: For testing, the extract pH value should be neutral (pH 6 to 8)).

6. Precautions

Use the test only for its intended purpose.

Only trained laboratory personnel may perform the test. Follow the instructions for use strictly.

Test strips are single-use items. Use a separate test strip for each sample.

If stored in a refrigerator, bring the test kit to room temperature (20 to 25 °C / 68 to 77 °F) before use. Do not open the test strip container while it is colder than its surroundings, since moisture or condensation can damage the strips.

Do not return removed test strips to the test strip container.

Test strips are sensitive to moisture. Moisture can negatively affect test results, so protect the strips from humidity at all times. This is particularly important when opening the test strip container. Remove test strips from the packaging only immediately before use and only after they have reached room temperature. Close the container tightly immediately after removing a test strip.

Do not use the test kit after the expiration date (See test kit outer label).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

Responsibly dispose all reagents and materials after use, with due regard for the protection of human health and the environment. Follow all applicable national disposal regulations (e.g. applicable waste management and hazardous substance regulations).

7. Storage instructions

Store the test at 2 - 25 °C (36 - 77 °F).

8. Indication of instability or deterioration of reagents

- A missing control band

9. Sample preparation

Clean all surfaces and equipment (For example mixer) thoroughly before and after each sample to prevent allergen contamination caused by airborne allergens or unclean laboratory equipment.

Bring all reagents to room temperature (20 to 25 °C / 68 to 77 °F) before use.


Perform sample preparations and test procedure at room temperature (20 to 25 °C / 68 to 77 °F).

Positive control (if needed):

Add 1 mL of distilled or deionized water to the positive control and incubate for 5 minutes. Mix thoroughly. Then use 0.1 mL as a sample into the test. The positive control can be stored for 12 months or more when stored at -20 °C (-4 °F). Avoid repeated freezing and thawing. Dissolved or thawed controls must be used within two days if stored at room temperature (20 to 25 °C / 68 to 77 °F).

9.1 Solid food samples

1. Thoroughly mix, chop or grind a representative amount of sample (At least 5 to 50 g) and dilute the sample 1:10 (1+9) in distilled or deionized water (e.g. 1 g sample plus 9 mL distilled or deionized water). Mix until a homogenous suspension is obtained (e.g. using a vortexer).
2. Centrifuge 2 min at 2,000 x g or filter. The supernatant or extract should be as particle-free as possible.

9.2 Liquid food samples

- Mix a representative amount of the sample (At least 5 to 50 mL) thoroughly and dilute 1:10 (1+9) in distilled or deionized water (e.g. 1 mL sample plus

9 mL distilled or deionized water). Mix until a homogenous suspension forms (e.g. using a vortexer).


- If the extract contains particles, centrifuge for 2 min at 2,000 x g or filter it. The supernatant or extract should be as particle-free as possible.

9.3 (CIP) water

- Use 0.1 mL of the water directly in the test.

9.4 Surfaces (swab samples)

1. Use one PBS pre-filled vial per sample. For dry surfaces, moisten a clean swab with the PBS.
2. Swab the area of interest using a crosshatch technique in all directions (See graphic on the right). A surface area of approximately 10 x 10 cm is recommended. Sampling a specific problem area (e.g. corners of equipment) is also possible.
3. Transfer the swab sample from Step 1 to the vial and rinse it out: Press the swab against the vial wall, gently rotate it, and move it up and down.






10. Test implementation

Write the sample ID on a reaction vial and open it.

1. Add 0.1 mL (three drops) of running buffer to the reaction vial.
2. Add 0.1 mL of the prepared sample (See Chapter 9) to the reaction vial (Single-use plastic pipettes with a 0.1 mL scale are included in the kit) Mix the content by gently pipetting up and down about 3-times. Alternatively, close the reaction vial and mix thoroughly by tapping the vial several times on the top of table.


After the dried antibody is dissolved only in running buffer, the solution should take on a faint red color. When the dried antibody is mixed with running buffer and sample, this faint red coloration may no longer be visible, depending on the coloring of the sample.

3. Incubate for 5 minutes.
4. Immediately before use, take an allergen test strip out of the test strip container and directly close the test strip container again. Insert the test strip into the reaction vial with arrows pointing downward.
5. Read the result exactly at 5 minutes. After that time, results are not valid.

11. Result interpretation

Allergen test strips include three lines in the reaction area: A control line (C, upper line), the hook line (H, middle line) and a test line (T, lower line). One to three lines may appear in the test strip reaction area. The control line (C) must be present for any valid test run, since it indicates that the test was performed correctly. If the control line is not present, the result is invalid, and tests must be repeated. The test line (T) indicates the presence of the target residue in the sample. Very high concentrations of allergen residues may reduce the intensity of the test line or suppress its formation completely. The hook line (H) serves as an “alarm line”. If it is missing, a high content of allergen is present in the sample. In this case, a missing test line must not be interpreted as a negative result.

12. Limits of the method

The instructions for use, in particular as to sample preparation, procedure, and storage, must be followed to obtain valid results within the stated limits of detection.

The limit of detection is dependent on sample type and extraction efficiency, surface characteristics of swabbed surfaces, and the type of contamination.

Because of the wide variety of food products, test kit manufacturers cannot validate all foods for the test kit. Users are ultimately responsible for verifying suitability of a specific food product for use with the test. Test results may vary depending on sample matrix.

Some solid samples containing cereal flour or thickening agents may impact flow behavior at dilutions of 1:10. Very viscous samples should preferably be tested at a dilution of 1:20. Reduced sensitivity must be expected due to higher sample dilution.

The bioavid Absorbent Buffer (Art. No. BS810-15) is recommended for testing samples containing polyphenols (e.g. coffee, cocoa, spices, wine) or samples containing oil.

Acidic or alkaline samples may impair the test performance and must be neutralized before or after extraction (pH 6 to 8).

Pure nuts may interfere with the hook line development and must be diluted 1:20 (Instead of 1:10; see Section 5.1.1) with distilled or deionized water. Reduced sensitivity must be expected due to higher sample dilution.

The bioavid lateral flow tests are qualitative tests. They do not allow quantitation of the test results based on test line intensity.

Improper weighing affects measurement results (For example, a deviation of +10 % leads to a measured concentration that is 10 percent higher). Reliable measurement results are generally achieved only when weighing deviations do not exceed ± 1 %.

13. Notes

Refer to the current validation report for detailed results and additional information. Data from comparative laboratory tests and interlaboratory comparisons may also be available.

For verification of unknown food samples, Standard Operating Procedures (SOP) are available upon request free of charge.

14. Recommendations

To ensure high analytical reliability, we recommend:

- Follow general laboratory quality assurance requirements set out in standards EN 15633-1 and 15842.
- Acidic or alkaline samples may impair the test performance and must be neutralized before or after extraction (pH 6 to 8).


- Perform spiking experiments to verify correct execution of the determination and the absence of interference.
- Use a RIDASCREEN® ELISA to quantify positive results.
- To investigate unexpected results, analyze the same using a RIDASCREEN® ELISA or run a PCR (e.g. SureFood®).
- Use the free “RIDA®SMART APP Allergen” documentation app available in the Google Play Store (Android devices only).

For additional product information and applications, please contact:
sales@r-biopharm.de.


Version overview

Version number	Chapter and title
2023-02-01	Release Version
2023-07-18	Revision
2023-10-20	Revision
2024-02-27	Revision
2024-07-30	Revision
2024-10-14	Revision
2024-11-06	Revision
2025-04-29	Revision
2026-04-27	Current Version Changes made: <ul style="list-style-type: none"> – Alignment with R-Biopharm design. – General linguistic revision. – Additions to food sample specifications; Revised specifications: changed result reporting to protein (instead of allergens). – Chapter 3: Additional note on test procedure. – Chapter 4: linguistic revision of provided reagents – Added chapter 5: Materials required, but not provided. – Chapter 6: Additional information on test procedure, use of separate test strips and hazardous substances. – Added chapter 8. Indication of instability or deterioration of reagents. – Chapter 9: Additional information on processing temperature. – Chapter 10: Additional note on validity of results. – Chapter 12: additional information on following the instructions, thickening agents and improper weighing. – Addition of an Icon Legend. – Updated Disclaimer Section.

1 Sample preparation

2 Test procedure

Icon Legend

- General Icons:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM-DD)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM-DD)

Manufacturer + address

Disclaimer

1. R-Biopharm AG provides a limited statutory warranty (“Gewährleistung”) as mandated and pursuant to German law against defects in title and material. For products with a defined shelf life, the warranty applies for the stated shelf life. For use cycle limited products, it applies until the stated cycle limit is reached, but not beyond the general warranty period. For all other products, the general warranty period is twelve (12) months from the date the risk transfers.
Such warranty is subject to timely notice of the alleged defect, giving sufficient detail to assess the claim. A product’s defectiveness is assessed based on the specifically agreed quality, its suitability for the contractually stated purpose, and whether all agreed-upon accessories, documentation, and instructions were provided.
2. R-Biopharm AG is not liable for issues arising from:
 - a) The failure to read, understand, or follow product instructions and information for use or safety
 - b) The use of unqualified personnel
 - c) The failure to follow industry standards such as Good Laboratory Practices
 - d) The use of unsuitable controls, sample matrices, or procedures
 - e) Other misuse
 - f) unilateral product modifications
 - g) Improper storage by others
 - h) External chemical, electromagnetic, mechanical, or electrolytic influences that exceed the defined parameters.
 - i) Harm caused by external events beyond R-Biopharm's control, including but not limited to burglary, theft, lightning, fire, water, or acts of nature.
3. R-Biopharm AG’s liability is limited to foreseeable, contract-typical losses in case of negligent breaches of obligations that are essential to achieving the agreement’s purpose and that the other party could reasonably rely on (material contract obligations). R-Biopharm AG is not liable for losses resulting from non-material contract obligations.
4. ALL OTHER WARRANTIES OR GUARANTEES OF ANY KIND, EXPRESS OR IMPLIED, ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, COURSE OF DEALINGS BETWEEN THE PARTIES OR OTHERWISE. R-Biopharm AG assumes no liability for consequential damages, in particular loss of profit, interruptions, or other indirect damages.
5. The foregoing is not intended to limit liability to the extent such liability cannot be excluded or limited under applicable German law or convention, including cases of intent, willful misconduct, or gross negligence, personal injury or death, express warranties or assumptions of procurement risks, strict liability, or other non-waivable legal mandates.
6. The foregoing liability provisions apply equally to R-Biopharm AG, its directors, officers, employees, and agents.

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© 2026 R-Biopharm AG
are registered trademarks of R-Biopharm AG.