bioavid Lateral Flow Mustard inkl. Hook-Linie

Art. Nr. BLH703-15

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

Art. No. BLH703-15

INTENDED USE

The Lateral Flow Mustard (Art. No. BLH703-15) with included hook line is an immunochromatographic test for the sensitive and qualitative detection of mustard residues on surfaces (e.g. swab test for the hygiene control in food production lines), in cleansing and process water (CIP water) and food samples. The integrated hook line prevents false negative interpretation of high positive samples.

Time requirement: Sample preparation (swab, CIP water) approx. 1 min
 Sample preparation (solid samples) approx. 5 min
 Test implementation (incubation time) 10 min

Limit of Detection: Mustard seed 1 ppm
 Spiked zero matrices 10 - 20 ppm
 Swab samples $\leq 0.08 \mu\text{g}$ mustard / cm^2
 CIP water 1 ppm

SPECIFICITY

The test detects mustard residues down to 1 parts per million (ppm = mg/kg). Very high concentrations of mustard residues may reduce the intensity of the test line (> 10,000 ppm) as well as the hook line (> 1,000 ppm), indicating high levels of contamination (high positive samples).

- Hook effect starting at 1,000 ppm mustard, detected by the disappearance of the hook line.
- 0.1 % cross-reactivity to rucola.
- 100 % cross-reactivity to rapeseed.

VALIDATED MATRICES

Spiked zero matrices	Further
Baby food	Mustard seed
Curry spices	Surface (swab samples)
Canned chicken-pasta	CIP water
Canned meat-pasta	
Canned creamy mushroom sauce-pasta	
Canned spaghetti sauce-pasta	
Canned mixed beans	
Mixed flour	

1. BACKGROUND INFORMATION

People allergic to mustard will react to any ingredient of black, Indian brown or yellow / white mustard. Mustard consists of oil, glycoside and more than 20 % protein, some of which have an allergenic effect. Even small traces of allergens in a food can cause allergic reactions in sensitive people. Mustard is used as a condiment in meat, to prepare pastes or as an ingredient for a large variety of sauces, soups, sausages, salad dressings, and marinades. The allergen can be present as an ingredient or as a contamination in raw and cooked products. Thus, it is frequently encountered in food as a hidden allergen. Considerable deficiencies in the labelling and quality assurance of allergens can be observed, particularly for mustard. Many products, which are supposed to be mustard free, are already contaminated on the field or during processing. Cross contamination of spices with the allergen likely occurred during processing (dust formation, during crushing) or packaging. In bread available in bakeries, a mustard contamination is often observed.

Food production lines and the food itself have to be free of residues to prevent allergic reactions in sensitive persons. The effectiveness of cleaning can be verified by using a lateral flow method for mustard detection. According to the regulation (EU) No. 1169/2011, mustard and products thereof must be declared on food labels.

2. BRIEF INFORMATION

The bioavid Lateral Flow Allergen Devices (LFD; dip stick) with included hook line are effective tools for the fast, sensitive, easy and reliable detection of allergen contamination in production lines, rinse water (CIP water), surfaces and in food products. Summarized information about the usage of LFDs within an allergen management can be found on our infographic, which can be downloaded on www.r-biopharm.com.

Most important challenges with the safe and fast LFD analysis are false negative results due to very high allergen concentrations in the test sample (so-called High-Dose-Hook-Effect). The hook effect (or overload effect; Standard EN 15633-1. 2019) is observed when a very high amount of an analyte is present in the sample. In this case, the development of the test line can be inhibited and there is a risk that the test result will be falsely evaluated as negative (or analyte undetectable). The hook line on the test strip makes this effect visible. A faint or missing hook line indicates a high allergen content in the sample and can no longer be evaluated as negative (see chapter 9).

3. TEST PRINCIPLE

The principle of the test is an antigen-antibody reaction. The reaction vial contains mustard-specific antibodies. After addition of running buffer, the sample (e.g. supernatant, swab sample, CIP water) is added to the reaction vial. If the sample contains the target antigen, an antigen-antibody complex will form. After incubation, a dip stick is placed into the reaction vial. Capillary forces cause the liquid to flow up the dip stick. The dip stick contains three areas, where reactions take place and lines for visual evaluation can form. First, the control line in order to verify the test run and which should appear every time. Second, the test line, where antigen-specific antibodies are immobilized. If target antigen is present, the previously formed antigen-antibody-complexes react with the immobilized antibodies on the test line resulting in the formation of an antibody-antigen-antibody complex. Third, the hook line, where target analyte is immobilized and will appear, as long as target-analyte-specific antibodies are present unbound. The result is read visually. The test is very rapid and reliable and does not require laboratory equipment for most applications.

4. STORAGE

Store the test kit at 2 - 25 °C (35 - 77 °F).

5. REAGENTS PROVIDED

Each kit contains sufficient materials for 15 measurements.

BLH703-15		
Content	Description	Amount
Reaction vials	Containing labeled antibodies in stabilized dried form, needs to be dissolved with running buffer	15 x
Allergen test strips	1 tube, reclosable	15 x
Dropper bottle with running buffer	8 mL	1 x
Positive control	Needs to be dissolved with 1 mL distilled or deionized water	1 x
Plastic pipettes	With 100 µL mark	15 x
Plastic swab	Used to swab surfaces	16 x
Pre-filled PBS vials	Containing 1.2 mL PBS; used for swabbing (moisten the swab, swab and flush swab)	15 x
Evaluation card		1 x

6. WARNINGS AND PRECAUTIONS FOR THE USERS

Bring the test kit to room temperature (20 - 25 °C / 68 - 77 °F) before use. Do not open the tube with test strips when it is colder than the environment, since strip deterioration may result from moisture.

Once removed, test strips should not be returned to the test strip tube.

The test strips are sensitive to moisture. Therefore, remove the test strips from the test strip tube immediately before use only.

Do not use the test kit after the expiration date (see test kit label).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).


7. PREPARATION OF SAMPLES

Positive control (if needed):

Add 1 mL distilled or deionized water to the positive control and leave for 5 minutes. Then apply 0.1 mL as a sample into the assay. The positive control can be stored at -20 °C (-4 °F) for at least one year. Do not repeatedly freeze and thaw. Dissolved or thawed controls are stable for 2 days at room temperature (20 - 25 °C / 68 - 77 °F).

7.1. Food samples:

1. Grind/mix a representative amount of sample (e.g. 5 - 50 g) thoroughly and solve the sample 1:10 (1+9) in distilled or deionized water (e.g. 1 g sample + 9 mL distilled or deionized water).
2. Centrifuge 2 min at 2,000 x g or filter. The supernatant or extract should be as particle-free as possible.


7.2. CIP water:

1. 0.1 mL of cleansing or process water (CIP water) can directly be used in the assay.

7.3. Swab samples (surfaces):





1. Use one vial (pre-filled with PBS) per sample. Moisten a clean cotton swab with the contained PBS.
2. Swab the area of interest with a crosshatch technique in all directions (see graphic on the right).
3. Release the sample from the swab into the vial from step 1: squeeze the swab against the vial, carefully rotate and move up and down.

8. ASSAY


Note the sample ID on a reaction vial and open it.

1. Add 0.1 mL (3 drops) of running buffer into the reaction vial.
2. Add 0.1 mL of sample (supernatant, CIP water or swab sample) into the reaction vial (single-use plastic pipettes with 0.1 mL scale are contained in the kit) and mix the content by carefully pipetting up and down (approx. 3-times) or close the reaction vial and mix intensively by carefully tapping the vial on the desk a few times.
3. Incubate for **5 minutes**.
4. Immediately before use, take an allergen test strip out of the test strip tube. Immediately close the test strip tube again. Insert the test strip with the arrows pointing down into the reaction vial.
5. Read the result at **5 minutes**.

9. INTERPRETATION OF THE RESULT

The allergen test strip includes three lines in the reaction field: the control line (C, top line), the hook line (H, middle line) and the test line (T, lowest line). In total, 1 to 3 purple lines may appear in the reaction field on the test strip. The control line (C) indicates that the test has been performed correctly and consequently must be present on every valid test strip. The test line (T) indicates the presence of the target residue in the sample. Very high concentrations of mustard residues may reduce the intensity of the test line or suppress its formation completely. The hook line (H) serves as an “alarm line”. If it is faint or missing, a high content of allergen is present in the sample. In this case, a missing test line must not be interpreted as a negative result.

- A C-line and H-line present, no T-line: **valid negative** result.
- B All 3 lines present: **valid positive** result.
- C C-line and T-line present, H-line faint or absent: **valid high positive** result.
- D Only C-line present: **suspected high positive**, repeat test with higher sample dilution.

No C-line: test is **invalid**.

10. LIMITATION

Due to the large amount of different foods, it is not possible for the manufacturer of a test kit to validate all foods for the test. The suitability of the respective food for use in the test must be verified by the user himself.

For verification of unknown food samples, an SOP is available on the product-specific website www.r-biopharm.com.

Some solid samples containing cereal flour and/or thickeners may interfere with a correct development of the hook line if applied in a 1:10 dilution. Thus, very viscous samples should be tested in a 1:20 dilution. In this case, lower sensitivity should be noted due to higher sample dilution.

Pure nuts interfere with the development of the hook line and must be diluted 1:20 (instead of 1:10) with distilled or deionized water.

To ensure high analytical reliability, it is recommended to set the pH value to neutral for extremely acidic ($\text{pH} \leq 5$) or basic ($\text{pH} \geq 11$) samples.

Bioavid lateral flow tests are qualitative tests. They do not allow quantitative evaluation even in the case of variably developed test lines (T).

11. NOTES

If a sample after extraction was too viscous to allow the control line to appear within 2 minutes, a higher dilution is to be tested.

Further information can be found in the validation report, which is available upon request.

The QC certificates are available upon request.

12. DISCLAIMER


The user assumes all risk in using bioavid Diagnostics GmbH & Co. KG products.

Bioavid Diagnostics GmbH & Co. KG will warrant that its products meet all quality control standards set by bioavid Diagnostics GmbH & Co. KG, and bioavid Diagnostics GmbH & Co. KG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. Bioavid Diagnostics GmbH & Co. KG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of bioavid Diagnostics GmbH & Co. KG.

1 Sample preparation

2 Test procedure