bioavid Swabbing kit

Art. Nr. BS803-30

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Art. No. BS803-30

1. INTENDED USE

The bioavid Swabbing kit (Art. No. BS803-30) is used for the sampling of allergen residues on surfaces (e.g. in food production lines).

2. STORAGE

Store the test kit at 2 - 25 °C (35 - 77 °F).

3. REAGENTS PROVIDED

Each kit contains sufficient materials for 30 samplings.

BS803-30		
Content	Description	Amount
Sample tube	Containing 1.2 mL PBS	30 x
Plastic swab	Used to swab surfaces	30 x
Plastic pipettes	With 100 µL mark	30 x

4. WARNINGS AND PRECAUTIONS FOR THE USERS


Bring the buffer to room temperature (20 - 25 °C / 68 - 77 °F) before use.

Do not use the test kit after the expiration date (see test kit label).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

5. PREPARATION OF SAMPLES

1. Use one vial (pre-filled with PBS) per sample. Moisten a clean cotton swab with the contained PBS.
2. Swab the area of interest with a crosshatch technique in all directions (see graphic below).
3. Release the sample from the swab into the vial from step 1: squeeze the swab against the vial, carefully rotate and move up and down.
4. 0.1 mL of the supernatant is used in the test for testing for allergens (for further procedure, see product-specific instruction for use chapter "Test procedure").

6. DISCLAIMER

The user assumes all risk in using bioavid Diagnostics GmbH & Co. KG products.

Bioavid Diagnostics GmbH & Co. KG will warrant that its products meet all quality control standards set by bioavid Diagnostics GmbH & Co. KG, and bioavid Diagnostics GmbH & Co. KG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. Bioavid Diagnostics GmbH & Co. KG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of bioavid Diagnostics GmbH & Co. KG.

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Der Anwender trägt das alleinige Risiko bei der Verwendung der Produkte der bioavid Diagnostics GmbH & Co. KG.