bioavid Absorbent Buffer

Art. Nr. BS810-15

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Art. No. BS810-15

INTENDED USE

The bioavid Absorbent Buffer (Art. No. BS810-15) is particularly suitable for preparing samples containing polyphenols (e.g. coffee, cocoa, red wine) for the bioavid rapid allergen tests. It is also suitable for neutralizing liquid samples with a strongly acidic or alkaline pH (orange juice, wine).

1. STORAGE

Store the bioavid Absorbent Buffer at 2 - 25 °C (35 - 77 °F).

2. REAGENTS PROVIDED

Each kit contains sufficient materials for 15 measurements.

BS810-15		
Content	Description	Amount
Sample tube	Containing 9 mL absorption buffer each	15 x

3. WARNINGS AND PRECAUTIONS FOR THE USERS

Bring the bioavid Absorbent Buffer to room temperature (20 - 25 °C / 68 - 77 °F) before use.

Do not use the test kit after the expiration date (see test kit label).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

4. PREPARATION OF SAMPLES

Shake the test tube before use to suspend the insoluble components.

1. Place 1 mL liquid sample or 1 g water-soluble or homogenized sample into a sample tube.
2. Mix for 10 - 15 minutes (overhead shaker, rotator or similar).
3. Centrifuge at 2,000 - 3,000 x g for 3 - 5 minutes or filter.
4. After centrifugation or filtration, a sediment of approx. 0.5 - 1 mL volume is visible.
5. 0.1 mL of the supernatant is used in the test for testing for allergens (for further procedure, see product-specific instruction for use chapter "Test procedure").

Note: As a use to neutralize liquid samples with strongly acidic or alkaline pH (orange juice, wine), no incubation (step 2) is required. Immediately after addition and mixing of 1 mL sample, centrifugation or filtration can be performed.

5. LIMITATION

Due to the large amount of different foods, it is not possible for the manufacturer of a test kit to validate all foods for the test. The suitability of the respective food for use in the test must be verified by the user himself.

The bioavid Absorption Buffer has not yet been fully tested. It is therefore recommended to analyze a spiked food sample in parallel. For this purpose, add 10 µg to 100 µg of the allergen to be tested to the sample in the sample tube.

6. DISCLAIMER

The user assumes all risk in using bioavid Diagnostics GmbH & Co. KG products.

Bioavid Diagnostics GmbH & Co. KG will warrant that its products meet all quality control standards set by bioavid Diagnostics GmbH & Co. KG, and bioavid Diagnostics GmbH & Co. KG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. Bioavid Diagnostics GmbH & Co. KG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of bioavid Diagnostics GmbH & Co. KG.