
CONGEN

SureFast® PREP DNA/RNA Virus

Art. No. F1051
100 extractions

User Manual

Efficient DNA preparation of viral or bacterial DNA/RNA

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Content

1	General Information	10
1.1	Description	10
1.2	Kit components (per box) and storage	10
1.3	Additionally required equipment and materials	11
1.4	Precautions for users	11
2	Protocol	12
2.1	Principle	12
2.2	Preparations	12
2.2.1	General	12
2.2.2	Before each preparation	12
2.3	Preparation of the basic material	12
2.3.1	Liquid samples	12
2.3.2	Stool samples	12
2.3.3	Swabs and filters	12
2.3.4	Sample preparation of solid and semi-solid samples	12
2.4	Sample preparation for surface swabbing for the detection of SARS-CoV-2	13
2.4.1	Surface swabbing	13
2.4.2	Sample preparation for surface swabbing for direct analysis without transport of the swab	13
2.4.3	Sample preparation for surface swabbing with transport medium (UTM) after the external transport (48 hours maximal transport duration)	13
2.5	Extraction protocol	13
3	Further Information	14
3.1	Product Information	14
3.2	Technical Support	14
4	Safety Information	15

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Kit Code	Reagenz /Material	Menge (per box)
B	Binding Buffer *	1 x 30 ml
E	Elution Buffer	1 x 5 ml
P	Pre-Wash Buffer **	1 x 40 ml
W	Wash Buffer **	1 x 60 ml
ET	Extraction Tubes	1 x 50 Tubes
R	Receiver Tubes 2,0 ml	1 x 50 Tubes
T	Receiver Tubes 1,5 ml	1 x 50 Tubes
S	Spin Filter Set	1 x 50 Filter

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1 General Information

1.1 Description

The kit is intended to be used for the isolation of DNA and RNA viruses nucleic acids from serum, urine, plasma, cell culture supernatants, foods (for example wash up fluids from fruits, salads etc.), filters from water samples as well as stool samples and swabs (from work surfaces etc.).

If an Internal Control RNA/DNA (ICR/ICD, not supplied in the kit) is used as an extraction control for the sample preparation procedure and as inhibition control, it should be added during extraction procedure. The ICR/ICD should always be pipetted to the already mixed sample inside the Extraction Tube (**Code ET**) (see 2.5.1 Lysis of the basic material as well as information in the corresponding SureFast® PCR User Manual).

1.2 Kit components (per box) and storage

Kit Code	Reagent /Material	Amount (per box)
B	Binding Buffer *	30 ml
E	Elution Buffer	1 x 5 ml
P	Pre-Wash Buffer **	1 x 40 ml
W	Wash Buffer **	1 x 60 ml
ET	Extraction Tubes	1 x 50 Tubes
R	Receiver Tubes 2,0 ml	1 x 50 Tubes
T	Receiver Tubes 1,5 ml	1 x 50 Tubes
S	Spin Filter Set	1 x 50 Filter

* After adding isopropanol (not supplied with the kit, see 1.3 Additionally required equipment and materials and 2.2 Preparations)

** After adding ethanol (not supplied with the kit, see 1.3 Additionally required equipment and materials and 2.2 Preparations)

The components are equally divided into two boxes.

All reagents of the kit should be stored dry and at room temperature (14-25°C). After opening the stability of the kit is reduced to 12 months.

1.3 Additionally required equipment and materials

- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes free from DNA and DNase 1.5 ml; 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to 100°C)
- micro centrifuge (up to 12,000 rpm)
- Ethanol for preparation of Wash Buffer and Pre-Wash Buffer (puriss, $\geq 96\%$)
- isopropanol for preparation of Binding Buffer (e.g. Carl Roth 2-Propanol Rotipur® >99.7%; Applichem 2-Propanol für Mikrobiologie; Sigma-Aldrich 2-Propanol)
- RNase free PCR grade H₂O
- disposal bags or waste bin
- autoclave
- biological safety cabinet

1.4 Precautions for users

- This test must only be performed by qualified and authorized laboratory personnel.
- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Samples must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.
- For more details see Material Safety Data Sheet, at www.congen.de/en/eifu/.

2 Protocol

2.1 Principle

1. Preparation of the basic material
2. Lysis at 65°C and 95°C
3. Setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

2.2 Preparations

2.2.1 General

Add 30 ml isopropanol to the Binding Buffer (**Code B**).

Add 20 ml ethanol to the Pre-Wash Buffer (**Code P**) and mix thoroughly.

Add 48 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

2.2.2 Before each preparation

Preheating the Elution Buffer (**Code E**) – Transfer the needed amount under calculation of a reserve volume of Elution Buffer into a reaction tube (not provided with the kit).

Equilibrate the Elution Buffer to 65°C. The Elution Buffer is necessary for step 6.

2.3 Preparation of the basic material

2.3.1 Liquid samples

Mix 200 µl samples fluids with 200 µl RNase free PCR grade H₂O in a reaction tube (not supplied with the kit). For samples which have a smaller volume than 200 µl please fill up to a total volume of 400 µl with RNase free PCR grade H₂O.

By small yields of nucleic acids use directly 400 µl sample fluid. Transfer the complete sample in an Extraction Tube (**Code ET**).

2.3.2 Stool samples

Dilute stool samples 10 to 20 fold in RNase free PCR grade H₂O (for a lower viscosity). Centrifuge the dilution for 2 min at maximum speed. Mix 200 µl supernatant with 200 µl RNase free PCR grade H₂O. Transfer the complete sample in an Extraction Tube (**Code ET**).

2.3.3 Swabs and filters

Place swabs and filters complete in the Extraction Tube (**Code ET**) (Cut off the shaft of the swab, so that the cap of the extraction tube can be closed). Add 400 µl RNase free PCR grade H₂O.

2.3.4 Sample preparation of solid and semi-solid samples

Solid and semi solid samples as dried fruits, fruits, salads, bivalvia (mussels, oysters) will be treated for sampling procedure according to special application notes.

For further questions please send an e-mail to marketing@congen.de

2.4 Sample preparation for surface swabbing for the detection of SARS-CoV-2

2.4.1 Surface swabbing

The swab should be pre-moistened with 100 µl UTM prior of swabbing the surface. Swab a surface of approximately 25 cm² with the wet swab.

For transport to an external laboratory, put the swab back into the Copan tube and send it to the designated laboratory. After transport follow with 2.3.3.

2.4.2 Sample preparation for surface swabbing for direct analysis without transport of the swab

Place the swab with the shaft in an Extraction Tube (**Code ET**) and cut the shaft, so that the lid of the Extraction Tube can be closed. Add 400 µl RNase free PCR grade H₂O.

2.4.3 Sample preparation for surface swabbing with transport medium (UTM) after the external transport (48 hours maximal transport duration)

Place the swab with the shaft in an Extraction Tube (**Code ET**) and cut the shaft, so that the lid of the Extraction Tube can be closed. Add 400 µl of the transported UTM. Transportation of the swab samples with UTM should take place at 2-8°C.

Note: After lysis time carefully squeeze out the swab on the wall of the tube and discard the swab.

2.5 Extraction protocol

1. Lysis of the basic material

Place the Extraction Tubes (**Code ET**) into a Thermomixer and incubate for 15 minutes at 65°C.

Incubation for 10 minutes at 95°C under continuously shaking.

Note: If an internal control RNA/DNA (not supplied in the kit) is used as an extraction control for the sample preparation procedure and as inhibition control, 20 µl of the ICR/ICD should be added during extraction procedure. The ICR/ICD should always be pipetted to the already mixed sample inside the Extraction Tube (**Code ET**) and must not be added directly to the raw sample.

2. Setting of optimal binding conditions

Add directly 400 µl Binding Buffer (**Code B**) to fluids without particles and mix the sample by vortexing.

Note: Does the sample lysate show particles it has to be centrifuged for 1 min at maximum speed. The liquid supernatant has to be transferred into a new reaction tube (not supplied with the kit). Add 400 µl Binding Buffer (**Code B**) to the supernatant and mix the sample by vortexing.

3. Binding of the nucleic acids on a Spin Filter

Transfer the complete sample to a Spin Filter Set (**Code S**). Close the cap and centrifuge for 1 min at 12,000 rpm. Discard the Receiver Tube with the filtrate and place the Spin Filter in a new 2.0 ml Receiver Tube (**Code R**).

4. Purification of the bound nucleic acids

Add 500 µl Pre-Wash Buffer (**Code P**) to the Spin Filter and centrifuge for 1 min at 10,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Add 700 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 10,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

5. Drying of the Spin Filter

Remove the residual ethanol by final centrifugation for 4 min at maximum speed.

6. Elution of nucleic acids from the Spin Filter

Place the Spin Filter into a clear 1.5 ml Receiver Tube (**Code T**) and add 60 µl of the preheated to 65°C Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate for 3 min and centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

The eluted nucleic acid is ready-to-use for the PCR or RT-PCR. For longer storage, the nucleic acids should be stored at -20°C or -80°C.

Note: The nucleic acids can also be eluted with a lower (but not lower than 40 µl) or a higher volume of Elution Buffer (depends on the expected yield or needed concentration of the nucleic acids).

3 Further Information

3.1 Product Information

- Flow chart (Download: www.congen.de/en/downloads)
- Validation data upon request
- Product-related documents (Download: www.congen.de/en/eifu/)

3.2 Technical Support

For further questions please send an e-mail to marketing@congen.de.

4 Safety Information

Pre-Wash Buffer

Warning

H302-H314-H412
 P260-P264-P270-P273-P280-P301+P312-
 P301+P330+P331-P303+P361+P353-P304+P340-
 P501

Extraction Tubes

Danger

H302-H315-H318-H373-H410
 P260-P264-P270-P273-P280-P301+P312-
 P302+P352-P305+P351+P338-P310-P314-P321-
 P330-P332+P313-P362+P364-P501-EUH208

H302:	Harmful if swallowed.
H314:	Causes severe skin burns and eye damage.
H315:	Causes skin irritation.
H318:	Causes serious eye damage.
H373:	May cause damage to organs through prolonged or repeated exposure.
H410:	Very toxic to aquatic life with long lasting effects.
H412:	Harmful to aquatic life with long lasting effects.
P260:	Do not breathe dust/fume/mist/vapours/spray.
P264:	Wash hands, forearms and face thoroughly after handling.
P270:	Do not eat, drink or smoke when using this product.
P273:	Avoid release to the environment.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P301+P312:	IF SWALLOWED: Call a POISON CENTRE or doctor if you feel unwell.
P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P302+P352:	IF ON SKIN: Wash with plenty of water.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P305+P351+P338:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310:	Immediately call POISON CENTER or doctor.
P314:	Get medical advice/attention if you feel unwell
P321:	Specific treatment (see supplemental first aid instruction on this label).
P330:	Rinse mouth.
P332+P313:	If skin irritation occurs: Get medical advice/attention.
P362+P364:	Take off contaminated clothing and wash it before reuse.
P501:	Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.
EUH208:	Contains Proteinase K, Tritirachium album serine (39450-01-6). May produce an allergic reaction.

For further information we offer a Material Safety Data Sheet, see at www.congen.de/en/eifu/. Alternatively please contact your distributor or send an e-mail to marketing@congen.de.