**CONGEN**

# **SureCycle®**

Art. No. F4001  
260 rxn

## **User Manual** **System Verification**

**April 2026**

 **Inhalt**

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## **Content**

1	General Information .....	14
1.1	Intended Use.....	14
1.2	Description .....	14
1.3	Kit components and storage .....	15
1.4	Additionally required equipment and materials .....	15
1.5	Precautions for users .....	15
1.6	Setup.....	16
1.7	Detection channel Set-up .....	16
2	Protocol .....	17
2.1	Recommended pipetting scheme for system verification .....	17
2.2	Preparation of the master-mix .....	18
2.3	Preparation of the real-time PCR mix.....	18
3	Interpretation of results .....	18
3.1	SureCycle® profile 1 Denaturation .....	18
3.1.1	Denaturation system (VIC channel).....	18
3.1.2	Annealing system (FAM-channel) .....	20
3.2	SureCycle® Profile 2 Annealing.....	21
3.2.1	Annealing system (FAM-channel) .....	21
3.2.2	Denaturation system (VIC-channel).....	23
4	Further Information .....	23
4.1	Product Information.....	23
4.2	Technical Support .....	23
4.3	Explanation of Symbols .....	23

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## **1 General Information**

### **1.1 Intended Use**

The SureCycle®-Kit is intended to use for real-time PCR routine control in all laboratories. Especially laboratories which need a constant control concerning their testing devices (DIN EN ISO 15189 and DIN EN ISO 9001, GMP and GLP) may use the SureCycle®-Kit as a reliable and fast verification method.

The kit is designed to control compatibility of the selected plastic material and to verify the correct temperature control of the thermocycler.

### **1.2 Description**

The real-time Polymerase chain reaction (PCR) is an established method for the specific and sensitive detection of DNA. The area of use spans the clinical infection diagnostics to the analysis of pathogens in food as well as forensics. Each PCR-reaction basically is defined by the primer sequence and the PCR protocol with specific temperature and time intervals. The detection is therefore dependent to the correct temperature system of the thermocycler in use. To prevent false positive or false negative results, the control and verification of the correct temperature system is obligatory.

The SureCycle®-System is a thermosensitive duplex real-time PCR. The temperature sensitivity of the system results from the instability of the primer binding. Binding of the primer to the complementary chain is being destabilized and thereby the annealing- and denaturation reaction become temperature sensitive. Such a system is reacting very sensitively to 'overshooting' of thermocyclers. The test is equally sensitive regarding the inadequate denaturation of double-stranded DNA (dsDNA) due to a too-low denaturation temperature ('undershooting'). Both, 'overshooting' and 'undershooting', of denaturation temperature reduces the efficiency of the PCR.

To control the correct temperature during the annealing- and denaturation reactions, two thermal profiles are used: SureCycle® Profile 1 Denaturation and SureCycle® Profile 2 Annealing.

With Profile 1 an optimal denaturation temperature is programmed, as well as a suboptimal annealing temperature.

Using Profile 2 the opposite situation is evaluated, here the optimal annealing temperature and the suboptimal denaturation temperature is performed.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for the detection of two fluorescence emissions at the channels FAM and VIC/HEX at the same time.

The internal technical verification of instruments was performed on Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96, Bio Molecular Systems MIC, Agilent AriaDx and Agilent MX3005P.

### 1.3 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	4 x 1400 µl	Yellow
2	Taq Polymerase	1 x 200 µl	Dark Red

**Store all reagents at –28°C to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple use on the same day.**

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

### 1.4 Additionally required equipment and materials

- real-time PCR instrument with two detection channels (510 nm and 580 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

### 1.5 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

**1.6 Setup**

	<b>SureCycle® Profile 1 Denaturation</b>	<b>SureCycle® Profile 2 Annealing</b>
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, <b>95°C</b>	15 sec, <b>92°C</b>
Annealing/Extension	30 sec, <b>63°C</b>	30 sec, <b>60°C</b>
Cooling	10 sec, 40 °C	10 sec, 40 °C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

**1.7 Detection channel Set-up**

**For additional details see cyler operating instructions for real-time PCR device.**

<b>Real-time PCR device</b>	<b>Detection channel</b>	<b>Quencher</b>	<b>Note</b>
<b>Agilent Mx3005P</b>	FAM	+	
	HEX	+	
<b>Agilent AriaDx / Mx</b>	FAM	+	
	HEX	+	
<b>Applied Biosystems 7500</b>	FAM	None	Check the passive reference option ROX is none.
	VIC	None	
<b>Bio-Rad CFX96 /Dx/Opus</b>	FAM	+	Baseline Settings: <ul style="list-style-type: none"> <li>• Baseline subtracted curve fit</li> <li>• Apply fluorescence drift correction</li> </ul>
	VIC / HEX	+	
<b>R-Biopharm RIDA®CYCLER</b>	green	+	<b>Ignore cycles before</b> , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cyler operating instructions, section
	yellow	+	
<b>Qiagen Rotor-Gene Q</b>	green	+	<b>Note:</b> Please use only 0.1 ml reaction tube. The gain settings must be set to 5 (factory default) for all channels.
	yellow	+	
<b>Roche LightCycler® 480 II</b>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	533-580	+	
<b>Roche cobas® z 480 Analyzer</b>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	540-580	+	

## 2 Protocol

### 2.1 Recommended pipetting scheme for system verification

**Blockcyler:**

For verification of the real-time PCR instrument it is recommended to pipette 4 x 4 wells in the corner positions (highlighted in red) and 4 wells in the center (highlighted in red) per each SureCycle® profile (see figure 1).


	1	2	3	4	5	6	7	8	9	10	11	12
A			/	/	/	/	/	/	/	/		
B			/	/	/	/	/	/	/	/		
C	/	/	/	/	/	/	/	/	/	/		
D	/	/	/	/	/			/	/	/		
E	/	/	/	/	/			/	/	/		
F	/	/	/	/	/			/	/	/		
G			/	/	/	/	/	/	/	/		
H			/	/	/	/	/	/	/	/		

**Figure 1: Blockcyler pipetting scheme in 96-well format**

**Rotorcyler:**

It is recommended to pipette 20 reactions per SureCycle®. 5 reactions should be used as reference reactions (see figure 2, position 1-5, highlighted in red).

In the figure below the Qiagen Rotor-Gene Q (72-well disc) is represented as an example. The positioning for use with a different rotor size (Bio Molecular Systems MIC, RIDA®CYCLER) can be adjusted according to the reactions.

**Figure 2: Rotorcyler pipetting scheme for e.g. 72-well disc, Qiagen Rotor-Gene Q**

## 2.2 Preparation of the master-mix

The total number of reactions needed for the PCR should be calculated according to the recommended plate set-up with 20 reactions per SureCycle® profile (see 2.1).

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

### Example for the calculation and preparation of 20 reactions:

Components of the master-mix	Amount per reaction	20 reactions (with 10 % excess)
Reaction Mix	19.3 µl	424.6 µl
Taq Polymerase	0.7 µl	15.4 µl
<b>Total volume</b>	<b>20 µl</b>	<b>440 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

## 2.3 Preparation of the real-time PCR mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start run according to the setup.

# 3 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

A reaction is evaluated as **positive**, if an amplification is shown in the respective channel.

A reaction is evaluated as **negative**, if no amplification is shown in the respective channel.

## 3.1 SureCycle® profile 1 Denaturation

### 3.1.1 Denaturation system (VIC-channel)

The SureCycle® Profile 1 Denaturation takes place with an optimal denaturation temperature at 95°C and suboptimal Annealing/Extension temperature of 63°C. In this case positive results are to be expected for the Denaturation system (VIC-channel).

The FAM-channel controls the correct Annealing/Extension. Since the temperature setting is suboptimal, no amplification should occur and hence negative results are to be expected.

**Blockcycler**

The target value is calculated with the mean value of the reference reactions (Mean value D6/D7/E6/E7, highlighted in red). The target range of control reactions is defined as  $\pm 2$  Cp of the reference reactions. The real-time PCR instrument passed the verification, if Cp-values of all control reactions are within the target range.

The real-time PCR instrument did not pass the verification, if one or more Cp-values of control reactions are out of the target range.

In case there is a Cp-value deviation of a control reaction from the target range, a complete inspection of all wells should be performed.

	1	2	3	4	5	6	7	8	9	10	11	12
A	26.8	26.6									26.8	26.4
B	26.8	26.8									26.9	26.7
C												
D						26.9	27.0					
E						27.0	27.0					
F												
G	26.7	26.8									26.8	26.7
H	26.6	26.8									26.6	26.7

**Figure 3: Example calculation SureCycle® Profile 1 Denaturation, VIC-channel**

Mean value D6/D7/E6/E7: **27.0** Target range  $\pm 2$  Cp: **25.0 – 29.0**


**Rotorcycler**

The target value is calculated with the mean value of the reference reactions (Mean value Position 1-5, highlighted in red). The target range of control reactions is defined as  $\pm 2$  Cp of the reference reactions.

The real-time PCR instrument passed the verification, if Cp-values of all control reactions are within the target range.

The real-time PCR instrument did not pass the verification, if one or more Cp-values of control reactions are out of the target range.

In the figure below the Qiagen Rotor-Gene Q (72-well disc) is represented as an example. The positioning for use with a different rotor size (Bio Molecular Systems MIC, RIDA®CYCLER) can be adjusted according to the reactions.

**Figure 4: Example calculation SureCycle® Profile 1 Denaturation, VIC-channel**

Mean value Pos. 1-5: **27.0** target-range  $\pm 2$  Cp: **25.0-29.0**

**Table 1: Results Interpretation of SureCycle® Profile 1 Denaturation (VIC-channel)**

<b>VIC-channel Denaturation-System</b>			
<b>Reference reaction</b>	<b>Control reaction</b>	<b>Result</b>	<b>Consequence for further analysis</b>
positive	Cp within target range	correct	Real-time PCR instrument temperature system works correctly
positive	Cp out of target range	incorrect	
negative	negative	incorrect	Denaturation temperature too low → potential false PCR-results

**3.1.2 Annealing system (FAM-channel)**

All results should be negative.

The real-time PCR instrument passed the verification, if all reactions are negative.

In case control or reference reactions show positive Cp-values, the results have to be interpreted according to table 2.

**Table 2: Results Interpretation of SureCycle® Profile 1 Denaturation (FAM-channel)**

<b>FAM-channel Annealing-System</b>			
<b>Reference reaction</b>	<b>Control reaction</b>	<b>Result</b>	<b>Consequence for further analysis</b>
negative	negative	correct	Real-time PCR instrument temperature system works correctly
positive	negative	incorrect	
positive	positive	incorrect	
negative	positive	incorrect	

**3.2 SureCycle® Profile 2 Annealing**

**3.2.1 Annealing system (FAM-channel)**

The SureCycle® Profile 2 Annealing takes place with a suboptimal denaturation temperature at 92°C and the optimal Annealing/Extension temperature at 60°C. In this case negative results are to be expected for the Denaturation system (VIC-channel).

The FAM-channel controls the correct Annealing/Extension. Since the temperature setting is optimal, amplification should occur and hence positive results are to be expected.

**Blockcycler**

The target value is calculated with the mean value of the reference reactions (Mean value D6/D7/E6/E7, highlighted in red). The target range of control reactions is defined as  $\pm 2 C_p$  of the reference reactions.

The real-time PCR instrument passed the verification, if  $C_p$ -values of all control reactions are within the target range.

The real-time PCR instrument did not pass the verification, if one or more  $C_p$ -values of control reactions are out of the target range.

In case there is a  $C_p$ -value deviation of a control reaction from the target range, a complete inspection of all wells should be performed.

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	26.8	26.6									26.8	26.4
<b>B</b>	26.8	26.8									26.9	26.7
<b>C</b>												
<b>D</b>						26.9	27.0					
<b>E</b>						27.0	27.0					
<b>F</b>												
<b>G</b>	26.7	26.8									26.8	26.7
<b>H</b>	26.6	26.8									26.6	26.7

**Figure 5: Example calculation SureCycle® Profile 2 Annealing, FAM-channel**

Mean value D6/D7/E6/E7: **27.0** Target-range  $\pm 2 C_p$ : **25.0-29.0**


**Rotorcycler**

The target value is calculated with the mean value of the reference reaction (Mean value Position 1-5, highlighted in red). The target range of control reactions is defined as  $\pm 2 C_p$  of the reference reactions.

The real-time PCR instrument passed the verification, if  $C_p$ -values of all control reactions are within the target range.

The real-time PCR instrument did not pass the verification, if one or more  $C_p$ -values of control reactions are out of the target range.

In the figure below the Qiagen Rotor-Gene Q (72-well disc) is represented as an example. The positioning for use with a different rotor size (Bio Molecular Systems MIC, RIDA®CYCLER) can be adjusted according to the reactions.

**Figure 6: Example calculation SureCycle® Profile 2 Annealing, FAM-channel**

Mean value Pos. 1-5: **27.0** target-range  $\pm 2 C_p$ : **25.0-29.0**

**Table 3: Result Interpretation of SureCycle® Profile 2 Annealing (FAM-channel)**

FAM-channel Annealing-System			
Reference reaction	Control reaction	Result	Consequence for further analysis
positive	$C_p$ within target range	correct	Real-time PCR instrument temperature system works correctly
positive	$C_p$ out of target range	incorrect	
negative	negative	incorrect	Annealing temperature too low → potential false PCR-results

**3.2.2 Denaturation system (VIC-channel)**

All results should be negative.

The real-time instrument passed the verification, if all reactions are negative.

In case control or reference reactions show positive Cp-values, the results have to be interpreted according to table 4.

**Table 4: Result Interpretation of SureCycle® Profile 2 Annealing (VIC-channel)**

FAM-channel Denaturation-System		Result	Consequence for further analysis
Reference reaction	Control reaction		
negative	negative	correct	Real-time PCR instrument temperature system works correctly  Denaturation temperature too high → uncritical, false PCR-results are not to be expected
positive	positive/negative	incorrect	
negative	positive	incorrect	Denaturation temperature too high → uncritical, false PCR-results are not to be expected

**4 Further Information**

**4.1 Product Information**

- Detailed information about setup of several real-time PCR devices (Download: [www.congen.de/en/company/downloads](http://www.congen.de/en/company/downloads))
- Product-related documents (Download: [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/)) (Download:

**4.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).

**4.3 Explanation of Symbols**

Lot number

Article number

Expiry date

Number of Tests

Storage temperature

Manufacturer