CONGEN

SureFast[®] VEGAN

Art. No. F4055
100 rxn

User Manual

April 2023

 **Inhalt**

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Content

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	DNA-Preparation	9
1.4	Kit components and storage	10
1.5	Additionally required equipment and materials	10
1.6	Setup.....	11
1.7	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Preparation of the master-mix	12
2.1.2	Preparation of the real-time PCR-mix	12
2.2	Interpretation of results	13
3	Further Information	13
3.1	Product Information.....	13
3.2	Technical Support	13

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1 General Information

1.1 Description

The SureFast® VEGAN real time PCR is a highly sensitive system for the detection of traces of vertebrate-DNA and differentiation between DNA of vertebrates and plants in food. The Positive Control corresponds to 0.1% vertebrate-DNA.

Insects, arachnids, molluscs and crustaceans are not detected by this system.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of fluorescence emissions at the channels **FAM** and **Cy5** at the same time. The technical verification of instruments was performed on Roche LightCycler® 480 II, Qiagen Rotor-Gene Q, Bio-Rad CFX96Dx, Bio-Rad CFX96 Opus, R-Biopharm RIDA®CYCLER, Agilent AriaDx.

1.2 Limit of Detection

The SureFast® VEGAN real-time PCR has a limit of detection of 0.01%.

The limit of detection of the complete method (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA-preparation and DNA-content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 DNA-Preparation

For DNA-preparation the use of SureFood® PREP Basic (Art. No. S1052), SureFast® Mag PREP Food (Art. No. F1060) and for highly processed food and feed the use of SureFood® PREP Advanced (Art. No. S1053) is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at –20°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA-Extraction kit
(e.g. SureFood® PREP Basic Art. No. S1052 / SureFood® PREP Advanced Art. No. S1053)/
SureFast® Mag PREP Food (Art. No. F1060)
- real-time PCR instrument with **four** detection channels (**510 nm**, 580 nm, 610 nm and **660 nm**)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler
Initial Denaturation (HOLD) Cycles	5 min, 95°C 30	1 min, 95°C 30
Denaturation Annealing/Extension (CYCLE)	15 sec, 95°C 30 sec, 60°C	10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.7 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent AriaMx / Dx	Plant	FAM	+	
	Vertebrates	Cy5	+	
Bio-Rad CFX96 / Dx/ Opus	Plant	FAM	+	
	Vertebrates	Cy5	+	
R-Biopharm RIDA®CYCLER	Plant	green	+	
	Vertebrates	red	+	
Qiagen Rotor-Gene Q	Plant	green	+	Note: Please use only 0.1 ml reaction tube. The gain settings must be set to 5 (factory default) for all channels.
	Vertebrates	red	+	
Roche LightCycler® 480 II	Plant	465-510	+	
	Vertebrates	618-660	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay as well as for the inhibition control.

Recommended control reactions for the specific PCR assay: negative control, negative extraction control, positive control.

Reactions needed for the qualitative detection of plant and vertebrates:

3 reactions for controls (1x no-template control, 1x extraction control, 1 positive control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

With this sensitive real-time PCR kit it is possible to detect minimal contaminations by the target DNA (for example airborne pollen, human DNA from the operator). These can lead to positive signals of the negative controls.

Plant DNA is detected in the FAM-channel and vertebrate DNA is detected in the Cy5-channel.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel.

If the sample DNA is detected before the signal of the Positive Control in the Cy5 channel, the concentration of the sample is over 0.1 % vertebrate DNA. If the sample-DNA in the Cy5 channel is detected after the signal of the Positive Control, the concentration of the sample is below 0.1 % vertebrate DNA.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel.

result in the respective channel		Interpretation
FAM channel plant	Cy5 channel vertebrates	
positive	negative	plant DNA detected
negative	positive	vertebrate DNA detected
positive	positive	plant and vertebrate DNA detected
negative	negative	invalid

3 Further Information

3.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/company/downloads)
- Verification Report upon request

3.2 Technical Support

For further questions please send an e-mail to info@congen.de.