
CONGEN

SureFast® Vibrio Virulence 4plex

Art. No. F5181
100 rxn

User Manual

July 2026

Inhalt

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Content

1	General Information	10
1.1	Description	10
1.2	Limit of Detection	10
1.3	DNA-preparation	11
1.4	Kit components and storage	11
1.5	Additionally required equipment and materials	11
1.6	Precautions for users	11
1.7	Setup.....	12
1.8	Detection channel Set-up	12
2	Qualitative Analysis	13
2.1	Protocol	13
2.1.1	Preparation of the master-mix	13
2.1.2	Preparation of the real-time PCR-mix	14
2.2	Interpretation of results	14
3	Limitations of the method	16
4	Further Information	16
4.1	Product Information	16
4.2	Technical Support	16

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1 General Information

1.1 Description

SureFast® Vibrio Virulence 4plex is a real-time PCR assay for the direct qualitative detection and differentiation of specific DNA sequences of the virulence genes *rtxA*, *vvh*, *tdh*, *toxR*, *trh* and *rtxA* in enriched food.

Detection of the cytotoxin gene *vvh* is performed as a marker for the presence of potentially pathogenic *Vibrio vulnificus*. The genes *tdh*, *toxR*, and *trh* are detected as relevant virulence markers associated with gastrointestinal disease. Their detection serves as an indication of the presence of *Vibrio parahaemolyticus*. Detection of the virulence factor *rtxA* is performed as a marker for the presence of *Vibrio cholerae*.

The selection of virulence genes is based on the requirements of DIN EN ISO 21872-1.

Each reaction contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The internal technical verification of instruments was performed on Roche LightCycler® 480 II, Bio-Rad CFX96 Dx, Bio-Rad CFX96 Opus, R-Biopharm RIDA®CYCLER, Qiagen Rotor-Gene Q and Agilent AriaDx.

1.2 Limit of Detection

The SureFast® Vibrio Virulence 4plex real-time PCR has a limit of detection of ≤ 5 DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP Bacteria (Art. No. F1021) is recommended.

To assess the process of bacterial growth, it is suggested to compare the samples at the beginning and at the end of the culturing (bacterial growth at Cp difference > 3).

For cultivation the method of the DIN EN ISO 21872-1 is recommended.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.

Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.

1.5 Additionally required equipment and materials

- DNA-Extraction kit (e.g. SureFast® PREP Bacteria Art. No. F1021)
- real- time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- Incubator
- Enrichment broth
- microcentrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples and enrichment cultures must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled and disposed according to the national safety regulations.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

SureFast® Vibrio Virulence 4plex (100 rxn)

Art. No. F5181

July 2026

1.7 Setup

	Blockcycler & R-Biopharm RIDA®CYCLER	Rotorcycler & LightCycler® 480 II
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel Set-up

For additional details see cycler operating instructions for real-time PCR device.

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent AriaDx /Mx	<i>V. parahaemolyticus</i> <i>tdh</i> , <i>toxR</i> , <i>trh</i>	FAM	+	
	IAC	HEX	+	
	<i>V. vulnificus</i> <i>vvh</i>	ROX	+	
	<i>V. cholerae</i> <i>rtxA</i>	Cy5	+	
Bio-Rad CFX96/Dx/Opus	<i>V. parahaemolyticus</i> <i>tdh</i> , <i>toxR</i> , <i>trh</i>	FAM	+	Baseline Settings: • Baseline subtracted curve fit • Apply fluorescence drift correction
	IAC	VIC/HEX	+	
	<i>V. vulnificus</i> <i>vvh</i>	ROX	+	
	<i>V. cholerae</i> <i>rtxA</i>	Cy5	+	
R-Biopharm RIDA®CYCLER	<i>V. parahaemolyticus</i> <i>tdh</i> , <i>toxR</i> , <i>trh</i>	green	+	Ignore cycles before , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter.
	IAC	yellow	+	
	<i>V. vulnificus</i> <i>vvh</i>	orange	+	
	<i>V. cholerae</i> <i>rtxA</i>	red	+	
Qiagen Rotor- Gene Q	<i>V. parahaemolyticus</i> <i>tdh</i> , <i>toxR</i> , <i>trh</i>	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	IAC	yellow	+	
	<i>V. vulnificus</i> <i>vvh</i>	orange	+	
	<i>V. cholerae</i> <i>rtxA</i>	red	+	
Roche LightCycler® 480 II	<i>V. parahaemolyticus</i> <i>tdh</i> , <i>toxR</i> , <i>trh</i>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	533-580	+	
	<i>V. vulnificus</i> <i>vvh</i>	533-610	+	
	<i>V. cholerae</i> <i>rtxA</i>	618-660	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control and Positive control.

For the analysis of enrichments additional controls are recommended: zero control and medium control.

The reaction mix contains an internal amplification control (Inhibition control) per reaction.

Reactions needed for the qualitative *V. parahaemolyticus* *tdh*, *toxR*, *trh*, *V. vulnificus* *vvh* and *V. cholerae* *rtxA* detection:

3 reactions for controls* (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

Reactions needed for the qualitative *V. parahaemolyticus* *tdh*, *toxR*, *trh*, *V. vulnificus* *vvh* and *V. cholerae* *rtxA* detection:

5 reactions for controls* (1x negative control, 1x extraction control, 1x Positive Control, 1x zero control, 1x medium control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

* Description of the controls

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – components from used Prep Kit
- Positive Control: master-mix and within the kit's provided Positive Control
- Zero control: master-mix and sample prior enrichment
- Medium control: only medium – no sample

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

In the FAM channel, the *Vibrio parahaemolyticus* virulence factors *tdh*, *toxR*, and *trh* are detected; in the ROX channel, the *Vibrio vulnificus* virulence factor *vvh* is detected, and in the Cy5 channel, the *Vibrio cholerae* virulence factor *rtxA* is detected (see table). An internal amplification control (IAC) is detected in the VIC/HEX channel.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the internal control (VIC/HEX-channel) of the sample is **positive**.

If the sample DNA shows no amplification or sigmoidal amplification curve in the VIC/HEX channel, this indicates that inhibitors are present in the sample DNA that inhibit the PCR. Under these circumstances DNA isolation and purification of the sample need to be improved. It is recommended to repeat the extraction process. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific PCR assay.

Result in the respective channel				Interpretation
FAM channel <i>V. parahaemolyticus</i> <i>tdh, toxR, trh</i>	ROX channel <i>V. vulnificus</i> <i>vvh</i>	Cy5 channel <i>V. cholerae</i> <i>rtxA</i>	VIC/HEX channel IAC	
positive	negative	negative	positive/negative	<i>V. parahaemolyticus</i> <i>tdh, toxR, trh</i> DNA detected
negative	positive	negative	positive/negative	<i>V. vulnificus</i> <i>vvh</i> DNA detected
negative	negative	positive	positive/negative	<i>V. cholerae</i> <i>rtxA</i> DNA detected
negative	negative	negative	positive	Negative, target DNA is not detected
negative	negative	negative	negative	invalid

Note: The results displayed in the table above represent merely an example. Additional combinations are also possible.

The following table shows the specification ranges of the kit controls

	Specification range	Internal amplification control (IAC)
negative control	negative	$28 \leq C_p \leq 36$
Positive Control (FAM – <i>V. parahaemolyticus</i> <i>tdh, toxR, trh</i>)	$25 \leq C_p \leq 33$	not relevant
Positive Control (FAM – <i>V. vulnificus</i> <i>vvh</i>)	$25 \leq C_p \leq 33$	not relevant
Positive Control (Cy5 – <i>V. cholerae</i> <i>rtxA</i>)	$25 \leq C_p \leq 33$	not relevant

3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.
- A positive test result does not necessarily indicate the presence of viable organism. It is indicative for the presence of the target DNA.
- The *in silico* analysis of the *toxR* system of *Vibrio vulnificus* revealed occasional hits for *Vibrio alginolyticus*, as approximately 1 % of the examined *V. alginolyticus* strains may also contain this gene. This indicates that PCR amplification of *V. alginolyticus* is theoretically possible. However, due to the very low number of hits, a relevant or frequent cross-reaction is considered unlikely based on the current data. Occasional amplification cannot be completely ruled out.

4 Further Information

4.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

4.2 Technical Support

For further questions please send an e-mail to info@congen.de.