
**CONGEN**

# **SureFast® Parasitic Water Panel 4plex**

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## **Content**

1	General Information .....	10
1.1	Description .....	10
1.2	Limit of Detection .....	10
1.3	DNA-preparation .....	11
1.4	Kit components and storage .....	11
1.5	Additionally required equipment and materials .....	11
1.6	Precautions for users .....	11
1.7	Setup.....	12
1.8	Detection channel Set-up .....	12
2	Qualitative Analysis .....	13
2.1	Protocol .....	13
2.1.1	Preparation of the master-mix .....	13
2.1.2	Preparation of the real-time PCR-mix .....	13
2.2	Interpretation of results .....	14
3	Limitations of the method .....	15
4	Further Information .....	16
4.1	Product Information .....	16
4.2	Technical Support .....	16

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## 1 General Information

### 1.1 Description

The SureFast® Parasitic Water Panel 4plex is a real-time PCR for the direct, qualitative detection and differentiation of specific DNA sequences of *Giardia intestinalis*, *Entamoeba histolytica* and *Cryptosporidium* spp. in water.

The test contains an internal amplification control (IAC). If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are algae, alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride).

The real-time PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of four fluorescence emissions at the channels FAM, VIC/HEX, ROX and Cy5 at the same time. The internal technical verification of instruments was performed on Agilent Mx3005P, Bio-Rad CFX96, Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Applied Biosystems 7500 and R-Biopharm RIDA®CYCLER.

### 1.2 Limit of Detection

The SureFast® Parasitic Water Panel 4plex real-time PCR has a limit of detection of  $\leq 5$  DNA copies.

The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target DNA is sufficient for a successful analysis. The concentration of total DNA in the sample does not allow a conclusion on the quantity and quality of the target DNA.

**Note:** Inconsistent mixing ratios may cause a loss of sensitivity in the low concentration channel in mixed samples\* especially with high amplicon concentrations (Cp value < 20).

\* e.g. 99.9 % *Giardia intestinalis* and 0.1 % *Cryptosporidium* spp.

### 1.3 DNA-preparation

For DNA-preparation the use of SureFast® PREP DNA/RNA Virus (Art. No. F1051) or SureFast® PREP Aqua Art. (No. F1023) is recommended.

### 1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 1050 µl	Yellow
2	Taq Polymerase	1 x 80 µl	Dark Red
3	Positive Control	1 x 200 µl	Light Blue

**Store all reagents at -28 to -16°C and protected from light. The Taq Polymerase can be stored at +2 to +8°C for multiple uses on the same day.**

**Note: The Taq Polymerase may be in a frozen or unfrozen state. This does not affect the quality of the Taq Polymerase or the performance of the real-time PCR.**

### 1.5 Additionally required equipment and materials

- DNA-Extraction kit  
(e.g. SureFast® PREP Aqua Art. No. F1023 / SureFast® PREP DNA/RNA Virus Art. No. F1051)
- real-time PCR instrument with four detection channels (510 nm, 580 nm, 610 nm and 660 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- microcentrifuge with a rotor for the reaction tubes

### 1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Samples shall be considered as potentially infectious and must be handled and disposed of in accordance with applicable regulations for potentially infectious materials. All reagents and materials that come into contact with such samples shall be treated as biohazardous waste and disposed of accordingly.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

# SureFast® Parasitic Water Panel 4plex (100 rxn)

Art. No. F5506

April 2026

## 1.7 Setup

	Blockcyler & R-Biopharm RIDA®CYCLER	Rotorcyler & LightCycler® 480 II
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

## 1.8 Detection channel Set-up

For additional details see cyler operating instructions for real-time PCR device.

Real-time PCR device	Detection	Detection channel	Quencher	Note
<b>Agilent Mx3005P</b>	<i>Giardia intestinalis</i>	FAM	+	
	IAC	HEX	+	
	<i>Entamoeba histolytica</i>	ROX	+	
	<i>Cryptosporidium</i> spp.	Cy5	+	
<b>Applied Biosystems 7500</b>	<i>Giardia intestinalis</i>	FAM	None	Check the passive reference option ROX is none.
	IAC	VIC	None	
	<i>Entamoeba histolytica</i>	ROX	None	
	<i>Cryptosporidium</i> spp.	Cy5	None	
<b>Bio-Rad CFX96/Dx/Opus</b>	<i>Giardia intestinalis</i>	FAM	+	Baseline Settings: <ul style="list-style-type: none"> <li>Baseline subtracted curve fit</li> <li>Apply fluorescence drift correction</li> </ul>
	IAC	VIC/HEX	+	
	<i>Entamoeba histolytica</i>	ROX	+	
	<i>Cryptosporidium</i> spp.	Cy5	+	
<b>R-Biopharm RIDA®CYCLER</b>	<i>Giardia intestinalis</i>	green	+	<b>Ignore cycles before</b> , if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cyler operating instructions, section 12.1.2 Cycling analysis parameter.
	IAC	yellow	+	
	<i>Entamoeba histolytica</i>	orange	+	
	<i>Cryptosporidium</i> spp.	red	+	
<b>Roche LightCycler® 480 II</b>	<i>Giardia intestinalis</i>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	533-580	+	
	<i>Entamoeba histolytica</i>	533-610	+	
	<i>Cryptosporidium</i> spp.	618-660	+	
<b>Roche cobas® z 480 Analyzer</b>	<i>Giardia intestinalis</i>	465-510	+	The SureCC Color Compensation Kit I (Art. No. F4009) is required.
	IAC	540-580	+	
	<i>Entamoeba histolytica</i>	540-610	+	
	<i>Cryptosporidium</i> spp.	610-670	+	

## 2 Qualitative Analysis

### 2.1 Protocol

#### 2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, extraction control and Positive Control.

The reaction mix contains an internal amplification control (Inhibition control) per reaction.

**Reactions needed for the qualitative *Giardia intestinalis*, *Entamoeba histolytica* and *Cryptosporidium* spp. detection:**

3 reactions for controls\* (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

**Example for the calculation and preparation of 10 reactions:**

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19.3 µl	212.3 µl
Taq Polymerase	0.7 µl	7.7 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

**Mix each master-mix well and centrifuge shortly before use.**

#### \* Description of the controls

- Negative control: only master-mix
- Extraction control: the extraction is performed without the sample – use components from used Prep Kit
- Positive Control: master-mix and within the kit's provided Positive Control

## 2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control.
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

## 2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

*Giardia intestinalis* DNA is detected in the FAM-channel, *Entamoeba histolytica* DNA is detected in the ROX-channel and *Cryptosporidium* spp. DNA is detected in the Cy5-channel (see table). In the VIC/HEX-channel the amplification control is detected.

A sample is stated **positive** for the respective parameter, if the sample DNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control (IAC).

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative** for the respective parameter, if the sample DNA shows no amplification in the respective channel and if the internal control (VIC/HEX-channel) of the sample is **positive**.

If the sample DNA shows no amplification or sigmoidal amplification curve in the VIC/HEX channel, this indicates that inhibitors are present in the sample DNA that inhibit the PCR. Under these circumstances DNA isolation and purification of the sample need to be improved. It is recommended to repeat the extraction process. Alternatively the DNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific PCR assay.

Result in the respective channel				Interpretation
FAM channel <i>Giardia intestinalis</i>	ROX channel <i>Entamoeba histolytica</i>	Cy5 channel <i>Cryptosporidium</i> spp.	VIC/HEX channel IAC	
<b>positive</b>	negative	negative	<b>positive/negative</b>	<i>Giardia intestinalis</i> DNA detected
negative	<b>positive</b>	negative	<b>positive/negative</b>	<i>Entamoeba histolytic</i> DNA detected
negative	negative	<b>positive</b>	<b>positive/negative</b>	<i>Cryptosporidium</i> spp. DNA detected
negative	negative	negative	<b>positive</b>	Negative, target DNA is not detected
negative	negative	negative	negative	invalid

**Note: The results displayed in the table above represent merely an example. Additional combinations are also possible.**

The following table shows the specification ranges of the kit controls

	Specification range	Internal amplification control (IAC)
negative control	negative	$25 \leq Cp \leq 34$
Positive Control (FAM – <i>Giardia intestinalis</i> )	$25 \leq Cp \leq 33$	not relevant
Positive Control (ROX – <i>Entamoeba histolytica</i> )	$25 \leq Cp \leq 33$	not relevant
Positive Control (Cy5 – <i>Cryptosporidium</i> spp.)	$25 \leq Cp \leq 33$	not relevant

### **3 Limitations of the method**

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.
- A positive test result does not necessary indicate the presence of viable organism. It is indicative for the presence of the target DNA (*Giardia intestinalis* DNA).

### **4 Further Information**

#### **4.1 Product Information**

- Detailed information about setup of several real-time PCR devices and additional information for sample preparations (Download: [www.congen.de/en/downloads](http://www.congen.de/en/downloads))
- Product-related documents (Download: [www.congen.de/en/eifu/](http://www.congen.de/en/eifu/))
- Validation Report upon request

#### **4.2 Technical Support**

For further questions please send an e-mail to [info@congen.de](mailto:info@congen.de).