CONGEN

SureFast® Norovirus/Hepatitis A 3plex

Art. No. F7124
100 rxn

User Manual

August 2025

 Inhalt

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Content

1	General Information	10
1.1	Description	10
1.2	Limit of Detection	10
1.3	RNA Preparation	11
1.4	Kit components and storage	11
1.5	Additionally required equipment and materials	11
1.6	Precautions for users	11
1.7	Setup.....	12
1.8	Detection channel Set-up	12
2	Qualitative Analysis	13
2.1	Protocol	13
2.1.1	Preparation of the master-mix	13
2.1.2	Preparation of the real-time PCR-mix	14
2.2	Interpretation of results	14
3	Limitations of the method	16
4	Further Information	16
4.1	Product Information	16
4.2	Technical Support	16

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1 General Information

1.1 Description

The SureFast® Norovirus/Hepatitis A 3plex is a multiplex real-time RT-PCR for the direct, qualitative detection of specific RNA sequences of norovirus (genogroup I and II) and hepatitis A in food.

Each reaction contains an Internal Control RNA (ICR, consisting of MS2-bacteriophage) as an internal control of sample preparation procedure and to monitor possible PCR-inhibition.

If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time RT-PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of three fluorescence emissions at the channels FAM, VIC/HEX and Cy5 at the same time. The internal technical verification of instruments was performed on Roche LightCycler® 480 II, Qiagen Rotor-Gene Q, R-Biopharm RIDA®CYCLER, Bio-Rad CFX96 Dx, Bio-Rad CFX96 OPUS and Agilent Aria Dx.

1.2 Limit of Detection

The SureFast® Norovirus/Hepatitis A 3plex real-time RT-PCR has a limit of detection of ≤ 25 RNA copies.

The assay limit of detection depends on sample matrix, processing grade, RNA preparation and RNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target RNA is sufficient for a successful analysis. The concentration of total RNA in the sample does not allow a conclusion on the quantity and quality of the target RNA.

1.3 RNA Preparation

For RNA-preparation the use of SureFast® PREP DNA/RNA Virus (Art. No. F1051) is recommended.

The test assay contains an Internal Control RNA (ICR), which can either be used as internal amplification control or as positive extraction control for the sample preparation procedure and as a PCR inhibition control.

If the ICR is used only as a PCR internal amplification control (IAC), 1 µl per reaction of the ICR should be added to the master-mix (see point 2.1.1).

If the ICR is used as an extraction control for the sample preparation procedure and as inhibition control, 20 µl of the ICR should be added during extraction procedure. The ICR should always be pipetted to the sample with added Lysis Buffer and must not be added directly to the raw sample.

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 700 µl	Yellow
2	PP-Mix	1 x 770 µl	Light Green
3	Enzyme Mix	1 x 80 µl	Red
R	Internal Control RNA	2 x 1700 µl	Brown
N	PCR Water	1 x 450 µl	White
P	Positive Control	1 x 90 µl	Light Blue

Store all reagents at -28 to -16°C and protected from light.

1.5 Additionally required equipment and materials

- RNA-Extraction kit (e.g. SureFast® PREP DNA/RNA Virus Art. No. F1051)
- real- time PCR instrument with three detection channels (510 nm, 580 nm and 660 nm)
- real-time PCR consumable (plates, tubes, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.
- All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulation for disposal.

1.7 Setup

	Blockcycler & Rotor-Gene & R-Biopharm RIDA®CYCLER	Rotorcycler & LightCycler® 480 II
Reverse Transcription	10 min, 58 °C	10 min, 58 °C
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent AriaDx /Mx	norovirus	FAM	+	
	ICR	HEX	+	
	hepatitis A	Cy5	+	
Bio-Rad CFX96/Dx/Opus	norovirus	FAM	+	Baseline Settings: - Baseline subtracted curve fit - Apply fluorescence drift correction
	ICR	VIC/HEX	+	
	hepatitis A	Cy5	+	
R-Biopharm RIDA®CYCLER	norovirus	green	+	Ignore cycles before, if there is a significant deviation in the baseline at the start of the run. Please see page 45 of the cycler operating instructions, section 12.1.2 Cycling analysis parameter.
	ICR	yellow	+	
	hepatitis A	red	+	
Qiagen Rotor-Gene Q	norovirus	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	ICR	yellow	+	
	hepatitis A	red	+	
Roche LightCycler® 480 II	norovirus	465-510	+	
	ICR	533-580	+	
	hepatitis A	618-660	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay.

Recommended control reactions for the specific PCR assay: negative control, negative extraction control and Positive Control.

The test assay contains an Internal Control RNA (ICR), which can either be used as internal amplification control or as positive extraction control for the sample preparation procedure and as a PCR inhibition control.

Reactions needed for the qualitative norovirus and hepatitis A detection:

3 reactions for controls* (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample RNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	12.5 µl	137.5 µl
PP-Mix	6.9 µl	75.9 µl
Enzyme Mix	0.7 µl	7.7 µl
Total volume	20.1 µl	221.1 µl

Mix each master-mix well and centrifuge shortly before use.

Example for the calculation and preparation of 10 reactions for ICR only as PCR inhibition control:

Components of the master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	12.5 µl	137.5 µl
PP-Mix	6.9 µl	75.9 µl
ICR	1.0 µl	11.0 µl
Enzyme Mix	0.7 µl	7.7 µl
Total volume	21.1 µl	232.1 µl

Mix each master-mix well and centrifuge shortly before use.

*** Description of the controls**

- Negative control: master-mix and 5 µl PCR Water (Kite Code N)
- Extraction control: the extraction is performed without the sample – components from used Prep Kit
- Positive Control: master-mix and within the kit's provided Positive Control

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Pipette 5 µl of PCR Water in the tubes/wells and close the negative control
- Pipette 5 µl of sample RNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Norovirus RNA is detected in the FAM-channel and hepatitis A RNA is detected in the Cy5-channel (see table). In the VIC/HEX-channel the internal amplification (IAC) or positive extraction control (ICR) is detected.

A sample is stated **positive** for the respective parameter, if the sample RNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the IAC.

A Cp value for the IAC is not needed to obtain a positive result of the Positive Control.

A sample is stated **negative** for the respective parameter, if the sample RNA shows no amplification in the respective channel and if the amplification or positive extraction control (VIC/HEX-channel) of the sample is **positive**.

If the sample RNA in the VIC/HEX-Channel shows **no amplification** or an irregular amplification curve, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances RNA isolation and purification of the sample need to be improved. Alternatively the RNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific norovirus or hepatitis A PCR assay.

Result in the respective channel			Interpretation
FAM channel norovirus	Cy5 channel hepatitis A	VIC/HEX channel ICR	
positive	negative	positive/negative	norovirus (GG I/GG II) RNA detected
negative	positive	positive/negative	hepatitis A RNA detected
positive	positive	positive/negative	norovirus (GG I/GG II) and hepatitis A RNA detected
negative	negative	positive	Negative, no target RNA detected
negative	negative	negative	invalid

The following table shows the specification ranges of the kit controls

	Specification range	Internal amplification control (IAC)
negative control	negative	$25 \leq C_p \leq 34$
Positive Control (FAM –norovirus)	$22 \leq C_p \leq 30$	not relevant
Positive Control (Cy5 – hepatitis A)	$22 \leq C_p \leq 30$	not relevant

3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- In highly processed samples, the limit of detection may be shifted. Factors such as high pressures, mechanical stresses, chemical treatment, extreme temperatures and/or extreme pH values during manufacturing process – such as in canning production – can damage or degrade nucleic acids. This means that the sensitivity of the test kit may be reduced and not all original components may be detected.
- A positive test result does not necessarily indicate the presence of a viable organism. It is indicative for the presence of the target RNA (norovirus RNA and hepatitis A RNA).

4 Further Information

4.1 Product Information

- Detailed information about setup of several real-time PCR devices (Download: www.congen.de/en/downloads)
- Product-related documents (Download: www.congen.de/en/eifu/)
- Validation Report upon request

4.2 Technical Support

For further questions please send an e-mail to info@congen.de.