CONGEN

SureFast® **Norovirus/Hepatitis A 3plex**

Art. No. F7124
100 rxn

User Manual

June 2024

 **Inhalt**

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Content

1	General Information	9
1.1	Description	9
1.2	Limit of Detection	9
1.3	RNA Preparation	10
1.4	Kit components and storage	10
1.5	Additionally required equipment and materials	10
1.6	Precautions for users	10
1.7	Setup.....	11
1.8	Detection channel Set-up	11
2	Qualitative Analysis	12
2.1	Protocol	12
2.1.1	Preparation of the master-mix	12
2.1.2	Preparation of the real-time PCR-mix	13
2.2	Interpretation of results	13
3	Limitations of the method	14
4	Further Information	14
4.1	Product Information	14
4.2	Technical Support	14

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1 General Information

1.1 Description

The SureFast® Norovirus/Hepatitis A 3plex is a multiplex real-time RT-PCR for the direct, qualitative detection of specific RNA sequences of norovirus (genogroup I and II) and hepatitis A in food.

The test assay contains an external Internal Control RNA (ICR) as an internal control of sample preparation procedure and/or to determine possible PCR inhibition.

If the DNA contains PCR inhibiting substances, the signal of the amplification control will be affected or the amplification will be suppressed. Examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). In addition spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

The real-time RT-PCR assay can be performed with commonly used real-time PCR instruments, equipped for detection of three fluorescence emissions at the channels FAM, VIC/HEX and Cy5 at the same time. The technical verification of instruments was performed on Roche LightCycler® 480 II, Roche cobas® z 480 Analyzer, Qiagen Rotor-Gene Q, Applied Biosystems 7500, Bio-Rad CFX96, Bio Molecular Systems MIC, Agilent AriaDx and Agilent Mx3005P.

1.2 Limit of Detection

The SureFast® Norovirus/Hepatitis A 3plex real-time RT-PCR has a limit of detection of ≤ 25 RNA copies.

The assay limit of detection depends on sample matrix, processing grade, RNA preparation and RNA content.

The SureFast® PCR systems are very sensitive and therefore even a small amount of target RNA is sufficient for a successful analysis. The concentration of total RNA in the sample does not allow a conclusion on the quantity and quality of the target RNA.

1.3 RNA Preparation

For RNA-preparation the use of SureFast® PREP DNA/RNA Virus (Art. No. F1051) is recommended. The test assay contains an Internal Control RNA (ICR), which can either be used as internal amplification control or as positive extraction control for the sample preparation procedure and as a PCR inhibition control. If the ICR is used only as a PCR internal amplification control (IAC), 1 µl per reaction of the ICR should be added to the master-mix (see point 2.1.1).

If the ICR is used as an extraction control for the sample preparation procedure and as inhibition control, 20 µl of the ICR should be added during extraction procedure. The ICR should always be pipetted to the sample with added Lysis Buffer and must not be added directly to the raw sample

1.4 Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	2 x 700 µl	Yellow
2	PP-Mix	1 x 770 µl	Light Green
3	Enzyme Mix	1 x 80 µl	Red
R	Internal Control RNA	2 x 1700 µl	Brown
N	PCR Water	1 x 450 µl	White
P	Positive Control	1 x 90 µl	Light Blue

Store all reagents at -28 to -16°C and protected from light.

1.5 Additionally required equipment and materials

- RNA-Extraction kit (e.g. SureFast® PREP DNA/RNA Virus Art. No. F1051)
- real-time PCR instrument with three detection channels (510 nm, 580 nm and 660 nm)
- real-time PCR consumable (plates, tubes, capillaries, foils, caps)
- pipettes with filter tips
- powder-free disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

1.6 Precautions for users

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Food samples must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

1.7 Setup

	Blockcycler & Rotor-Gene & R-Biopharm RIDA®CYCLER	LightCycler® 480 II
Reverse Transcription	10 min, 58 °C	10 min, 58 °C
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum

1.8 Detection channel Set-up

Real-time PCR device	Detection	Detection channel	Quencher	Note
Agilent Mx3005P	norovirus	FAM	+	
	ICR	HEX	+	
	hepatitis A	Cy5	+	
Agilent AriaDx /Mx	norovirus	FAM	+	
	ICR	HEX	+	
	hepatitis A	Cy5	+	
Applied Biosystems 7500	norovirus	FAM	None	Check the passive reference option ROX is none.
	ICR	VIC	None	
	hepatitis A	Cy5	None	
Bio-Rad CFX96/Dx/Opus	norovirus	FAM	+	
	ICR	VIC/HEX	+	
	hepatitis A	Cy5	+	
R-Biopharm RIDA®CYCLER	norovirus	green	+	
	ICR	yellow	+	
	hepatitis A	red	+	
Qiagen Rotor-Gene Q	norovirus	green	+	Note: Please use only 0.1 ml reaction tubes. The gain settings must be set to 5 (factory default) for all channels.
	ICR	yellow	+	
	hepatitis A	red	+	
Roche LightCycler® 480 II	norovirus	465-510	+	
	ICR	533-580	+	
	hepatitis A	618-660	+	
Roche cobas® z 480 Analyzer	norovirus	465-510	+	
	ICR	540-580	+	
	hepatitis A	610-670	+	

2 Qualitative Analysis

2.1 Protocol

2.1.1 Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions) for the specific PCR assay. Recommended control reactions for the specific PCR assay: negative control, negative extraction control, Positive Control.

The test assay contains an Internal Control RNA (ICR), which can either be used as internal amplification control or as positive extraction control for the sample preparation procedure and as a PCR inhibition control. We also recommend to add 1 µl of the ICD to the negative control and positive control PCR Mix.

Reactions needed for the qualitative norovirus and hepatitis A detection:

3 reactions for controls (1x negative control, 1x extraction control, 1x Positive Control)

For each sample: at least 1 reaction for each sample RNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	12.5 µl	137.5 µl
PP-Mix	6.9 µl	75.9 µl
Enzyme Mix	0.7 µl	7.7 µl
Total volume	20.1 µl	221.1 µl

Mix each master-mix well and centrifuge shortly before use.

Example for the calculation and preparation of 10 reactions for ICR only as PCR inhibition control:

Components of the master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	12.5 µl	137.5 µl
PP-Mix	6.9 µl	75.9 µl
ICR	1.0 µl	11.0 µl
Enzyme Mix	0.7 µl	7.7 µl
Total volume	21.1 µl	232.1 µl

Mix each master-mix well and centrifuge shortly before use.

2.1.2 Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample RNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

2.2 Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions have to show the correct results.

Norovirus RNA is detected in the FAM-channel and hepatitis A RNA is detected in the Cy5-channel (see table). In the VIC/HEX-channel the internal amplification (IAC) or positive extraction control (ICR) is detected.

A sample is stated **positive** for the respective parameter, if the sample RNA shows amplification in the respective channel. High amplicon concentrations can result in a weak or absent signal of the internal amplification control.

A Cp value for the internal amplification control is not needed to obtain a positive result of the positive control.

A sample is stated **negative** for the respective parameter, if the sample RNA shows no amplification in the respective channel and if the amplification or positive extraction control (VIC/HEX-channel) of the sample is **positive**.

If the sample RNA in the VIC/HEX-Channel shows no amplification or an irregular amplification curve, it contains PCR inhibiting substances. A significant decrease in the fluorescence signal can also show the presence of PCR inhibiting substances. Under these circumstances RNA isolation and purification of the sample need to be improved. Alternatively the RNA can be diluted (recommendation 1:2 in PCR-water) and analysed again for inhibition. Please note that the dilution factor also affects the detection limit of the specific norovirus or hepatitis A PCR assay.

Result in the respective channel			Interpretation
FAM channel norovirus	Cy5 channel hepatitis A	VIC/HEX channel ICR	
positive	negative	positive/negative	norovirus (GG I/GG II) RNA detected
negative	positive	positive/negative	hepatitis A RNA detected
positive	positive	positive/negative	norovirus (GG I/GG II) and hepatitis A RNA detected
negative	negative	positive	Negative, no target RNA detected
negative	negative	negative	invalid

3 Limitations of the method

- The presence of PCR inhibitors may cause invalid results.
- Extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
- A positive test result does not necessarily indicate the presence of a viable organism. It is indicative for the presence of the target DNA (norovirus RNA and hepatitis A RNA).

4 Further Information

4.1 Product Information

- Detailed information about setup of several real-time PCR devices
(Download: www.congen.de/en/company/downloads)
- Validation Report upon request

4.2 Technical Support

For further questions please send an e-mail to info@congen.de.