

## Flow Chart for SureFast® PREP Bacteria

Art. No. F1021

June 2023

**(1) Preparation of the basic material**

Transfer 1.0 ml of the enrichment into a 1.5 ml reaction tube  
 Centrifuge for 5 min at 12.000 rpm.  
 Discard the fluid

**(2) Lysis of basic material**




Add **400 µl Lysis Buffer** (Code L)  
 Mixing  
 Incubation **10 min at 99°C**.  
 Centrifugation **1 min 12000 rpm**

**Setting of optimal binding conditions**



Transfer **ca. 300 µL supernatant** into a new **1.5 ml reaction tube**

**(3) Binding of nucleic acids on a spin filter**



Add **200 µl Binding Buffer** (Code B) to sample  
 Mixing  
 Place a **Spin Filter** (Code S) into a new **clear Receiver Tube** (Code R)  
 Transfer the **complete solution** onto the Spin Filter (Code S)  
 Incubation **1 min at room temperature**  
 Centrifugation **1 min 12000 rpm**  
 Discard the filtrate

## Flow Chart for SureFast® PREP Bacteria

Art. No. F1021

June 2023

Information

(4) **Purification of bound nucleic acids &**

(5) **Drying of the Spin Filter**



Place the **clear Spin Filter** back into the **clear Receiver Tube** (Code R)

Add **550 µL Wash Buffer** (Code W)

Centrifugation **1 min 12000 rpm**

Discard filtrate  
Place Spin Filter back into the Receiver Tube (Code R)



Add **550 µL Wash Buffer** (Code W)


Centrifugation **1 min 12000 rpm**

Discard filtrate  
Place Spin Filter back into the Receiver Tube (Code R)



Centrifugation **2 min 12000 rpm**

(6) **Elution of nucleic acids**



Place the **Spin Filter** into a **clear 1.5 ml Receiver Tube** (Code T)

Add **100 µL preheated Elution Buffer** (Code E)

Incubation **3 min room temperature**

Centrifugation **1 min 10000 rpm**

Discard Spin Filter

The eluted DNA is ready-to-use for the PCR