

GENTAMICIN ELISA

5111GEN[18]02.24

A competitive enzyme immunoassay for
screening and quantitative analysis of
Gentamicin in various matrices

EUROPROXIMA GENTAMICIN ELISA

A competitive enzyme immunoassay for screening and quantitative analysis of Gentamicin in various matrices

TABLE OF CONTENTS

	PAGE:
Brief Information	2
1. Introduction	2
2. Principle of the Gentamicin-ELISA	3
3. Specificity and sensitivity	4
4. Handling and storage	5
5. Kit contents	6
6. Equipment and materials required but not provided.....	7
7. Safety precautions	7
8. Sample preparations	8
9. Preparations of reagents	10
10. Assay Procedure	11
11. Interpretation of results.....	12
12. Literature.....	14
13. Ordering information	14
14. Revision history	14


The quality management system of R-Biopharm Nederland B.V. is
ISO 9001:2015 certified

BRIEF INFORMATION

The Gentamicin ELISA is a competitive enzyme immunoassay for measurement of the concentration of gentamicin in various samples. With this ELISA kit 96 analyses can be performed. Samples and standards are measured in duplicate which means that a total 40 samples can be analyzed.

The ELISA kit contains all reagents to perform the assay. Reagents for sample preparation are not included in the kit.

1. INTRODUCTION



Gentamicin belongs to a group of carbohydrate containing antibiotics called aminoglycosides. All the aminoglycosides are potentially toxic compounds causing significant damage in vestibular and auditory functions in human as well as in animals. Nevertheless, they are used in practice because of their antibacterial and antifungal activities. These compounds have been found to be useful for the treatment of serious infections due to Gram negative micro-organisms. However, the range between therapeutic effectiveness and toxicity is narrow, therefore, dosage must be monitored. Aminoglycoside residues may occur in products of animal origin for several reasons such as deliberate feeding, inadvertent feeding to prevent infections in cows or to avoid outbreak of diseases of digestive and respiratory tracts of poultry.

Within the European Union, provisional Maximum Residue Limits for aminoglycosides have been fixed (see Table I).

Table I: Provisional Maximum Residue Limits (mg/kg) for aminoglycosides.

Aminoglycosides	Kidney	Liver	Muscle	Milk	Fat	Eggs
Streptomycin	1.0	0.5	0.5	0.2	0.5	-
Dihydrostreptomycin	1.0	0.5	0.5	0.2	0.5	-
Gentamicin	0.75	0.2	0.05	0.1	0.05	-
Neomycin	5.0	0.5	0.5	1.5	0.5	0.5

EuroProxima has also available a Neomycin- and (Dihydro) Streptomycin-ELISA.

2. PRINCIPLE OF THE GENTAMICIN-ELISA

The microtiter plate based gentamicin ELISA consists of one precoated plate (12 strips, 8 wells each). Antibody, horseradish peroxidase (-HRP) labelled gentamicin and standard solution or sample are added to the wells. Free gentamicin from the samples or standards and gentamicin-HRP conjugate compete for the specific antibody binding sites (competitive enzyme immunoassay).

After an incubation time of 1 hour, the non-bound reagents are removed in a washing step. The amount of bound gentamicin-HRP conjugate is visualized by the addition of a substrate/chromogen solution (H_2O_2 /TMB). Bound gentamicin-HRP conjugate transforms the colourless chromogen into a coloured product.

The substrate reaction is stopped by the addition of sulfuric acid. The colour intensity is measured photometrically at 450 nm. The optical density is inversely proportional to the gentamicin concentration in the sample.

3. SPECIFICITY AND SENSITIVITY

The Gentamicin ELISA utilizes antibodies raised in rabbits against protein conjugated gentamicin. The reactivity pattern of the antibody is:

Cross- reactivity:	Gentamicin	100%
	Sisomycin	25%
	Neomycin	< 0.1%
	Kanamycin	< 0.1%
	Tobramycin	< 0.1%
	Lincomycin	< 0.1%

The cross-reactivities are determined in a buffer system. The reported values may be different in samples due to matrix effects.

The test cannot discriminate between analytes and cross-reactive substances

The Limit of detection (LOD) and the detection capability (CC β) are determined under optimal conditions. Cut-off criteria need critical consideration.

Matrix	Procedure	LOD (ppb)	CC β (ppb)
Milk	8.1	0.4	10
Tissue	8.2	12.5	-
Serum	8.4	2	-
Honey	8.5	5	-
Egg	8.6	1	-
Urine	8.7	4	-
Feed	8.8	10	-

If the sample is found to be non-compliant, the results shall be verified by re-analysis of the sample using a confirmatory method.

4. HANDLING AND STORAGE

- Kit and kit components should be stored at 2°C to 8°C in a dark place. For repeated use store kit components as specified under chapter 9.
- After the expiry date of the kit and/or components has passed, no further quality guarantee is valid.
- Bring all kit components including the microtiter plate to ambient (room) temperature before use.
- Dilute the kit components immediately before use, but after the components are brought to ambient temperature.
- Avoid condensation in the wells of the plate. Bring the sealed plate to ambient temperature before opening the plate sealing.
- The substrate chromogen solution can be stored in a refrigerator (2°C to 8°C) until the expiry date stated on the label.
- Exposure of the chromogen solution to light should be avoided.

Degeneration of the reagents may have occurred when the following phenomena are observed:


- A blue colouring of the chromogen solution before transferring it into the wells.
- A weak or no colour reaction in the zero standard wells ($E_{450nm} < 0.8$).

5. KIT CONTENTS

Manual

One sealed (96-wells) microtiter plate (12 strips, 8 wells each), coated with antibody. Plate is ready-to-use.

Position of the reagents in the kit. For preparation of the reagents see Chapter 9.



1. **Dilution buffer** (20 ml, ready-to-use)
2. **Rinsing buffer** (30 ml, 20x concentrated)
3. **Substrate solution** (12 ml, ready-to-use)
4. **Stop solution** (15 ml, ready-to-use)
5. **Conjugate** (lyophilized, blue cap)
6. **Antibody** (lyophilized, yellow cap)
7. not in use
8. **Standard solution 1000 ng/ml** (1ml, ready-to-use)
9. not in use
10. not in use
11. not in use
12. not in use
13. **Zero Standard solution** (2ml, ready-to-use)
14. **Standard solution 1** (1ml, ready-to-use) **0.25 ng/ml**
15. **Standard solution 2** (1ml, ready-to-use) **0.5 ng/ml**
16. **Standard solution 3** (1ml, ready-to-use) **1 ng/ml**
17. **Standard solution 4** (1ml, ready-to-use) **5 ng/ml**
18. **Standard solution 5** (1ml, ready-to-use) **10 ng/ml**
19. **Standard solution 6** (1ml, ready-to-use) **50 ng/ml**

6. EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- Scales and weighing vessels
- Gloves
- Fume hood
- Homogeniser (blender, Ultra Turrax, mixer)
- Centrifuge (2000 x g)
- Vortex
- Automated microplate washer or 8 channel micropipette 100 - 300 μ l
- Microtiter plate shaker
- Microtiter plate reader with 450 nm filter
- Siliconised glass test tubes or plastic tubes
- Micropipettes 20 - 200 μ l, 100 - 1000 μ l
- Multipipette with 2.5 ml combitips
- Trichloroacetic acid
- Disodium hydrogen phosphate
- Potassium dihydrogen phosphate
- Potassium chloride
- Sodium chloride
- Tween 80

7. SAFETY PRECAUTIONS

- This kit may contain hazardous substances. For hazard notes please refer to the appropriate safety data sheets (SDS).
- Avoid contact of all biological materials with skin and mucous membranes.
- Do not pipette by mouth.
- Do not eat, drink, smoke, store or prepare foods, or apply cosmetics within the designated work area.
- Do not use components past expiration date and do not use components from different lots.
- Each well is ultimately used as an optical cuvette. Therefore, do not touch the under surface of the wells, prevent damage and dirt.
- All components should be completely dissolved before use. Take special attention to the substrate and rinsing buffer, which crystallize at +4°C.
- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting and washing throughout this procedure are necessary to maintain good precision and accuracy.

R-Biopharm Nederland makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm Nederland shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

8. SAMPLE PREPARATIONS

8.1 Milk samples

- Defat milk by centrifugation for 10 minutes at 4°C and 2000 x g.
- Eliminate the upper fat layer.
- Dilute and homogenize the defatted milk sample 10 times in sample dilution buffer (see Chapter 9).
- Check pH 7.4 ± 0.4 .
- Pipette 50 μ l of this solution into the ELISA.

8.2 Tissue samples

- Weigh 5 gram finely cut subsequently homogenized tissue in a plastic tube
- Add 20 ml sample dilution buffer (see Chapter 9)
- Homogenise (for instance using an Ultra Turrax or head over head mixer) for 30 minutes
- Centrifuge a part of the mixture 10 minutes at 4000 x g at 4°C
- Remove the upper fat layer
- Pipette 50 μ l supernatant into a plastic tube, add 450 μ l sample dilution buffer, vortex
- Use 50 μ l of this solution into the ELISA

8.3 Fat samples

- Weigh 1 g of homogenized (melted) fat in a plastic tube.
- Add 10 ml of sample dilution buffer (see Chapter 9) and heat in a water bath at 70°C for 30 minutes.
- Centrifuge 15 minutes, 2000 x g, at 4°C.
- Remove the upper fat layer and pipette 1 ml of sample extract in a plastic tube.
- Add 4 ml of sample dilution buffer mix well.
- Pipette 50 μ l of this solution into the ELISA.

8.4 Serum samples

- Dilute serum samples 10 times in sample dilution buffer (see Chapter 9).
- E.g. pipette 50 μ l of serum sample into a clean siliconised glass tube or plastic vial.
- Add 450 μ l of sample dilution buffer.
- Mix well using a vortex
- Use 50 μ l of this solution into the ELISA.

8.5 Honey samples

- Weigh 1 g of homogenized honey in a plastic tube.
- Add 4 ml of sample dilution buffer (see Chapter 9) and mix well using a vortex
- Wait for a minute to obtain a separation between the solid part and the liquid part in the sample.
- Pipette 1 ml of the clear upper liquid into a clean tube and add 4 ml of sample dilution buffer.
- Mix well using a vortex
- Use 50 μ l of the dilution into the ELISA.

8.6 Egg samples

- Homogenise an egg, both egg-protein and yolk.
- Pipette 1 ml of the homogenised sample into a clean tube.
- Add 4 ml of sample dilution buffer (see Chapter 9).
- Mix well using a vortex.
- Use 50 μ l of the mixture in the ELISA.

8.7 Urine samples

- Dilute urine samples 10 times in sample dilution buffer (see Chapter 9).
- Check pH 7.4 ± 0.4 .
- Use 50 μ l of the dilution in the ELISA.

8.8 Feed samples

- Weigh 1 g finely cut and subsequently homogenised feed sample in the plastic tube.
- Add 5 ml of a trichloroacetic acid solution (3%)
- Homogenise (for instance using an Ultra Turrax) for 1 minute.
- Mix head over head for 30 minutes.
- Centrifuge 10 minutes at 2000 x g and 4°C.
- Remove the upper fat layer eventually present.
- Pipette 100 μ l of the supernatant into a plastic tube.
- Add 900 μ l of sample dilution buffer (see Chapter 9) and mix.
- Adjust pH to 7.4 ± 0.4 .
- Pipette 50 μ l of this solution into the microtiter plate.

9. PREPARATION OF REAGENTS

Before beginning the test, the reagents should be brought up to ambient temperature. Any reagents not used should be put back into storage immediately at 2°C to 8 °C. Prepare reagents fresh before use.

Microtiter plate

Return unused strip into resealable bag with desiccant and store at 2°C to 8°C for use in subsequent assays. Retain also the strip holder.

Rinsing buffer

The rinsing buffer is delivered 20 times concentrated. Prepare dilutions freshly before use. Per strip 20 ml of diluted rinsing buffer is used (1 ml concentrated rinsing buffer + 19 ml distilled water).

Substrate solution

The substrate solution (ready-to-use) precipitates at 4°C. Take care that this vial is at room temperature (keep in the dark) and mix the content before pipetting into the wells.

Standard solution 1000 ng/ml

To prepare standards in the appropriate matrix or to prepare spikes use the standard solution containing 1000 ng gentamicin per ml. Dilute the standard solution in the appropriate matrix to make a dilution range of 50, 10, 5, 1, 0.5, 0.25 ng/ml. Also the zero standard should be of the same matrix.

Conjugate solution

Reconstitute the vial of lyophilised conjugate (gentamicin-HRP) with 4 ml of dilution buffer, mix thoroughly and keep in the dark until use.

Antibody solution

Reconstitute the vial of lyophilised antibodies with 4 ml of dilution buffer, mix thoroughly and keep in the dark until use.

Sample dilution buffer

Sample dilution buffer is not provided in the kit. Prepare this buffer as follows:

Dissolve in 1 litre distilled water:

Na ₂ HPO ₄	1.15 g
KH ₂ PO ₄	0.2 g
KCl	0.2 g
NaCl	30 g
Tween 80	0.5 ml
pH	7.4 (7.3-7.5)

10. ASSAY PROCEDURE

Rinsing protocol

In ELISA's, between each immunological incubation step, unbound components have to be removed efficiently. This is reached by appropriate rinsing. It should be clear that each rinsing procedure must be carried out with care to guarantee good inter- and intra-assay results.

Basically, manual rinsing or rinsing with automatic plate wash equipment can be performed as follows:

Manual rinsing

1. Empty the contents of each well by turning the microtiter plate upside down and remove residual liquid by striking the plate against a paper towel.
2. Fill all the wells to the rims (300 μ l) with rinsing solution.
3. This rinsing cycle (1 and 2) should be carried out 3 times.
4. Turn the plate upside down and empty the wells by a firm short vertical movement.
5. Place the inverted plate on absorbent paper towels and tap the plate firmly to remove residual washing solution from the wells.
6. Take care that none of the wells dry out before the next reagent is dispensed.

Rinsing with automatic microtiter plate wash equipment

When using automatic plate wash equipment, check that all wells can be aspirated completely, that the rinsing solution is nicely dispensed reaching the rim of each well during each rinsing cycle. The washer should be programmed to execute three rinsing cycles.

Assay Protocol

1. Prepare samples according to Chapter 8 and prepare reagents according to Chapter 9.
2. Pipette 100 μ l of zero standard in duplicate (wells H1, H2, blank).
Pipette 50 μ l of zero standard in duplicate (wells A1, A2).
Pipette 50 μ l of each of the standard dilutions in duplicate (wells B1,2 to G1,2 i.e. 50,10, 5, 1, 0.5 and 0.25 ng gentamicin/ml).
3. Pipette 50 μ l of each sample solution in duplicate into the remaining wells of the microtiter plate.
4. Pipette 25 μ l of conjugate (gentamicin-HRP) to all wells, except H1 and H2.
5. Add 25 μ l of antibody solution into all wells, except wells H1 and H2.
6. Seal the microtiter plate and shake the plate for a few seconds on a microtiter plate shaker.

7. Incubate for 1 hour in the dark in a refrigerator (2°C to 8°C).
8. Discard the solution from the microtiter plate and wash 3 times with rinsing buffer.
9. Pipette 100 µl of substrate solution into each well.
10. Incubate 30 minutes in the dark at room temperature (20°C to 25°C).
11. Pipette 100 µl of stop solution to each well.
12. Read the absorbance values immediately at 450 nm.

11. INTERPRETATION OF RESULTS

In general, the O.D. values of the samples, that are measured in this ELISA test, are converted into corresponding concentrations in the following way:

The mean optical density (O.D.) value of the blank wells is subtracted from the individual O.D. of the wells containing the standards and the samples. Then these O.D. values of the standards and the samples (mean values of the duplicates) are divided by the mean O.D. value of the zero standard and multiplied by 100. The zero standard is thus made equal to 100% (maximal absorbance) and the other O.D. values are quoted in percentages of the maximal absorbance. The values calculated for the standards are entered in a semi-logarithmic plot against the concentration of the analyte.

Alternative for calibration curve:

The values of absorption (logit) calculation of the standards are plotted on Y-axis versus the analyte equivalent concentration on a logarithmic X-axis.

As customer-friendly options R-Biopharm offers 2 software programs that can be used for the interpretation of the ELISA test:

- RIDASOFT@Win.NET Food & Feed (art. nr. Z9996FF). In this software program all R-Biopharm ELISA tests (including EuroProxima) are preprogrammed in the database. A list of compatible readers is available on request.
- Simplefit. This Excel based program is developed to calculate the results of all EuroProxima tests.

8.1 Milk samples

The amount of gentamicin in the milk samples is expressed as gentamicin equivalents (ng/ml). The gentamicin equivalents in the milk (ng/ml) corresponding to the % maximal absorbance of each sample can be read from the calibration curve. These calculated gentamicin equivalents have to be multiplied by 10 to obtain the gentamicin equivalents (ng/ml) in the undiluted milk sample.

8.2 Tissue samples

The gentamicin equivalents in the tissue extract corresponding to the % maximal absorbance of each sample can be read from the calibration curve. These equivalents have to be multiplied by 50 to obtain the gentamicin equivalents (ng/g) in the tissue samples.

8.3 Fat samples

Multiply the calculated gentamicin equivalents by 50 to obtain the gentamicin equivalents (ng/g) in fat samples

8.4 Serum samples

Multiply the calculated gentamicin equivalents by 10 to obtain the gentamicin equivalents (ng/ml) in undiluted serum samples.

8.5.Honey samples

Multiply the calculated gentamicin equivalents by 25 to obtain the gentamicin equivalents (ng/g) in honey samples.

8.6 Egg samples

Multiply the calculated gentamicin equivalents by 5 to obtain the gentamicin equivalents (ng/ml) in egg samples.

8.7 Urine samples:

Multiply the calculated gentamicin equivalents by 10 to obtain the gentamicin equivalents (ng/ml) in undiluted urine samples.

8.8 Feed samples:

Multiply the calculated gentamicin equivalents by 50 to obtain the gentamicin equivalents (ng/g) in feed samples.

12. LITERATURE

Berdy J., Aszalos A., Bostian M and McNitt K.L. Handbook of Antibiotic Compounds, Vol. 1, CRC Press Inc. Boca Raton, Florida U.S.A. 1980.
Standefor J.C. and Saunderson G.C. Enzyme immunoassay for Gentamicin, Clin. Chem.(1978), 24, 1903.
Analysis of Antibiotic Drug Residues in Food Products of Animal Origin., Ed. Vipin K. Agarwal U.S.A. 1992.
Haasnoot W., Stouten P., Cazemier G., Lommen A., Nouws J.F.M. and Keukens H.J. Immunochemical detection of Aminoglycosides in milk, and kidney. The Analyst, 1999, 124, 301-305.

13. ORDERING INFORMATION

For ordering the Gentamicin ELISA kit, please use cat. code 5111GEN.

14. REVISION HISTORY

Chapter 11 has been updated.