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Simplex[®] Easy DNA Kit

1. Intended use

Fast and optimal DNA-extraction from bacteria and yeast

2. Content

The kit contains reagents for 100 or 10 extractions:

1 x Easy-Reagent (10 mL- 100 rxn; 1 mL- 10 rxn)

3. Materials required but not provided

3.1 Instruments:

Heating block or waterbath 37 °C – 95 °C (98 ° - 203 °F)

Microcentrifuge for 1.5 – 2.0 mL reaction tubes

Pipettes

“Vortex”

3.2 Reagents and plastic ware:

for yeast DNA-extraction: Lyticase (Sigma without DNA and RNA) in **PCR-Kits with freeze-dried stripes the PCR-tubes contain lyticase (s. PCR Kit manual).**

sterile safe-lock reaction tubes 1.5 – 2.0 mL

suitable filter tips

single use gloves

4. Warnings

Careful use with personal protection according to good laboratory practice is recommended.

Do not incorporate. Avoid skin and eye contact with all solutions.

5. Storage

Store solution at room temperature 15 – 30 °C (59 – 86 °F)

6. Indications of deterioration of reagents

In case of accurate handling deterioration unknown

7. Preliminary preparations

Prepare Lyticase solution:

5U / μ L in lyticase buffer (50 % 1xTE+50 % Glycerol, pH-Wert 7.5 – 8.0). **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).**

Store solution at -20 °C (- 4 °F)

8. DNA-Extraction

8.1 DNA-Extraction from bacteria


Hint:

Samples that contain a lot of yeast should be centrifuged as follows:

Centrifuge 1.5 mL sample for 5 minutes at 400 x g. Transfer 1 mL of the supernatant to a 1.5 mL safe-seal reaction tube and follow the protocol 8.1


Protocol:

1. mix the sample and transfer 1.0 mL in a 1.5 mL safe-lock reaction tube
2. place the tube always in the same manner into the centrifuge (lid of the tube), so that even when no pellet will be visible, it will be clear where the cells/pellet will be

3. centrifuge 5 minutes at 15500 x g

4. discard carefully and completely the supernatant. Aspirate from the opposite side of the pellet. The bacterial pellet should not be greater than a pinhead (e.g. 10^7 cfu / mL), otherwise take less sample volume (100 μ L of the sample)

5. resuspend the pellet completely in **100 μ L Easy-Reagent** (vortex or pipetting)
*For DNA-extraction from colonies: pick a colony or a part of it with a sterile inoculating loop and resuspend it in **100 μ L Easy-Reagent***
6. incubate 10 minutes at 95 °C (203 °F)
7. centrifuge 5 minutes at 15500 x g
8. transfer the supernatant (contains DNA) in a new 1.5 mL safe-lock reaction tube
Attention: transfer no pellet
9. use 2.5 μ L DNA for PCR

8.2 DNA-Extraction from yeast

1. mix the yeast sample and transfer 1 mL in a 1.5 mL safe-seal reaction tube
2. place the tubes always in the same manner into the centrifuge (lid of the tube), so that even when no pellet will be visible, it will be clear where the cells/pellet will be
3. centrifuge 5 minutes at 15500 x g
4. discard carefully and completely the supernatant. Aspirate from the opposite side of the pellet. The yeast pellet should not be greater than a pinhead (e.g. 10^6 cfu / mL), otherwise take less sample volume (50-100 μ L)
5. resuspend the pellet completely in **100 μ L Easy-Reagent** (e.g. vortex)
*For DNA-Extraction from yeast colonies: pick colony or a part of it with a sterile inoculating loop and resuspend it in **100 μ L Easy-Reagent***
6. add **15 U Lyticase** and mix the sample. **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).** Continue with step 8.
7. incubate 15 minutes at 37 °C (98 °F)
8. incubate 10 minutes at 95 °C (203 °F)
9. centrifuge 5 minutes at 15500 x g
10. transfer the supernatant (contains DNA) in a new 1.5 mL safe-lock reaction tube
Attention: transfer no pellet
11. use 2.5 μ L DNA for PCR

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.