---

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# QuickGEN Sample Preparation Centrifugation

## 1. Intended use

Fast and optimal DNA-extraction from bacteria and yeast out of beverages with low inhibitory components. For samples with and without preenrichment, colonies from agar plates and swaps by centrifugation.

## 2. Content

The kit contains reagents for 50 - 100 extractions, depending on the protocol:

1 x QuickGEN buffer

1 x Washing solution

## 3. Materials required, but not provided

### 3.1 Instruments

- Centrifuge
- Pipettes
- "Vortex"

### 3.2 Reagents and plastic ware

- for yeast DNA-extraction: Lyticase (Sigma L2524) **in PCR-Kits with freeze-dried stripes the PCR-tubes contain lyticase (s. PCR Kit manual).**
- safe-lock reaction tubes 1.5 – 2.0 mL
- suitable filter tips
- single use gloves

## 4. Warnings

Careful use with personal protection according to good laboratory practice is recommended. Do not incorporate. Avoid skin and eye contact with all solutions

## 5. Storage

Store all solutions at room temperature 15 – 30 °C (59 – 86 °F)

## 6. Indications of deterioration of reagents

In case of accurate handling deterioration unknown

## 7. Preliminary preparations

Prepare Lyticase solution:

5U /  $\mu$ L in lyticase buffer (50 % 1xTE+50 % Glycerol, pH-Wert 7.5 – 8.0). **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).**

Store solution at -20 °C (- 4 °F)

## 8. DNA-Extraction

**Hint:** For the detection of bacteria or spores in samples with high concentrations of yeast cells the following centrifugation protocol should be used:

Centrifuge 1.5 mL for 5 min. at 400 x g. Carry over 1 mL / 100  $\mu$ L or 10  $\mu$ L of the supernatant to a 1.5 mL reaction tube according to the hints for sampling of preenriched cultures

For samples with high yeast cell counts e.g. yeast tanks, propagation etc. the QuickGEN yeast Sample Preparation Centrifugation Kit (Q005) should be used.

### Hints for sampling of preenriched cultures:

Clear culture:

Weak turbidity:

Bacteria:  $\sim 10^7$  cfu/mL  
Yeast:  $\sim 10^6$  cfu/mL

Strong turbidity:

Bacteria:  $\sim 10^9$  cfu/mL  
Yeast:  $\sim 10^{8-9}$  cfu/mL


Very strong turbidity:

Bacteria:  $> 10^9$  cfu/mL  
Yeast:  $> 10^8$  cfu/mL

Prepare 1 mL

Protocol A

Prepare 100  $\mu$ L

Protocol A

Prepare 10  $\mu$ L without  
centrifugation

Protocol D

Dilute 1:10  
Prepare 10  $\mu$ L without  
centrifugation

Protocol D

### **Protocol A (for 1 mL / 100 µL sample e.g. preenriched sample)**

1. Mix the sample and transfer 1 mL/ 100 µL in a 1.5 mL reaction tube (1 mL for clear samples, 100 µL for weak turbid samples)
2. Centrifuge the sample for 5 min. at 15.500 x g
3. Remove the liquid  
*Optional* for samples with potential inhibitors (e.g. dark beer, red wine): add 100 µL Washing solution, resuspend the pellet, centrifuge the sample for 5 min. at 15.500 x g and remove the liquid
4. Add 100 µL QuickGEN buffer, resuspend the pellet (e.g. vortex or pipette tip)
5. ***Optional for yeast:* Add 15 U Lyticase** and mix the sample. **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).**
6. Use 2.5 - 5 µL of the DNA for PCR depending on the used PCR-Kit (if PCR is inhibited dilute the sample 1:10 in QuickGEN buffer)

### **Protocol B (for 30 mL beverage sample)**

1. Centrifuge 30 mL sample for 20 min. at 5.500 x g in a 50 mL falcon tube
2. Remove the liquid and turn the falcon tube upside down on a waterleaf paper  
*Optional* for samples with potential inhibitors (e.g. dark beers, red wine): add 500 µL washing solution, resuspend the pellet, centrifuge the sample for 15 min. at 5.500 x g and remove the liquid
3. Add 100 µL QuickGEN buffer to the pellet and resuspend (e.g. vortex or pipette tip)
4. Fill the volume into a 1.5 mL reaction tube to avoid cross contamination later on
5. ***Optional for yeast:* Add 15 U Lyticase** and mix the sample. **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).**
6. Use 2.5 - 5 µL of the DNA for PCR depending on the used PCR-Kit (if PCR is inhibited dilute the sample 1:10 in QuickGEN buffer)

### Protocol C (for bacteria or yeast colony)

1. Add QuickGEN buffer in a 1.5 mL reaction tube according to the size of the colony (see step 2)
2. Pick a colony ( $\leq 1$ mm) with an inoculating loop (size 1  $\mu$ L) and solve it in 100  $\mu$ L of the QuickGEN buffer. For colonies  $> 1$ mm, take only a part of the colony and solve it in 100  $\mu$ L QuickGEN buffer. Very small colonies have to be solved in 50  $\mu$ L QuickGEN buffer.

3. **Optional for yeast: Add 15 U Lyticase** and mix the sample. **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).**
4. Use 2.5 - 5  $\mu$ L of the DNA for PCR depending on the used PCR-Kit (if PCR is inhibited dilute the sample 1:10 in QuickGEN buffer)

### Protocol D (e.g. for culture yeast)

1. Transfer 10  $\mu$ L sample in a 1.5 mL reaction tube
2. Add 100  $\mu$ L QuickGEN buffer and vortex
3. **Optional for yeast: Add 15 U Lyticase** and mix the sample. **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).**
4. Use 2.5 - 5  $\mu$ L of the DNA for PCR depending on the used PCR-Kit (if PCR is inhibited dilute the sample 1:10 in QuickGEN buffer)

### Protocol E (for swabs)

1. Place the swab into a 1.5 mL reaction tube
2. Add 200  $\mu$ L QuickGEN buffer (the swab should be completely be covered) and vortex 30 sec.
3. **Optional for yeast: Add 15 U Lyticase** and mix the sample. **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).**
4. Use 2.5 - 5  $\mu$ L of the DNA for PCR (if PCR is inhibited dilute the sample 1:10 in QuickGEN buffer)