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QuickGEN Yeast Sample Preparation Centrifugation

1. Intended use

Microbial DNA-Extraction from samples with high yeast cell counts e.g. yeast tanks, propagation etc., that are older and may contain dead yeast cells.

2. Content

The kit contains reagents for 100 extractions:

1 x QuickGEN yeast buffer

2 x Enzyme

3. Materials required, but not provided

3.1 Instruments

Heating block or waterbath 65 °C – 95 °C (149 °F - 203 °F)

Centrifuge

Pipettes

“Vortex”

3.2 Reagents and plastic ware

for yeast DNA-extraction: Lyticase (Sigma L2524), **in PCR-Kits with freeze-dried stripes the PCR-tubes contain lyticase (s. PCR Kit manual).**

sodium chlorid or Ringer solution

steril Falcon tubes 15 mL, 50 mL

safe-lock reaction tubes 1.5 – 2.0 mL

suitable filter tips

single use gloves

4. Warnings

Careful use with personal protection according to good laboratory practice is recommended.

Do not incorporate. Avoid skin and eye contact with all solutions.

5. Storage

Store QuickGEN yeast buffer at room temperature 15 °C – 30 °C (59 °F – 86 °F) and Enzyme after opening at -20 °C (- 4 °F).

6. Indications of deterioration of reagents

In case of accurate handling deterioration unknown

7. Preliminary preparations

Prepare Lyticase solution:

5 U / μ L in lyticase buffer (50 % 1 x TE + 50 % Glycerol, pH-Wert 7.5 – 8.0). **In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual).**

Store solution at -20 °C (- 4 °F)

8. DNA-Extraction

A. Protocol for yeast samples, yeast tanks ($>10^9$ cfu/mL)

1. Transfer 1 mL sample in a sterile 50 mL Falcon tube
2. Dilute the sample **1:20** (1 mL sample + 19 mL dilution buffer*)
3. Transfer 1 mL of the dilution to a 1.5 mL reaction tube
4. Add 20 μ L enzyme
5. Incubate 30 min. at 65 °C
6. Centrifuge the sample for 5 min. at 15500 x g
7. Remove the supernatant
8. Add 250 μ L QuickGEN yeast buffer and resuspend the pellet
9. Incubate 5 min. at 95 °C

Optional for yeast: After incubation add 15 U Lyticase and mix the sample. In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual)

10. Use 2.5 μ L of the DNA for PCR

* e.g. sterile Ringer Solution or NaCl, do not use ddH₂O

B. Protocol for yeast samples, propagation (10^5 - 10^9 cfu/mL)

1. Transfer 1 mL sample in a sterile 15 mL Falcon tube
2. Dilute the sample **1:10** (1 mL sample + 9 mL dilution buffer*)
3. Transfer 1 mL of the dilution to a 1.5 mL reaction tube
4. Add 20 μ L enzyme
5. Incubate 30 min. at 65 °C
6. Centrifuge the sample for 5 min. at 15500 x g
7. Remove the supernatant
8. Add 250 μ L QuickGEN yeast buffer and resuspend the pellet
9. Incubate 5 min. at 95 °C

Optional for yeast: After incubation add 15 U Lyticase and mix the sample. In PCR-Kits with freeze-dried PCR stripes the tubes contain lyticase (s. PCR Kit manual)

10. Use 2.5 μ L of the DNA for PCR

* e.g. sterile Ringer Solution or NaCl, do not use ddH₂O