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Advanced Settings MyGo Pro	Integration time (s)	0.7
Advanced Settings MyGo Pro ESR	Integration time (s)	1.0

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
QuickGEN PCR Kit

P1 and *S. cerevisiae* var. *diastaticus* Screening

1. Examples of successful detected microorganisms

- *Lactobacillus* spp.
 - L. acetotolerans*
 - L. acidophilus*
 - L. amylolyticus*
 - L. amylovorus*
 - L. backii*
 - L. brevis*
 - L. brevisimilis*
 - L. buchneri*
 - L. casei* / *paracasei*
 - L. collinoides* / *paracollinoides*
 - L. coryniformis* ssp. *coryniformis*
 - L. coryniformis* ssp. *torquens*
 - L. curvatus*
 - L. delbrückii* ssp. *delbrückii* / *lactis*
 - L. harbinensis*
 - L. helveticus*
 - L. johnsonii*
 - L. lindneri*
 - L. malefermentans*
 - L. parabrevis*
 - L. parabuchneri* (*frigidus*)
 - L. paucivorans*
 - L. pentosus*
 - L. perolens*
 - L. plantarum*
 - L. paraplantarum*
 - L. rossiae*
- *Lactococcus lactis*
- *Megasphaera* spp.
- *Micrococcus kristiniae*
 - *Pectinatus* spp.
 - Pectinatus cerevisiiphilus*
 - Pectinatus frisingensis*
 - Pectinatus haikarae*
 - Pectinatus* ssp.
- *Pediococcus* spp.
 - Pediococcus acidilactici*
 - Pediococcus damnosus*
 - Pediococcus dextrinicus*
 - Pediococcus inopinatus*
 - Pediococcus pentosaceus*
 - Pediococcus parvulus*
 - Pediococcus claussenii*
- *Saccharomyces cerevisiae* var. *diastaticus*

2. Intended use

Qualitative detection of bacterial and *S. cerevisiae* var. *diastaticus* contaminations in filtrated beer, beer mixes and yeast containing samples.

3. Test principle

The TaqMan® real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ; ATTO/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument. To avoid false negative PCR-results an inhibition control (ROX/DQ) is amplified together in one reaction vessel with the specific sequence. **The tubes contain lyticase.**

4. Kit contents

The kit contains sufficient reagents for 48 reactions:

1	x Premix	white cap
6	x Dye Strips (freeze-dried, incl. IC-DNA and lyticase)	tube strips
6	x Cap Strips	
1	x ddH ₂ O	colourless cap
1	x Control-DNA (freeze-dried)	yellow cap

5. Storage conditions

The Control-DNA is freeze-dried and has to be solved in ddH₂O prior to use (see 7.1).

Do **not** freeze the lyophilized Dye Strips and lyophilized Control-DNA.

The PCR-reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) **after arrival**. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

6. Materials required, but not provided

6.1. Instruments

Real-time PCR machine for low profile PCR tubes clear

Centrifuge for 1.5 – 2.0 mL reaction tubes

Centrifuge for strips

Pipettes

“Vortex”

6.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

7. PCR

7.1. PCR-Setup

When using the kit for the first time, the freeze-dried Control-DNA has to be shortly centrifuged and carefully resolved:

- add 55 µL sterile ddH₂O to the freeze-dried Control-DNA
- after 15 minutes mix well

Before every use thoroughly mix all PCR-components and centrifuge briefly.

PCR-reaction Setup:

Remove the film from the needed tube strips and pipette the PCR-Components. After pipetting close the tube strips with the provided cap strips.

PCR-components	volume (µL)
Premix	15.0
Sample-DNA	5.0*
Total volume	20.0

* if using the Simplex easy DNA- or QuickGEN Yeast Sample Preparation Kit:
add 2.5 µL DNA and 2.5 µL sterile ddH₂O

1. Pipette 15 µL of the Premix into each PCR-tube, making sure that, prior to the first filling, the tip has been moistened.
2. Add 5 µL sample-DNA, add 5 µL of the Control-DNA for the positive control reaction, add 5 µL of the extraction control and 5 µL of ddH₂O for the negative control* reaction. Use a fresh tip with each DNA filling.
3. Close the tubes immediately and centrifuge them shortly.
4. Place the tubes in the PCR-machine and start run.

Very important: * Please fill up the negative control with 5 µL ddH₂O to avoid unspecific amplification.

Work swiftly to avoid warming up and keep away from light.

7.2 PCR-Program

Software MyGo Pro/ ESR 3.5 / 3.6

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30* - 35 x
Cycling Annealing/ Elongation	30 sec	63 °C	

*** for preenriched samples program 30 cycles**

Advanced Settings MyGo Pro	Integration time (s)	0.7
-----------------------------------	----------------------	-----

Advanced Settings MyGo Pro ESR	Integration time (s)	1.0
---------------------------------------	----------------------	-----

8. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

Lactobacillus/ Pediococcus-DNA: FAM-channel

Megasphaera/ Pectinatus-DNA: HEX-channel

Inhibition Control-DNA: ROX-channel

***S. cerevisiae* var. *diastaticus*-DNA:** ATTO-channel


A sample is ***Lactobacillus/ Pediococcus* positive**, if there is a detectable fluorescence increase in the **FAM-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the ROX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is ***Megasphaera/ Pectinatus* positive**, if there is a detectable fluorescence increase in the **HEX-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the ROX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is ***S. cerevisiae* var. *diastaticus* positive**, if there is a detectable fluorescence increase in the **ATTO-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the ROX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is **negative**, if there is no detectable fluorescence increase in the different channels and the positive controls have a positive fluorescence signal. The negative controls show no amplification. The Inhibition Control in the ROX-channel has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

analysis flowchart

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc..

The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.