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QuickGEN PCR Kit

Yeast and bacteria differentiation

1. Intended use

Differentiation of beer spoilage bacteria and yeast in beer and beer mixing drinks.

The following bacteria and yeasts are differentiated:

Lactobacilli, *Pediococci*, Acetic acid bacteria, *Enterobacteriaceae*, Bottom-fermented yeast, Top-fermented yeast and wild yeast.

2. Test principle

The TaqMan® real-time PCR is based on hot-start PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ; ROX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

To avoid false negative PCR-results an Inhibition Control is amplified together in one tube with the specific sequence. **Tubes 4, 5, 7 and 8 contain lyticase.**

3. Kit contents

The kit contains sufficient reagents for 96 PCR reactions (24 samples):

1 x Premix	white cap
12 x Dye Strips (freeze-dried, tubes 4, 5, 7 und 8 contain lyticase)	tube strips
12 x Cap Strips	
1 x ddH ₂ O (for PCR negative controls)	colourless cap

4. Storage conditions

Do **not** freeze the lyophilized Dye Strips.

The PCR reagents should be stored at 2– 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) **after arrival**. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR instrument for white profile tubes

Centrifuge for strips

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

Colour Compensation Set for LC480 (GEN-IAL[®], REF Q800)

LC480: Use a Colour compensation prior to first use of the kit

6. PCR

6.1. PCR-Setup

Before every use thoroughly mix all PCR-components and centrifuge briefly.

Remove the film from the needed tube strips and pipette the PCR-Components. After pipetting close the tube strips with the provided cap strips.

PCR-reaction setup:

PCR-Components	Amount (µL)
Premix	15.0
Sample-DNA	5.0*
Total volume	20.0

* if using the Simplex easy DNA- or QuickGEN Sample Preparation Kit:

add 2.5 µL DNA and 2.5 µL sterile ddH₂O

Pipetting scheme:

Tube 1*	Premix 15.0 µL	ddH ₂ O 5.0 µL, Attention: pipette no sample-DNA (negative control)
Tube 2*	Premix 15.0 µL	ddH ₂ O 5.0 µL, Attention: pipette no sample-DNA (positive control)
Tube 3*	Premix 15.0 µL	sample-DNA of sample 1 5.0 µL
Tube 4*	Premix 15.0 µL	sample-DNA of sample 1 5.0 µL
Tube 5*	Premix 15.0 µL	sample-DNA of sample 1 5.0 µL
Tube 6*	Premix 15.0 µL	sample-DNA of sample 2 5.0 µL
Tube 7*	Premix 15.0 µL	sample-DNA of sample 2 5.0 µL
Tube 8*	Premix 15.0 µL	sample-DNA of sample 2 5.0 µL

* please notice the orientation of the tubes. They are marked on the strip at Tube 1 with an A and at Tube 8 with an H.

6.2 PCR-Program

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30 ¹ -35 x
Cycling Annealing/ Elongation	20-30 ² sec	62 °C	

¹ for preenriched samples program 30 cycles

² please program 20 sec, if the real-time-PCR machine allows it

6.2.1 PCR Program LC480

1. Click the button **tool** on the right side in the window **LightCycler 480 Software release 1.5.0. SP1**
2. Click the button **Detection formats** at the left side of the menu bar
3. Click **New** in the window **Detection formats** and name the experiment
4. Open the window **Filter Combination Selection** and choose the following filter combinations: 465-510 / 533-580 / 533-610
5. Open the window **Selected Filter Combination** list and add the following amounts

Excitation Filter	Emission Filter	Name	Melt factor	Quant factor	Max. Integ. Time
465	510	465-510 FAM	1	10	2
533	580	533-580 HEX	1	10	2
533	610	533-610 610	1	10	2

6. **Close** the window
7. Click the button **New Experiment** on the right side of the menu bar
8. From the pull-down menu **Detection formats** choose the defined experiment, click the button **Customize** and check the detection formats. All of them have to be activated
9. Click the button **ok**
10. Define the following program:

1. Programm Name: **Heat**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
37	None	00:15:00	4.40		0	0	0
95	None	00:15:00	4.40		0	0	0

2. Programm Name: **Ampli**

Cycles **30*-35** Analysis Mode **Quantification**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
95	None	00:00:15	4.40		0	0	0
62	Single	00:00:20	2.20		0	0	0

* for preenriched samples program 30 cycles

3. Programm Name: **Cool**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
40	None	00:00:20	2.20		0	0	0

Optional: Saving the program as **run template**. The template button allows to select and apply a template to the currently open object and to save the currently open object as a template. Click the clamp beside the button **Apply Template** . Click **Save as template** and save the file in **templates**.

11. Click the button **Subset editor** on the left side

12. Click the button **+** and **New Subset 1** appears

13. Mark the wells in the **New Subset 1 Settings** window

14. Click the button **Apply**

15. Click the button **Sample editor** on the left side

16. **Very important:** Activate In the window **Step 1 Select Workflow: Abs.Quant**

17. Choose the subset **New Subset** in the window **Step 2 Select Samples**

18. Define your probes in the **Sample table**

19. Switch to the button **Experiment** and start the run with the button **start run**

7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

LC480: For analysis activate Colour Compensation

The negative control (NTC, tube 1) has to be negative in the FAM-, and HEX-channel and positive in the ROX-channel (Inhibition Control). Is the Inhibition Control in the ROX-channel negative and/or the results in the FAM-/ HEX-channel positive, PCR has to be repeated.

The Inhibition Controls in the samples (tubes 4/7) have to be positive in the ROX-channel. Is the Inhibition Control in the ROX-channel negative or weak inhibitors or high amount of DNA are in the sample reaction. The sample should be repeated with a dilution.

The Positive Control (tube 2) has to be positive in the FAM-channel, if it is negative PCR, has to be repeated.

*** Wild yeast 1 detects:**

<i>Dekkera anomala</i>	<i>Dekkera bruxellensis</i>	<i>Dekkera custersiana</i>	<i>Dekkera naardenensis</i>	<i>Debaromyces hansenii</i>
<i>Hanseniaspora guilliermondii</i>	<i>Hanseniaspora osmophila</i>	<i>Hanseniaspora uvarum</i>	<i>Issotchenkia orientalis</i>	<i>Kazachstania exigua</i>
<i>Kluyveromyces marxianus</i>	<i>Metschnikowia pulcherrima</i>	<i>Pichia anomala</i>	<i>Pichia fermentans</i>	<i>Pichia membranaefaciens</i>
<i>S. cerevisiae</i> var. <i>diastaticus</i>	<i>Saccharomycodes ludwigii</i>	<i>Torulaspora delbrückii</i>		

*** Wild yeast 2 detects:**

<i>Candida glabrata</i>	<i>Candida albicans</i>	<i>Candida kefyr</i>	<i>Candida intermedia</i>	<i>Candida parapsilosis</i>
<i>Candida sake</i>	<i>Candida tropicalis</i>	<i>Naumovozyma dairenensis</i>	<i>Pichia guilliermondii</i>	<i>Zygosaccharomyces bailii</i>
<i>Zygosaccharomyces rouxii</i>				

PCR-Analysis:

please notice the orientation of the tubes. They are marked on the strip at Tube 1 with an A and at Tube 8 with an H.

Tube	Sample	FAM	HEX	ROX
1	NTC	-	-	Inhibition Control
2	PTC	Positive Control	-	-
3	1	Enterobacteriaceae	<i>Lactobacillus / Pediococcus</i>	<i>Pediococcus</i>
4	1	Wild yeast 1*	Bottom fermented yeast	Inhibition Control
5	1	Wild yeast 2*	Top fermented yeast	Acetic acid bacteria
6	2	Enterobacteriaceae	<i>Lactobacillus / Pediococcus</i>	<i>Pediococcus</i>
7	2	Wild yeast 1*	Bottom fermented yeast	Inhibition Control
8	2	Wild yeast 2*	Top fermented yeast	Acetic acid bacteria

NTC: negative Control PTC: positive Control

Example for interpretation:

If sample 1 shows a positive FAM-signal in tube 4, wild yeast 1* is identified. If the sample contains mixtures, you have positive signals in different channels.

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.