---

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Pipettierschema:

Tube 1*	Premix 15,0 µL	<b>ddH<sub>2</sub>O 5,0 µL,</b> <b><u>Achtung:</u> keine Proben-DNA pipettieren</b>
Tube 2*	Premix 15,0 µL	Proben-DNA 5,0 µL
Tube 3*	Premix 15,0 µL	Proben-DNA 5,0 µL
Tube 4*	Premix 15,0 µL	Proben-DNA 5,0 µL
Tube 5*	Premix 15,0 µL	Proben-DNA 5,0 µL
Tube 6*	Premix 15,0 µL	Proben-DNA 5,0 µL
Tube 7*	Premix 15,0 µL	Proben-DNA 5,0 µL
Tube 8*	Premix 15,0 µL	Proben-DNA 5,0 µL

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**Advanced Settings MyGo Pro**    Integration time (s)    0.7

**Advanced Settings MyGo Pro ESR**    Integration time (s)    1.0

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# QuickGEN PCR Kit

## Screening and differentiation of 14 wine spoilage bacteria and yeast

### 1. Intended use

Screening and differentiation of 14 bacteria and yeast in wine. The following yeast are differentiated:

*Lactobacillus spp.*

*Pediococcus spp.*

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*Leuconostoc mesenteroides*

*Oenococcus oeni*

*Lactococcus lactis*

*Alicyclobacillus spp.*

*Diastatische Saccharomyces cerevisiae*

*Schizosaccharomyces pombe*

*Dekkera spp. :*

*D. anomala*

*D. custersiana*

*D. bruxellensis*

*D. naardenensis*

*D. nanus*

*Candida spp. :*

*C. albicans*

*C. parapsilosis*

*C. glabrata*

*C. tropicalis*

*C. sake*

*C. intermedia*

*C. pulcherrima*

*Pichia spp. :*

*W. anomala*

*P. fermentans*

*P. membranaefaciens*

*P. guilliermondii*

## 2. Test principle

The TaqMan® real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ; ROX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument. To avoid false negative PCR-results an Inhibition Control is amplified together in one tube with the specific sequence. **Tubes 5, 6, 7 und 8 contain lyticase.**

## 3. Kit contents

The kit contains sufficient reagents for 96 PCR reactions (12 samples):

1 x Premix	white cap
12 x Dye Strips (freeze-dried, tubes 5, 6, 7 und 8 contain lyticase)	tube strips
12 x Cap Strips	
1 x ddH <sub>2</sub> O (for PCR negative controls)	colourless cap

## 4. Storage conditions

Do **not** freeze the lyophilized Dye Strips.

The PCR reagents should be stored at 2 - 8 °C (35 – 46 °F). Keep Premix for storage at - 20 °C (- 4 °F) **after arrival**. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

## 5. Materials required, but not provided

### 5.1. Instruments

Real-time PCR instrument for low profile tubes

Centrifuge for tube strips

Pipettes

“Vortex”

### 5.2. Reagents and plastic ware

sterile ddH<sub>2</sub>O

sterile filter tips

## 6. PCR

### 6.1. PCR-Setup

Before every use thoroughly mix Premix and centrifuge briefly. Remove the film from the needed tube strips and pipette the PCR-Components. After pipetting close the tube strips with the provided cap strips.

#### Example for pipetting an 8-well strip (1 sample):

Mix the premix well before use and centrifuge briefly if necessary. Prepare a sterile 1.5 mL reaction tube and fill it with 127.5  $\mu$ L Premix (8.5 x 15  $\mu$ L). Remove the foil from the required tube strips and add 5  $\mu$ L ddH<sub>2</sub>O to tube 1 of the strip and pipette 15  $\mu$ L premix from the reaction tube (NTC). Then pipette 37.5  $\mu$ L of the corresponding sample into the PCR tubes, vortex briefly and pipette 20  $\mu$ L of the sample premix mixture into each of the tubes 2 - 8. After pipetting, close the tube strips with the cap strips provided and centrifuge briefly.

#### PCR-reaction setup per sample:

PCR-Components	Volume ( $\mu$ L)
Premix	15.0
Sample-DNA	5.0
Total volume	<b>20.0</b>

Pipetting schema per sample:

Tube 1*	Premix 15.0 µL	<b>ddH<sub>2</sub>O 5.0 µL</b> <b><u>Attention:</u> pipette no sample-DNA</b>
Tube 2*	Premix 15.0 µL	sample-DNA 5.0 µL
Tube 3*	Premix 15.0 µL	sample-DNA 5.0 µL
Tube 4*	Premix 15.0 µL	sample-DNA 5.0 µL
Tube 5*	Premix 15.0 µL	sample-DNA 5.0 µL
Tube 6*	Premix 15.0 µL	sample-DNA 5.0 µL
Tube 7*	Premix 15.0 µL	sample-DNA 5.0 µL
Tube 8*	Premix 15.0 µL	sample-DNA 5.0 µL

\* please notice the orientation of the tubes: they are marked with A - H

## 6.2 PCR-Program

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 35x
Cycling Annealing/ Elongation	30 sec	60 °C	

## MyGo User

### Software MyGo Pro/ ESR 3.5 / 3.6

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 35x
Cycling Annealing/ Elongation	30 sec	60 °C	

**Advanced Settings MyGo Pro**                      Integration time (s)                      0.7

**Advanced Settings MyGo Pro ESR**                      Integration time (s)                      1.0

## **7. Evaluation**

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

The negative control (NTC, Tube 1) has to be negative in the FAM channel and positive in the HEX-channel (Inhibition Control). Is the Inhibition Control in the HEX-channel negative and/or the results in the FAM-channel positive, PCR has to be repeated.

The Inhibition Control (Tube 5) has to be positive in the HEX-channel. Is the Inhibition Control in the HEX-channel negative inhibitors or high amount of DNA are in the sample reaction.

The Positive Control (Tube 8) has to be positive in the HEX-channel. If it is negative, PCR has to be repeated.

## PCR-Analysis:

please notice the orientation of the tubes: they are marked with A - H

	Kanäle	Spezies	Signal
<b>Tube 1</b>	FAM	NTC	negativ
	HEX	Inhibitionskontrolle	positiv
<b>Tube 2</b>	FAM	<i>Lactobacillus spp./ Pediococcus spp.</i>	
	HEX	Essigsäurebakterien	
<b>Tube 3</b>	FAM	<i>Leuconostoc mesenteroides</i>	
	HEX	<i>Oenococcus oeni</i>	
<b>Tube 4</b>	FAM	<i>Lactococcus lactis</i>	
	HEX	<i>Alicyclobacillus spp.</i>	
<b>Tube 5</b>	FAM	<i>Candida spp.</i>	
	HEX	Inhibitionskontrolle	positiv
<b>Tube 6</b>	FAM	<i>Dekkera bruxellensis</i>	
	HEX	<i>Saccharomyces cerevisiae</i>	
<b>Tube 7</b>	FAM	diastatic <i>Saccharomyces cerevisiae</i>	
	HEX	<i>Schizosaccharomyces pombe</i>	
	ROX	<i>Zygosaccharomyces spp.</i>	
<b>Tube 8</b>	FAM	<i>Pichia spp.</i>	
	HEX	PTC	positiv

### Example for interpretation:

If sample shows a positive FAM-signal in Tube 5, *Candida spp.* is identified. If the sample contains mixtures, you have positive signals in different channels.

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.