---

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Cycle 35 x

|                          |                 |                      |     |
|--------------------------|-----------------|----------------------|-----|
| <b>Advanced Settings</b> | <b>MyGo Pro</b> | Integration time (s) | 0.7 |
|--------------------------|-----------------|----------------------|-----|

|                          |                     |                      |     |
|--------------------------|---------------------|----------------------|-----|
| <b>Advanced Settings</b> | <b>MyGo Pro ESR</b> | Integration time (s) | 1.0 |
|--------------------------|---------------------|----------------------|-----|

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# QuickGEN PCR Kit

## Screening and differentiation of 8 wine spoilage bacteria and yeast

### 1. Intended use

Screening and differentiation of 8 bacteria and yeast in wine. The following bacteria and yeasts are differentiated:

*Lactobacilli / Pediococci*, *Oenococcus oeni*, Acetic acid bacteria, *Zygosaccharomyces bailii*, *Zygosaccharomyces rouxii*, *Dekkera bruxellensis*, *Saccharomyces cerevisiae* and yeast.

### 2. Test principle

The TaqMan<sup>®</sup> real-time PCR is based on hot-start PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ; ROX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

To avoid false negative PCR-results an Inhibition Control (ROX/DQ) is amplified together in one tube with the specific sequence. **Tubes 4, 5, 7 and 8 contain lyticase.**

### 3. Kit contents

The kit contains sufficient reagents for 96 PCR reactions (24 samples):

|                                                                      |                |
|----------------------------------------------------------------------|----------------|
| 1 x Premix                                                           | white cap      |
| 12 x Dye Strips (freeze-dried, tubes 4, 5, 7 und 8 contain lyticase) | tube strips    |
| 12 x Cap Strips                                                      |                |
| 1 x ddH <sub>2</sub> O (for PCR negative controls)                   | colourless cap |

#### **4. Storage conditions**

Do **not** freeze the lyophilized Dye Strips.

The PCR reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) **after arrival**. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

#### **5. Materials required, but not provided**

##### **5.1. Instruments**

Real-time PCR instrument for low profile tubes

Centrifuge for strips

Pipettes

“Vortex”

##### **5.2. Reagents and plastic ware**

sterile ddH<sub>2</sub>O

sterile filter tips

## 6. PCR

### 6.1. PCR-Setup

**Before every use thoroughly mix all PCR-components and centrifuge briefly.**

Remove the film from the needed tube strips and pipette the PCR-Components. After pipetting close the tube strips with the provided cap strips.

PCR-reaction setup:

| PCR-Components | volume (µL) |
|----------------|-------------|
| Premix         | 15.0        |
| Sample-DNA     | 5.0*        |
| Total volume   | <b>20.0</b> |

\* if using the Simplex Easy Wine or QuickGEN Yeast Sample Preparation Kit:  
add 2.5 µL DNA and 2.5 µL sterile ddH<sub>2</sub>O

Pipetting scheme:

|         |                |                                                                                       |
|---------|----------------|---------------------------------------------------------------------------------------|
| Tube 1* | Premix 15.0 µL | ddH <sub>2</sub> O 5.0 µL, <b>Attention:</b> pipette no sample-DNA (negative control) |
| Tube 2* | Premix 15.0 µL | ddH <sub>2</sub> O 5.0 µL, <b>Attention:</b> pipette no sample-DNA (positive control) |
| Tube 3* | Premix 15.0 µL | sample-DNA of sample 1 5.0 µL                                                         |
| Tube 4* | Premix 15.0 µL | sample-DNA of sample 1 5.0 µL                                                         |
| Tube 5* | Premix 15.0 µL | sample-DNA of sample 1 5.0 µL                                                         |
| Tube 6* | Premix 15.0 µL | sample-DNA of sample 2 5.0 µL                                                         |
| Tube 7* | Premix 15.0 µL | sample-DNA of sample 2 5.0 µL                                                         |
| Tube 8* | Premix 15.0 µL | sample-DNA of sample 2 5.0 µL                                                         |

\* Please notice the orientation of the tubes: they are marked with 1,2,3.....

## 6.2 PCR-Program

| Step                          | Time   | Temp. |            |
|-------------------------------|--------|-------|------------|
| Lyticase treatment            | 15 min | 37 °C |            |
| Initial denaturation of DNA   | 15 min | 95 °C |            |
| Cycling Denaturation          | 15 sec | 95 °C | Cycle 35 x |
| Cycling Annealing/ Elongation | 30 sec | 60 °C |            |

## MyGo Pro User

### Software MyGo Pro/ ESR 3.5 / 3.6

| Step                          | Time   | Temp. |                  |
|-------------------------------|--------|-------|------------------|
| Lyticase treatment            | 15 min | 37 °C |                  |
| Initial denaturation          | 15 min | 95 °C |                  |
| Cycling Denaturation          | 15 sec | 95 °C | Cycle 30* - 35 x |
| Cycling Annealing/ Elongation | 30 sec | 62 °C |                  |

**Advanced Settings MyGo Pro**                      Integration time (s)                      0.7

**Advanced Settings MyGo Pro ESR**                      Integration time (s)                      1.0

## 7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

The Negative Control (NTC, tube 1) has to be negative in the FAM-, and HEX-channel and positive in the ROX-channel (Inhibition Control). Is the Inhibition Control in the ROX-channel negative and/or the results in the FAM-/ HEX-channel positive, PCR has to be repeated.

The Inhibition Controls in the samples (tubes 4/7) have to be positive in the ROX-channel. Is the Inhibition Control in the ROX-channel negative or weak inhibitors or high amount of DNA are in the sample reaction. The sample should be repeated with a dilution.

The Positive Control (tube 2) has to be positive in the FAM-channel, if it is negative PCR, has to be repeated.

## PCR-Analysis:

please notice the orientation of the tubes: they are marked with 1,2,3.....

| Tube | Probe | FAM                                | HEX                             | ROX                             |
|------|-------|------------------------------------|---------------------------------|---------------------------------|
| 1    | NTC   | -                                  | -                               | Inhibition control              |
| 2    | PTC   | Positive control                   | -                               | -                               |
| 3    | 1     | <i>Lactobacillus / Pediococcus</i> | <i>Oenococcus oeni</i>          | Acetic acid bacteria            |
| 4    | 1     | yeast                              | <i>Zygosaccharomyces bailii</i> | Inhibition control              |
| 5    | 1     | <i>Dekkera bruxellensis</i>        | <i>Saccharomyces cerevisiae</i> | <i>Zygosaccharomyces rouxii</i> |
| 6    | 2     | <i>Lactobacillus / Pediococcus</i> | <i>Oenococcus oeni</i>          | Acetic acid bacteria            |
| 7    | 2     | yeast                              | <i>Zygosaccharomyces bailii</i> | Inhibition control              |
| 8    | 2     | <i>Dekkera bruxellensis</i>        | <i>Saccharomyces cerevisiae</i> | <i>Zygosaccharomyces rouxii</i> |

NTC: Negative Control    PTC: Positive Control

### Example for interpretation:

If sample 1 shows a positive FAM-signal in tube 5, *Dekkera bruxellensis* is identified. If the sample contains mixtures, you have positive signals in different channels

#### Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.