---

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## Wine Screening

### 1. Intended use

Qualitative detection and differentiation of bacterial and yeast contaminations in wine and grape must.

### 2. Test principle

The (TaqMan®) real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ; ROX/DQ; CY5/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

To avoid false negative PCR-results an Inhibition Control (HEX/DQ) is amplified together in one reaction vessel with the specific sequence. The system contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit). **The tubes contain lyticase.**

### 3. Kit contents

The kit contains sufficient reagents for 48 PCR reactions:

|                                                           |                |
|-----------------------------------------------------------|----------------|
| 1 x Premix                                                | white cap      |
| 48 x Dye Strips (freeze-dried, incl. IC-DNA and lyticase) | tube strips    |
| 1 x ddH <sub>2</sub> O                                    | colourless cap |
| 1 x Solution S                                            | green cap      |
| 1 x Control-DNA (freeze-dried)                            | yellow cap     |

## 4. Storage

**The Control-DNA is freeze-dried and has to be solved in Solution S prior to use (see 6.1).**

Do **not** freeze the lyophilized Dye Strips and lyophilized Control-DNA.

The PCR-reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) after opening. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

## 5. Materials required, but not provided

### 5.1. Instruments:

Real-time PCR instrument for low profile tubes

Centrifuge for 1.5 – 2.0 mL reaction vessels

Centrifuge for strips

Pipettes

“Vortex”

### 5.2. Reagents and plastic ware:

sterile ddH<sub>2</sub>O

sterile filter tips

**optional:** Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction)

## 6. PCR

### 6.1. PCR-Setup

*When using the kit for the first time, the freeze-dried Control-DNA has to be shortly centrifuged and carefully resolved:*

- add 55 µL sterile Solution S to the freeze-dried Control-DNA
- after 15 minutes mix well

**Before every use thoroughly mix all PCR-components and centrifuge briefly.**

PCR-reaction Setup:

| PCR-Components | volume (µL) |
|----------------|-------------|
| Premix         | 15.0        |
| Sample-DNA     | 5.0*        |
| Total volume   | <b>20.0</b> |

\* if using the Simplex Easy Wine Kit: add 2.5 µL DNA and 2.5 µL sterile ddH<sub>2</sub>O

1. Pipette 15 µL of the Premix into each PCR-tube, making sure that, prior to the first filling, the tip has been moistened.
2. Add 5 µL sample DNA, add 5 µL of the Control-DNA for the positive control reaction, add 5 µL of the extraction control and 5 µL of ddH<sub>2</sub>O for the negative control\* reaction. Use a fresh tip with each DNA filling.
3. Close the tubes immediately and centrifuge them shortly.
4. Place the tubes in the PCR-machine and start run.

**Very important: \* Please fill up the negative control with 5 µL ddH<sub>2</sub>O to avoid unspecific amplification.**

**Work swiftly to avoid warming up and keep away from light.**

## 6.2 PCR-Program

### Note:

For the use of UNG the thermal cycler programm has to be changed according to manufacturers` instructions.

| Step                          | Time   | Temp. |            |
|-------------------------------|--------|-------|------------|
| Lyticase treatment            | 15 min | 95 °C |            |
| Initial denaturation          | 15 min | 95 °C |            |
| Cycling Denaturation          | 15 sec | 95 °C | Cycle 40 x |
| Cycling Annealing/ Elongation | 25 sec | 58 °C |            |

## 7. Evaluation

The evaluation has to be made according to the data analysis programme recommended by the real-time instrument manufacturer.

**Lactobacillus/Pediococcus/Oenococcus oeni-DNA:** FAM-channel

**Acetic acid bacteria-DNA:** ROX-channel

**Yeast specific-DNA:** CY5-channel

**Inhibition Control-DNA:** HEX-channel


A sample is ***Lactobacillus/Pediococcus/Oenococcus oeni* positive**, if there is a detectable fluorescence increase in the **FAM-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the HEX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is **Acetic acid bacteria positive**, if there is a detectable fluorescence increase in the **ROX-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the HEX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is **yeast positive**, if there is a detectable fluorescence increase in the **CY5-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition Control in the HEX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). In the negative controls it has to be positive.

A sample is **negative**, if there is no detectable fluorescence increase in the different channels and the positive controls have a positive fluorescence signal. The negative controls show no amplification. The Inhibition Control in the HEX-channel has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

## analysis flowchart

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method. GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.