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Wild Yeast 1

1. Intended use

Detection of wild yeast in beverages..

The following yeasts are detected:

<i>Dekkera anomala</i>	<i>Dekkera bruxellensis</i>	<i>Dekkera custersiana</i>	<i>Dekkera naardenensis</i>	<i>Debaromyces hansenii</i>
<i>Hanseniaspora guillermondii</i>	<i>Hanseniaspora osmophila</i>	<i>Hanseniaspora uvarum</i>	<i>Issatchenkia orientalis</i>	<i>Kazachstania exigua</i>
<i>Kluyveromyces marxianus</i>	<i>Metschnikowia pulcherrima</i>	<i>Wickerhamomyces anomalus</i>	<i>Pichia fermentans</i>	<i>Pichia membranaefaciens</i>
<i>Saccharomyces diastaticus</i>	<i>Saccharomycodes ludwigii</i>	<i>Torulasporea delbrueckii</i>		

2. Test principle

The TaqMan® real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument. To avoid false negative PCR-results an Inhibition Control (HEX/DQ) is amplified together in one reaction vessel with the specific sequence. **The tubes contain Lyticase.**

3. Kit contents

The kit contains sufficient reagents for 48 reactions:

1 x Premix	white cap
6 x Dye Strips (lyophilisiert, tubes incl. IC-DNA)	tube strips
6 x Cap Strips	
1 x ddH ₂ O	colourless cap
1 x Control-DNA (lyophilisiert)	yellow cap

4. Storage conditions

The Control-DNA is lyophilized, it has to be solved in ddH₂O prior to use (see 6.1).

Do **not** freeze the lyophilized Dye Strips and lyophilized Control-DNA.

The PCR reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) **after arrival**. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR machine for low profile tubes

Centrifuge for 1.5 – 2.0 mL vials

Centrifuge for strips

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

6. PCR

6.1. PCR-Setup

When using the kit for the first time, the freeze-dried Control-DNA has to be shortly centrifuged and carefully resolved in ddH₂O:

- add 55 µL ddH₂O to the freeze-dried Control-DNA
- after 15 minutes mix well

Before every use thoroughly mix all PCR-components and centrifuge briefly.

Remove the film from the needed tube strips and pipette the PCR-Components. After pipetting close the tube strips with the provided cap strips.

PCR-reaction Setup:

PCR-Components	amount (µL)
Premix	17.5
Sample DNA	2.5
total volume	20.0

1. Pipette 17.5 µL of the Premix into each PCR-tube, making sure that, prior to the first filling, the tip has been moistened.
2. Add 2.5 µL sample DNA, for the PCR positive control, add 2.5 µL of the Control-DNA, add 2.5 µL of the extraction control and 2.5 µL of ddH₂O for the negative control* reaction. Use a fresh tip with each DNA filling.
3. Close the tubes immediately and centrifuge them shortly.
4. Place the tubes in the PCR-machine and start run.

Very important: * Please fill up the negative control with 2.5 µL ddH₂O to avoid unspecific amplification.

Work swiftly to avoid warming up and keep away from light.

6.2 PCR-Program

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	30 sec	62 °C	

MyGo Pro User

Software MyGo Pro 3.5 / 3.6 / ESR

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	30 sec	62 °C	

Advanced Settings MyGo Pro Integration time (s) 0.7

Advanced Settings MyGo Pro ESR Integration time (s) 1.0

7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.


Wild yeast 1-DNA: FAM-channel

Inhibition Control-DNA: HEX-channel

A sample is **Wild Yeast 1** positive, if there is a detectable fluorescence increase in the **FAM-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition-Control in the HEX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). For negative controls it has to be positive.

A sample is negative, if there is no detectable fluorescence increase in the FAM-channel and the positive controls have a positive fluorescence signal. The negative controls have to show no amplification in the FAM-channel. The Inhibition-Control in the HEX-channel has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

analysis flowchart

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method. GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.