---

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|                                       |                      |     |
|---------------------------------------|----------------------|-----|
| <b>Advanced Settings MyGo Pro</b>     | Integration time (s) | 0.7 |
| <b>Advanced Settings MyGo Pro ESR</b> | Integration time (s) | 1.0 |

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
## Wild Yeast 1

### 1. Intended use

Detection of wild yeasts in beverages.

The following yeasts are detected:

|                                                            |                                  |                               |                                |                                |
|------------------------------------------------------------|----------------------------------|-------------------------------|--------------------------------|--------------------------------|
| <i>Dekkera anomala</i>                                     | <i>Dekkera bruxellensis</i>      | <i>Dekkera custersiana</i>    | <i>Dekkera naardenensis</i>    | <i>Debaromyces hansenii</i>    |
| <i>Hanseniaspora guillermundii</i>                         | <i>Hanseniaspora osmophila</i>   | <i>Hanseniaspora uvarum</i>   | <i>Issotchenkia orientalis</i> | <i>Kazachstania exigua</i>     |
| <i>Kluyveromyces marxianus</i>                             | <i>Metschnikowia pulcherrima</i> | <i>Pichia anomala</i>         | <i>Pichia fermentans</i>       | <i>Pichia membranaefaciens</i> |
| <i>Saccharomyces cerevisiae</i><br><i>var. diastaticus</i> | <i>Saccharomycodes ludwigii</i>  | <i>Torulaspora delbrückii</i> |                                |                                |

### 2. Test principle

The TaqMan® real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument. To avoid false negative PCR-results an Inhibition Control (HEX/DQ) is amplified together in one reaction vessel with the specific sequence.

### 3. Kit contents

The kit contains sufficient reagents for 50 reactions:

|                                           |                    |
|-------------------------------------------|--------------------|
| 1 x Premix                                | white cap          |
| 1 x Dye Mix (lyophilisiert, incl. IC-DNA) | dark vial, red cap |
| 1 x ddH <sub>2</sub> O                    | colourless cap     |
| 1 x Control-DNA (lyophilisiert)           | yellow cap         |

#### 4. Storage conditions

**The Dye Mix and the Control-DNA are lyophilized, they have to be solved in ddH<sub>2</sub>O prior to use (see 6.1).**

Do **not** freeze the lyophilized Dye Mix and lyophilized Control-DNA.

The PCR reagents should be stored at 2 – 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) **after arrival**. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Mix contains the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

#### 5. Materials required, but not provided

##### 5.1. Instruments

Real-time PCR instrument with channels FAM and HEX

Centrifuge for 1.5 – 2.0 mL reaction tubes

Centrifuge for multiwell plates or stripes

Pipettes

“Vortex”

##### 5.2. Reagents and plastic ware

sterile ddH<sub>2</sub>O

sterile reaction vessels 1.5 – 2.0 mL

sterile optical tubes, plates, capillaries

sterile filter tips

## 6. PCR

### 6.1. PCR-Setup

*When using the kit for the first time, the freeze-dried kit components have to be shortly centrifuged and carefully resolved:*

- add 80 µL sterile ddH<sub>2</sub>O to freeze-dried Dye Mix
- add 55 µL sterile ddH<sub>2</sub>O to the freeze-dried Control-DNA
- after 15 minutes mix well

**Before every use thoroughly mix all PCR components and centrifuge briefly.**

PCR-reaction Setup:

| PCR-components | amount (µL) |
|----------------|-------------|
| Premix         | 13.5        |
| Dye Mix        | 1.5         |
| Sample-DNA     | 5.0*        |
| Total volume   | <b>20.0</b> |

\* if using the Simplex Easy DNA- or QuickGEN Yeast Sample Preparation kit: add 2.5 µL DNA and 2.5 µL sterile ddH<sub>2</sub>O

1. Prepare a mastermix by mixing Premix and Dye Mix.
2. Multiply said volumes with the number of PCR preparations including controls (positive control, negative control, extraction control), taking into account pipette reserves of approximately 5-10 %.
3. Divide 15 µL of the PCR-mastermix among the individual reaction vessels, making sure that, prior to the first filling, the tip of the pipette has been moistened.
4. Add 5 µL sample DNA, add 5 µL of the Control-DNA for the PCR positive control, add 5 µL of the extraction control and 5 µL of ddH<sub>2</sub>O for the negative control\* reaction. Use a fresh tip with each DNA filling.
5. Close the tubes immediately and centrifuge them shortly.
6. Place the tubes in the PCR-machine and start run.

**Very important: \* Please fill up the negative control with 5 µL ddH<sub>2</sub>O to avoid unspecific amplification.**

**Work swiftly to avoid warming up and keep away from light.**

## 6.2 PCR-Program

### 6.2.1 PCR-Program LC480

1. Click the button **tool** on the right side in the window **LightCycler 480 Software release 1.5.0. SP1**.
2. Click the button **Detection formats** at the left side of the menu bar.
3. Click **New** in the window **Detection formats** and name the experiment
4. Open the window **Filter Combination Selection** and choose the following filter combinations: 465-510 / 533-580.
5. Open the window **Selected Filter Combination** list and add the following amounts.

| Excitation Filter | Emission Filter | Name        | Melt factor | Quant factor | Max. Integ. Time |
|-------------------|-----------------|-------------|-------------|--------------|------------------|
| 465               | 510             | 465-510 Fam | 1           | 10           | 2                |
| 533               | 580             | 533-580 Hex | 1           | 10           | 2                |

6. **Close** the window.
7. Click the button **New Experiment** on the right side of the menu bar.
8. From the pull-down menu **Detection formats** choose the defined experiment, click the button **Customize** and check the detection formats. All of them have to be activated.
9. Click the button **ok**
10. Define the following program:

1. Programm Name: **Heat**

Cycles **1** Analysis Mode **None**

| Target<br>(°C) | Acquisition<br>Mode | Hold<br>(hh:mm:ss) | Ramp<br>Rate (°C/s) | Acquisitions<br>(per °C) | Sec Target<br>(°C) | Step size<br>(°C) | Step Delay<br>(cycles) |
|----------------|---------------------|--------------------|---------------------|--------------------------|--------------------|-------------------|------------------------|
| 37             | None                | 00:15:00           | 4.40                |                          | 0                  | 0                 | 0                      |
| 95             | None                | 00:15:00           | 4.40                |                          | 0                  | 0                 | 0                      |

2. Programm Name: **Ampli**

Cycles **40** Analysis Mode **Quantification**

| Target<br>(°C) | Acquisition<br>Mode | Hold<br>(hh:mm:ss) | Ramp<br>Rate (°C/s) | Acquisitions<br>(per °C) | Sec Target<br>(°C) | Step size<br>(°C) | Step Delay<br>(cycles) |
|----------------|---------------------|--------------------|---------------------|--------------------------|--------------------|-------------------|------------------------|
| 95             | None                | 00:00:10           | 4.40                |                          | 0                  | 0                 | 0                      |
| 62             | Single              | 00:00:30           | 2.20                |                          | 0                  | 0                 | 0                      |

3. Programm Name : **Cool**

Cycles **1** Analysis Mode **None**

| Target<br>(°C) | Acquisition<br>Mode | Hold<br>(hh:mm:ss) | Ramp<br>Rate (°C/s) | Acquisitions<br>(per °C) | Sec Target<br>(°C) | Step size<br>(°C) | Step Delay<br>(cycles) |
|----------------|---------------------|--------------------|---------------------|--------------------------|--------------------|-------------------|------------------------|
| 40             | None                | 00:00:20           | 2.20                |                          | 0                  | 0                 | 0                      |

**Optional:** Saving the program as **run template**. The template button allows to select and apply a template to the currently open object and to save the currently open object as a template. Click the clamp beside the button **Apply Template**. Click **Save as template** and save the file in **templates**.

11. Click the button **Subset editor** on the left side

12. Click the button **+** and **New Subset 1** appears

13. Mark the wells in the **New Subset 1 Settings** window

14. Click the button **Apply**

15. Click the button **Sample editor** on the left side

16. **Very important:** Activate In the window **Step 1 Select Workflow: Abs.Quant**

17. Choose the subset **New Subset** in the window **Step 2 Select Samples**

18. Define your probes in the **Sample table**

19. Click the button **Experiment** and start the run with **Start run**

## 6.2.2 PCR-Program for other real-time machines

The Premix contains no ROX. This has to be done according to the settings of the real-time machine. **ABI 7500:** Under “Assign Targets and Samples” in ‘Select the dye to use as the passive reference’ choose „none“.

| Step                          | Time   | Temp. |            |
|-------------------------------|--------|-------|------------|
| Lyticase treatment            | 15 min | 37 °C |            |
| Initial denaturation          | 15 min | 95 °C |            |
| Cycling Denaturation          | 10 sec | 95 °C | Cycle 40 x |
| Cycling Annealing/ Elongation | 30 sec | 62 °C |            |

## MyGo Pro User

### Software MyGo Pro 3.5 / 3.6 / ESR

| Step                          | Time   | Temp. |            |
|-------------------------------|--------|-------|------------|
| Lyticase treatment            | 15 min | 37 °C |            |
| Initial denaturation          | 15 min | 95 °C |            |
| Cycling Denaturation          | 10 sec | 95 °C | Cycle 40 x |
| Cycling Annealing/ Elongation | 30 sec | 62 °C |            |

|                                   |                      |     |
|-----------------------------------|----------------------|-----|
| <b>Advanced Settings MyGo Pro</b> | Integration time (s) | 0.7 |
|-----------------------------------|----------------------|-----|

|                                       |                      |     |
|---------------------------------------|----------------------|-----|
| <b>Advanced Settings MyGo Pro ESR</b> | Integration time (s) | 1.0 |
|---------------------------------------|----------------------|-----|

## 7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

For LC480: Activate Colour Compensation

**Wild yeast 1-DNA:** FAM-channel (LC480: 465-510 nm)


**Inhibition Control-DNA:** HEX-channel (LC480: 533-580 nm)

A sample is **Wild Yeast 1** positive, if there is a detectable fluorescence increase in the **FAM-channel** and the negative controls show no amplification. The positive controls should have a positive fluorescence signal. The Inhibition-Control in the HEX-channel may be positive or negative (depending on the amount of DNA or inhibitors in the sample reaction). For negative controls it has to be positive.

A sample is negative, if there is no detectable fluorescence increase in the FAM- channel and the positive controls have a positive fluorescence signal. The negative controls have to show no amplification in the FAM-channel. The Inhibition-Control in the HEX-channel has to be positive in the sample and in the negative controls, a false negative result due to inhibitory effects is then excluded.

# analysis flowchart

## LC480: analysis after activating Colour Compensation

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.