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Yeast Differentiation

1. Intended use

Differentiation of beer spoilage yeasts in beer and beer mixing drinks.
The following yeast are differentiated:

Saccharomyces exiguus

Saccharomyces cerevisiae var. *diastaticus*

Saccharomyces bayanus / *pastorianus*

Saccharomycodes ludwigii

Torulaspora delbrückii

Klyveromyces marxianus

Debaromyces hansenii

<i>Dekkera</i> spp.:	<i>D. anomala</i>
	<i>D. custersiana</i>
	<i>D. bruxellensis</i>
	<i>D. naardenensis</i>
	<i>D. nanus</i>

<i>Rhodotorula</i> spp.:	<i>R. mucilaginosa</i>
	<i>R. glutinis</i>
	<i>R. graminis</i>
	<i>R. toruloides</i>
	<i>R. bacarum</i>

<i>Candida</i> spp.:	<i>C. albicans</i>
	<i>C. parapsilosis</i>
	<i>C. glabrata</i>
	<i>C. tropicalis</i>
	<i>C. sake</i>
	<i>C. intermedia</i>
	<i>C. pulcherrima</i>
	<i>Klyveromyces marxianus</i>

<i>Pichia</i> spp.:	<i>P. anomala</i>
	<i>P. fermentans</i>
	<i>P. membranaefaciens</i>

<i>Hanseniaspora</i> spp.:	<i>H. guilliermondii</i>
	<i>H. uvarum</i>
	<i>H. osmophila</i>

2. Test principle

The TaqMan® real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument. To avoid false negative PCR-results an Inhibition Control is amplified together in one tube with the specific sequence. **The tubes contain lyticase.**

3. Kit contents

The kit contains sufficient reagents for 96 PCR reactions (12 samples):

1 x Premix	white cap
12 x Dye Strips (freeze-dried, tubes contain lyticase)	tube strips
12 x Cap Strips	
1 x ddH ₂ O (for PCR negative controls)	colourless cap

4. Storage conditions

Do **not** freeze the lyophilized Dye Strips.

The PCR reagents should be stored at 2 - 8 °C (35 – 46 °F).

Keep Premix for storage at - 20 °C (- 4 °F) **after arrival**. Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The Dye Strips contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 12 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR instrument for low profile tubes

Centrifuge for strips

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

6. PCR

6.1. PCR-Setup

Before every use thoroughly mix Premix and centrifuge briefly.

Remove the film from the needed tube strips and pipette the PCR-Components. After pipetting close the tube strips with the provided cap strips.

Pipetting example for an 8-well strip (1 sample):

Mix the premix well before use and centrifuge briefly if necessary. Prepare a sterile 1.5 mL reaction tube and fill it with 149 μL Premix (8.5 x 17.5 μL). Remove the foil from the required tube strips and add 2.5 μL ddH₂O to tube 1 of the strip and pipette 17.5 μL premix from the reaction tube (NTC). Then pipette 19 μL of the corresponding sample into the PCR tubes, vortex briefly and pipette 20 μL of the sample premix mixture into each of the tubes 2 - 8. After pipetting, close the tube strips with the cap strips provided and centrifuge briefly.

PCR-reaction setup:

PCR-Components	Volume (μL)
Premix	17.5
Sample-DNA	2.5
Total volume	20.0

Pipetting schema per sample:

Tube 1*	Premix 17.5 μL	ddH₂O 2.5 μL <u>Attention:</u> pipette no sample-DNA
Tube 2*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 3*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 4*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 5*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 6*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 7*	Premix 17.5 μL	sample-DNA 2.5 μL
Tube 8*	Premix 17.5 μL	sample-DNA 2.5 μL

* please notice the orientation of the tubes: they are marked with A - H

6.2 PCR-Program

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30 ¹ -35x
Cycling Annealing/ Elongation	20-30 sec ²	63 °C	

¹ for preenriched samples program 30 cycles

² please program 20 sec, if the real-time PCR machine allows it

MyGo User

Software MyGo Pro/ ESR 3.5 / 3.6

Step	Time	Temp.	
Lyticase treatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	15 sec	95 °C	Cycle 30*-35x
Cycling Annealing/ Elongation	30 sec	63 °C	

* for preenriched samples program 30 cycles

Advanced Settings MyGo Pro Integration time (s) 0.7

Advanced Settings MyGo Pro ESR Integration time (s) 1.0

7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer.

The negative control (NTC, Tube 1) has to be negative in the FAM channel and positive in the HEX-channel (Inhibition Control). Is the Inhibition Control in the HEX-channel negative and/or the results in the FAM-channel positive, PCR has to be repeated.

The Inhibition Control (Tube 7) has to be positive in the HEX-channel. Is the Inhibition Control in the HEX-channel negative inhibitors or high amount of DNA are in the sample reaction.

The Positive Control (Tube 8) has to be positive in the HEX-channel. If it is negative, PCR has to be repeated.

PCR-Analysis:

please notice the orientation of the tubes: they are marked with A - H

	Channels	Species	Signal
Tube 1	FAM	NTC	negative
	HEX	Inhibition Control	positive
Tube 2	FAM	<i>Rhodotorula spp.</i>	
	HEX	<i>Saccharomyces exiguus</i>	
Tube 3	FAM	<i>Candida spp.</i>	
	HEX	<i>Saccharomyces cerevisiae var. diastaticus</i>	
Tube 4	FAM	<i>Saccharomycodes ludwigii</i>	
	HEX	<i>Debaromyces hansenii</i>	
Tube 5	FAM	<i>Torulaspora delbrückii</i>	
	HEX	<i>Saccharomyces bayanus / pastorianus</i>	
Tube 6	FAM	<i>Kluyveromyces marxianus</i>	
	HEX	<i>Hanseniaspora spp.</i>	
Tube 7	FAM	<i>Dekkera spp.</i>	
	HEX	Inhibition Control	positive
Tube 8	FAM	<i>Pichia spp.</i>	
	HEX	PTC	positive

Example for interpretation:

If sample shows a positive FAM-signal in Tube 3, *Candida spp.* is identified. If the sample contains mixtures, you have positive signals in different channels.

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.