Simplex[®] Easy Spin

DNA-Extraktion aus *Alicyclobacillus* spp.

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Simplex[®] Easy Spin

1. Intended use

Fast *Alicyclobacillus* spp. DNA extraction from e.g. juices, concentrates and tomato products.

2. Content

The kit contains sufficient reagents for 50 preparations:

- 1 x Lysis buffer SES
- 1 x Binding buffer
- 3 x Elution buffer
- 50 x Columns incl. Collection tubes

3. Materials required, but not provided

3.1. Instruments

Heating block or water bath 37 °C – 95 °C (98.6 °F – 203 °F)

Centrifuge for 1.5 – 2.0 mL reaction tubes

Pipettes

3.2. Reagents and plastic ware

sterile safe-seal reaction tubes 1.5 – 2.0 mL

suitable filter tips

single use gloves

Ethanol (96-100 %)

4. Storage

All reagents can be stored at room temperature (18-25 °C, 64.4 °F – 77 °F) and are stable up to one year. If there is any precipitate present in the buffers, warm the buffer to 37 °C (98.6 °F) to dissolve the precipitate before use.

Attention: The binding buffer contains Guanidin hydrochloride.

5. Risk and safety phrases (R & S)

Component	Hazard components	Hazard symbol	Risk phrases	Safety phrases
Binding buffer	Guanidin hydrochloride 50 – 66%	 Warnung	R 22-36/38	S26-37/39

Risk phrases

R 22 Harmful if swallowed
R 36/38 Irritating to eyes and skin

Safety phrases

S 26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S 37/39 Wear suitable gloves and eye/face protection

6. Indications of deterioration of reagents

In case of accurate handling deterioration unknown.

7. DNA-Extraction from *Alicyclobacillus* spp.

7.1 Sample preparation according to IFU-method 12

For non filterable solutions:

- Add 90 mL BAT-Bouillon to 10 mL product
- For inactivation of vegetative cells and activation of spores incubate the sample for 10 min at 80 °C and cool down to 45 °C (113 °F)
- Incubate the sample (3-7 days, aerobe at 45 °C +/- 1 °, 113 °F)

For filterable solutions:

- For inactivation of vegetative cells and activation of spores incubate 100 mL sample for 10 min at 80 °C (176 °F) and cool down to 45 °C (113 °F)
- Membrane filtration of the sample (0.45 µm filter)
- Incubate the filter in BAT-Bouillon (3-7 days, aerobe at 45 °C +/- 1 °, 113 °F)

7.2 DNA-Extraction protocol for non filterable solutions (e.g. concentrates)

1. add 2 mL of the pre-enriched sample in a 2.0 mL safe-lock reaction tube
2. centrifuge 5 min at 400 x g
3. transfer 1 mL supernatant in a 1.5 mL safe-seal reaction tube
4. centrifuge 3 min at 15500 x g
5. discard the supernatant
6. resuspend the pellet in 250 µL Lysis buffer SES
7. incubate 30 min at 95 °C (203 °F)
8. add 240 µL Binding buffer and 160 µL Ethanol
9. vortex 10 sec
10. transfer the lysate to the column placed in a collection tube
11. centrifuge 2 min at 11000 x g
12. transfer the column in a new 1.5 mL safe-seal reaction tube
13. add 100 µL Elution buffer (preheated to 95 °C, 203 °F) onto the column
14. incubate for 15 min at room temperature
15. centrifuge 1 min at 11000 x g
16. use 5 µL of the eluted DNA for PCR

7.3 DNA-Extraction protocol for filterable solutions

1. add 2 mL of the pre-enriched sample in a 2.0 mL safe-lock reaction tube
2. centrifuge 5 min at 15500 x g
3. discard the supernatant
4. resuspend the pellet in 250 µL Lysis buffer SES
5. incubate 30 min at 95 °C (203 °F)
6. add 240 µL Binding buffer and 160 µL Ethanol
7. vortex 10 sec
8. transfer the lysate to the column placed in a collection tube
9. centrifuge 2 min at 11000 x g
10. transfer the column in a new 1.5 mL safe-seal reaction tube
11. add 100 µL Elution buffer (preheated to 95 °C, 203 °F) onto the column
12. incubate for 5 min at room temperature
13. centrifuge 1 min at 11000 x g
14. use 5 µL of the eluted DNA for PCR

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.