Simplex[®] Easy Spin Plus

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Simplex[®] Easy Spin Plus

1. Intended use

Alicyclobacillus spp. DNA-extraction from e.g. juices, concentrates and tomato products without preenrichment.

2. Content

The kit contains sufficient reagents for 50 preparations:

- 1 x Lysis buffer SESP
- 1 x Binding buffer
- 1 x Washing buffer
- 2 x Elution buffer
- 50 x Columns incl. Collection tubes

3. Materials required, but not provided

3.1. Instruments

Heating block or water bath 37 °C – 95 °C (98.6 °F – 203 °F)

Centrifuge for 1.5 – 2.0 mL reaction tubes

Pipettes

3.2. Reagents and plastic ware

Polycarbonate-filter, 50 mm Ø, 0.4 µm

sterile safe-lock reaction tubes 1.5 – 2.0 mL

suitable filter tips

single use gloves

Ethanol (96-100 %)

sterile water

4. Storage

All reagents can be stored at room temperature (18-25 °C, 64.4 °F – 77 °F) and are stable up to one year. If there is any precipitate present in the buffers, warm the buffer to 37 °C (98.6 °F) to dissolve the precipitate before use.

Attention: The binding buffer contains Guanidin hydrochloride.

5. Risk and safety phrases (R & S)

Component	Hazard components	Hazard symbol	Risk phrases	Safety phrases
Binding buffer	Guanidin hydrochloride 50 – 66%	 Warnung	R 22-36/38	S26-37/39

Risk phrases

R 22 Harmful if swallowed
R 36/38 Irritating to eyes and skin

Safety phrases

S 26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S 37/39 Wear suitable gloves and eye/face protection

6. Indications of deterioration of reagents

In case of accurate handling deterioration unknown.

7. Preliminary preparations

Add 18 mL 100% ETOH to the washing buffer and mix

8. DNA-Extraction from *Alicyclobacillus spp.* without preenrichment

DNA-Extraction protocol for filterable solutions

1. filter 0.33 L - 0.5 L of liquid through a 0.4 µm polycarbonate filter
 - a. The filterable quantity depends on the sample matrix
2. add 50 mL filtration buffer (Ref. No. Q704) and filter off
3. add 50 mL sterile water and filter off
4. place the polycarbonate filter in a 1.5 mL sterile safe-lock reaction tube
5. add 300 µL Lysis buffer
 - a. **Very important:** the filter must be completely covered
6. vortex the sample 30 sec.
7. incubate the sample for 20 min. min at 95 °C (203 °F)
8. remove the sample from the heating block, add 240 µL Binding buffer and 160 µL Ethanol
(*alternative: prepare a mixture of Binding buffer and 100% EtoH (ratio 1,5 Vol : 1,0 Vol) and pipette 400 µl of this mixture to the sample*)
9. vortex the sample 10 sec.
10. place the column in an empty 1.5 mL reaction vessel and transfer the complete sample volume (approx. 700 µL) onto the column
11. incubate the sample for 5 min. min at room temperature
12. centrifuge 2 min at 10000 x g
13. transfer the column in a new 1.5 mL reaction vessel
14. pipette 400 µL Washing buffer onto the column
15. centrifuge for 2 min at 10000 x g
16. place the column in a new 1.5 mL safe-lock reaction tube
17. add 50 µL Elution buffer (preheated to 95 °C, 203 °F) onto the column
18. incubate for 5 min at room temperature
19. centrifuge 1 min at 10000 x g
20. discard the column
21. use 5 µL of the eluted DNA for PCR