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Simplex[®] Easy Spin Plus

1. Intended use

Alicyclobacillus spp. DNA-extraction from e.g. juices, concentrates and tomato products without preenrichment..

2. Content

The kit contains sufficient reagents for 50 preparations:

- 1 x Lysis buffer SESP
- 1 x Binding buffer
- 1 x Washing buffer
- 3 x Elution buffer
- 50 x Columns incl. Collection tubes

3. Materials required, but not provided

3.1. Instruments

- Heating block or water bath 37 °C – 95 °C (98.6 °F – 203 °F)
- Microcentrifuge for 1.5 – 2.0 mL reaction tubes
- Pipettes: 10 – 100 µL, 100 – 1000 µL, 0.5 – 10 µL

3.2. Reagents and plastic ware

- sterile safe-seal reaction tubes 1.5 – 2.0 mL
- suitable filter tips
- single use gloves
- optional: lysozyme
- Ethanol (96-100 %)

4. Storage

All reagents can be stored at room temperature (18-25 °C, 64.4 °F – 77 °F) and are stable up to one year. If there is any precipitate present in the buffers, warm the buffer to 37 °C (98.6 °F) to dissolve the precipitate before use.

Attention: *The binding buffer contains Guanidin hydrochloride.*

5. Risk and safety phrases (R & S)

Component	Hazard components	Hazard symbol	Risk phrases	Safety phrases
Binding buffer	Guanidin hydrochloride 50 – 66%	 Warnung	R 22-36/38	S26-37/39

Risk phrases

R 22 Harmful if swallowed
R 36/38 Irritating to eyes and skin

Safety phrases

S 26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
S 37/39 Wear suitable gloves and eye/face protection

6. Indications of deterioration of reagents

In case of accurate handling deterioration unknown.

7. DNA-Extraction from *Alicyclobacillus spp.*

Test preparations: Add 18 mL 100% ETOH to the washing buffer and mix

DNA-Extraction protocol for filterable solutions (e.g. concentrates)

- filter 0.33 L - 0.5 L of liquid through a 0.4 µm polycarbonate filter
The filterable quantity depends on the sample matrix
- add 50 mL filtration buffer (Ref. No. Q704) and filter off
- add 50 mL sterile water and filter off
- place the polycarbonate filter in a 1.5 mL sterile safe-seal reaction vial
- add 300 µL Lysis buffer
Very important: the filter must be completely covered
- vortex the sample 30 sec.
- incubate the sample for 20 min. min at 95 °C (203 °F)
- remove the sample from the heating block, add 240 µL Binding buffer and 160 µL Ethanol
(*alternative: prepare a mixture of Binding buffer and 100% EtoH (ratio 1,5 : 1,0 V) and pipette 400 µl of this mixture to the sample*)
- vortex the sample 10 sec.
- place the column in an empty 1.5 mL reaction vessel and transfer the complete sample volume (approx. 700 µL) onto the column
- incubate the sample for 5 min. min at room temperature
- centrifuge 2 min at 10000 x g
- transfer the column in a new 1.5 mL safe-seal reaction tube
- pipette 400 µL Washing buffer onto the column
- centrifuge for 2 min at 10000 x g
- place the column in a new 1.5 mL reaction tube
- add 50 µL Elution buffer (preheated to 95 °C, 203 °F) onto the column
- incubate for 5 min at room temperature
- centrifuge 1 min at 10000 x g
- discard the column
- use 5 µL of the eluted DNA for PCR

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.