Vitamin B7 (Biotin)

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VitaFast® Vitamin B7 (Biotin)

Brief information

Easy to use microbiological microtiter plate test for the quantification of total biotin (enriched and natural biotin) in food and pharmaceutical products. All required reagents and standard are contained in the test kit. The kit is sufficient for 96 determinations incl. standard. Evaluation requires a microtiter plate reader (610 - 630 nm, alternatively at 540 - 550 nm).

Sample preparation

Liquid samples (added biotin):
sterile filtration and dilution

Solid samples (added biotin):
sample homogenization, extraction, centrifugation
and dilution

Liquid and solid samples (native and added biotin):
sample homogenization, enzymatic treatment,
extraction, centrifugation and dilution

Time requirement: Test conductionapprox. 60 min
Evaluation.....2 min

Incubation: 44 - 48 h in the dark at 37 °C (98.6 °F)

AOAC-RI Certified matrices:
infant formula, cereals, pills, powders, juice and milk

Standard range: 0.08 - 0.72 µg / 100 g (ml)

Limit of Quantitation: 0.08 µg / 100 g (ml)

Limit of Detection: 0.013 µg / 100 g (ml)

Recovery: 90 - 105 %

Intra-assay CV for standards: < 10 %

Inter-assay CV for standards: < 10 %

1. Principle of the test

The VitaFast® Biotin microtiter plate test is a microbiological method for the quantitative determination of total biotin (added and natural biotin) in food and pharmaceutical products. The microbiological test system is in accordance with international norms.

Biotin is extracted from the sample and the extract is diluted. The biotin assay-medium and the diluted extract are pipetted into the wells of a microtiter plate which is coated with *Lactobacillus plantarum*. The growth of *Lactobacillus plantarum* is dependent on the supply of biotin. Following the addition of biotin as a standard or as a compound of the sample, the bacteria grow until the vitamin is consumed. The incubation is done in the dark at 37 °C (98.6 °F) for 44 - 48 h.

The intensity of metabolism or growth of *Lactobacillus plantarum* in relation to the extracted biotin is measured as turbidity and compared to a standard curve. The measurement is done using a microtiter plate reader at 610 - 630 nm (alternatively at 540 - 550 nm).

2. Reagents provided

Each test kit contains sufficient materials for 96 determinations incl. standards. Each test kit contains:

- 1 x microtiter plate with 96 wells, coated with *Lactobacillus plantarum*
- 3 x redist., sterile water (30 ml) for the preparation of assay medium, standards and dilution of extracted samples
- 3 x biotin assay medium (solid)
- 3 x biotin standard (solid)
- 3 x adhesive foil (1 complete and 2 halves of adhesive foils, sufficient for 3 test runs)
- 1 x additional holder for microtiter strips

Remark:

No quality guarantee can be assumed after expiration of the shelf life.

3. Required reagents and instruments, not provided

- **sterile bench** (sterile working for recommended)
- **microtiter plate reader** 610 - 630 nm (540 - 550 nm)
- incubator with dark chamber, 37 °C (98.6 °F)
- water bath heatable to 95 °C (203 °F)
- pH meter
- centrifuge > 8000 x g (if the sample cannot be filtrated)
- sterile tips for graduated micropipette 20 - 200 µl and 100 - 1000 µl,
- sterile graduated centrifuge vials with screw cap (15 and 50 ml) and sterile reaction vials 1.5 or 2.0 ml
- 500 ml screw glass jar and volumetric flasks (100 and 1000 ml), 100 ml beaker
- sterile filters polyethersulfon 0.2 µm with syringe
- redist. or deionized water for sample extraction
- HCl 1.0 mol / l and 0.1 mol / l
- NaOH 1.0 mol / l und 0,1 mol / l
- NaOH concentrated (40 g NaOH ad 100 ml redist. or deionized water)
- NaOH 2 mol / l (8 g NaOH ad 100 ml redist. or deionized water)

Reagents for the determination of native biotin

- taka diastase, aspergillus oryzae (e. g. Fluka 86250)
- H₂SO₄ 1 mol / l and 0.5 mol / l
- citrate buffer pH 4.5 (weigh 1.5 g citric acid monohydrate (e.g. Roth 5110.3) in a 100 ml beaker with magnetic stirrer; solve the citric acid with about 50 ml redist. or. deionised water under stirring; thereafter add 12 ml NaOH 1 mol / l (or 0.48 g NaOH); the pH should be 4.5 (correct with HCl 0.1 mol / l); transfer the solution quantitatively with redist. or deionised water in a 100 ml volumetric flask and fill up with redist. or deionised water to the mark; the buffer can be stored 3 days at 2 - 8 °C (35.6 - 46.4 °F)

4. Warning and precautions for the user

- the assay medium could evoke irritations of mucosa, eyes and skin
- after running the test the strips used must be disposed of according to regulations (e.g. autoclaved)

5. Storage instructions

Store the kit / reagents at 2 - 8 °C (35.6 - 46.4 °F).

Use the prepared reagents (standard, medium) immediately and reject them after the assay.

6. Sample preparation

For the determination of **added biotin** in nutrient enriched solid samples, a hot, alkali extraction is usually sufficient. Liquid samples can be used directly after sterile filtration and dilution with sterile water from the test kit.

For the determination of the **total biotin (native and added)** content of most foods (with the exception of cereals), the sample has to be treated with enzyme. For cereal samples with native biotin content, the sample must be extracted by autoclaving.

Samples should be stored protected from light at 4 °C (39.2 °F). Standard and samples should be run in triplicate. **Unknown sample matrices** should be analyzed **with two dilutions** of the sample extract. Sample extracts have to be used within one day and should be stored in the dark until analyzing.

Sample extraction is carried out with 1 g (ml) homogenized sample in 40 ml redist. or deionized water or extraction solution. This equals a **sample extraction dilution factor of 40**. This factor is already included in the standard curve (see Quality Assurance Certificate). **For low biotin concentrations a sample weight of up to 5 g (ml) can be used** (this has to be considered in the evaluation).

Note: when an enzymatic sample preparation is done, it is recommended to carry out a reagent blank to ensure that no biotin is present in the used reagents (enzymes).

Spike recovery studies may also be performed in order to identify the presence of any inhibitory compounds, which may negatively affect the result. An aliquot of spiking solution containing a known amount of biotin should be added to the sample prior to extracting. Spiking levels should generally be approximately 100 % of the expected vitamin content of the sample, in order to obtain an accurate result. When low spike recoveries are obtained, this may indicate the presence of an inhibitory compound. The sample should be diluted further where possible to eliminate these inhibitory effects. Spiking standards are available from R – Biopharm.

Only sterile sample extracts or sterile dilutions thereof should be pipetted onto the microtiter plate. Dilutions have to be prepared with sterile water from the test kit. **Therefore after the sample extraction sterile working conditions and sterile consumables are necessary.** A sterile filtration of the sample or the sample extract is **always** necessary for:

- samples like fruit juices and fitness drinks, which are not heated during sample extraction (except when the sample is heated 30 min at 95 °C (203 °F) in a water bath
- samples containing herbs and spices as well as honey and tea
- vitamin mixes, premixes, tablets (highly enriched samples, see 6.3) (except when the sample is heated 30 min at 95 °C (203 °F) in a water bath
- samples with a low vitamin concentration that are highly coloured; (the filtration step eliminates the colouring)
- if filtration is not possible due to solid particles or due to cloudiness, centrifugation should be carried out before the sterile filtration step (greater than 8000 x g for 5 min)

The sterile filtration of the sample is not necessary if the sample extraction is carried out at 95 °C (203 °F) for 30 min. Nevertheless the extracted samples have to be diluted with sterile water from the test kit (**the assay medium has to be always filtered**).

Example for the dilution factors of the sample extract

Solid sample with a labeled concentration of 80 µg / 100 g

The dilution of the extract should be in the middle of the standard curve. Therefore, the concentration is divided by standard 2.

Calculation:

$80 \mu\text{g} / 0.24 \mu\text{g} = 333$
→ dilution factor 300 (1 : 300)

Dilution steps:

- a) 1 : 10 → 0.1 ml sample extraction solution + 0.9 ml sterile water from the test kit
- b) 1 : 10 → 0.1 ml from a) + 0.9 ml sterile water from the test kit
- c) 1 : 3 → 0.2 ml from b) + 0.4 ml sterile water from the test kit

6.1. Added biotin in liquid samples (multivitamin juices, fitness drinks)

Add 1 ml sample into a 50 ml sterile centrifuge vial and fill up exactly to 40 ml with redist. or deionized water. Thereafter shake, sterile filter (or heat the sample 30 min at 95 °C in a water bath, thereafter chill down quickly below 30 °C (86 °F)) and, depending on the concentration range, further dilute in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

6.2. Added biotin in fruit gums and candies

Weigh about 15 - 20 g fruit gums or candies in a 50 ml sterile centrifuge vial, add about 40 ml redist. or deionized water and solve the sample at 95 °C (203 °F) in water bath. Chill down quickly to below 30 °C (86 °F), transfer the extraction solution with redist. or deionized water quantitatively into a 100 ml volumetric flask and fill up to the mark with redist. or deionized water. Transfer the volume corresponding to 1 g of test material into a 50 ml centrifuge vial and fill up to exactly 40 ml with redist. or deionized water, shake, sterile filter (or heat the sample 30 min at 95 °C in a waterbath, thereafter chill down quickly below 30 °C (86 °F)) and, depending on the concentration range, further dilute in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

Example: sample weight 17 g fruit gums

calculation: $100 \text{ ml} / 17 \text{ g sample weight} = 5.88 \text{ ml} / \text{g}$

1 g sample is contained in 5.88 ml. Transfer 5.88 ml into a 50 ml centrifuge vial and continue as described above.

6.3. Added biotin in capsules, pills and vitamin mixes

First the capsule or pill weight has to be determined (average of 5 - 10 capsules or pills). Then pills are homogenized in a mortar or mixer. Capsules are open and extracted.

6.3.1 Sample preparation with 1 g sample size and pre-extraction

Weigh 1 g of pills, vitamin mix, premix or cut open capsule (take weight into consideration) into a 500 ml screw glass jar, add about 400 ml redist. or deionized water, shake well. Adjust pH to 8.0 +/- 0.2 with HCl or NaOH. Extract for 30 minutes at 95 °C (203 °F) in a water bath. During extraction the glass jar has to be shaken well at least five times. Chill down quickly to below 30 °C (86 °F), transfer the extraction solution with redist. or deionized water quantitatively into a 1000 ml volumetric flask and fill up to the mark with redist. or deionized water. Transfer 1 ml into a 50 ml sterile centrifuge vial and fill up exactly to 40 ml with redist. or deionized water. Thereafter shake, sterile filter (or heat the sample 30 min at 95 °C in a water bath, thereafter chill down quickly below 30 °C (86 °F) and, depending on the concentration range, further dilute in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

Attention: For the calculation of the result, the pre-dilution factor of 1000 has to be considered (1 g up to 1000 ml). The dilution step 1 ml up to 40 ml is already included in the standard curve.

6.3.2 Sample preparation with 0.2 sample size

Weigh 0.2 g of pills, vitamin mix, premix or cut open capsule (take weight into consideration) into a 50 ml centrifuge vial, add about 20 ml redist. or deionized water, shake well, adjust pH 8.0 +/- 0.2 with HCl or NaOH and fill up to exactly 40 ml with redist. or deionized water. Extract for 30 minutes at 95 °C (203 °F) in a water bath. During extraction the vial has to be shaken well at least five times. It is important to make sure that the centrifuge vial is tightly closed. Chill down quickly to below 30 °C (86 °F). Thereafter centrifuge and, depending on the concentration range, further dilute the clear supernatant in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

Note: For the evaluation the weight of sample has to be considered.

6.4. Added biotin in cereals, baby food and flour

Weigh 1 g (ml) homogenized sample into a 50 ml centrifuge vial, add about 30 ml redist. or deionized water, adjust to pH 8.0 +/- 0.2 with HCl or NaOH, shake well and fill up to exactly 40 ml with redist. or deionized water. Extract for 30 minutes at 95 °C (203 °F) in a water bath. During extraction the vial has to be shaken well at least five times. It is important to make sure that the centrifuge vial is tightly closed. Chill down quickly to below 30 °C (86 °F). Thereafter centrifuge and, depending on the concentration range, further dilute the clear supernatant in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

6.5. Total vitamin content (native and added biotin) in milk products and infant formula

To extract the bound, native biotin as well or to determine it in non fortified samples, the sample has to be extracted with enzyme. Procedures with additional treatments using enzymes are described in the literature. The following extraction procedure has proved itself:

Weigh 1 g (ml) homogenized sample into a 50 ml sterile centrifuge vial, add 20 ml redist. or deionized water, shake and adjust pH 4.5 with HCl.

Alternatively instead of water a citrate buffer can be used for the extraction (no pH adjustment is necessary); to 1 g sample add 20 ml citrate buffer pH 4.5 and shake.

Add 300 mg taka diastase, shake well and incubate 1 h in the dark at 37 °C (98.6 °F) (incubator or water bath, shake at times). Fill up exactly to 40 ml with redist. or deionized water. Thereafter, heat the extract for 30 min in a water bath at 95 °C (203 °F). During extraction the vial has to be shaken well at least five times. It is important to make sure that the centrifuge vial is tightly closed. Chill down quickly to below 30 °C (86 °F). Thereafter centrifuge and, depending on the concentration range, further dilute the clear supernatant in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

6.6. Total vitamin content (native and added biotin) in cereals

For the extraction of native biotin from cereal matrix an extraction with acid in an autoclave is necessary.

Weigh 1 g homogenized sample into a 50 ml sterile centrifuge vial, add 20 ml 0.5 mol / l H₂SO₄ and shake well. Autoclave at 121 - 123°C for 30 minutes. Chill down quickly to below 30 °C (86 °F). Adjust the pH to 6.0 - 7.0 with NaOH. Fill up exactly to 40 ml with redist. or deionized water. Thereafter centrifuge and, depending on the concentration range, further dilute the clear supernatant in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

6.7. Total vitamin content (native and added biotin) in yeasts and yeast products

Weigh 1 g of sample into a 50 ml sterile centrifuge vial, fill up to 30 ml with 1 mol / l H₂SO₄ and shake well. Extract for 2 h at 95 °C (203 °F) in a water bath. During extraction the vials have to be shaken well at least five times. It is important to make sure that the centrifuge vials are tightly closed. Chill down quickly to below 30 °C (86 °F) and adjust to pH 6.0 to 7.0 with concentrated and 2 mol / l NaOH. Fill up exactly to 40 ml with redist. or deionized water. Thereafter shake, sterile filter (or heat the sample 30 min at 95 °C in a water bath, thereafter chill down quickly below 30 °C (86 °F) and centrifuge) and, depending on the concentration range, further dilute the clear supernatant in 1.5 ml (or 2.0 ml) sterile reaction vials with sterile water from the test kit.

7. Test implementation

7.1. Test preparation

The **bottle with sterile water**: push the coloured lid up, pull off right up to the glass rim and turn entire lid to remove it.

Biotin standards should be dissolved and diluted freshly. Each dilution is sufficient for three wells.

- open the biotin standard bottle, place the lid down with the opening facing upwards
- add **x ml (x = see quality assurance certificate and label standard bottle)** sterile water (from the test kit) to the standard bottle. Close the standard bottle with the lid and dissolve the standard by shaking = **standard concentrate**
- take 6 sterile vials (1.5 or 2.0 ml) and prepare from the dissolved standard concentrate a **standard curve** according to the following scheme:

standard curve* in µg / 100 g (ml)	sterile water in µl		standard concentrate in µl		total volume in µl
blank: 0	900	+	0	=	900
standard 1: 0.08	900	+	100	=	1000
standard 2: 0.24	350	+	150	=	500
standard 3: 0.40	250	+	250	=	500
standard 4: 0.56	150	+	350	=	500
standard 5: 0.72	50	+	450	=	500

* the sample extraction dilution of 1 : 40 is already included in the standard curve

The **biotin assay medium** is sufficient for 6 microtiter strips. Open the assay medium bottle and remove the desiccant using tweezers (discard the desiccant).

- add 10 ml sterile water from the test kit to the biotin assay medium bottle
- close the assay medium bottle carefully and shake well
- heat the bottle in a water bath at 95 °C (203 °F) for 5 min, shaking at least twice; always make sure that the bottle is tightly closed
- quickly chill down to room temperature below 30 °C (86 °F)
- filter the medium through a 0.2 µm filter into a sterile 15 ml centrifuge vial

7.2. Test procedure

Only sterile samples which are diluted with sterile water from the test kit should be pipetted onto the microtiter plate.

- remove the **required strips** of the microtiter plate and place them into the additional holder; return the unused strips together with the desiccant to the foil bag and seal it well, store at 2 - 8 °C (35.6 - 46.4 °F)

Pipette first the assay – medium and then standard or diluted samples, as followed:

- pipette 150 µl biotin assay medium into the wells
- pipette 150 µl standard or diluted sample into the assigned wells (flush the pipette tip with standard or sample solution)
- cover the strips / cavities with adhesive foil: pull off the protective layer of the foil, place the foil flat onto the strips, smoothing it down by hand, press the foil firmly onto the strips
important: make sure the cavities are sealed airtight by smoothing down the foil over the cavities, take special care with the wells around the edges
- incubate at **37 °C (98.6 °F)** in the dark for **44 - 48 h** in an incubator

7.3. Measurement

- press down the adhesive foil once more, place the microtiter plate upside down on a table and dissolve the microorganisms thoroughly by shaking the plate on the surface of the desk
- invert the plate to the regular position and remove the adhesive foil; **the foil is strongly adhesive, so pulling it off the microtiter plate must be done with great care: hold the strips firmly in the frame with one hand while you pull off the foil**
- destroy any bubbles on the surface of liquid in the wells (by means of a pipette tip or a needle)
- measure the turbidity with a microtiter plate reader at 610 - 630 nm (alternatively at 540 - 550 nm)

Note:

- after 44 - 48 h of incubation, the microtiter plate can be stored for max. 48 h in the refrigerator, thereafter the turbidity should be measured
- to avoid any time losses due to weekends or bank holidays, the microtiter plate can be evaluated after 60 h. It is recommended to use a timer to turn off the incubator after 44 - 48 h.

8. Evaluation

A **4-parameter** evaluation is recommended, e. g. the RIDA[®]SOFT Win from R-Biopharm.

The test evaluation is correct on condition that

- OD blank < OD standard 1
- OD standard 5 > 0.6 OD

The sample dilution factor of 40 is already included in the standard curve. In the below formula merely the further dilution factor of the extract and a differing sample weight need to be taken into consideration.

$$\begin{array}{l} \text{biotin} = \\ \text{(in } \mu\text{g / 100 g)} \\ \text{(in } \mu\text{g / 100 ml)} \end{array} = \frac{\text{conc. standard curve} \times \text{dilution factor}}{\text{sample weight in g (ml)}}$$

Example:

Sample weight:	1 g
Sample extraction dilution:	1 : 40 (must not be considered)
Dilution of the sample extraction:	1 : 300 (has to be considered)
Measured concentration from the standard curve:	0.24 μg biotin / 100 g (ml)

$$0.24 * 300 / 1 = 72 \mu\text{g Biotin} / 100 \text{ g (ml)}$$

Evaluation for capsules, pills, vitamin mixes:

biotin in μg / pill or capsule =

$$\frac{\text{conc. standard curve} \times \text{dilution factor} \times \text{pill or. capsule weight in g}}{\text{sample weight in g} \times 100}$$

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

9. Literature

- Official Methods of Analysis (AOAC) 960.46
- Schweizer Lebensmittelbuch SLMB Kapitel 62 / 10.2.1
- Review Analytical techniques for determining biotin, Livaniou et al., Journal for Chromatography A, 881 (2000) 331-343
- Vitamin-Bestimmungen, Verlag Chemie GmbH Weinheim / Bergstr.; Strohecker und Henning
- The Vitamins, 2. Auflage, Vol. VII. Academic Press, New York / London (1967): V. Herbert und J.R.Bertino

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is a registered trademark of the ifp Institut für Produktqualität GmbH.