r-biopharm®

RIDA®CHECK

REF R1091

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1. Intended use

RIDA®CHECK is a rapid swab test for monitoring the efficiency of surface cleaning procedures in the production environment. The test is intended to be used as a support for hygiene monitoring as well as general protein screening in allergen management.

2. General information

Within the scope of a reasonable HACCP concept it is not sufficient to only carry out microbiological status within a production environment after the food is manufactured, as this control would only monitor the success of sanitation steps within the hygiene program of the food factory. Residues of foods left remaining on surfaces may additionally provide the ideal nutrient base for airborne microbes which can re-contaminate disinfected surfaces. Hence an efficient cleaning of surfaces and equipment is needed and therefore it is expedient to regularly monitor the effectiveness of these procedures.

3. Test principle

During production processes protein residues from raw materials or food stuffs may remain on surfaces. These should be removed during the normal cleaning process. However cleaning actions are not always successful and consequently contamination may still remain and possibly lead to the cause of future health issues.

RIDA®CHECK's methodology is based on the detection of protein residues and is independent of ATP. During the assay the swabbed protein residues react with the colour indicator. This interaction induces a pH decrease which leads to a rapid colorimetric reaction from yellow through to green. According to the intensity of the colour change (yellow → light green → green → dark green) the guideline determination of the contamination level is possible. The operator then determines whether the level of contamination in the tested area is still acceptable (according to routinely observed and recorded colour levels) to decide if corrective action has to be initiated.

4. Reagents provided

Each package of RIDA[®]CHECK (Art. No. R1091) contains 50 swabs for the performance of 100 tests / determinations.

Component	Format	Volume
Swabs	Ready to use	50 pieces

5. Reagents required but not provided

For the testing with RIDA[®]CHECK indicators no additional reagents are needed. Laboratory equipment is not required.

6. Warnings and precautions for the users

Used RIDA[®]CHECK indicators can be disposed with the regular household refuse. Even though the chemicals of the reaction liquid are not harmful to humans, animals or the environment it is recommended to rinse the sampled surface with alcohol or water subsequently to the application of the indicator.

For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the test kit at 4 - 24 °C (39.2 - 75.2 °F). Do not freeze the swabs.

Do not use the test kit after the expiration date (see test kit label).

8. Indication of instability or deterioration of reagents

Do not use swabs which are actually dried out or where the colour already changed to green.

9. Sample preparation

For the swab test on surfaces with RIDA[®]CHECK indicators no sample preparation is needed.

10. Test procedure

RIDA[®]CHECK indicators are ready-to-use and applicable without any preparation of reagents. Touching the cotton heads which are moisturized with reagent liquid should be avoided.

Performance:

1. Rip the packaging of the swab at its perforation.
2. Hold the swab on one side with one hand and remove the other side of the packaging with the other hand.
3. Firmly wipe the surface with the unpackaged part of the RIDA[®]CHECK indicator (test 1).
4. Return the used part of the swab into its packaging.
5. Remove the packaging of the second indicator and firmly wipe the surface of another sampling area (test 2)*.

* **Note:** After opening the indicator package, both indicators have to be used within 1 hour.

The size of the sampling area should equal to approximately 20 cm² (4.5 cm x 4.5 cm). The developing time of the colour reaction depends on the level of contamination with protein residues on the sampled surface. In case of high pollution the colour change will be visible after only a few seconds. In case of only low pollution the time until the reaction ends could last 1 - 2 minutes.

11. Evaluation

The interpretation of the results has to be done after the colour development is ended. In case of the contamination of the surface with protein residues the colour of the indicator changes from yellow to green. The relative degree of surface contamination can be determined according to the following colour scale:

Colour	Contamination value
	Clean
	Borderline
	Lightly contaminated
	Highly contaminated

A yellow “clean” result indicates that the amount of protein is below the detection limit of 20 µg. The threshold has been determined with serum albumin. Due to the fact that food borne proteins are normally mixed combinations of different protein residues, a direct quantitative correlation between the intensity of indicator dye and amount of proteins isn’t possible. The result of each test is only relative to albumin and for this reason the interpretation of the degree of contamination is individual and cannot be compared directly with any individual protein standard. By routinely comparing test results to historical cleaning data the person in charge for the decontamination has to decide personally whether the cleansing of surfaces is acceptable or not and if corrective actions have to be undertaken.

Note: Possibly the RIDA®CHECK indicator may also detect residues of disinfectants. The strength of the colorimetric reaction depends on the amount and the kind of the disinfectant residue. For further information on the detection of disinfectants please refer to R-Biopharm AG.

For further product information, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
2009-04-28	Release version
2023-06-27	Current version Changes made: <ul style="list-style-type: none">– Layout adaption– 4. Reagents provided,– 6. Warnings & precautions for the users– 7. Storage instructions– Version overview– Explanation of symbols– Disclaimer

Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.