r-biopharm®

## RIDA® Ochratoxin A column

**REF** R1303

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# RIDA® Ochratoxin A column

## Brief information

RIDA® Ochratoxin A column (Art. No. R1303) are immunoaffinity columns for sample clean up prior to analysis of ochratoxin A in food and feed (see chapter 1. Intended use).

R1303 contains 10 immunoaffinity columns for single use.

Gel suspension:	sepharose
Antibody:	monoclonal antibody conjugated to sepharose
Max. applied volume:	60 mL
Flow rate:	1 drop/s
Capacity of the column:	approx. 200 ng ochratoxin A

## 1. Intended use

RIDA® Ochratoxin A columns are usable in combination with enzyme immunoassays (e.g. Art. No. R1312 RIDASCREEN® Ochratoxin A 30/15 and R5402 RIDASCREEN®FAST Ochratoxin A) for the quantitative determination of ochratoxin A.

The sample preparation with the immunoaffinity columns simplifies and enhances the sample clean up procedure. Pure extracts are obtained, which can be analyzed by different analytical methods.

## 2. Purification by immunoaffinity column

The basis is the antigen-antibody reaction.

The column contains a gel suspension to which monoclonal antibodies are attached covalently. The antibodies are specific for the ochratoxin A.

- A: The sample is applied and is passed through the column. If ochratoxin A is present in the sample, it will be bound by the monoclonal antibodies.
- B: All other substances will not be retained on the column.
- C: Using methanol as eluent the ochratoxin A will be released from the antigen-antibody-complex. Methanol causes a denaturation of the antibodies. Therefore, the antigen (ochratoxin A) is set free and can be eluted.


Fig. 1: Principle of an immunoaffinity column

### 3. Reagents provided

With the immunoaffinity columns 10 sample preparations respectively could be done (one sample preparation per column). Each test kit contains:

Component	Format	Amount
Immunoaffinity column	Ready to use	10 columns

### 4. Equipment / Reagents required but not provided

#### 4.1 Equipment

- Gloves
- Chemical hood
- Stand / holder for columns (e.g. Immunoaffinity Column Rack, Art. No. RBRCR1);
- Immunoaffinity Column Accessory Pack (contains glass barrels, syringes and adapters, Art. No. RBRAP01);  
alternatively one way plastic pipette tips (5 or 10 mL, as sample reservoir)
- Graduated pipettes
- Evaporator; alternatively water bath (65 °C)
- Variable 20 - 200 µL and 200 - 1000 µL micropipettes

#### 4.2 Reagents

- Methanol p.a.
- Wash buffer Tween
- Distilled water (dist. water) or deionized water

### 5. Warnings and precautions for the users

This kit should only be used by trained laboratory personnel. The instruction for use must be strictly followed.

Methanol is harmful to health. It should be worked under a chemical hood. Ochratoxin A is a toxic and carcinogenic substance. Particular care should be taken. Avoid contact with the skin (use gloves).

Decontamination of the glassware and toxin-content solutions is best carried out using a sodium hypochlorite (bleach) solution (10 % (v/v)) overnight (adjust solution with HCl to pH 7).

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 6. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze columns.

Return any unused columns to their original foil bag and further store at 2 - 8 °C (35 - 46 °F).

Do not use the kit after the expiration date (see kit label).

## 7. Sample preparation

The samples should be stored in a cool place, protected against light.

Ochratoxin A is light sensitive, therefore, avoid exposure of the samples and sample extracts to direct light.

R-Biopharm can provide you different applications for RIDASCREEN® Ochratoxin A 30/15 (Art. No. R1312) and RIDASCREEN®FAST Ochratoxin A (Art. No. R5402) in combination with RIDA® Ochratoxin A column upon request.

**Please contact your local distributor or [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## 8. Separation with RIDA® Ochratoxin A column

### 8.1 Preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

Unused columns should be stored in the original foil bag at 2 - 8 °C (35 - 46 °F) when no longer required.

### Note

The columns are delivered without an adapter for sample reservoir (see chapter 4. Equipment /Reagents required but not provided).

The columns are plugged with caps on top and tip, which have to be removed prior to use.

The columns must not dry up during usage.

Excessive pressure on the column may cause compression of the gel and consequently low recoveries.

## 8.2 Procedure (clean up)

Do not let the columns dry up, especially during sample filling and washing.

### Note

Follow the instructions for sample preparation and clean-up (instruction for use/application note/recommendation) for the specific test kits and matrices in combination with RIDA® Ochratoxin A column (not included in the kit). In case of deviations from the following general procedure, matrix-related instructions must be followed.

1. Rinse the storage buffer of the column or remove it by pressing air.
2. Fill the column with prepared particle free extract and run through without pressure (flow rate: approx. 1 drop/s) (see instructions). Discard passed solution.
3. Wash the column twice with 5 mL wash buffer and remove buffer residues by pressing air. Discard passed solution.
4. Place a clean and closable vial directly below the column and elute the column with 1 mL of 100 % methanol. Press air thoroughly through the column to collect all eluate residues.

(Remark: backflushing of methanol is recommended. If the eluent passed too fast (faster than 10 s), collect the eluate and pass it again through the column.)

5. Process the eluate containing the toxin according to the instructions:
  - Evaporate the eluate to dryness at 65 °C.
  - Dissolve the dried residue (see instructions).
6. Use 50 µL of the reconstituted sample per well in the assay.

## 9. Application notes and recommendations

Application notes and recommendations are available on request.

**For further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

## Version overview

Version number	Chapter and title
1999-01-14	Release version
2009-07-16	General revision
2022-09-16	Current version General revision Changes made: <ul style="list-style-type: none"><li>– Chapter titles adapted to product</li><li>– Updating chapter content:<ul style="list-style-type: none"><li>– 4. Equipment / Reagents required but not provided</li><li>– 5. Warnings and precautions for the users</li><li>– 7. Sample preparation</li><li>– 8. Separation with RIDA® Ochratoxin A column</li></ul></li><li>– Chapter 9. Detection and quantification removed and replaced by 9. Application notes and recommendations</li></ul>

## Explanation of symbols

General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.