# **RIDASCREEN<sup>®</sup> Histamine (enzymatic)**

**Art. No. R1605**

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- Wine – More sensitive analysis with RIDASCREEN<sup>®</sup> Histamine (enzymatic)
- Reduced sample preparation time for fresh fish and canned fish

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# RIDASCREEN<sup>®</sup> Histamine (enzymatic)

## Brief information

RIDASCREEN<sup>®</sup> Histamine (enzymatic) (Art. No. 1605) is an enzymatic test in microtiter plate format for the quantitative determination of histamine in fresh fish, canned fish, fish meal, wine (use RIDA<sup>®</sup> Sample Decolorant, Art. No. R1699, for sample preparation), and cheese.

All reagents required for the enzyme assay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

|                          |  |
|--------------------------|--|
| Sample preparation:      | homogenization and extraction  |
| Time requirement:        | sample preparation (10 fish samples) approx. 35 min<br>test implementation (incubation time) ..... 15 min  |
| Standard:                | histamine  |
| Limit of detection:      | fresh fish, canned fish ..... 0.75 mg/kg<br>fish meal ..... 3.75 mg/kg<br>wine..... 0.54 mg/kg<br>cheese.... 0.75 mg/kg  |
| Limit of quantification: | fresh fish, canned fish ..... 2 mg/kg<br>fresh fish incl. ascorbic acid ..... 3 mg/kg<br>fish meal ..... 10 mg/kg<br>wine..... 1.4 mg/kg<br>cheese... 2 mg/kg  |
| Recovery:                | approx. 85 - 105 %   |
| Specificity:             | The enzyme is specific for histamine.<br>There are no side reactions to L-histidine,<br>L-tryptophan, serotonin, L-tyrosine, tyramine, L-lysine,<br>trimethylamine, putrescine and cadaverine. Further<br>information on side reactions are contained in the<br>updated validation report. |

The specificity of RIDASCREEN<sup>®</sup> Histamine (enzymatic) has been evaluated by the determination of the side activities of the respective substances in a buffer system (see validation report).

## Related products

RIDA<sup>®</sup> Sample Decolorant (Art. No. R1699) for wine analysis

## 1. Intended use

RIDASCREEN® Histamine (enzymatic) is used for the quantitative determination of histamine in fresh fish, canned fish, fish meal, wine and cheese. The test kit will be approved by the AOAC *Performance Tested Methods Program*<sup>SM</sup>.

## 2. General

Histamine is a product of the decomposition of histidine caused by the growth of certain bacteria in protein rich food like fish, cheese, and wine.

The amount of histamine formed depends on the bacterial species, the temperature, and the time of exposure and may exceed in fish 1000 mg/kg. Fish of good quality contains less than 10 mg/kg histamine. The limit value for fish and products thereof is between 50 and 200 mg/kg depending on the fish species and the country. Similar high concentrations can be found in cheese, while histamine in wine usually does not exceed 15 mg/l.

## 3. Test principle

The enzymatic determination is based on the histamine dehydrogenase which catalyzes the oxidative deamidation of histamine in the presence of an electron carrier that converts a dye to a color product. The color intensity is directly proportional to histamine concentration and is measured at 450 nm. The electron carrier and dye are coated on the microtiter plate.

## 4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). All reagents are ready to use. Each test kit contains:

| Component   | Cap color   | Format       | Histamine conc. | Volume   |
|---|-------------|--------------|-----------------|----------|
| <b>Microtiter plate</b>                             | -           | ready to use |                 | 96 wells |
| <b>Buffer</b>                                       | transparent | ready to use |                 | 15 ml    |
| <b>Standard 1</b>                                   | white       | ready to use | 0 mg/l          | 1.3 ml   |
| <b>Standard 2</b>                                   | white       | ready to use | 1 mg/l          | 1.3 ml   |
| <b>Standard 3</b>                                   | white       | ready to use | 5 mg/l          | 1.3 ml   |
| <b>Standard 4</b>                                   | white       | ready to use | 10 mg/l         | 1.3 ml   |
| <b>Standard 5</b>                                   | white       | ready to use | 15 mg/l         | 1.3 ml   |
| <b>Standard 6</b>                                   | white       | ready to use | 20 mg/l         | 1.3 ml   |
| <b>Enzyme solution</b>                              | red         | ready to use |                 | 1 ml     |
| <b>Spiking solution</b>                             | blue        | ready to use | 500 mg/l        | 3 ml     |
| <b>Catalase</b><br>Catalase to remove ascorbic acid | black       | ready to use |                 | 1 ml     |

## 5. Reagents and equipment required but not provided

### 5.1. Equipment and reagents

**For the sample preparation only plastic vials (no glass) e.g. polypropylene should be used.**

- distilled or deionized water
- vortex
- mixer, blender (e.g. Mr. Magic)
- centrifuge, centrifuge vials
- variable 100 - 200 µl micropipettes
- multistep-pipette with 10 µl, 100 µl, 150 µl, 200 µl pipetting volume (e.g. Eppendorf combitips)
- microtiter plate spectrophotometer (450 nm)
- RIDASOFT® Win.NET (Z9996); use **version 1.104** or later for calculation

### **Fresh fish analysis**

- 50 ml polypropylene (PP) screw-cap vials (e.g. Greiner; Product No. 227261)
- homogenizer: laboratory mincer / grinder, Ultra-Turrax or mixer
- water bath (100 °C / 212 °F)

### **Fresh fish analysis: As required for degradation of ascorbic acid**

- 1 N NaOH (e.g. Roth, Art. No. CN 58.1)
- 1 N HCl (e.g. Roth, Art. No. CN 63.1)
- H<sub>2</sub>O<sub>2</sub> solution, 30 % in water (e.g. Sigma, Art. No. 31642)
- Paper filter (e.g. Whatman no.1 or similar quality)

### **Wine analysis**

- 1.5 - 2 ml polypropylene vials (e.g. Eppendorf tube)
- RIDA® Sample Decolorant (R1699) for sample preparation

### **Cheese analysis**

- 1.5 - 2 ml polypropylene vials (e.g. Eppendorf tube)
- 1 N perchloric acid (e.g. Bernd Kraft, Art. No. 20126.2700)
- 1 N KOH (e.g. Carl Roth, Art. No. K017.1)

## 6. Warnings and precautions for the users

This test should be carried out only by trained laboratory technicians. The instructions for use must be strictly followed. Use clean equipment and wear gloves before starting and during the assay in order to avoid cross-contamination. Take precautions when using the boiling water bath to avoid being burnt or scalded.

The analytical procedure might include the use of hazardous chemicals. Please refer to the appropriate safety data sheets (SDS) available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## 7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

No quality guarantee is given after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Indication of instability or deterioration of reagents

If the difference between absorption of  $A_2 - A_1$  of standard 6 (e.g.  $\Delta A < 1.7$ ) is differing strongly from the certificate of analysis, then a deterioration of reagents has likely occurred. The slope of the standard curve should lie in the range of 0.0847 – 0.1025 (mean 0.0936).

## 9. Preparation of samples

The histamine concentration increases also in refrigerated food, therefore it is important to carry out the laboratory analysis directly after sampling. **Use only plastic vials (e.g. polypropylene, PP) for the sample preparation. Histamine adsorbs on glass surfaces which causes reduced recoveries.**

## 9.1. Fresh fish and canned fish

- weigh 5 g of homogenized sample into a 50 ml PP-screw cap vial
- pour 20 ml of dist. water over the sample
- close the vial and shake or vortex until the sample is **evenly suspended**
- incubate the sample for 20 min in a boiling water bath, take vials out of the water bath while wearing safety gloves after 10 min and shake immediately for 3 s so that samples are well suspended
- allow the sample to cool to room temperature in an ice bath for at least 2 min
- centrifuge at least 2,500 g for 10 min  
Alternatively transfer 2 ml of the extract into centrifuge vials and centrifuge at high speed (> 10,000 g) for 2 min in a micro-centrifuge. The extract should be as clear as possible. If particles are present, then filtrate the extract.  
(If a layer of fat has formed on top of the extract after centrifugation then hold the vial angular and take the extract from the side. Alternatively, carefully withdraw with a pipette the lower layer and deliver to a new vial, filter / centrifuge again)
- pipette 100 µl of the undiluted **clear extract** per well. (The extracted samples are stable at room temperature for 2 h, or at 2 - 8 °C (36 - 46 °F) for 7 days.)

### 9.1.1 Degradation of ascorbic acid (> 250 mg/kg) in fish

If fresh fish has been treated with ascorbic acid (analysis using test kits Art. No. E1267 or 10409677035), an intense yellow color may occur (due to a creep reaction) during the subsequent enzymatic analysis within seconds ( $A_1$  measurement e.g. 3.0).

The ascorbic acid is degraded by  $H_2O_2$  / catalase as follows:

- Transfer 1 ml of the clear supernatant (see 9.1) in a 50 ml PP screw cap vial
- Add 100 µl 1 N NaOH
- Add 100 µl 30 %  $H_2O_2$
- Mix softly and incubate 10 min at room temperature
- Add 100 µl 1 N HCl
- Mix gently
- Add 10 µl Catalase
- Mix softly and incubate 10 min at room temperature
- Put the mixture for 5 min into a boiling water bath
- Let cool down and pipette 100 µl of the **clear extract** per well (The extracted samples are stable at room temperature for 2 h, or at 2 - 8 °C (36 - 46 °F) for 7 days.)

For validation spike the sample and determine the recovery.

## 9.2. Fish meal

- weigh 1 g of homogenized fish meal into a 50 ml PP-screw cap vial
- pour 25 ml of dist. water over the sample
- close the vial and shake or vortex until the sample is evenly suspended
- incubate the sample for 20 min in a boiling water bath, take vials out of the water bath while wearing safety gloves after 10 min and shake immediately for 3 s so that samples are well suspended
- allow the sample to cool to room temperature in an ice bath for at least 2 min
- centrifuge at least 2,500 g for 10 min

Alternatively transfer 2 ml of the extract into centrifugal vials and centrifuge at high speed (> 10,000 g) for 2 min in a micro-centrifuge. The extract should be as clear as possible. If particles are present, then filtrate the extract.

(If a layer of fat has formed on top of the extract after centrifugation then hold the vial angular and take the extract from the side. Alternatively, carefully withdraw with a pipette the lower layer and deliver to a new vial, filter / centrifuge again)

Pipette 100 µl of the undiluted **clear extract** per well. (The sample extracts can be stored at room temperature for 2 h, at 2 - 8 °C (36 - 46 °F) for 3 days or at -20 °C (- 4 °F) for 3 months).

## 9.3. Wine

Wine (red wine, white wine, rosé) contains anthocyanins and polyphenols which have to be removed prior to the enzymatic analysis.

- Add 200 µl Reagent 1 of RIDA<sup>®</sup> Sample Decolorant to 200 µl wine (rinse pipette tip)
- Add 200 µl Reagent 2 of RIDA<sup>®</sup> Sample Decolorant, vortex, and incubate 5 min at room temperature, then centrifuge 2 min at a minimum of 14,000 g (room temperature, 20 - 25 °C, 68 - 77 °F)
- Transfer 500 µl of the supernatant in a new vial
- Add 100 µl Reagent 3 of RIDA<sup>®</sup> Sample Decolorant, vortex and incubate 5 min at room temperature
- centrifuge at 14,000 g (room temperature) for 2 min
- Use 100 µl of the **clear extract** per well (the sample extracts can be stored at room temperature for 2 h or at 2 - 8 °C (36 - 46 °F) for 1 day)

## 9.4. Cheese

If the cheese rind is edible then it should also be part of the representative sample.

- Weigh 60 g of cheese into the plastic cup for blending and add 140 ml of water. Mix for 15 s to suspend all particles. Repeat mixing two more times to suspend all particles
- Remove the fat layer with a spatula. Transfer 50 ml of the unclear solution to a plastic vial
- Centrifuge at 3,500 g for 10 min (at 4 °C if possible otherwise cool sample down in an ice bath)
- Remove the fat layer with a spatula and transfer 0.8 ml of the clear solution to a plastic vial
- Add 200 µl 1 N perchloric acid, mix well
- Centrifuge at 20,000 g for 2 min at room temperature
- Transfer 500 µl of the clear supernatant to a new centrifugal plastic vial and add 100 µl 1 N KOH, mix
- Centrifuge 2 min at of least 20,000 g at room temperature (20-25 °C, 68-77 °F)
- Use 100 µl of **clear extract** per well. (The sample extracts can be stored at room temperature for 2 h. or at 2 - 8 °C (36 - 46 °F) for 18 h)

## 10. Test implementation

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) and mix carefully before use.

### 10.2. Test implementation

Insert a sufficient number of wells into the microwell holder for all standard and samples to be run in duplicate. Record standard and sample positions.

1. Using a multi-stepper, add 150 µl of buffer to the wells and carefully shake the plate manually for 3 s
2. Add 100 µl of standards, controls or prepared samples to separate wells in duplicate. Flush the pipette tips. Thereafter, carefully shake the plate manually for 3 s
3. After 3 min\* (for wine samples at 10 min) the absorption A1 is measured at 450 nm **without** a reference filter (e.g. 620 nm) (avoid any bubbles in the wells)
4. Add 10 µl of the blue dyed enzyme solution with a multi-stepper to each well. Thereafter, carefully shake the plate manually for 3 s

5. After 10 min the absorption A2 is measured at 450 nm **without** a reference filter (e.g. 620 nm). For evaluation use the RIDASOFT® Win.NET, select A1 and A2 and calculate.

\*For **fish** an **increase of yellow color** may occur (creep reaction)

- If A1 is e.g. 0.3 then the incubation time is increased to 10 min.

- If A1 is e.g. 3.0 then the fish sample should be prepared with catalase (see 9.1.1)

## 11. Results

At R-Biopharm AG a special software, the RIDASOFT® Win.NET (Art. Nr. Z9996, use version 1.104 or later), is available for evaluation of the RIDASCREEN® assays. The calculation is done by using a linear regression, the dilution factor of the samples has to be taken into account. (From software version 1.105 onwards the default factor is 5 which needs to be modified if required).

| Sample                                  | Dilution factor | Measuring range  |
|---|-----------------|------------------|
| Fresh/canned fish                       | 5               | 2 - 100 mg/kg    |
| Fish incl. degradation of ascorbic acid | 6.55            | 2.62 - 131 mg/kg |
| Fish meal                               | 25              | 10 - 500 mg/kg   |
| Wine                                    | 3.6             | 1.4 - 72 mg/kg   |
| Cheese                                  | 5               | 2 - 100 mg/kg    |

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit. The slope of the standard curve should lie in the range of 0.0847 – 0.1025 (mean 0.0936).

The residuals for  $\Delta A$  should lie between - 0,05 und 0,05 (refer to the residual plot, last page of the RIDASOFT® Win report). Larger deviation of single wells are indicating (almost) always an insufficient precision of pipetting, therefore the pipetting technique needs to be improved.

The coefficient of variation should not exceed 10 % for samples and standards (except for CVs calculated from standard 1 and samples close to standard 1). The recovery of check samples is usually between 75 and 125 %.

Higher values of the absorbance ( $A_{450\text{ nm}}$ ) for the calibration curve, especially for the zero standard, may be a result of histamine contamination. If absorbance values ( $A_{450\text{ nm}}$ ) > standard 6, dilute the final test extract with water and perform another assay.

### **Recommendation:**

In order to ensure a high analytical reliability:

- Each final test extract should be assayed in duplicate, **ensure accurate pipetting**
- Employ also histamine-free and histamine-containing (spiked) samples as test controls (e.g. FAPAS<sup>®</sup> reference samples)

To verify the correct and interference-free implementation of the test procedure the spike solution contained in the kit should be used for spiking.

### **10 mg/l spike control**

- 1:50 dilution of the Spiking Solution Histamine (500 mg/l) e.g. 100  $\mu$ l + 4.9 ml water
- Use 100  $\mu$ l in the assay.

### **50 mg/kg histamine spiked fish sample**

- pipette 500  $\mu$ l of the Spiking Solution Histamine (500 mg/l) to 5 g homogenized fish
- extract the sample with 19.5 ml of water (see 9.2.)

### **250 mg/kg histamine spike fish meal sample**

- pipette 500  $\mu$ l of the Spiking Solution Histamine (500 mg/l) to 1 g homogenized fish meal sample
- extract the sample with 24.5 ml of water (see 9.3.)

### **10 mg/l histamine spiked wine sample**

- pipette 200  $\mu$ l of the Spiking Solution Histamine (500 mg/l) to 9.8 ml wine
- extract the sample with RIDA<sup>®</sup> Sample Decolorant (see 9.4)

**Further application notes:**

- Wine – More sensitive analysis with RIDASCREEN® Histamine (enzymatic)
- Reduced sample preparation time for fresh fish and canned fish.

For additional applications (e.g. milk, yeast extract, pet food) and product information (e.g. automation, training videos) please contact your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG