RIDASCREEN®FAST Biotin

REF R2202

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RIDASCREEN®FAST Biotin

Brief information

RIDASCREEN®FAST Biotin (Art. No. R2202) is a competitive enzyme immunoassay for the quantitative determination of biotin in milk, milk powders, food for special medical purpose, grain and cereals, vitamin tablets, vitamin mixtures, vitamin powders, vitamin gummies, vitamin waters and vitamin juices.

All reagents required to perform the enzyme-linked immunosorbent assay (ELISA), including standards, are included in the test kit. One test kit is sufficient for up to 48 measurements (including standard measurements). A microplate reader is required for data analysis.

Sample preparation:	liquid milk: dilution and neutralisation vitamin waters and juices: dilution, pH adjustment vitamin tablets/gummies: extraction, dilution milk powders and cereals: extraction, neutralisation, dilution
Time required:	Sample preparation (10 samples).....10 - 60 min Test execution (Incubation time)25 min
Limit of detection:	1.3 µg/L (ppb) (in standard solution*) (corresponding to standard substance)
Specificity:	Biotin (Vitamin B7)..... 100 % Biocytin...approx. 250 % Biotin-PABA.....approx. 130 % Biotin (S)-sulfoxide.....approx. 35 % Biotin-sulfone.....approx. 14 % Desthiobiotin.....approx. 7 %

The specificity of the RIDASCREEN®FAST Biotin test was determined by assessing cross-reactivity with the corresponding substances in the buffer system.

In samples, specificity may differ from the results determined in the buffer system due to matrix effects. Before analyzing cross-reactive substances, the user must determine their limit of detection and recovery rate in the respective matrix. The test cannot distinguish between analytes and cross-reactive substances.

Note: The test exhibits cross-reactivity with biocytin, a biotin metabolite. Biocytin is a biotin-lysine conjugate that is primarily bound to proteins. Natural levels of biotin in vitamin juices and water are low or negligible; therefore, the presence of its metabolite, biocytin is also minimal. For infant milk powder, liquid milk, and cereal samples, an extensive validation study confirmed that interference from biocytin is negligible. This enables accurate quantification of biotin in these samples using the Biotin ELISA.

This conclusion aligns with the findings of Joseph G. et al. (2018), *Determination of Total Biotin by Liquid Chromatography Coupled with Immunoaffinity Column Cleanup Extraction: Multilaboratory Testing, Final Action 2016.02. J AOAC Int*, 101(3):831–842. In that study, no samples of infant milk powder were found to contain free biocytin.

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It contains minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website

<https://food.r-biopharm.com/media/technical-guides/>.

1. Intended use

RIDASCREEN®FAST Biotin is a competitive enzyme immunoassay for the quantitative determination of biotin in milk, milk powders, food for special medical purpose, grain and cereals, vitamin tablets, vitamin mixtures, vitamin powders, vitamin gummies, vitamin waters and vitamin juices.

2. General

Biotin belongs to the group of water-soluble vitamins. It cannot be produced by the human body and must, therefore, be obtained through one's diet. While biotin is required in only small amounts, it plays a crucial role in metabolism. Serving as a coenzyme, biotin assists in the conversion of carbohydrates and proteins into energy and facilitates the synthesis of fatty acids. Additionally, biotin is vital for the growth and repair of skin, hair, and nails. Although biotin deficiency rarely occurs in healthy adults with a balanced diet, genetic conditions can lead to serious health problems when there is a lack of biotin.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-biotin antibodies. Biotin standards or sample solutions, biotin enzyme conjugate and anti-biotin antibodies are added. Free biotin and biotin enzyme conjugate compete for the biotin antibody binding sites (competitive enzyme immunoassay). At the same time, the anti-biotin antibodies are also bound by the immobilized capture antibodies. Any unbound enzyme conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is inversely proportional to the biotin concentration in the sample.

4. Reagents provided

The reagents in one package are sufficient for up to 48 tests (including standard tests). Each test kit contains:

Component	Lid color	Condition		Volume
Microtiter plate	-	ready to use		48 wells
Sample buffer Probenpuffer	transparent	2x concentrate	2x	125 mL
Standard 1	white	ready to use	0 µg/L (ppb)	1.3 mL
Standard 2	white	ready to use	3 µg/L (ppb)	1.3 mL
Standard 3	white	ready to use	10 µg/L (ppb)	1.3 mL
Standard 4	white	ready to use	30 µg/L (ppb)	1.3 mL
Standard 5	white	ready to use	100 µg/L (ppb)	1.3 mL
Standard 6	white	ready to use	300 µg/L (ppb)	1.3 mL
Wash buffer salt Tween	-	salt for dissolving		
Conjugate	red	ready to use		3 mL
Substrate/Chromogen Red Chromogen Pro	brown	ready to use		10 mL
Stop solution	yellow	ready to use		14 mL
Neutralisation buffer	blue	ready to use		3 mL

5. Materials required but not provided

5.1. Equipment:

- Microtiter plate spectrophotometer (450 nm)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge + centrifugal vials
- Shaker
- Optional: paper filter and funnel
- Graduated pipettes
- Variable 20 - 200 µL and 200 - 1000 µL micropipette
- Optional: multi-dispenser or 8-channel pipette
- Optional: microtiter plate washer

5.2. Reagents:

- Distilled or demineralized water
- 1 M NaOH

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at www.r-biopharm.com.

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant bag provided and further store at 2 - 8 °C (36 - 46 °F).

The standards, the conjugate and the substrate/chromogen are light sensitive; therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers (exception: sample buffer, wash buffer (salt), stop solution).

8. Indication of instability or deterioration of reagents

- Bluish coloration of the red stained substrate/chromogen prior to test implementation
- An absorption smaller than 0.6 ($A_{450\text{ nm}} < 0.6$) for the zero standard

9. Preparation of samples

The samples should be stored in a cool place, protected from light. The samples should be brought to room temperature before analysis (20 - 25 °C / 68 - 77 °F). For quality control, it is recommended to analyse control samples.

Remark: The samples should be measured at the day of extraction.

Note: The sample buffer can also be used for the sample preparation of folic acid or vitamin B12 containing samples in combination with the RIDASCREEN®FAST Folic Acid (Art. No. R3203) or RIDASCREEN®FAST Vitamin B12 (Art. No. R2103), respectively.

9.1. Milk, milk powder and food for special medical purpose

Milk:

- Collect 10 mL of the sample in a tube.
- Centrifuge for (or alternatively filter) in order to remove undissolved components.
- Collect 200 µL of the supernatant. Add 50 µL of the neutralisation buffer and 350 µL of the sample buffer. Mix well.
- The sample can be now analysed directly. Add 50 µL per well in the assay. The calculation factor is 3.
- The sample can be also diluted further with the sample buffer to fit within the testing range of the assay. Note that the calculation factor will change!

Milk powder, food for special medical purpose:

- Suspend 1 g sample in approx. 5 ml distilled water, fill up to 10 ml with distilled water.
- Shake / mix for 10 min.
- Heat the diluted sample for 3 min at 100 °C (212 °F) in a water bath, cool down quickly (ice bath).
- Centrifuge (or alternatively filter) in order to remove undissolved components.
- Collect 200 µL of the supernatant. Add 50 µL of the neutralisation buffer and 350 µL of the sample buffer. Mix well.
- The sample can be now analysed directly. Add 50 µL per well in the assay. The calculation factor is 30.
- The sample can be also diluted further with the sample dilution buffer to fit within the testing range of the assay. Note that the calculation factor will change!

9.2. Grain and cereals

- Weigh 1 g of ground sample into a suitable vial, suspend with approx. 5 mL of sample buffer and fill up to 10 mL with sample buffer (dilution factor = 10).
- Shake / mix for 10 min.
- Centrifuge (or alternatively filter) in order to remove undissolved components.
- Take supernatant, if supernatant is still turbid, filtrate.
- Collect 200 µL of the clear supernatant and add 50 µL of the neutralisation buffer and 350 µL of the sample buffer.
- The sample can be now analysed directly. Add 50 µL per well in the assay. The calculation factor is 30.
- The sample can be also diluted further with sample buffer to fit within the testing range of the assay. Note that the calculation factor will change!

Remark: If the sample has a gel-like consistency after centrifugation it is difficult to pipette the upper layer. The sample should be re-extracted with up to 15 mL or 20 mL of the sample buffer allowing for the better separation of the phases. Note that the calculation factor will change!

9.3. Vitamin powders, vitamin mixtures, vitamin tablets, and vitamin gummies

Vitamin powders / vitamin mixtures:

- Dissolve 1 g of vitamin powder / vitamin mixture in 5 mL sample buffer and fill up to 10 mL with sample buffer (dilution factor = 10).

Vitamin tablets:

- Crush the tablets and dissolve 1 g in 5 mL sample buffer and fill up to 10 mL with sample buffer (dilution factor = 10).

Vitamin gummies:

- Dissolve 1 gummy in the sample buffer to obtain 1 g in 10 mL solution. For instance, if the gummy weighs 2.2 g, add up to 22 mL (dilution factor = 10).
- Heat the mixture in a water bath at 50°C or on a hot plate stirrer (50°C) with controlled temperature until fully dissolved

Continue with these samples as described below:

- Shake / mix for 10 min.
- Centrifuge (or alternatively filter) in order to remove undissolved components.
- Dilute the supernatant or filtrate if necessary with the sample buffer into the testing range (new calculation of the dilution factor is necessary!).
- Add 50 µL per well in the assay.

9.4. Vitamin juices

- Centrifuge (or alternatively filter) in order to remove undissolved components.
- Dilute at least 5 times in the sample buffer, for example 2 mL of juice with 8 mL of the sample buffer (dilution factor = 5).
- Adjust the pH to 8.0 using 1 M NaOH.
- Dilute the supernatant or filtrate if necessary with the sample buffer into the testing range (new calculation of the dilution factor is necessary!).
- Add 50 µL per well in the assay.

9.5. Vitamin waters

- Centrifuge (or alternatively filter) in order to remove undissolved components.
- Dilute at least 5 times with the sample buffer, for example 2 mL of sample with 8 mL of the sample buffer (dilution factor = 5).
- Add 50 µL per well in the assay.

Remark: acidic vitamin water samples (such as carbonated energy drinks) require the same procedure (pH adjustment to 8.0) as vitamin juices.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **sample buffer** is provided as a 2-fold concentrate. The concentrate is diluted 1:2 with distilled water (e.g. 50 mL sample buffer concentrate + 50 mL distilled water). The diluted sample buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

A PBS-Tween buffer is required as **wash buffer**, please use the wash buffer salt (pouch) contained in the kit (see chapter 4. Reagents provided). Dissolve the total content of the pouch in one liter of distilled water. The ready to use wash buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternatively: Dissolve the content of the pouch in 100 mL of distilled water to obtain a 10-fold concentrated wash buffer. This solution expires after approx. 8 - 12 weeks, stored at room temperature (20 - 25 °C / 68 - 77 °F). Use one part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.

10.2. Test procedure

Standards and samples must be tested on one microtiter plate at one time. Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µL of each standard or prepared sample in duplicate to the wells; use a new pipette tip for each standard or sample.
3. Add 50 µL of conjugate to each well (multi-dispenser or 8-channel pipette) and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Using an 8-channel pipette, fill the wells with ready to use wash buffer (250 µL per well) (see chapter 10.1). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.

5. Add 100 µL of the substrate/chromogen to each well (multi-dispenser or 8-channel pipette) and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
6. Add 100 µL of stop solution to each well (multi-dispenser or 8-channel pipette) and measure the absorbance at 450 nm within 5 min after addition of stop solution.

11. Results

A special software, the **RIDASOFT® Win.NET** (Art. No. Z9996FF), is available for evaluation of the RIDASCREEN® enzyme immunoassays. The results are given in µg/kg respectively µg/L and µg/100 g respectively µg/100 mL.

The course of the standard curve is shown in the Certificate of analysis (CoA) enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages (%). The values calculated for the standards are entered in a system of coordinates on semi-logarithmic graph paper against the biotin concentration (µg/kg or µg/L).

In order to obtain the biotin concentration actually contained in a sample, the concentration read from the standard curve must be further multiplied by the corresponding dilution factor.

12. Limitations of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

13. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

For further product information and applications, please contact your local distributor or R-Biopharm at this address: sales@r-biopharm.de.

Version overview

Version number	Chapter and title
1.0	Release version

Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM-DD)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM-DD)
	Manufacturer + address

Disclaimer

1. R-Biopharm AG provides a limited statutory warranty ("Gewährleistung") as mandated and pursuant to German law against defects in title and material. For products with a defined shelf life, the warranty applies for the stated shelf life. For use cycle limited products, it applies until the stated cycle limit is reached, but not beyond the general warranty period. For all other products, the general warranty period is twelve (12) months from the date the risk transfers.
Such warranty is subject to timely notice of the alleged defect, giving sufficient detail to assess the claim. A product's defectiveness is assessed based on the specifically agreed quality, its suitability for the contractually stated purpose, and whether all agreed-upon accessories, documentation, and instructions were provided.
2. R-Biopharm AG is not liable for issues arising from:
 - a) The failure to read, understand, or follow product instructions and information for use or safety
 - b) The use of unqualified personnel
 - c) The failure to follow industry standards such as Good Laboratory Practices
 - d) The use of unsuitable controls, sample matrices, or procedures
 - e) Other misuse
 - f) unilateral product modifications
 - g) Improper storage by others
 - h) External chemical, electromagnetic, mechanical, or electrolytic influences that exceed the defined parameters.
 - i) Harm caused by external events beyond R-Biopharm's control, including but not limited to burglary, theft, lightning, fire, water, or acts of nature.
3. R-Biopharm AG's liability is limited to foreseeable, contract-typical losses in case of negligent breaches of obligations that are essential to achieving the agreement's purpose and that the other party could reasonably rely on (material contract obligations). R-Biopharm AG is not liable for losses resulting from non-material contract obligations.
4. ALL OTHER WARRANTIES OR GUARANTEES OF ANY KIND, EXPRESS OR IMPLIED, ARE EXCLUDED, WHETHER IMPLIED BY CUSTOM, PRACTICE, COURSE OF DEALINGS BETWEEN THE PARTIES OR OTHERWISE. R-Biopharm AG assumes no liability for consequential damages, in particular loss of profit, interruptions, or other indirect damages.
5. The foregoing is not intended to limit liability to the extent such liability cannot be excluded or limited under applicable German law or convention, including cases of intent, willful misconduct, or gross negligence, personal injury or death, express warranties or assumptions of procurement risks, strict liability, or other non-waivable legal mandates.
6. The foregoing liability provisions apply equally to R-Biopharm AG, its directors, officers, employees, and agents.

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are registered trademarks of R-Biopharm AG.