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## Brief information

RIDASCREEN®FAST Folsäure (Folic Acid) (Art. No. R3203) is a competitive enzyme immunoassay for the quantitative determination of added folic acid in milk, milk powders, food for special medical purpose\*, grain and cereals, fortified flour, vitamin powders, vitamin mixtures, vitamin tablets, and vitamin juices.

All reagents required for the enzyme immunoassay, including standards, are contained in the test kit. The test kit is sufficient for a maximum of 48 determinations (including standards). A microtiter plate spectrophotometer is required for quantification.

Sample preparation:           milk: can be applied directly, dilution if necessary  
  juices: can be applied directly, dilution if necessary  
  other matrices: extraction, dilution if necessary

Time requirement:           sample preparation (10 samples).....10 - 60 min  
  test implementation (incubation time).....25 min

Limit of detection:           0.6 µg/L (ppb) (in measurement solution\*\*)  
(corresponding to standard substance)

Specificity:

Folic acid.....	100 %
Dihydrofolic acid.....	approx. 14 %
Tetrahydrofolic acid.....	< 0.5 %
5-Methyltetrahydrofolic acid.....	< 0.5 %
Formyltetrahydrofolic acid.....	< 0.5 %
Tetrahydrobiopterin.....	< 0.5 %
6-Biopterin.....	< 0.5 %

The specificity of the RIDASCREEN®FAST Folic Acid test was determined by analyzing the cross-reactivities to corresponding substances in a buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

\* Commission Directive 1999/21/EC on dietary foods for special medical purposes

\*\* see product information

In order to increase the quality of assessment when performing ELISA procedures, we refer additionally to our Good ELISA Practice brochure. It contains minimum standards and conditions that are required when using test kits of R-Biopharm AG to perform ELISA analysis. The brochure can be retrieved, printed and downloaded from the website <https://food.r-biopharm.com/media/technical-guides/>.

## **1. Intended use**

RIDASCREEN®FAST Folsäure (Folic Acid) is a competitive enzyme immunoassay for the quantitative determination of added folic acid in milk, milk powder, food for special medical purpose, grain and cereals, fortified flour, vitamin powders, vitamin mixtures, vitamin tablets, and vitamin juices.

## **2. General**

Folic acid also known as vitamin B9 or vitamin B11 belongs to the group of water-soluble vitamins. It cannot be produced by the human body and therefore must be included in the diet in order to meet the daily requirements. Folic acid plays a key role in all growth and cell division processes. Adequate folic acid intake during pregnancy is important to avoid neural tube defects, for example. Some countries have agreed to fortify either bread, cereals, flour, pasta, rice or other grain products with folic acid in different concentrations as a preventive approach to improve micronutrient status of the population.

## **3. Test principle**

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with specific antibodies against folic acid. Folic acid standards, respectively sample and enzyme labelled folic acid (conjugate) are added. Free folic acid and enzyme-conjugated folic acid compete for the antibody binding sites (competitive enzyme immunoassay). Any unbound conjugate is then removed in a washing step. Substrate/chromogen is added to the wells, bound conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The resulting absorbance values are inversely proportional to the folic acid concentration of the sample.

## 4. Reagents provided

Each kit contains sufficient materials for a maximum of 48 measurements (including standard analyses). Each test kit contains:

Component	Cap Color	Format		Volume
<b>Microtiter plate</b>	-	ready to use		48 wells
<b>Sample buffer</b>	transparent	2-fold concentrate	2x	125 mL
<b>Standard 1</b>	white	ready to use	0 µg/L (ppb)	1.3 mL
<b>Standard 2</b>	white	ready to use	0.75 µg/L (ppb)	1.3 mL
<b>Standard 3</b>	white	ready to use	1.5 µg/L (ppb)	1.3 mL
<b>Standard 4</b>	white	ready to use	3 µg/L (ppb)	1.3 mL
<b>Standard 5</b>	white	ready to use	6 µg/L (ppb)	1.3 mL
<b>Standard 6</b>	white	ready to use	12 µg/L (ppb)	1.3 mL
<b>Wash buffer salt Tween</b>	-	salt for dissolving		
<b>Conjugate</b>	red	ready to use		3 mL
<b>Substrate/Chromogen Red Chromogen Pro</b>	brown	ready to use		10 mL
<b>Stop solution</b>	yellow	ready to use		14 mL

## 5. Materials required but not provided

### 5.1. Equipment:

- Microtiter plate spectrophotometer (450 nm)
- Laboratory mincer / grinder, mortar, ultra-turrax or homogenizer
- Centrifuge + centrifugal vials
- Water bath (100 °C / 212 °F)
- Ice bath (0 °C / 32 °F)
- Shaker
- Optional: paper filter and funnel
- Graduated pipettes
- Variable 20 - 200 µL and 200 - 1000 µL micropipettes
- Optional: multi-dispenser or 8-channel pipette
- Optional: microtiter plate washer

### 5.2. Reagents:

- Distilled water or demineralized water

## 6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (SDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

All reagents and materials must be recovered or disposed after use at customers own responsibility according to the protection of human health and the environment. Please observe the applicable national regulations concerning waste disposal (e.g. Waste Management Act, Regulations on Dangerous Chemicals, etc.).

## 7. Storage instructions

Store the kit at 2 - 8 °C (36 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant bag provided and further store at 2 - 8 °C (36 - 46 °F).

The standards, the conjugate and the substrate/chromogen are light sensitive; therefore, avoid exposure to direct light.

Do not use the test kit after the expiration date (see test kit label).

Do not interchange individual reagents between kits of different lot numbers (exception: sample buffer, wash buffer (salt), stop solution).

## 8. Indication of instability or deterioration of reagents

- Bluish coloration of the red stained substrate/chromogen prior to test implementation
- An absorption smaller than 0.6 ( $A_{450\text{ nm}} < 0.6$ ) for the zero standard

## 9. Preparation of samples

The samples should be stored in a cool place, protected from light. The samples should be brought to room temperature before analysis (20 - 25 °C / 68 - 77 °F). For quality control, it is recommended to analyze control samples.

**Remark:** The samples should be measured at the day of extraction.

**Note:** The sample buffer can also be used for the sample preparation of vitamin B12 containing samples in combination with the RIDASCREEN®FAST Vitamin B12 (Art. No. R2103).

### 9.1. Milk, milk powder and food for special medical purpose

#### Milk:

- Homogenized milk can be directly applied in the test, dilute if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- Add 50 µL per well in the assay

#### Milk powder, food for special medical purpose:

- Suspend 1 g sample in approx. 5 mL distilled water, fill up to 10 mL with distilled water (dilution factor = 10)
- Shake / mix for 10 min
- Heat the diluted samples for 3 min at 100 °C (212 °F) in a water bath, cool down quickly (ice bath)
- Centrifuge (or alternatively filter) in order to remove undissolved components
- Dilute the supernatant or filtrate if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- Add 50 µL per well in the assay

### 9.2. Grain and cereals

A representative sample should be ground and thoroughly mixed before proceeding with the extraction procedure (particle size < 250 µm).

- Weigh 1 g of ground sample into a suitable vial, suspend with approx. 5 mL sample buffer and fill up to 10 mL with sample buffer (dilution factor = 10)
- Shake / mix for 10 min
- Centrifuge (or alternatively filter) in order to remove undissolved components
- Take supernatant, if supernatant is still turbid, filtrate
- Dilute the supernatant or filtrate if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- Add 50 µL per well in the assay

### 9.3. Fortified flour

A representative sample should be ground and thoroughly mixed before proceeding with the extraction procedure (particle size < 250 µm).

- Weigh 1 g of the homogenized flour sample into a suitable vial, suspend with approx. 5 mL sample buffer and fill up to 10 mL with sample buffer (dilution factor = 10)
- Shake / mix for 10 min
- Centrifuge (alternative: filter) in order to remove undissolved components
- Take supernatant, if supernatant is still turbid, filtrate
- Dilute the supernatant or filtrate if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- Add 50 µL per well in the assay

### 9.4. Vitamin powders, vitamin mixtures, and vitamin tablets

#### Vitamin powders / vitamin mixtures:

- Dissolve 1 g of vitamin powder / vitamin mixture in 5 mL sample buffer and fill up to 10 mL with sample buffer (dilution factor = 10)

#### Vitamin tablets:

- Crush the tablets and dissolve 1 g in 5 mL sample buffer and fill up to 10 mL with sample buffer (dilution factor = 10)

Continue with these samples as is described below:

- Shake / mix for 10 min
- Centrifuge (alternative: filter) in order to remove undissolved components
- Dilute the supernatant or filtrate if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- Add 50 µL per well in the assay

### 9.5. Vitamin juices

- The homogenized juice can be used directly in the test
- Centrifuge (alternative: filter) in order to remove undissolved components
- Dilute if necessary with sample buffer into the testing range (new calculation of the dilution factor is necessary!)
- Add 50 µL per well in the assay

## 10. Test implementation

### 10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **sample buffer** is provided as a 2-fold concentrate. The concentrate is diluted 1:2 with distilled water (e.g. 50 mL sample buffer concentrate + 50 mL distilled water). The diluted sample buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

A PBS-Tween buffer is required as **wash buffer**, please use the wash buffer salt (pouch) contained in the kit (see chapter 4. Reagents provided). Dissolve the total content of the pouch in one liter of distilled water. The ready to use wash buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternatively: Dissolve the content of the pouch in 100 mL of distilled water to obtain a 10-fold concentrated wash buffer. This solution expires after approx. 8 - 12 weeks, stored at room temperature (20 - 25 °C / 68 - 77 °F). Use one part of this concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.

### 10.2. Test procedure

Standards and samples must be tested on one microtiter plate at one time. Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µL of each standard or prepared sample in duplicate to the wells; use a new pipette tip for each standard or sample.
3. Add 50 µL of conjugate to each well (multi-dispenser or 8-channel pipette) and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Using an 8-channel pipette, fill the wells with ready to use wash buffer (250 µL per well) (see chapter 10.1). Empty the wells again and remove all remaining liquid. Repeat the washing step two more times.

5. Add 100 µL of the substrate/chromogen to each well (multi-dispenser or 8-channel pipette) and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
6. Add 100 µL of stop solution to each well (multi-dispenser or 8-channel pipette) and measure the absorbance at 450 nm within 5 min after addition of stop solution.

## 11. Results

A special software, the **RIDASOFT® Win.NET** (Art. No. Z9996FF), is available for evaluation of the RIDASCREEN® enzyme immunoassays. The results are given in µg/kg respectively µg/L and µg/100 g respectively µg/100 mL.

The course of the standard curve is shown in the Certificate of analysis (CoA) enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages (%). The values calculated for the standards are entered in a system of coordinates on semi-logarithmic graph paper against the folic acid concentration (µg/kg or µg/L).

In order to obtain the folic acid concentration actually contained in a sample, the concentration read from the standard curve must be further multiplied by the corresponding dilution factor.

## 12. Limits of the method

Test results may vary depending on the sample matrix, the actual test procedure and the laboratory environment.

Detection and quantification limits depend on the respective sample matrix, the degree of processing and the extraction method.

## 13. Recommendation

In order to ensure a high analytical performance we recommend to analyze each sample material in duplicates. Each laboratory may decide to perform the test in single determinations after a qualified risk management analysis. This has no influence on the function of the test kit. However, it should be noted that this increases the risk of overlooking errors in the performance of the test (e.g. pipetting errors). Moreover, a higher result variation will occur when pipetting in single determinations.

**For further product information and applications, please contact your local distributor or R-Biopharm at this address: [sales@r-biopharm.de](mailto:sales@r-biopharm.de).**

### Version overview

Version number	Chapter and title
2016-05-23	Release version
2023-10-09	Current version General revision Changes made: <ul style="list-style-type: none"><li>– Brief information</li><li>– Limit of Detection</li><li>– 6. Warnings and precautions for the users</li><li>– 11. Results</li></ul>

### Explanation of symbols

General symbols:

	Follow the instructions for use
	Batch number
	Expiry date (YYYY-MM)
	Storage temperature
	Article number
	Number of test determinations
	Manufacturing date (YYYY-MM)
	Manufacturer + address

## Disclaimer

The user assumes all risk in using R-Biopharm AG's products and services.

R-Biopharm AG will warrant that its products and services meet all quality control standards set by R-Biopharm AG, and R-Biopharm AG will, at its option, replace or repair any components, product or repeat services which prove to be defective in workmanship or material within product specific warranty periods or expiration dates and which our examination shall disclose to our satisfaction to be defective as such.

This warranty is expressly in lieu of all other warranties, expressed or implied, as to quality, description, fitness for any particular purpose, merchantability, productiveness, or any other matter. R-Biopharm AG shall be in no way responsible for the proper use of its products and hereby disclaims all other remedies, warranties, guarantees or liabilities, expressed or implied, arising by law or otherwise, and it shall have no liability for any lost profits or damage, direct, indirect or otherwise, to person or property, in connection with the use of any of its products or services.

This warranty shall not be extended, altered or varied except by a written instrument signed by an authorized representative of R-Biopharm AG.

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are registered trademarks of R-Biopharm AG.